Bioclin

BIOLISA HIV 1/2/O

REF **K118**

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INSTRUCTIONS FOR USE

FUNCTION

Third-generation test for the qualitative detection of the presence of total antibodies (IgG, IgM and IgA) of HIV-1, HIV-2, and/or subtype O in biological samples (human serum and plasma) through enzyme-immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or enzyme immunoassay.
The BIOLISA HIV 1/2/O kit is a solid-phase enzyme immunoassay based on the "sandwich" principle for the detection of total antibodies (IgG, IgM and IgA) to HIV-1, HIV-2 and/or subtype O in samples of serum and plasma. The microplate is coated with recombinant HIV antigens, specific for HIV-1 (p24, gp41), HIV-2 (gp36) and subtype O. Antibodies against HIV-1, HIV-2 and/ or subtype O, present in the samples , bind to recombinant HIV antigens coated on the microplate, forming anti-HIV antigen-antibody complexes. After the initial incubation, the microplate is washed to remove unbound material. Recombinant HIV antigens conjugated to the enzyme peroxidase are added to the microplate and then incubated. Peroxidase-conjugated recombinant HIV antigens bind to immobilized anti-HIV antigen-antibody complexes present. A further wash is performed to remove unbound materials. After this step, the substrate is added and incubated, producing a blue color that indicates the detection of anti-HIV antibodies present in the samples. Stop Solution is added to stop the reaction and a color change to yellow, measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8 °C. Contains: Microplate coated with recombinant HIV antigens.
- 2- Conjugate** – Store between 2 and 8 °C. Contains: Recombinant HIV antigens linked to peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** – Store between 2 and 8 °C. Contains: Buffer solution (phosphate < 0.5 mol/L, potassium chloride < 100 mmol/L, sodium chloride < 5 mol/L), surfactant and preservative.
- 4- Substrate** – Store between 2 and 8 °C. Contains: Buffer solution containing urea peroxide, tetramethylbenzidine (TMB) and preservative.
- 5- Stop Solution** – Store between 2 and 8 °C. Contains: 1 M hydrochloric acid solution.
- 6- Negative Control** – Store between 2 and 8 °C. Non-reactive solution for HIV-1 and HIV-2 and preservative. **Potentially infectious.**
- 7- HIV-1 Positive Control** – Store between 2 and 8 °C. Reactive solution for HIV-1, stabilizer, dye and preservative. **Potentially infectious.**
- 8- HIV-2 Positive Control** – Store between 2 and 8 °C. HIV-2 reactive solution, stabilizer, dye and preservative. **Potentially infectious.**
- 9- Plate sealers.**

PRESENTATION

REAGENTS	1	2	3
	96 cavities	192 cavities	480 cavities
1- Sensitized Plate	1 Unit	2 Units	5 Units
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
5- Stop Solution	1 Vial x 8 mL	1 Vial x 16 mL	1 Vial x 40 mL
6- Negative Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 2 mL
7- HIV-1 Positive Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 2 mL
8- HIV-2 Positive Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 2 mL
9-Plate sealers	3 Units	1 Unit	1 Unit

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 500 µL with a coefficient of variation less than 1.5%.
- 2- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer for technique on filter paper.
- 4- ELISA reader with absorbance capacity at 450 and 630 nm of wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 ± 2 °C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate should only be opened after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free of contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw material of the product is tested and must be non-reactive for HbsAg and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product containing serum is potentially capable of transmitting diseases. Therefore, it is necessary to take the proper biosafety precautions when handling these products.
- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.
- 9- As a protection measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the assay.
- 11- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advice Service) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed and highly lipemic samples should not be used. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C. Plasma samples stored for periods longer than recommended should not be used.

PROCESS DESCRIPTION

Stability After Opening

The stability test results prove that the BIOLISA HIV 1/2/O kit is stable after opening up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial Nº 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, contact SAC (Customer Advice Service).

- Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Separate the wells to be used considering: HIV-1 Positive Control, HIV-2 Positive Control, Negative Control, and Samples (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.
- 2- Separate the first cavity for the blank (OPTIONAL).
- 3- Pipette 100 µL of HIV-1 Positive Control, HIV-2 Positive Control, Negative Control and samples into previously determined wells.
- 4- Gently homogenize for ± 10 seconds, cover the wells with plate sealer.
- 5- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 6- Remove the sealer from the cavities.
- 7- Discard the contents of the wells by aspiration in a Microplate Washer. Use approximately 300 µL of **previously prepared** Washing Solution and perform a total of five (5) 5 second shake wash cycles. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause poor results.
- 8- Pipette 100 µL of Conjugate into all wells, **including the Blank well.**
- 9- Gently homogenize for ± 10 seconds. Cover the wells with the plate sealer.
- 10- Incubate for 20 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 11- Remove the sealer from the cavities.
- 12- Repeat item 7.
- 13- Pipette 100 µL of Substrate in all wells, **including the Blank well.**
- 14- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 15- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 16- Remove the plate sealer from the wells.
- 17- Pipette 50 µL of Stop Solution into all wells.
- 18- Gently homogenize for ± 30 seconds.
- 19- Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.100
Negative Control (R6)	< 0.100
HIV-1 Positive Control (R7)	> 0.500
HIV-2 Positive Control (R8)	> 0.500

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:
Cut Off = (Average Absorbance of the Negative Control (R6)) + 0.180.

Example:

ITEM	ABSORBANCES
Negative Control (R6)	0.048
	0.050
Cut Off = (Average Absorbance of the Negative Control (R6)) + 0.180.	((0.048 + 0.050) / 2) + 0.180 = 0.229

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCES
Sample	3.102
Cut Off Value	0.229
Index = Sample / Cut Off Value	3.102 / 0.229 = 13.54

Note: Data shown in examples is for illustration only and cannot be used for calculation of results.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (Non-Reagent)	< 0.9
Undetermined	Entre 0.9 e 1.1
Positive (Reagent)	> 1.1

Non-Reactive: Sample with absorbance less than the cut-off is considered non-reactive for antibodies to HIV-1, HIV-2, and/or subtype O and may be considered negative.

Reagent: A sample with an absorbance greater than the cut-off is considered initially reactive for antibodies to HIV-1, HIV-2, and/or subtype O. The sample must be reanalyzed in duplicate before the end of the interpretation. The sample that is reactive in at least one of the reanalyses is presumed to be reactive and must be confirmed through another diagnostic method.

Notes: In case of indeterminate results, the sample must be reanalyzed in duplicate. Samples that give repeatedly indeterminate results should be retested using an alternate method. If results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected shall prevail. The interpretation of a diagnostic test should not be based on a single test. Additional confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, prior to making a descriptive diagnosis of the disease.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/ dL and Anti-Streptolysin O 1023 IU/mL. Plasma samples stored for periods longer than recommended (30 days) may show fibrin and fibronectin precipitation that may interfere with the test. Samples collected in a serum tube, EDTA, heparin and citrate do not interfere with the assay.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 71 serum and plasma samples negative for HIV but positive for other infections. Among them 5 positive samples for HBsAg, 6 positive samples for HTLV, 5 positive samples for HCV, 10 positive samples for COVID-19, 5 positive samples for Dengue, 10 positive samples for Rubella, 5 positive samples for CMV, 8 positive samples for Chagas disease, 8 positive samples for Syphilis and 9 positive samples for Toxoplasmosis. No cross-reactivity was observed with positive samples for HBsAg, HTLV, HCV, COVID-19, Dengue, Rubella, CMV, Chagas Disease, Syphilis and Toxoplasmosis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Amostra		
	1	2	3
Average	2.430	0.697	0.007
Standard Deviation	0.033	0.024	0.001
Coefficient of Variation (%)	1.35	3.44	14.29

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Amostra		
	1	2	3
Average	1.715	0.368	0.033
Standard Deviation	0.144	0.039	0.005
Coefficient of Variation (%)	8.41	10.58	15.13

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HIV 1/2/O kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA HIV 1/2/O kit is >99.9% and the clinical specificity is >99.9%

		Result Reference		
		Positive	Negative	Total
BIOLISA HIV 1/2/O	Positive	87	0	87
	Negative	0	398	398
	Total	87	398	485

Clinical Sensitivity: >99.99% (87/87) - IC 95% (95.8 – 100%)

Clinical Specificity: >99.9% (398/398) - IC 95% (99.1 – 100%)

DIAGNOSTIC SIGNIFICANCE

HIV is the etiological agent of Acquired Immunodeficiency Syndrome (AIDS). The main routes of transmission include exposure to blood and blood products, including sharing needles and syringes, sexual contact and mother-to-child transmission. The virion is surrounded by a lipid envelope, derived from the host cell membrane. Several viral glycoproteins are on the envelope. Each virus contains two copies of genomic RNA. HIV-1 has been isolated from AIDS patients, and from healthy people at high potential risk for developing AIDS. HIV-1 consists of subtype M and subtype O. Highly divergent strains of HIV-1 were recognized in 1990 and provisionally grouped as subtype O, as this variation was similar to HIV-1 glycoprotein markers, but with a slight variation for the protein marker. Although rarely compared to HIV-1 and HIV-2, infections caused by subtype O have so far been identified in Africa (Cameroon), France and Germany. HIV-2 has been isolated from AIDS patients and asymptomatic HIV-positive individuals in West Africa. HIV-1, HIV-2, and Subtype O induce an immune response. Detection of HIV antibodies in serum, plasma, or whole blood is the most efficient and common way to determine whether an individual has been exposed to HIV. Despite differences in their biological characteristics, serological activities and genome sequences, HIV-1, HIV-2 and subtype O show strong antigenic cross-reactivity. Most HIV-2 positive sera can be identified using HIV-1 serological tests. The test uses recombinant HIV antigens to selectively detect HIV-1, HIV-2 and/or subtype O antibodies.

BIBLIOGRAPHIC REFERENCES

- 1- Chang, SY, Bowman, BH, Weiss, JB, Garcia, RE and White, TJ. The origin of HIV-1 isolate HTLV-IIIB. Nature (1993) 3;363:466-9.
- 2- Arya, SK, Beaver, B, Jagodzinski, L, Ensoli, B, Kanki, PJ,Albert, J, Fenyo, EM, Biberfeld, G, Zagury, JF and Laure, F. New human and simian HIV related retroviruses possess functional transactivator (tat) gene. Nature (1987) 328:548-550.
- 3- Caetano JA Immunologic Aspects of HIV infection. Acta Med Port (1991) 4 Suppl 1:52S-58S.
- 4- Janssen, RS, Satten, GA, Stramer, SL, Rawal, BD, O'Brien, TR, Weiblen, BJ, Hecht, FM, Jack, N, Cleghorn, FR, Kahn, JO, Chesney, MA and Busch MP. New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes. JAMA (1998) 280(1):42-48.
- 5- Travers, K, Mboup, S, Marlink, R, Gueye-Nidaye, A, Siby, T, Thior, I, Traore, I, Dieng-Sarr, A, Sankale, JL and Mullins, C. Natural protection against HIV-1 infection provided by HIV-2. Science (1995) 268:1612-161.
- 6- Greenberg, AE, Wiktor, SZ, DeCock, KM, Smith, P, Jaffe HW and Dondero, TJ, Jr. HIV-2 and natural protection against HIV-1 infection. Science (1996) 272:1959-1960.
- 7- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR.
- 8- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 9- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the BIOLISA HIV 1/2/O kit at ANVISA: 10269360200

Review: September/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER