BIOLISA HTLV I/II DBS

REF **K257**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Enzyme-linked immunosorbent test (ELISA) for qualitative determination of total antibodies (IgG, IgA and IgM) to the HTLV I and HTLV II viruses in dried blood samples collected on filter paper (DBS). Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

**Methodology:** Enzyme immunoassay or enzyme immunoassay  
The BIOLISA HTLV I/II DBS kit is a solid phase enzyme immunoassay based on the sandwich principle of qualitative detection of total antibodies (IgG, IgA and IgM) to HTLV I/II in human dried blood samples collected on filter paper. Total antibodies (IgG, IgA and IgM) present in the dried blood sample are eluted and bind to the recombinant HTLV I/II antigens coated on the microplate, forming immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Recombinant antigens conjugated to Peroxidase are added to the microplate, which is then incubated. The enzyme-conjugated recombinant antigens bind to the Anti-HTLV I/II total antibodies (IgG, IgA and IgM) present, bound to the antigen-coated plate. New washing is carried out to remove excess. After this step, the Substrate is added and incubated, producing a blue color that indicates the amount of total Anti-HTLV I/II Antibodies (IgG, IgA and IgM) present in the sample. The Stop Solution is added to stop the reaction and there is a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Sensitized Plate with recombinant HTLV I/II antigen and preservative.
- 2- Conjugate** - Store between 2 and 8°C. Contains: Buffer solution containing recombinant HTLV I/II antigens linked to peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Substrate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 5- Stop Solution** - Store between 2 and 8°C. Contains: 1 M Hydrochloric Acid.
- 6- Elution Buffer** - Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 7- DBS Negative Control** - Store between 2 and 8 °C. Contains: Non-reactive sample for total antibodies (IgG, IgA and IgM) Anti-HTLV I/II impregnated on filter paper. **Potentially infectious.**
- 8- DBS Positive Control** - Store between 2 and 8 °C. Contains: Reactive sample for total antibodies (IgG, IgA and IgM) Anti - HTLV I/II impregnated on filter paper. **Potentially infectious.**

PRESENTATION

| REAGENTS                | 1              | 2               | 3               |
|-------------------------|----------------|-----------------|-----------------|
|                         | 96 cavities    | 192 cavities    | 480 cavities    |
| 1- Sensitized Plate     | 1 Unit         | 2 Units         | 5 Units         |
| 2- Conjugate            | 1 Vial x 12 mL | 1 Vial x 24 mL  | 1 Vial x 60 mL  |
| 3- Concentrated Washing | 1 Vial x 50 mL | 1 Vial x 100 mL | 1 Vial x 250 mL |
| 4- Substrate            | 1 Vial x 12 mL | 1 Vial x 24 mL  | 1 Vial x 60 mL  |
| 5- Stop Solution        | 1 Vial x 12 mL | 1 Vial x 24 mL  | 1 Vial x 60 mL  |
| 6- Elution Buffer       | 1 Vial x 24 mL | 1 Vial x 48 mL  | 1 Vial x 120 mL |
| 7- DBS Negative Control | 1 Unit         | 2 Units         | 3 Units         |
| 8- DBS Positive Control | 1 Unit         | 2 Units         | 3 Units         |

EQUIPMENTS AND OPERATIONAL INPUTS

- Reagents described in the previous table

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 50 to 500 µL with a coefficient of variation of less than 1.5%.
- 2- Repipettor for repetitive pipetting of 500 µL volumes with a coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader with absorbance capacity at 450 and 620 nm wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 °C ± 2 °C.
- 11- Paper shredder (diameter 3 + 0.2 mm).

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8 °C. Transport at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate must only be opened after reaching room temperature. Return the unused microcavity strips to the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free from contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw materials for the product are tested and are non-reactive for HBsAg, Anti-HIV 1&2 and Anti HCV. However, these tests do not offer complete assurance of the absence of infectious agents. The handling of any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take appropriate biosafety precautions when handling these products.

- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 9- As a protective measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
- 11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation stages.
- 12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through a request through SAC (Customer Advisory Service) from Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

**Dried blood collected on filter paper (DBS) – EDTA or Puncture**  
The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use.  
Blood samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, samples must remain refrigerated, between 2 and 8 °C. Storage for more than 2 years must be carried out at a temperature of -20 °C.

PROCESS DESCRIPTION

Stability After Opening

The stability test results prove that the BIOLISA HTLV I/II DBS kit is stable after opening for up to 30 days. This stability may vary depending on test and environmental conditions. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of bottle No. 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, consult SAC (Customer Advisory Service).

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.
- 1- Separate the cavities to be used considering: DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.
- 2- Separate the first cavity for the Blank (OPTIONAL).
- 3- Add a 3mm disc of DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper, previously perforated, into the determined cavities.
- 4- Pipette 200 µL of Elution Buffer into all wells.
- 5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 7- Remove the sealer from the cavities.

- 8- Discard the contents of the cavities by vacuuming (Washer for filter paper technique). Use approximately 500 µL of Washing Solution, **previously prepared**, and perform a total of five (5) wash and soak cycles of 30 seconds for each cycle. To ensure the plate dries, at the end of washing, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause inadequate results.
- 9- Pipette 100 µL of Conjugate into all wells, including the Blank well.
- 10- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the plate sealer from the cavities.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells.
- 15- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 17- Remove the plate sealer from the cavities.
- 18- Pipette 100 µL of Stop Solution into all wells.
- 19- Gently homogenize for ± 30 seconds.
- 20- Read using a double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE CHECK

Check whether the results obtained for reading the Blank and Controls are compatible with the values presented below:

| ITEM                 | ABSORBANCES |
|----------------------|-------------|
| Blank                | < 0.300     |
| DBS Negative Control | < 0.300     |
| DBS Positive Control | > 1.000     |

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = Mean Absorbance of DBS Negative Control + 0.200.

Example:

| ITEM                                                      | ABSORBANCES                          |
|-----------------------------------------------------------|--------------------------------------|
| DBS Negative Control                                      | 0.102                                |
|                                                           | 0.104                                |
| Cut Off = Mean Absorbance of DBS Negative Control + 0.200 | [(0.102 + 0.104) /2] + 0.200 = 0.303 |

Calculate the Index by dividing the sample absorbance by the Cut Off value. Example:

| ITEM                           | ABSORBANCES          |
|--------------------------------|----------------------|
| Sample                         | 0.953                |
| Cut Off Value                  | 0.303                |
| Index = Sample / Cut Off Value | 0.953 / 0.303 = 3.14 |

**Note:** The data presented in the examples are for illustration only and cannot be used to calculate results.  
Each laboratory must validate the cut-off according to the instrumentation used and the population researched.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

| RESULTS                 | QUALITATIVE |
|-------------------------|-------------|
|                         | INDEX       |
| Negative (Non-Reactive) | ≤ 0.8       |
| Undetermined            | 0.8 and 1.2 |
| Positive (Reagent)      | ≥ 1.2       |

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative and samples with an index between 0.8 and 1.2 are considered indeterminate. Therefore, the sample mentioned in the example, whose absorbance was 0.953 and index of 3.77, presents a positive result (index ≥ 1.2).

**Note:** In the case of an indeterminate result, the sample must be re-analyzed. Samples that yield repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected within two weeks. The result of the last sample collected must prevail. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient’s diagnosis and/or treatment.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established based on a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before definitive diagnosis of the disease.

INTERFERENCES

No interference was observed for concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A study was carried out with 95 samples that were negative for HTLV I/II, but positive for other infections. Among them 15 positive samples for Human Immunodeficiency Virus (HIV), 10 positive samples for HBsAg, 10 positive samples for Syphilis, 10 positive samples for Hepatitis C (HCV), 10 positive samples for Chagas Disease, 10 positive samples for Zika, 10 samples positive for Cytomegalovirus (CMV), 10 samples positive for Rubella, 10 samples positive for Toxoplasmosis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient’s clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

| Repeatability                | SAMPLE |        |         |
|------------------------------|--------|--------|---------|
|                              | 1      | 2      | 3       |
| Average                      | 0.4680 | 1.6336 | 0.1242  |
| Standard Deviation           | 0.0425 | 0.0553 | 0.0131  |
| Coefficient of Variation (%) | 9.0889 | 3.3824 | 10.5328 |

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

| Reproducibility              | SAMPLE |        |        |
|------------------------------|--------|--------|--------|
|                              | 1      | 2      | 3      |
| Average                      | 0.4653 | 1.6090 | 0.1210 |
| Standard Deviation           | 0.0307 | 0.0570 | 0.0073 |
| Coefficient of Variation (%) | 6.6015 | 3.5453 | 6.0205 |

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HTLV I/II DBS kit was used for the analysis of clinical samples previously characterized and confirmed using another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA HTLV I/II DBS kit is 96.0% and the clinical specificity is 98.73%.

|                          |          | Reference Result |          |       |
|--------------------------|----------|------------------|----------|-------|
|                          |          | Positive         | Negative | Total |
| BIOLISA<br>HTLV I/II DBS | Positive | 48               | 2        | 50    |
|                          | Negative | 1                | 78       | 79    |
|                          | Total    | 49               | 80       | 129   |

Clinical Sensitivity: 96.0% (48/50) (95% CI = 86.3 – 99.5%)

Clinical Specificity: 98.73% (78/79) (95% CI = 93.1 – 100%)

DIAGNOSTIC MEANING

Human T-Cell Lymphotropic Virus (HTLV) is a member of the Retroviridae family and is classified into two types: HTLV-1 and HTLV-2. These retroviruses infect human T cells that are important for the body’s defense system. Therefore, HTLV-1 has a preference for CD4+ cells and HTLV-2 for CD8 cells, both with the ability to integrate into the host cell’s genome and are now called proviruses. Its transmission can be both vertical and horizontal through sexual intercourse, breast milk, blood transfusion and/ or intravenous drug use. The main pathologies caused by these viruses are adult T-cell leukemia (ALLcTA), Uveitis, in addition to other manifestations such as infectious dermatitis.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Tel.: (31) 3439.5454 | E-mail:bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HTLV I/II DBS kit: 10269360448

**Review:** March/2024

UNIVERSAL SYMBOLOGY

|  |                                          |  |                                        |
|--|------------------------------------------|--|----------------------------------------|
|  | CATALOG NUMBER                           |  | MADE BY                                |
|  | LOT NUMBER                               |  | CONTROL                                |
|  | MANUFACTURING DATE                       |  | POSITIVE CONTROL                       |
|  | VALIDITY DATE<br>(last day of the month) |  | NEGATIVE CONTROL                       |
|  | TEMPERATURE LIMIT<br>(store)             |  | BIOLOGICAL RISK                        |
|  | CONTENT IS SUFFICIENT<br>FOR <N> TEST    |  | FLAMMABLE                              |
|  | SEE INSTRUCTIONS<br>FOR USE              |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC<br>PRODUCT           |  | TOXIC                                  |
|  | KEEP AWAY<br>FROM SUNLIGHT               |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |
|  | DO NOT REUSE                             |  | PRODUCT<br>STERILIZED                  |
|  | CAUTION                                  |  | DANGER                                 |