BIOLISA HTLV I/II

REF **K237**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for the qualitative determination of Total Antibodies (IgG, IgA and IgM) for the HTLV I and HTLV II viruses in biological samples of serum or plasma through enzyme-immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or enzyme immunoassay
The BIOLISA HTLV I/II kit is a solid phase immunoenzymatic assay based on the sandwich principle for the qualitative detection of Total Antibodies (IgG, IgA and IgM) to HTLV I/II in human serum and plasma. Total Antibodies (IgG, IgA and IgM) present in the sample bind to the Recombinant HTLV I/II Antigens coated on the microplate, forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated Recombinant Antigens are added to the microplate which is then incubated. The enzyme-conjugated Recombinant Antigens bind to the Anti-HTLV I/II Total Antibodies (IgG, IgA and IgM) present, bound to the Antigen-coated plate. A new wash is carried out to remove the surplus. After this step, the Substrate is added and incubated, producing a blue color that indicates the amount of Total Anti-HTLV I/II Antibodies (IgG, IgA and IgM) present in the sample. Stop Solution is added to stop the reaction and there is a color change from blue to yellow as measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with Recombinant HTLV I/II Antigens and preservative.
- 2- Conjugate** - Store between 2 and 8°C. Contains: Buffer Solution containing Recombinant HTLV I/II Antigens linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Sample Diluent** - Store between 2 and 8°C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
- 5- Substrate** - Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: Hydrochloric Acid 1 M.
- 7- Negative Control** - Store between 2 and 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: Buffer Solution, Total Antibodies (IgG, IgM and IgA) Anti-HTLV I/II, dye, stabilizers, surfactant and preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 – Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2 – Conjugate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
3 – Concentrated Washing	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
4 – Sample Diluent	1 Vial x 42 mL	2 Vials x 42 mL	5 Vials x 42 mL
5 - Substrate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
9 – Stop Solution	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
7 –Negative Control	1 Vial x 300 µL	2 Vials x 300 µL	5 Vials x 300 µL
8 –Positive Control	1 Vial x 300 µL	2 Vials x 300 µL	5 Vials x 300 µL
9 – Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes of 5 to 500 µL with lower coefficient of variation than 1.5%.
- 2- Re-pipettor for repetitive pipetting volumes of 500 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Stopwatch or watch.
- 7- Flask to store the washing solution after dilution.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.
- 7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.
- 8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9- As a safety measure, you should cover the plate during the reaction.
- 10- You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.

- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

PROCESS DESCRIPTION

Stability After Opening

The stability test results prove that the Biolisa HTLV I/II kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial N°3 (Concentrated Washing) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30°C for up to 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37°C until dissolution.

Substrate

The Substrate is ready to use.

TECHNIQUE

For use in automatic equipment, contact SAC (Customer Advice Service).

- Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Separate the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.
- 2- Separate the first cavity for the Blank (OPTIONAL).
- 3- Pipette 100 µL of Sample Diluent. **Including the cavity for Blank.**
- 4- Pipette 5 µL of Samples and Controls into previously determined wells. Observe the color change of the diluent when adding the sample. The color change indicates that the sample was added to the well correctly.
- 5- Gently homogenize for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- Remove the sealer from the cavities.
- 8- Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual). Use approximately 500 µL of Washing Solution, **previously prepared**, and perform a total of five (5) 30 second wash and soak cycles for each cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause poor results.
- 9- Pipette 100 µL of Conjugate into all wells.
- 10- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the plate sealer from the wells.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells.
- 15- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.

- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 17- Remove the plate sealer from the wells.
- 18- Pipette 100 µL of Stop Solution into all wells.
- 19- Gently homogenize for ± 30 seconds.
- 20- Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.250
Negative Control	< 0.250
Positive Control	> 1.000

If the values are out the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance of the Positive Control x 0.01) + 0.200

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	2.156
	2.194
(Average Absorbance of the Positive Control x 0.01) + 0.200	((2.156 + 2.194)/2 x 0.01) + 0.200= 0.222

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.564
Cut Off Value	0.222
Index = Sample / Cut Off Value	1.564 / 0.222 = 7.045

INTERPRETATION OF RESULTS

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE FOR SERUM SAMPLE AND PLASMA
	INDEX
Negative	< 0.90
Positive	> 1.1
Undetermined	0.9 – 1.1

Note: In the case of an undetermined result, the sample must be reanalyzed. Samples that obtain results repeatedly indeterminate conditions must be retested using an alternative method. If results remain indeterminate, a new sample in two weeks. The result of the last sample must prevail collected. The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for the determining the diagnosis and /or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established with basis of a single trial. Other confirmatory tests should be included, before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before the definitive diagnosis of the disease.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL. Plasma samples stored for periods longer than recommended (30 days) may show fibrin and fibronectin precipitation that may interfere with the test. These samples should not be used, as well as samples collected with other anticoagulants other than those mentioned (EDTA and heparin).

CROSS REACTIVITY

A study was performed with 90 samples negative for HTLV I/II but positive for other infections. Among them 10 positive samples for Zika, 10 positive samples for Dengue, 10 positive samples for Chikungunya, 10 positive samples for Toxoplasmosis, 10 positive samples for Cytomegalovirus, 10 positive samples for Rubella, 10 positive samples for HIV, 10 positive samples for HBV and 10 samples positive for HCV. No cross-reactivity was observed with positive samples for Zika, Dengue, Chikungunya, Toxoplasmosis, Cytomegalovirus, Rubella, HIV, HBV and HCV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

ACCURACY

Repetibility

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repetibility	Sample		
	1	2	3
Average	0.700	2.539	0.103
Standard Deviation	0.051	0.068	0.004
Coefficient of Variation (%)	7.372	2.669	4.297

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	0.764	2.672	0.118
Standard Deviation	0.049	0.056	0.008
Coefficient of Variation (%)	6.425	2.102	6.892

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HTLV I/II kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA HTLV I/II kit is >99% and the clinical specificity >99%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	BIOLISA HTLV I/II
Positive Sample	50	50
Negative Sample	183	183
Total Samples Tested	233	233

Clinic Sensitivity: >99% (50/50) (IC 95% = 92.9 – 100%)

Clinic Specificity: >99% (183/183) (IC 95% = 98.0 – 100%)

DIAGNOSTIC SIGNIFICANCE

Human T-cell lymphotropic virus (HTLV) is a member of the Retroviridae family and is classified into two types: HTLV-1 and HTLV-2. These retroviruses infect human T cells that are important for the body's defense system. Since, HTLV-1 with preference for CD4 + cells and HTLV-2 with CD8 cells, both with the ability to integrate into the host cell's genome, becoming known as a provirus. Its transmission can be either vertical or horizontal through sexual intercourse, breast milk, blood transfusion and / or the use of intravenous drugs. The main pathologies caused by these viruses are adult T cell leukemia (LLcTA), Uveitis, in addition to other manifestations such as infectious dermatitis.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HTLV I/II kit: 10269360361

Review: June/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER