BIOLISA LH

[REF] K260

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for the quantitative determination of Luteinizing Hormone (LH) in biological samples of serum or plasma, through an enzyme-immunoassay test. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA LH kit is a solid-phase enzyme immunoassay based on the “sandwich” principle for the quantitative detection of Luteinizing Hormone (LH) in human serum or plasma samples. The microplate is coated with specific Antibodies against an antigenic site present exclusively on LH molecules. At first, the sample is incubated with the Conjugate, also made up of specific Anti-LH Antibodies, conjugated to Peroxidase. The LH molecules present in the sample bind to the Anti-LH Antibodies immobilized on the plate and to the Anti-LH Antibodies present in the Conjugate, forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Then the Substrate is added and incubated, producing a blue color that indicates the detection of Luteinizing Hormone in the sample. Stop Solution is added to stop the reaction, with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Reference Standards (A – F)** - Store between 2 and 8 °C. Contains: Six (6) vials of Reference Standards (A - F), containing Luteinizing Hormone (LH) in different concentrations, Buffer Solution, stabilizers, preservatives, surfactant and dye. The concentrations vary for each batch (from 0 to 100 mIU/mL), see the vial labels. **Potentially infectious.**
- 2- Conjugate** - Store between 2 and 8 °C. Contains: Buffer Solution containing specific Anti-LH Antibodies conjugated to Peroxidase, stabilizers, preservatives, surfactant and dye.
- 3- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with specific Anti-LH Antibodies.
- 4- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1M Hydrochloric Acid.

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Reference Standards (A-F)	6 Vials x 1 mL	6 Vials x 1 mL	6 Vials x 2 mL
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Sensitized Plate	1 Unit	2 Units	5 Units
4- Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the table above.
- Instructions for use (manual).

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 500 µL with a coefficient of variation less than 1.5%.
- 2- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer.
- 4- ELISA reader with absorbance capacity at 450 and 630 nm of wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or watch.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- 37 °C ± 2 °C incubator.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- **For *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate should only be opened after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free of contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All product raw material is tested and must be non-reactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product containing serum is potentially capable of transmitting diseases. Therefore, it is necessary to take the proper biosafety precautions when handling these products.
- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.
- 9- As a protection measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the assay.
- 11- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Safety Service) Customer Advice) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA)

Hemolyzed and highly lipemic samples should not be used. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C. Plasma samples stored for periods longer than recommended and samples collected with other anticoagulants other than those mentioned above should not be used.

PROCESS DESCRIPTION

Stability After Opening

The stability test results prove that the BIOLISA LH kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute 50 mL of Concentrated Washing (vial N° 4) in 1000 mL of Distilled or Deionized Water. After preparation, the solution can be stored at 2 to 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, contact SAC (Customer Advice Service).

Before starting the assay, place all reagents, Reference Standards and Samples to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

- 1- Separate the cavities to be used considering: Reference Standards (from A to F) and Samples. It is recommended to test in duplicate. Return unused microplate strips to the original sealed package.
- 2- Separate the first cavity for the Blank (Optional).
- 3- Pipette 50 µL of each of the Reference Standards (A-F) and the Samples in the previously determined wells.
- 4- Pipette 100 µL of conjugate into all wells, **including the Blank well.**
- 5- Gently homogenize for ± 10 seconds. Cover the wells with the plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes at 37°C.
- 7- Remove the sealer from the cavities.
- 8- Discard the contents of the wells by aspiration in a Microplate Washer. Use approximately 300 µL of previously prepared Washing Solution and perform a total of five (5) 3-second shake wash cycles. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.

Note: Poor washing/drying may cause poor results.

- 9- Pipette 100 µL of Substrate in all wells, **including the Blank well.**
- 10- Gently homogenize for ± 10 seconds. Cover the wells with the plate sealer.
- 11- Incubate for 10 minutes ± 2 minutes at 37°C.
- 12- Remove the sealer from the cavities.
- 13- Pipette 100 µL of Stop Solution into all wells, **including the Blank well.**
- 14- Gently homogenize for ± 10 seconds.
- 15- Read using a double filter: 450 nm / 630 nm in up to 10 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.100
Standard Reference A	< 0.100
Standard Reference B	> A and < C
Standard Reference C	> B and < D
Standard Reference D	> C and < E
Standard Reference E	> D and < F
Standard Reference F	> E

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

A calibration curve is used to determine the LH concentration in unknown samples.

Preparation of the Calibration Curve

If they were performed in duplicate, calculate the absorbance averages obtained in the microplate reader of each of the Reference Standards (A-F), as in the example:

Reference Standard	Concentration (See Label)	Absorbance	Average Absorbance
A	0.0	0.018	0.018
		0.017	
B	7.5	0.138	0.136
		0.134	
C	15.0	0.283	0.283
		0.283	
D	30.0	0.636	0.640
		0.643	
E	60.0	1.260	1.271
		1.281	
F	100.0	2.096	2.110
		2.124	

Plot the mean absorbances of each Reference Standard (A-F) versus the corresponding concentration in mIU/mL on graph paper or using computer programs. Trace the point-to-point curve.

Calculation of the LH Concentration of the Samples

Calculate the concentration of each analyzed sample from the Calibration Curve made on the graph paper, interpolating the absorbance value obtained in the sample reading with the corresponding concentration value (mIU/mL). The calculation of the concentration of the tested samples can also be performed using adequate computer programs, through linear regression.

The measuring range of the kit is between 0.4 and 100.0 mIU/mL (analytical sensitivity and maximum point of the curve, respectively).

Results falling below the lower limit of detection should be expressed as < 0.4 mIU/mL. Results above the upper limit of detection should be expressed as > 100.0 mIU/mL. Samples that show values greater than the last point on the curve can be diluted and reanalyzed, as described below.

Sample Dilution

Samples with concentration >100.0 mIU/mL can be diluted and reanalyzed. To do so, add 10 µL of sample in 90 µL of 0.9% Sodium Chloride Solution (1:10 dilution). The concentration obtained from the diluted sample must be multiplied by the dilution factor (if diluted 1:10, multiply by 10). If the sample concentration is still greater than 100.0 mIU/mL, a higher dilution must be performed.

Note: The data shown are for illustration only and cannot be used as a substitute for the calibration curve, which must be constructed in the laboratory. The Reference Standards must always be tested with each new assay.

REFERENCE VALUES

Sex	Reference Value (IC 95%) mIU/mL
Men	0.4 to 13.0
Women	
Follicular Phase	0.4 to 14.3
Midcycle	10.2 to 100.0
Luteal Phase	0.4 to 14.3
Postmenopausal	12.2 to 68.3

These values should be used as guidelines, and each laboratory should create its range of reference values, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PROCESS LIMITATIONS

Follicular and luteal phase samples may present values lower than the analytical and functional sensitivity of the product. Midcycle samples or ovulatory phase may show results higher than the last curve calibration point. In this case, the sample can be diluted and the value found must be multiplied by the dilution factor to determine the correct concentration value.

LH is usually suppressed by estrogens and progestogens, but in women using hormonal contraceptives, the level can be low or normal. Excessive dieting and weight loss can reduce LH levels.

The LH concentration in a sample can vary due to differences in assay methods, calibration, and reagent specificity. The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis.

INTERFERENTS

No interference was observed by Triglycerides 1,500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1,000 mg/dL, Acid Oxalic 60 mg/dL, C-Reactive Protein 41.2 mg/dL, anti-Streptolysin O 1,023 IU/mL and Rheumatoid Factor 1,080 IU/mL. Plasma samples collected in EDTA can be used for LH quantification. Samples collected with other anticoagulants other than the aforementioned cannot be used. Plasma samples stored for periods longer than recommended (30 days) may show fibrin and fibronectin precipitation that may interfere with the test.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 50 serum and plasma samples free of LH but with known concentrations of other analytes. Among them, 10 samples with Human Chorionic Gonadotropin (hCG) up to 50,000 mIU/mL, 10 samples with Follicle Stimulating Hormone (FSH) up to 1,000 mIU/mL, 10 samples with Thyroid Stimulating Hormone (TSH) up to 1,000 mIU/mL, 10 samples with prolactin up to 1,000 ng/mL and 10 samples with Prostate Specific Antigen (PSA) up to 1,000 ng/mL. No false positive results were observed for the above analytes, up to the maximum concentration tested. Despite the results found, the possibility of cross-reactivity with other analytes or with higher concentrations cannot be completely ruled out. Due to elevated levels of hCG during pregnancy, use of enzyme-linked immunoassay for LH is not recommended during pregnancy or immediately after delivery. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	1.242	1.517	0.074
Standard Deviation	0.039	0.065	0.005
Coefficient of Variation (%)	3.15	4.26	7.32

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	1.426	1.659	0.087
Standard Deviation	0.147	0.182	0.012
Coefficient of Variation (%)	10.32	11.00	13.45

ANALYTICAL SENSITIVITY

The analytical sensitivity of the BIOLISA LH kit is 0.4 mIU/mL and the functional sensitivity is 0.49 mIU/mL.

LINEARITY

The reaction is linear up to the highest point concentration on the calibration curve. For samples with higher values, dilute the samples as described in "Sample Dilution", repeat the dosage and multiply the result obtained by the dilution factor.

HIGH DOSE HOOK EFFECT

No high dose hook effect was observed at up to 1,000 mIU/mL LH.

DIAGNOSTIC SIGNIFICANCE

Luteinizing Hormone (LH) is produced in both men and women by the anterior pituitary gland in response to luteinizing hormone-releasing hormone (LH-RH or Gn-RH), which is released by the hypothalamus. LH is a glycoprotein composed of two chains of amino acids, alpha and beta. The alpha chain is similar to that found in Thyroid Stimulating Hormone (TSH), Follicle Stimulating Hormone (FSH) and Human Chorionic Gonadotropin (hCG). The difference between these hormones lies in the amino acid composition of their beta subunits, which are responsible for their immunological differentiation. Basal LH secretion in men is episodic and its primary function is to stimulate interstitial cells to produce testosterone. Variation in LH concentrations in women is subject to the complex ovulatory cycle and depends on a sequence of hormonal events along the hypothalamic-pituitary-gonadal axis. LH is released as a result of the rapid increase in estradiol due to follicle development, with direct stimulation of the pituitary gland. Ovulation occurs approximately 12 to 18 hours after LH peaks. Patients suffering from hypogonadism have increased concentrations of serum LH. High concentrations of LH can also be found in primary testicular failure and Klinefelter syndrome. Increased LH concentrations are also present during renal failure, cirrhosis, hyperthyroidism, and severe starvation. Decreased LH concentrations can result in infertility in both men and women, and may indicate some pituitary or hypothalamic dysfunction. In the differential diagnosis of hypothalamic, pituitary, or gonadal dysfunction, LH concentration assays are routinely performed in conjunction with FSH assays, as their roles are closely related. In addition, hormone levels are used to determine menopause, identify ovulation, and monitor endocrine therapy.

BIBLIOGRAPHIC REFERENCES

- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Harris, G.W. and Naftolin., The Hypothalamus and Control of Ovulation, Brit. Med. Bullet., 6, 1-9 (1970).
- 3- Knobil, E., The Neuroendocrine Control of the Menstrual Cycle, Rec. Prog. Horm Res., 36, 52-88 (1980).
- 4- Jeffcoate, S.L., Clinics in Endocrinol. Metab., 4, 521-543 (1975).
- 5- Whitely, R.J., Keutmann, H.T. and Ryan, R.J., Endocrinology, 102, 1874 (1978).
- 6- Pierce, J.G. and Parsons, T.F., Glycoprotein hormones: Structure and Function, Annual Rev. Biochem., 50, 465-495 (1981).
- 7- Bardin, C.W. and Paulsen, C.A., "The Testes" in Textbook of Endocrinology, (ed.) R.H. Williams M.D., W.B. Saunders Co., (1981).
- 8- Shome, B. and Parlow, A.F., J. Clin. Endocrinol. Metabl., 39, 199-202 (1974).
- 9- Shome, B. and Parlow, A.F., J. Clin. Endocrinol. Metabl., 39, 203-205 (1974).
- 10- Ross, G.T., VandeWeile, R.L., and Frantz, A.G. Chapter 7 in "The Ovaries and the Breasts" in "Textbook of Endocrinology" (R. H. Williams, Ed.), W.B. Saunders Co. (1981).
- 11- Marshall, J.C., Clinics in Endocrinal Metab., 4, 545-567 (1975).
- 12- Acosta, A.A.M.D. and Wright, G.L., Journal of Clinical Immunoassays, 6, 41 (1983).
- 13- Brown, J. Wood (eds.), Churchill Livingston Co., New York, 7-42 (1980).
- 14- Yen, S.S.C., Vela, P. and Rankin, J., J. Clin. Endocrinol. Metab., 30, 434-442 (1970).
- 15- Cohen, K. L., Metabolism, 26, 1165-1177 (1977).
- 16- Engvall, E., Methods in Endocrinology, Volume 70, Van Vunakis, H. and Langone, J.J. (eds.), Academic Press, New York, 419-492 (1980).
- 17- Uotila, M. Ruoslahti,E. and Engvall, E., J. Immunol. Methods, 42, 11-15 (1981).
- 18- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the BIOLISA LH kit at ANVISA: 10269360353

Review: March/2023

UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER