Bioclin

BIOLISA RUBÉOLA IgM

REF **K320**

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FINALIDADE

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- Bioclin – Dados de arquivos.

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2. 2002:31.
- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.

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
BIOLISA RUBELLA IgM

REF K320 INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Capture enzyme-linked immunosorbent assay (ELISA) for the determination of IgM antibodies to Rubella in human serum or plasma samples. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic assay based on the principle of qualitative detection by capture of Rubella IgM antibodies in human serum or plasma samples. IgM antibodies present in the sample bind to the anti-IgM antibodies coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Rubella antigens conjugated to peroxidase are added to the microplate which is then incubated. The enzyme-conjugated antigens bind to the anti-Rubella IgM antibodies present, bound to the plate coated with anti-IgM antibodies. A new wash is performed to remove the excess. After this step, the substrate is added and incubated, producing a blue color that indicates the amount of anti-Rubella IgM antibodies present in the sample. Stop solution is added to stop the reaction and there is a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with anti-IgM antibodies.
- 2- Conjugate** - Store between 2 and 8 °C. Contains: Buffer solution, Rubella antigens conjugated to peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Wash** - Store between 2 and 8 °C. Contains: Buffer solution (phosphate < 0.5 mol/L, potassium chloride < 100 mmol/L, sodium chloride < 5 mol/L), surfactant and preservative.
- 4- Sample Diluent** - Store between 2 and 8 °C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1 M Hydrochloric Acid.
- 7- Negative Control** - Store between 2 and 8 °C. Contains: Buffer solution, stabilizer, surfactant and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8 °C. Contains: IgM anti-Rubella antibodies, buffer solution, dye, stabilizers, surfactant and preservative. **Potentially infectious.**

PRESENTATION

REACTIVES	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unidad	2 Unidades	5 Unidades
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Concentrated Wash	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4- Sample Diluent	1 Vial x 25 mL	1 Vial x 50 mL	1 Vial x 100 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
7- Negative Control	1 Vial x 300 µL	1 Vial x 300 µL	1 Vial x 1 mL
8- Positive Control	1 Vial x 300 µL	1 Vial x 300 µL	1 Vial x 1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- Pipettes capable of dispensing volumes from 5 to 300 µL with a coefficient of variation lower than 1.5%.
- Repipette for repetitive pipetting of volumes of 300 µL with a coefficient of variation lower than 1.5% or multichannel pipette (optional).
- Microplate washer (optional).
- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- Absorbent paper to dry the microwells.
- Stopwatch or clock.
- Vial to store the Wash Solution after dilution.
- Distilled or deionized water.
- Quality Control Tools.
- Incubator 37 °C ± 2 °C.

STORAGE AND TRANSPORTATION CONDITIONS

The storage temperature should be between 2 and 8 °C. Transportation at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- For professional *in vitro* diagnostic use only.
- Strictly follow the proposed methodology to obtain accurate results.
- The sachet containing the microplate should only be opened after it has reached room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- The water used to clean the material must be fresh and free of contaminants.
- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- The Stop Solution contains Hydrochloric Acid, which is a strong acid. Therefore, handle it with due care.
- All raw materials of the product are tested and found to be

non-reactive for HBsAg, Anti-HIV and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take the appropriate biosafety precautions when handling these products.

- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.
- As a protective measure, the plate must be covered during the reaction.
- Ensure that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during the storage or incubation stages.
- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or by requesting it through Quibasa's SAC (Customer Advisory Service).
- Do not use the product if the packaging is damaged.
- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed or highly lipemic samples should not be used. Samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C. Plasma samples stored for periods longer than the recommended period (30 days) and samples collected with anticoagulants other than those mentioned should not be used.

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test show that the BIOLISA RUBELLA IgM kit is stable after opening for up to 30 days. This stability may vary depending on the test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of Vial No. 4 (Concentrated Wash) in a ratio of 1:20 with distilled or deionized water; for example: 50 mL of Concentrated Washing Solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready for use.

TECHNIQUE

For use in automatic equipment, consult SAC (Customer Advisory Service).

Before starting the test, place reagents, Controls and Samples to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

- Separate the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed packaging.
- Separate the first well for the Blank (OPTIONAL).
- Pipette 100 µL of Sample Diluent, including in the well for the Blank.

4- Pipette 5 µL of Sample and Controls in the previously determined wells. Observe the color change of the diluent when adding the sample. The color change indicates that the sample was added to the well correctly.

5- Gently mix for ± 30 seconds, cover the wells with plate sealer.

6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

7- Remove the sealer from the wells.

8- Discard the contents of the wells by aspiration (Washer) or by decantation (Manual). Use approximately 300 µL of Washing Solution, **previously prepared**, and perform a total of five (5) washing cycles with shaking for 5 seconds. To ensure that the plate is dry, at the end of the washing, tap the plate for a few seconds on absorbent paper.

Note: Inadequate washing/drying may cause inadequate results.

9- Pipette 100 µL of Conjugate into all wells, including the Blank well.

10- Gently mix for ± 30 seconds. Cover the wells with the plate sealer.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

12- Remove the plate sealer from the wells.

13- Repeat item 8.

14- Pipette 100 µL of Substrate into all wells, including the Blank well.

15- Gently mix for ± 30 seconds. Cover the wells with the plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

17- Remove the plate sealer from the wells.

18- Pipette 100 µL of Stop Solution into all wells, including the Blank well.

19- Gently mix for ± 30 seconds.

20- Read using double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check whether the results obtained for reading the Blank and Controls are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.150
Negative Control	< 0.150
Positive Control	> 1.000

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance of Positive Control x 0.1) + 0.200.

Example:

ITEM	ABSORBANCIAS
Positive Control	2.125
	2.120
Cut Off = (Average Absorbance of Positive Control x 0.1) + 0.200	Cut Off = ((2.125 + 2.120) / 2 x 0.1) + 0.200 Cut Off = 0.412

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCIA
Sample	1.421
Cut Off Value	0.412
Index = Sample / Cut Off Value	1.421/ 0.412 = 3.45

Note: The data presented in the examples are for illustration purposes only and cannot be used to calculate results. Each laboratory must validate the cutoff according to the instrumentation used and the population studied.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (Non Reactive)	≤ 0,8
Indeterminate	0,8 – 1,2
Positive (Reactive)	≥ 1,2

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative and samples with an index between 0.8 and 1.2 are considered indeterminate. Thus, we have that the sample mentioned in the example, whose absorbance was 1.421 and index of 3.45, presents a positive result (index ≥ 1.2).

Notes: In the case of an indeterminate result, the sample must be reanalyzed in duplicate. Samples that repeatedly obtain indeterminate results must be retested using an alternative method. If the results remain indeterminate, a new sample must be collected in two weeks. The result of the last sample collected must prevail. The interpretation of a diagnostic test should not be established based on a single assay. Further confirmatory tests must be included before a sample is considered positive. A negative result does not rule out the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before a descriptive diagnosis of the disease is made. The results provided by this kit must be interpreted by the responsible medical professional and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

PROCESS LIMITATIONS

IgM antibodies are present in serum and plasma for a short period of time, disappearing up to six months after infection, while IgG antibodies remain present for a long period of time. However, the high sensitivity of new diagnostic methodologies, in some cases, allows very low levels of IgM antibodies, called residual IgM, to be found for a longer period of time. The detection of these residual IgM antibodies makes clinical interpretation of the period of infection difficult. In these cases, to confirm the result, it is recommended to perform an IgG avidity test. The interpretation of a diagnostic test should not be established based on a single assay. Other confirmatory tests should be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other available clinical information before definitive diagnosis of the disease. The results provided by this kit should be interpreted by the responsible medical professional and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS-REACTIVITY

A study was performed with 49 serum/plasma samples negative for Rubella IgM but positive for other infections in order to evaluate the possibility of cross-reactivity of these interferents in the result of the BIOLISA RUBELLA IgM kit. Among them, 5 samples were positive for HTLV, 6 samples were positive for Syphilis, 6 samples were positive for HBsAg, 5 samples were positive for HCV, 6 samples were positive for Chagas disease, 6 samples were positive for Zika, 5 samples were positive for Toxoplasmosis, 5 samples were positive for CMV and 5 samples were positive for HIV. No cross-reactivity was observed with samples positive for HTLV, Syphilis, HBsAg, HCV, Chagas disease, Zika, Toxoplasmosis, CMV and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results.

Repeatability	SAMPLE		
	1	2	3
Mean	1.791	0.085	0.053
Standard Deviation	0.049	0.005	0.003
Coefficient of Variation (%)	2.76	6.29	5.69

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Mean	0.069	0.053	1.582
Standard Deviation	0.004	0.002	0.060
Coefficient of Variation (%)	5.80	3.77	3.79

Clinical Sensitivity and Specificity

The BIOLISA RUBELLA IgM kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA RUBELLA IgM kit is 99.9% and the clinical specificity is 98.2%.

RESULT	REFERENCE			
	Positive	Negative	Total	
BIOLISA RUBELLA IgM	Positive	52	2	54
	Negative	0	112	112
	Total	52	114	166

Clinical Sensitivity: 99.9% (52/52) 95% CI = 93.2 to 99.9 %
Clinical Specificity: 98.2% (112/114) 95% CI = 93.8 to 99.8%

DIAGNOSTIC SIGNIFICANCE

Rubella is a small, enveloped, spherical RNA virus belonging to the family Togaviridae. It is commonly known as German or 3-day measles. Rubella virus infection is transmitted through salivary droplets, resulting in a mild contagious rash in children or young adults. In childhood, the infection is a self-limiting, benign disease characterized by low-grade fever, headache, lymphadenopathy, arthralgia, and conjunctivitis. However, infection during pregnancy, especially in the first trimester, can lead to spontaneous abortion, intrauterine infection causing fetal death, or congenital anomalies. Congenital rubella depends on the period at which the infection occurs and can result in serious complications, including deafness, eye problems including cataracts and glaucoma, congenital heart disease, and mental retardation. IgM antibodies against rubella are produced initially, may reach detectable levels within 2-3 days and peak 14-21 days after the onset of symptoms that remain detectable during the next 4 - 8 weeks. The diagnosis of active or recent infection can be obtained by the presence of IgM antibodies in an initial sample. After several days, IgG antibodies appear after IgM, peaking at 14 - 21 days, with varying levels persisting throughout life. The presence of IgG anti-Rubella antibodies is indicative of previous infection and/or immunity.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
4. Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
5. Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
6. Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
7. Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
8. Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
9. Bioclin – Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date indicated on the packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Rubella IgM kit:
10269360465

Review: June/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N>-TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER