Bioclin

BIOLISA RUBÉOLA IgM DBS

REF **K275**

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FINALIDADE

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- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB).
- Bioclin – Dados de arquivos.

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INSTRUCTIONS FOR USE

FUNCTION

Enzymatic immunoassay capture test (ELISA) for determining IgM antibodies to the Rubella virus in biological samples of **dried blood collected on filter paper (DBS)**, for screening in adults and newborns, using an enzyme immunoassay test. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic immunoassay based on the principle of qualitative detection by capture of IgM Antibodies to the Rubella virus in human dried blood samples collected on filter paper. IgM antibodies to the Rubella virus, present in the dried blood sample, are eluted and bind to the antibodies coated on the microplate, forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated Rubella antigens are added to the microplate, which is then incubated. The enzyme-conjugated Antigens bind to the Anti-Rubella IgM Antibodies present, attached to the plate coated with Anti-IgM Antibodies. New washing is carried out to remove excess. After this step, the Substrate is added and incubated, producing a blue color that indicates the amount of Anti-Rubella IgM Antibodies present in the sample. The Stop Solution is added to stop the reaction and there is a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8 °C. Contains: Plate sensitized with anti-IgM antibodies.
- 2- Conjugate** – Store between 2 and 8 °C. Contains: Buffer Solution, Rubella Antigen conjugated to Peroxidase surfactant, stabilizers, coloring and preservative.
- 3- Concentrated Washing** – Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Substrate** – Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 5- Stopping Solution** – Store between 2 and 8 °C. Contains: 1 M Hydrochloric Acid.
- 6- Elution Buffer** – Store between 2 and 8 °C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 7- DBS Negative Control** – Store between 2 and 8 °C. Contains: Non-reactive sample for anti-Rubella IgM Antibodies impregnated on Filter Paper. **Potentially infectious.**
- 8- DBS Positive Control** – Store between 2 and 8 °C. Contains: Reactive sample for IgM Anti-Rubella Antibodies impregnated on Filter Paper. **Potentially infectious.**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unit	2 Units	5 Units
2- Conjugate	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
3- Concentrated Wash	1 Vial x 50 mL	1 Vials x 100 mL	1 Vials x 500 mL
4- Substrate	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
5- Stop Solution	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
6- Elution Buffer	1 Vial x 20 mL	1 Vials x 40 mL	1 Vials x 100 mL
7- DBS Negative Control	1 Unit	2 Units	3 Units
8- DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

Required materials not contained in the kit:

- Pipettes capable of dispensing volumes of 5 to 500 µL with a coefficient of variation of less than 1.5%.
- Repipettor for repetitive pipetting of volumes of 300 µL with a coefficient of variation of less than 1.5% or multichannel pipette (optional).
- Microplate washer for filter paper technique.
- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- Absorbent paper to dry the microcavities.
- Stopwatch or clock.
- Bottle to store the washing solution after dilution.
- Distilled or deionized water.
- Quality control tools.
- Incubator 37 °C ± 2 °C.
- Paper mincer (diameter 3 ± 0.2 mm).

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8 °C. Transportation at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the proposed methodology in order to obtain accurate results.
- The sachet containing the microplate should only be opened after it has reached room temperature. Place the unused microcavity strips back in the sachet, seal and store between 2 and 8 °C.
- The water used to clean the material must be fresh and free of contaminants.
- Saturated deionizing columns release alkaline water, various ions, oxidizing and reducing agents which can significantly alter the results.
- The Stop Solution contains hydrochloric acid, which is a strong acid. Therefore, handle it with due care.

- All the product's raw materials are tested as non-reactive for HBsAg and Anti-HCV. However, these tests do not provide complete assurance of the absence of infectious agents. The handling of any product containing a biological sample is potentially capable of transmitting diseases. Biosafety precautions must therefore be taken when handling these products.
- Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.
- As a protective measure, cover the plate during the reaction.
- Ensure that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
- Do not expose reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- We recommend applying local, state and federal environmental protection regulations so that reagents and biological material are disposed of in accordance with current legislation.
- For information related to biosafety or in the event of accidents with the product, consult the SDS (Safety Data Sheet) available at www.bioclin.com.br or by contacting Quibasa's Customer Service Center.
- Do not use the product if the packaging is damaged.
- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Dried blood collected on filter paper (DBS) - EDTA, Puncture or Heel stick (for neonates)

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use. Stability studies, carried out by Bioclin's Research and Development department, show that dried blood samples on filter paper can be stored for up to 35 days at room temperature, 2 years at 2-8 °C or -20°C. The Manual "Blood Collection: Diagnosis and monitoring of STDs, AIDS and Viral Hepatitis: Brasília, Ministry of Health, Department of STD, AIDS and Viral Hepatitis. 2010, 98 p. (TELELAB Series)[®], reports that dried blood samples on filter paper can be stored at room temperature, as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, samples must remain refrigerated, between 2 and 8 °C. Storage for more than 2 years must be carried out at a temperature of -20 °C.

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test prove that the BIOLISA RUBELLA IgM DBS kit is stable, after opening, for 30 days. This stability may vary depending on test and environmental conditions. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of bottle No. 3 (Concentrated Wash) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Wash solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved

All other reagents are ready for use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

- Separate the cavities to be used considering: DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.
- Separate the first cavity for the White (OPTIONAL).
- Add a 3 mm disc of DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper, previously perforated, into the determined cavities.
- Pipette 100 µL of Elution Buffer into all wells, including the blank well.
- Homogenize gently for ± 30 seconds, cover the cavities with plate sealer.
- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- After incubation, discard the contents of the cavities by aspiration (Washer for filter paper technique). Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) washing cycles with 5-second shaking. To ensure the plate dries, at the end of washing, tap the plate for a few seconds on absorbent paper.
- Poor washing/drying may cause inadequate results.
- Pipette 100 µL of Conjugate into all wells, including the White well.
- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- Remove the plate sealer from the cavities.
- Repeat item 7.
- Pipette 100 µL of Substrate into all wells. Including in the White cavity.
- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- Remove the plate sealer from the cavities.
- Pipette 100 µL of Stopping Solution into all wells.
- Gently homogenize for ± 30 seconds.
- Read using double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check whether the results obtained for reading the Blank and controls are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.150
DBS Negative Control	< 0.150
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

QUALITATIVE

Calculate the Cut-Off according to the following formula:
Cut Off = Average Absorbance of Negative Control DBS + 0.200.

Example:

ITEM	ABSORBANCE
DBS Negative Control	0.020
	0.022
Cut Off = Average Absorbance of Negative Control DBS + 0.200	$(0.020 + 0.022) / 2 + 0.200 = 0.221$

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCIA
Sample	1.254
Cut-Off Value	0.221
Index = Sample / Cut-Off Value	$1.254 / 0.221 = 5.67$

Note: The data presented in the examples is for illustrative purposes only and cannot be used to calculate results.

Each laboratory must validate the cut-off according to the instrumentation used and the population surveyed.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	INDEX
Negative (No Reactive)	≤ 0.8
Indeterminate	0.8 - 1.2
Positive (Reactive)	≥ 1.2

Thus, samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative and samples with an index between 0.8 and 1.2 are considered indeterminate. Therefore, the sample mentioned in the example, whose absorbance was 1.254 and index of 5.67, presents a positive result (index > 1.2).

Each laboratory must validate the cut according to the instrumentation used and the population researched.

Observations: In case of an indeterminate result, the sample must be re-analyzed in duplicate. Samples that yield repeatedly indeterminate results should be retested using an alternative method. If results remain indeterminate, a new sample should be collected within two weeks. The result of the last sample collected must prevail. The interpretation of a diagnostic test should not be established based on a single assay. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before descriptive diagnosis of the disease. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

LIMITATIONS OF THE PROCESS

The interpretation of a diagnostic test should not be established on the basis of a single assay. All results must be interpreted in conjunction with other available clinical information before a definitive diagnosis can be made. The results provided by this kit must be interpreted by the medical professional responsible and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed for concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS-REACTIVITY

A study was carried out with 69 dried blood samples collected on filter paper that were negative for Rubella, but positive for other infections, in order to evaluate the possibility of cross-reactivity of these interferences in the BIOLISA RUBELLA IgM DBS result. Among them, 5 samples for HTLV, 10 positive samples for Syphilis, 5 positive samples for HBsAg, 5 positive samples for HCV, 10 positive samples for Chagas Disease, 4 positive samples for Zika, 10 positive samples for CMV, 10 positive samples for Toxoplasmosis and 10 positive samples for HIV. No cross-reactivity was observed with samples positive for Syphilis, HBsAg, HCV, Chagas Disease, Zika, CMV, Toxoplasmosis and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is advisable to use controls to assess the precision and accuracy of dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	1.5004	0.882	0.013
Standard Deviation	0.028	0.099	0.002
Coefficient of Variation (%)	1,856	11,209	15,008

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	1.513	0.907	0.013
Standard Deviation	0.014	0.032	0.001
Coefficient of Variation (%)	0.895	3.517	2.670

Considering the variation of all samples tested in the precision study, it was observed that the kit presented a coefficient of variation of less than 20%.

Clinical Sensitivity and Specificity

The BIOLISA RUBELLA IgM kit analyzed clinical samples and compared it to other EIA methods. The results show that the clinical sensitivity of the BIOLISA Rubella IgM kit is > 99.99% and the clinical specificity is > 99.99%.

RESULT	REFERENCE			
	Positive	Negative	Total	
BIOLISA RUBELLA IgM DBS	Positive	53	0	53
	Negative	0	67	67
	Total	53	67	120

Clinical Sensitivity: > 99.99% (53/53) - IC 95% = 95.60 a 100%

Clinical Specificity: > 99.99% (67/67) – IC 95% = 94.45 a 100 %

DIAGNOSTIC SIGNIFICANCE

Rubella is a spherical, small enveloped RNA virus belonging to the Togaviridae family. It is commonly known as German or 3-day measles. Rubella virus infection is transmitted through droplets of saliva, resulting in a mild contagious rash in children or young adults. In childhood, the infection is a self-limiting, benign disease, characterized by low fever, headache, lymphadenopathy, arthralgia and conjunctivitis. However, infection during pregnancy, especially in the first trimester, can lead to miscarriage, intrauterine infection causing fetal death, or congenital anomalies. Congenital Rubella depends on the period in which the infection occurs and can result in serious complications, including deafness, eye problems including cataracts and glaucoma, congenital heart disease and mental retardation. IgM antibodies against Rubella are produced initially, reaching detectable levels within 2 - 3 days and peaking 14 - 21 days after the onset of symptoms that remain detectable for the next 4 - 8 weeks. The diagnosis of active or recent infection can be obtained by the presence of IgM antibodies in the initial sample. After several days, IgG antibodies appear after IgM, peaking at 14 - 21 days, persisting at varying levels for life. The presence of anti-Rubella IgG antibodies is indicative of previous infection and/or immunity. Therefore, the Biolisa Rubella IgM DBS kit is a valuable tool in monitoring prenatal and neonatal screening.

BIBLIOGRAPHIC REFERENCES

- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE. A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB).
- Bioclin – Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for RUBELLA IgM DBS kit: 10269360459

Review: January/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER