BIOLISA RUBÉOLA IgG DBS

REF K276

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FINALIDADE

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- Hermann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- Cancer Epidemiol Biomarkers Prev; 26(8): 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology, 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Indirect enzyme immunoassay (ELISA) for the quantitative and qualitative determination of IgG Antibodies to the Rubella virus in **dried blood spot collected on filter paper (DBS)**. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or enzyme immunoassay. The BIOLISA RUBELLA IgG DBS kit is a solid phase immunoenzymatic assay based on the principle of indirect qualitative and quantitative detection of IgG Antibodies to Rubella in dried human blood samples collected on filter paper. Rubella IgG Antibodies, present in the dried blood sample are eluted and bind to the Antigens coated on the microplate forming immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Anti-IgG human Antibodies conjugated to Peroxidase are added to the microplate, which is then incubated. The enzyme-conjugated Anti-IgG human Antibodies bind to the IgG Antibodies present, attached to the Antigen coated plate. A new wash is carried out to remove the surplus. After this step, the Substrate is added and incubated producing a blue color, which indicates the amount of Anti-Rubella IgG Antibodies present in the sample. Stop Solution is added to stop the reaction with a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS

- 1- Reference Standards (A - E)** - Store between 2 and 8°C. Contains: Five (5) vials (A - E) of Reference Standards containing IgG Anti-Rubella Antibodies at different concentrations in Buffer Solution containing surfactant, stabilizers, dye and preservative. **Potentially infectious.** The concentrations of the Reference Standards (A-E) vary from lot to lot. See vial labels.
- 2- Conjugate** - Store between 2 and 8°C. Contains: Buffer Solution, Anti-IgG human Antibody linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Sensitized Plate** - Store between 2 and 8°C. Contains: Plate sensitized with purified Rubella Antigens.
- 4- Concentrated Washing** - Store between 2 and 8°C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 5- Substrate** - Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: Hydrochloric Acid 1 M.
- 7- Elution Buffer** - Store between 2 and 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 8- DBS Negative Control** - Store between 2 and 8°C. Contains: Non-reactive sample for Anti-Rubella IgG Antibodies impregnated on filter paper. **Potentially infectious.**
- 9- DBS Positive Control** - Store between 2 and 8°C. Contains: Reactive sample for Anti-Rubella IgG Antibodies impregnated on filter paper. **Potentially infectious.**

PRESENTATION

Reagents	1	2	3
	96 cavities	192 cavities	480 cavities
1- Reference Standards (A – E)	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 2 mL
2- Conjugated	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Sensitized Plate	1 Unit	2 Units	5 Units
4- Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
7- Elution Buffer	1 Vial x 20 mL	1 Vial x 40 mL	1 Vial x 100 mL
8- DBS Negative Control	1 Unit	2 Units	3 Units
9- DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 100 to 500 µL with a coefficient of variation less than 1.5%.
- 2- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer for technique on filter paper.
- 4- ELISA reader with absorbance capacity at 450 and 620 nm of wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or watch.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 ± 2 °C.
- 11- Paper shredder (diameter 3 ± 0.2 mm).

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be between 2 and 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

1- For *in vitro* diagnostic use only.

- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate should only be opened after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free of contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw material of the product is tested and non-reactive for HBsAg and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take the proper biosafety precautions when handling these products.
- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.
- 9- As a protection measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the assay.
- 11- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Dried Blood collected on filter paper (DBS) - EDTA or Puncture

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use.

Blood samples dried on filter paper should be stored away from direct sunlight and in low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8 °C. Storage for longer than 2 years must be carried out at a temperature of -20 °C.⁹

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test prove that the BIOLISA RUBELLA IgG DBS kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial Nº 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, consult the SAC (Customer Advisory Service).

- Before starting the assay, place all Reagents, Samples and Controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Separate the wells to be used considering: Reference Standards (A - E), DBS Positive Control, DBS Negative Control and Dried Blood Spot collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.
 - 2- Separate the first cavity for the Blank (OPTIONAL).
 - 3- Add a 3mm disk of DBS Positive Control, DBS Negative Control and Dried Blood Spot collected on filter paper, previously perforated, in the determined wells.
 - 4- Pipette 100 µL and Reference Standards (A - E) into previously determined wells.
 - 5- Pipette 100 µL of Elution Buffer, including in the well for the Blank. Do not pipette into the Reference Standards wells (A - E).
 - 6- Gently homogenize for ± 30 seconds, cover the wells with plate sealer.
 - 7- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
 - 8- Remove the sealer from the cavities.
 - 9- After incubation, discard the contents of the wells by aspiration (Washer for filter paper technique). Use approximately 300 µL of **previously prepared** Washing Solution and perform a total of five (5) 5 second shake wash cycles. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
 - Note:** Poor washing/drying may cause poor results.
 - 10- Pipette 100 µL of Conjugate into all wells, **including the Blank well**.
 - 11- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
 - 12- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
 - 13- Remove the plate sealer from the wells.
 - 14- Repeat item 9.
 - 15- Pipette 100 µL of Substrate into all wells, **including the Blank well**.
 - 16- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
 - 17- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
 - 18- Remove the plate sealer from the wells.
 - 19- Pipette 100 µL of Stop Solution into all wells.
 - 20- Gently homogenize for ± 30 seconds.
 - 21- Read using a double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.100
Reference Standards A	< 0.100
Reference Standards B	> 0.200 and < 0.600
Reference Standards C	> B and < D
Reference Standards D	> C and < E
Reference Standards E	> 1.200
DBS Negative Control	< 0.150
DBS Positive Control	> 0.500

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Consider as Cut Off the mean absorbance obtained with Reference Standard B.

Example:

ITEM	ABSORBANCES
Reference Standard B	0.417
	0.419
Cut Off = Average Absorbance of Reference Standard B	Cut Off = (0.417 + 0.419) / 2 Cut Off = 0.418

Calculate the Index by dividing the Sample absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCES
Sample	1.346
Cut Off Value	0.418
Index = Sample / Cut Off Value	1.346 / 0.418 = 3.22

Note: Data shown in examples is for illustration only and cannot be used for calculation of results. Each laboratory must validate the Cut Off according to the instrumentation used and the researched population.

QUANTITATIVE

A calibration curve is used to determine the concentration of Rubella IgG Antibodies in unknown samples.

Preparation of the Calibration Curve

Record the absorbances obtained in the microplate reader, as shown in example 1. Calculate the averages of the duplicates (if duplicates are performed). Plot the mean absorbances of each Reference Standard versus the Corresponding concentration in IU/mL (described on the label) on graph paper (before plotting them on the graph) trace the curve.

Example 1:

Reference Standard	Concentration in IU/mL, see label	Absorbance	Average of duplicates
A	0	0.029	0.030
		0.031	
B	15	0.417	0.418
		0.419	
C	30	0.822	0.823
		0.824	
D	60	1.420	1.414
		1.408	
E	120	2.142	2.143
		2.144	

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE	QUANTITATIVE
	INDEX	CONCENTRATION
Negative (Non-Reagent)	≤ 0.8	≤ 12
Undetermined	0.8 and 1.2	12.0 – 18.0
Positive (Reagent)	≥ 1.2	≥ 18.0

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative, and samples with an index between 0.8 and 1.2 are considered indeterminate. Thus, we have that the sample mentioned in the example, whose absorbance was 1.346 and index of 3.22, presents a positive result (index ≥ 1.2).

The calculation of the concentration of the tested samples can also be performed using adequate computer programs, through linear regression.

Sample	Absorbance	Calculated concentration (IU/mL)	Quantitative Result
1	0.054	0.0	< 1.7 IU/mL
2	2.562	Above the upper limit of detection	> 120.0 IU/mL
3	0.242	7.4	74.0 IU/mL

Notes: In case of indeterminate results, the sample must be reanalyzed in duplicate. Samples that give repeatedly indeterminate results should be retested using an alternate method. If results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected shall prevail. The interpretation of a diagnostic test should not be based on a single test. Additional confirmatory tests must be included before a specimen is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, prior to making a descriptive diagnosis of the disease. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient. The product is not indicated for neonatal screening, as the detection of IgG antibodies in neonates has no clinical significance. The product is indicated for the analysis of whole blood samples on filter paper, and the dosage of anti-rubella IgG antibodies is usually indicated for the analysis of samples from women in the prenatal period.

INTERFERING

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A study was carried out with 52 dried blood samples collected on filter paper negative for Rubella IgG, but positive for other infections, in order to assess the possibility of cross-reactivity of these interferents in the result of BIOLISA RUBELLA IgG DBS. Among them 5 positive samples for HTLV, 10 positive samples for Syphilis, 6 positive samples for HBsAg, 5 positive samples for HCV, 7 positive samples for Chagas Disease, 6 positive samples for Zika, 4 positive samples for CMV, 4 positive samples for Toxoplasmosis and 5 HIV positive samples. No cross-reactivity was observed with positive samples for HTLV, Syphilis, HBsAg, HCV, Chagas disease, Zika, CMV, Toxoplasmosis and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	1.7699	0.8366	0.0225
Standard Deviation	0.0244	0.0169	0.0027
Coefficient of Variation (%)	1.3810	2.0193	11.898

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	1.7696	0.8404	0.0300
Standard Deviation	0.0271	0.0235	0.0029
Coefficient of Variation (%)	1.5349	2.7983	10.502

ANALYTICAL SENSITIVITY

The analytical sensitivity of the BIOLISA RUBELLA IgG DBS kit is 1.7 IU/mL.

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA RUBELLA IgG DBS kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA RUBELLA IgG DBS kit is 99.9% and the clinical specificity is 99.9%.

		Result Reference		
		Positive	Negative	Total
BIOLISA RUBELLA IgG DBS	Positive	151	0	151
	Negative	0	41	41
	Total	151	41	192

Clinical Sensitivity: >99.9% (151/151) CI 95% = 97.59 to 100%

Clinical Specificity: >99.9% (41/41) CI 95% = 91.40 to 100%

DIAGNOSTIC SIGNIFICANCE

Rubella is a spherical, small-enveloped RNA virus belonging to the Togaviridae family. It is commonly known as German or 3-day measles. Rubella virus infection is spread through droplets of saliva, resulting in a mild contagious rash in children or young adults. In childhood, the infection is a self-limiting, benign illness characterized by low-grade fever, headache, lymphadenopathy, arthralgia, and conjunctivitis. However, infection during pregnancy, especially in the first trimester, can lead to miscarriage, intrauterine infection causing fetal death, or congenital anomalies. Congenital Rubella depends on the period in which the infection occurs and can result in serious complications including deafness, eye problems including cataracts and glaucoma, congenital heart disease and mental retardation. Rubella IgM antibodies are produced initially, can reach detectable levels within 2 - 3 days and peak 14 - 21 days after the onset of symptoms which remain detectable for the next 4 - 8 weeks. The diagnosis of active or recent infection can be obtained by the presence of IgM antibodies in the initial sample. After several days, IgG antibodies appear after IgM, peaking at 14 - 21 days, with varying levels persisting throughout life. The presence of anti-Rubella IgG antibodies is indicative of previous infection and/or immunity.

BIBLIOGRAPHIC REFERENCES

- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.
- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- Cancer Epidemiol Biomarkers Prev; 26(8); 1337-44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA RUBELLA IgG DBS kit: 10269360434

Review: October/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER