BIOLISA RUBÉOLA IgG

REF K125

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FINALIDADE

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- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.

- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Herrmann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
- Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996), 275-286.
- Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
- Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
- Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.

- Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative and qualitative determination of IgG antibodies to Rubella virus in human samples of serum, plasma and whole blood on microplate enzyme immunoassay filter paper. Only for *in vitro* diagnostic use.

Methodology: Enzyme immunoassay or immunoenzyme
The BIOLISA Rubella IgG kit is an immunoenzymatic assay in solid phase based on the principle of indirect qualitative and quantitative detection of IgG Antibodies to Rubella virus in serum, plasma and whole blood on filter paper. IgG Anti-Rubella virus Antibodies present in the sample bind to the coated Antigens on the microplate, forming IgG Antigen-Antibody complexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated Anti-Human IgG Antibodies are added to the microplate and then incubated. Enzyme Antibodies-Human Anti-IgG Conjugate bind to IgG Antibodies immobilized by Antigens present. A new Wash is performed to remove the surplus. After this stage, the Substrate is added and incubated producing a blue color that indicates the amount of IgG Antibodies present in the samples. The Stop Solution is added to stop the reaction with a color change to yellow, measured in a microplate reader.

REAGENTS

1- Reference Standards (A - E) - Store between 2 and 8°C. Contains: Five (5) vials (A - E) of Reference Standards containing IgG Antibodies for Anti-Rubella in different concentrations in Buffer Solution, surfactant, stabilizers, dye and preservative.
Potentially infective.
The concentrations of the Reference Standards (A - E) vary with each batch. See vial label.
2- Conjugate - Store between 2 and 8°C. Contains: Buffer Solution, Human Anti-IgG Antibody linked to Peroxidase, surfactant, stabilizers, dye and preservative.
3- Sensitized Plate - Store at 2 and 8°C.
4- Concentrated Washing - Store at 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.
5- Sample Diluent - Store at 2 and 8°C. Contains: Buffer Solution, stabilizers, surfactant, chelator and preservative.
6- Substrate - Store at 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
7- Negative Control - Store at 2 and 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. Potentially infective.
8- Positive Control - Store between 2 and 8°C. Contains: Buffer Solution, IgG Anti-Rubella Antibodies, dye, stabilizers, surfactant and preservative. **Potentially infective.**
9- Stop Solution - Store at 2 and 8°C. Contains: Hydrochloric Acid 1 M.
10- Plate Sealers.
11- Filter Paper Elution Buffer - Store at 2 and 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
12- Negative Extraction Control - Store at 2 and 8°C. Contains: Non-reactive sample for Anti-Rubella IgG antibodies impregnated on Filter Paper. **Potentially infective.**
13- Positive Extraction Control - Store at 2 and 8°C. Contains: Reactive sample for Anti-Rubella IgG antibodies impregnated in Filter Paper. **Potentially infective.**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1 - Reference Standards (A - E)	1 Vial (A - E) x 300 µL	2 Vials (A - E) x 300 µL	5 Vials (A - E) x 300 µL
2 - Conjugate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
3 - Sensitized Plate	1 Unit (96 Cavities)	2 Units (192 Cavities)	5 Units (480 Cavities)
4 - Concentrated Washing	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
5 - Sample Diluent	1 Vial x 42 mL	2 Vials x 42 mL	5 Vials x 42 mL
6 - Substrate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
7 - Negative Control	1 Vial x 300 µL	2 Vials x 300 µL	5 Vials x 300 µL
8 - Positive Control	1 Vial x 300 µL	2 Vials x 300 µL	5 Vials x 300 µL
9 - Stop Solution	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
10 - Plate Sealers	3 Units	6 Units	15 Units
11- Filter Paper Elution Buffer	1 Vial x 20 mL	2 Vials x 20 mL	5 Vials x 20 mL
12- Negative Extration Control	1 Unit	2 Units	5 Units
13- Positive Extration Control	1 Unit	2 Units	5 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table.
- Operating instructions (manual).

Required materials not contained in the kit:

1- Pipette capable of dispensing volumes of 5 to 300 µL with lower coefficient of variation than 1.5%.
2- Re-pipettor for repetitive pipetting volumes of 300 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
3- Microplate washer (optional).
4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
5- Paper towel to dry cavities
6- Stopwatch or watch.
7- Vial to store the washing solution after diluted.
8- Distilled water.
9- Tools of Quality Control.
10- Incubator 37°C ± 2°C.
11- Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For in vitro diagnostic use only.
2- Strictly follow the methodology proposed to obtain accurate results.
3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
4- The water used in material cleaning must to be recent and free of contaminants.
5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.

7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.
8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
9- As a safety measure, you should cover the plate during the reaction.
10- You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
14- Do not use the product in case of damaged packaging.
15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

Whole Blood on Filter Paper

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Samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8°C. Storage for more than 2 years must be carried out the temperature of -20°C.

PROCESS DESCRIPTION

Stability After Open

The results of the stability test show that the BIOLISA Rubella IgG kit is stable after being opened for 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Reagent N°4 (Concentrated Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

Sample of Serum or Plasma

Before starting the assay, place all reagents, Reference Standards (A - E), Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.
1- Select the wells to be used considering: Reference Standards (A - E), Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate that will not be used for the original sealed packaging.
2- Select the first cavity for Blank (OPTIONAL).
3- Pipette 100 µL of Sample Diluent. Even in the cavity for Blank.
4- Pipette 5 µL of Reference Standards (A – E), Controls and Samples into previously determined wells. Observe the color change of the diluent at the time of sample addition. The color change indicates that the sample was added to the well correctly.
5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

7- Remove the sealing of wells.

8- Discard the contents of the wells by aspiration (Washer) or by decanting (Manual).

Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.

Note: Poor washing and drying can cause inadequate results.

9- Pipette 100 µL of Conjugate into each well, included in the Blank cavity.

10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

12- Remove the sealer from plate cavities.

13- Repeat item 8.

14- Pipette 100 µL of Substrate into all wells, even in the cavity for Blank.

15- Shake gently for ± 30 seconds. Cover wells with plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

17- Remove the sealer from plate cavities.

18- Pipette 100 µL of Stop Solution into all wells, even in the cavity for Blank.

19- Mix gently for ± 30 seconds.

20- Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

Whole Blood Sample on Filter Paper

Before starting the assay, place all reagents, Reference Standards (A - E), Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

1- Select the wells to be used considering: Reference Standards (A – E), Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
2- Select the first cavity for Blank (OPTIONAL).

3- Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper, previously perforated, in the determined wells.

4- Pipette 100 µL of Filter Paper Elution Buffer. Even in the cavity for Blank.

5- Pipette 5 µL of Reference Standards (A – E) and Controls into previously determined wells.

6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

7- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

8- Remove the sealing of wells.

9- Discard the contents of the wells by aspiration (Washer) or by decanting (Manual).

Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.

Note: Poor washing and drying can cause inadequate results.

10- Pipette 100 µL of Conjugate into each well, included in the Blank cavity.

11- Mix gently for ± 30 seconds. Cover cavities with plate sealer.

12- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

13- Remove the sealer from plate cavities.

14- Repeat item 9.

15- Pipette 100 µL of Substrate into all wells, even in the cavity for Blank.

16- Shake gently for ± 30 seconds. Cover wells with plate sealer.

17- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

18- Remove the sealer from plate cavities.

19- Pipette 100 µL of Stop Solution into all wells, even in the cavity for Blank.

20- Mix gently for ± 30 seconds.

21- Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION
Samples of Serum, Plasma and Whole Blood Sample on Filter Paper
Verify if the results obtained by the reading of the Blank, Reference Standards (A - E) and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.100
Reference Standard A	< 0.100
Reference Standard B	> 0.250 e < 0.600
Reference Standard C	> B e < D
Reference Standard D	> C e < E
Reference Standard E	> 1.400
Negative Control	< 0.100
Positive Control	> 1.000
Negative Extraction Control	< 0.300
Positive Extraction Control	> 1.000

If the values are out of the expected values. you must repeat the technique.

CALCULATIONS

QUALITATIVE
Serum, Plasma and Whole Blood Samples on Filter Paper
Cut Off regarded as the mean absorbance obtained with the Reference Standard B.

Example:

ITEM	ABSORBANCE
Reference Standard B	0.451
	0.449
Cut Off = Average Absorbance from Reference Standard B	Cut Off = (0.451 + 0.449) / 2 Cut Off = 0.450

Calculate the index by dividing the Sample absorbance by the value of the Cut Off.

Example:

ITEM	ABSORBANCE
Sample	1.346
Cut off Value	0.450
Index = Sample / Cut Off Value	Index = 1.346 / 0.450 Index = 2.99

QUANTITATIVE
A calibration curve is used to determine the concentration of Anti-Rubella Virus Antibodies in unknown samples.

Preparation of Calibration Curve
Serum and Plasma Technique
Record the absorbances obtained in the microplate reader, as shown in example 1. Calculate the averages of duplicates (if duplicates are performed). Plot the average absorbances of each Reference Standard versus the corresponding concentration in IU/mL (described on the label) on graph paper (before plotting them on the graph) plot the curve.

Whole Blood Technique on Filter Paper
Record the absorbances obtained in the microplate reader, as shown in example 1. Calculate the averages of duplicates (if duplicates are performed). Define the new concentration values for the filter paper technique by multiplying the concentration of each Reference Standard (described on the label) by a factor of 2.5.
Plot the average absorbances of each Reference Standard versus the corresponding concentration in IU/mL on graph paper (before plotting them on the graph) plot the curve.

Example 1

STANDARDS	ABSORBANCE	AVERAGE ABSORBANCE	CONCENTRACIÓN FOR SERUM AND PLASMA TECHNIQUE	CONCENTRACIÓN FOR TOTAL BLOOD TECHNIQUE ONFILTER PAPER
A	0.007	0.009	0	0 x 2.5 = 0
	0.011			
B	0.451	0.450	15	15 x 2.5 = 37.5
	0.449			
C	0.848	0.850	30	30 x 2.5 = 75
	0.852			
D	1.490	1.491	60	60 x 2.5 = 150
	1.493			
E	2.006	2.100	120	120 x 2.5 = 300
	2.195			

Samples that have an absorbance above the Reference Standard E, must be pre-diluted using Sample Diluent and must be retested. The concentration must be multiplied by the dilution factor. Automatic reading and calculation can be performed using the linear regression function in suitable computer programs.

Note: The data presented in the examples are for illustration only and cannot be used to replace the calibration curve, which must be constructed in the laboratory.

INTERPRETATION OF RESULTS
Serum, Plasma and Whole Blood Samples on Filter Paper

RESULTS	QUALITATIVE RESULTS FOR SERUM, PLASMA AND WHOLE BLOOD SAMPLES ON FILTER PAPER	QUANTITATIVE FOR SAMPLES OF SERUM OR PLASMA	QUANTITATIVE FOR SAMPLES OF WHOLE BLOOD ON FILTER PAPER
	INDEX	CONCENTRATION	CONCENTRATION
Negative	≤ 0.80	≤ 12	≤ 30
Positive	≥ 1.20	≥ 18	≥ 45
Indeterminate	0.81 – 1.19	12.1 – 17.9	30.1 – 44.9

Note: In case of indeterminate results, the sample must be retested. Samples that results have been obtained repeatedly indeterminate should be retested using an alternative method. If the results remain uncertain, one must collect a new sample in two weeks. Should prevail the result of the last sample collected. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/ or treatment of the patient.

PROCEDURE LIMITATIONS
The interpretation of a diagnostic test there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE
PRECISION
Repeatability
The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPETIBILITY	SAMPLE		
	1	2	3
Average	1.048	2.092	0.148
Standard Deviation	0.034	0.086	0.019
Coefficient of Variation (%)	3.215	4.106	12.580

Reproducibility
The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	1.049	2.092	0.148
Standard Deviation	0.007	0.010	0.003
Coefficient of Variation (%)	0.625	0.454	1.697

ANALYTICAL SENSITIVITY
The analytical sensitivity was calculated from the dosage of two samples diluted successively until the last concentration detected. The analytical sensitivity found was 0.161 IU/mL..

CLINICAL SENSITIVITY AND SPECIFICITY
Serum and Plasma Sample
BIOLISA Rubella IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA Rubella IgG kit is > 99.99%and clinical specificity is > 99.99%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	Biolisa Rubella IgG
Positive Sample	195	195
Negative Sample	50	50
Total Tested Sample	245	

Clinical Sensitivity: > 99.99% (195/195) IC 95% = 98.48 a 100 %
Clinical Specificity: > 99.99% (50/50) IC 95% = 94.12 a 100 %

Whole Blood Samples on Filter Paper
BIOLISA Rubella IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA Rubella IgG kit is > 99.99% and clinical specificity is > 99.99%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	Biolisa Rubella IgG
Positive Sample	95	95
Negative Sample	51	51
Total Tested Sample	146	

Clinical Sensitivity: > 99.99% (95/95) IC 95% = 96.90 a 100 %
Clinical Specificity: >99.99% (51/51) IC 95% = 94.23 a 100 %

LINEARITY
The reaction is able to detect concentrations up to the highest point on the calibration curve. For samples with higher values, dilute the same with Sample Diluent, repeat the dosage and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE
Rubella is an RNA virus, spherical, small envelope, belonging to the Togaviridae family. It is commonly known as German or 3 day measles. Rubella virus infection is transmitted through droplets of saliva, resulting in a mild contagious rash in children or young adults.
In childhood, infection is a benign, self-limiting disease, characterized by low fever, headache, lymphadenopathy, arthralgia and conjunctivitis. However, infection during pregnancy, especially in the first trimester, can lead to spontaneous abortion, intrauterine infection causing fetal death or congenital anomalies.
Congenital Rubella depends on the period in which the infection occurs and can result in serious complications, including deafness, eye problems, including cataracts and glaucoma, congenital heart disease and mental retardation. IgM Antibodies against Rubella are initially produced, reaching detectable levels within 2 - 3 days and peaking in 14 – 21 days after the onset of symptoms that remain detectable for the next 4 - 8 weeks.
The diagnosis of active or recent infection can be obtained by the presence of IgM Antibodies in an initial sample. After several days, IgG Antibodies appear after IgM, peaking at 14 - 21 days, with varying levels persisting throughout life. The presence of Anti-Rubella IgG Antibodies is indicative of previous infection and/or immunity.

BIBLIOGRAPHIC REFERENCES
1. Hermann, KL. Rubella Virus. In: Manual of Clinical Microbiology. American Society for Microbiology 4th Edition (1985) 779-784.
2. Turgeon, ML. Rubella Infection. In: Immunology and Serology in Laboratory Medicine. 2nd Edition (1996). 275-286.
3. Chernesky, MA, Mahony JB. Rubella Virus. In: Manual of Clinical Microbiology. 6th Edition (1995) 968-973.
4. Voller, A, Bidwell, DE, A simple Method for Detecting Antibodies to Rubella. Brit. J. Exp. Pathol. (1975) 56:338-339.
5. Rawls WE, Chernesky MA. Rubella Virus. Manual Clinical Immunology (1976) 452-455.

6. Millian, SJ, Wegman D. Rubella Serology: Applications, Limitations, and Interpretations. Amer. J. Pub. Health (1972) 170-176.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Rubella IgG kit: 10269360191

Review: September/2024

UNIVERSAL SYMBOLOGY

 REF	CATALOG NUMBER		MADE BY
 LOT	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
 VALIDITY DATE (last day of the month)			NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER