BIOLISA TOXOPLASMOSE IgG DBS

REF **K271**

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FINALIDADE

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- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
- Cancer Epidemiol Biomarkers Prev; 26(8): 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Indirect enzyme-linked immunosorbent assay (ELISA) for quantitative and qualitative determination of IgG antibodies to *Toxoplasma gondii* in dried blood samples collected on filter paper (DBS). Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzymeimmunoassay or enzyme immunoassay.

The BIOLISA TOXOPLASMOSIS IgG DBS kit is a solid phase enzyme immunoassay based on the principle of indirect qualitative and quantitative detection of IgG Antibodies to *Toxoplasma gondii* in human dried blood samples collected on filter paper. IgG antibodies to *Toxoplasma gondii*, present in the dried blood sample, are eluted and bind to the antigens coated on the microplate, forming IgG Antigen-Antibody immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Anti-human IgG antibodies conjugated to Peroxidase are added to the microplate, which is then incubated. The enzyme-conjugated Anti-Human IgG Antibodies bind to the IgG Antibodies present, bound to the plate coated with Antigens. New washing is carried out to remove excess. After this step, the Substrate is added and incubated, producing a blue color, which indicates the amount of Anti-*Toxoplasma gondii* IgG Antibodies present in the sample. The Stop Solution is added to stop the reaction and there is a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS

- 1- Reference Standards (A - E)** - Store between 2 and 8°C. Contains: Five (5) vials (A - E) of Reference Standards containing IgG Anti-*Toxoplasma gondii* Antibodies in different concentrations in buffer solution containing surfactant, stabilizers, dye and preservative. Potentially infectious. The concentrations of Reference Standards (A-E) vary with each batch. **See bottle labels.**
- 2- Conjugate** - Store between 2 and 8°C. Contains: Buffer solution, Anti-human IgG Antibody linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Sensitized Plate** - Store between 2 and 8°C. Contains: Plate sensitized with purified *Toxoplasma gondii* antigens.
- 4- Concentrated Washing** - Store between 2 and 8°C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 5- Substrate** – Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** – Store between 2 and 8°C. Contains: 1 M Hydrochloric Acid.
- 7- Elution Buffer** – Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 8- DBS Negative Control** – Store between 2 and 8°C. Contains: Non-reactive sample for IgG Anti-*Toxoplasma gondii* antibodies impregnated on filter paper. **Potentially infectious.**
- 9- DBS Positive Control** – Store between 2 and 8°C. Contains: Reactive sample for IgG Anti-*Toxoplasma gondii* antibodies impregnated on filter paper. **Potentially infectious.**

PRESENTATION

Reagents	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Reference Standards (A – E)	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 2 mL
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Sensitized Plate	1 Unit	2 Units	5 Units
4- Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
7- Elution Buffer	1 Vial x 20 mL	1 Vial x 40 mL	1 Vial x 100 mL
8- DBS Negative Control	1 Unit	2 Units	3 Units
9- DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 100 to 500 µL with a coefficient of variation of less than 1.5%.
- 2- Repipettor for repetitive pipetting of 300 µL volumes with a coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer for filter paper technique.
- 4- ELISA reader with absorbance capacity at 450 and 620 nm wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 °C ± 2 °C.
- 11- Paper shredder (diameter 3 ± 0.2 mm).

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8 °C. Transport at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate must only be opened after reaching room temperature. Return the unused microcavity strips to the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free from contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw materials for the product are tested and are non-reactive for HBsAg and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. The handling of any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take appropriate biosafety precautions when handling these products.
- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 9- As a protective measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.

- 11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation stages.
- 12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 13- To obtain information related to biosafety or in the event of accidents with the product, consult the SDS (Chemical Product Safety Information Sheet) available on the website www.bioclin.com.br or through a request through the SAC (Service Customer Advisory) from Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Dried blood collected on filter paper (DBS) - EDTA or Puncture

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use. Stability studies, carried out by Bioclin's Research and Development department, show that dried blood samples on filter paper can be stored for up to 35 days at room temperature, 2 years at 2-8 °C or 2 years at -20° W. The Manual "Blood Collection: Diagnosis and monitoring of STDs, AIDS and Viral Hepatitis: Brasília, Ministry of Health, Department of STD, AIDS and Viral Hepatitis. 2010, 98 p. (TELELAB Series)" ^[9], reports that dried blood samples on filter paper can be stored at room temperature, as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, samples must remain refrigerated, between 2 and 8 °C. Storage for more than 2 years must be carried out at a temperature of -20 °C.

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test prove that the BIOLISA TOXOPLASMOSIS IgG DBS kit is stable after being opened for up to 30 days. This stability may vary depending on test and environmental conditions. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of bottle No. 4 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, consult SAC (Customer Advisory Service).

- Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Separate the cavities to be used considering: Reference Standards (A - E), DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.
- 2- Separate the first cavity for the Blank (OPTIONAL).
- 3- Add a 3mm disc of DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper, previously perforated, into the determined cavities.
- 4- Pipette 100 µL and Reference Standards (A - E) into the previously determined wells.
- 5- Pipette 100 µL of Elution Buffer, **including** into the Blank cavity. Do not pipet into the Reference Standards wells (A - E).
- 6- Homogenize gently for ± 30 seconds, cover the cavities with plate sealer.
- 7- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 8- Remove the sealer from the cavities.

- 9- After incubation, discard the contents of the cavities by aspiration (Washer for filter paper technique). Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) washing cycles with 5-second shaking. To ensure the plate dries, at the end of washing, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause inadequate results.
- 10- Pipette 100 µL of Conjugate into all wells, **including** the Blank well.
- 11- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
- 12- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 13- Remove the plate sealer from the cavities.
- 14- Repeat item 9.
- 15- Pipette 100 µL of Substrate into all wells, **including** the Blank well.
- 16- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
- 17- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 18- Remove the plate sealer from the cavities.
- 19- Pipette 100 µL of Stop Solution into all wells.
- 20- Gently homogenize for ± 30 seconds.
- 21- Read using a double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE CHECK

Check whether the results obtained for reading the Blank and Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.100
Reference Standard A	< 0.100
Reference Standard B	> 0.200 and < 0.600
Reference Standard C	> B and < D
Reference Standard D	> C and < E
Reference Standard E	> 1.500
DBS Negative Control	< 0.150
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique must be repeated. The absorbances for the Controls and Reference Standards were obtained without decreasing the absorbance of the Blank.

CALCULATIONS

QUALITATIVE

Consider the average absorbance obtained with Reference Standard B as Cut Off.

Example:

ITEM	ABSORBANCE
Reference Standard B	0.356
	0.352
Cut Off = Average Absorbance of the Reference Standard B	Cut Off = (0.356 + 0.352) / 2 Cut Off = 0.354

Calculate the Index by dividing the sample absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	0.933
Cut Off Value	0.354
Index = Sample / Cut Off Value	0.933 / 0.354 = 2.635

Note: The data presented in the examples are for illustration only and cannot be used to calculate results. Each laboratory must validate the Cut Off according to the instrumentation used and the population researched.

QUANTITATIVE
A calibration curve is used to determine the concentration of IgG antibodies to Toxoplasmosis in unknown samples.

Preparation of the Calibration Curve
Record the absorbances obtained in the microplate reader, as shown in example 1. Calculate the averages of the duplicates (if duplicates are performed).
Plot the average absorbances of each Reference Standard (A-E) versus the corresponding Concentration in IU/mL (described on the label) on graph paper (before plotting them on the graph) and plot the curve.

Example 1

Standard Reference	Concentration in IU/mL, see label	Absorbance	Average duplicates
A	0	0.006	0.006
		0.007	
B	32	0.356	0.354
		0.352	
C	75	0.872	0.872
		0.873	
D	150	1.521	1.527
		1.534	
E	300	1.851	1.837
		1.824	

INTERPRETATION OF RESULTS
After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE	QUANTITATIVE
	INDEX	CONCENTRATION
Negative (Non-Reactive)	≤ 0.8	≤ 12
Undetermined	Between 0.8 and 1.2	Between 12.0 and 18.0
Positive (Reagent)	≥ 1.2	≥ 18.0

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative and samples with an index between 0.8 and 1.2 are considered indeterminate. Therefore, the sample mentioned in the example, whose absorbance was 0.933 and index of 2.635, presents a positive result (index ≥ 1.2).

The calculation of the concentration of the tested samples can also be carried out using suitable computer programs, through linear regression.

Sample	Absorbance	Calculated concentration (UI/mL)	Quantitative Result
1	0.117	Below the upper limit of detection	<13.5 UI/mL
2	3.256	Above the upper limit of detection	>300.0 UI/mL
3	1.528	150.5	150.5 UI/mL

Observations: In case of an indeterminate result, the sample must be re-analyzed in duplicate. Samples that yield repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected within two weeks. The result of the last sample collected must prevail. All results must be interpreted in conjunction with other available clinical information, before descriptive diagnosis of the disease. A negative result does not exclude the possibility of exposure.

The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PROCESS LIMITATIONS
The interpretation of a diagnostic test should not be established based on a single assay. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before descriptive diagnosis of the disease.
The product is not recommended for neonatal samples. The product is indicated for measuring IgG antibodies to Toxoplasma gondii in dried blood samples collected on filter paper from women in the prenatal period.

Traceability
The Biolisa Toxoplasmosis IgG DBS kit is traceable to the reference material NIBSC 01/600 (International Standard Anti-Toxoplasma IgG, Human – WHO International Standard.

INTERFERENCES
No interference was observed for concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY
A study was carried out with 78 dried blood samples collected on filter paper that were negative for Toxoplasmosis IgG, but positive for other infections, in order to evaluate the possibility of cross-reactivity of these interferences in the BIOLISA TOXOPLASMOSE IgG DBS result. Among them 7 positive samples for HTLV, 10 positive samples for Syphilis, 5 positive samples for HBsAg, 6 positive samples for HCV, 10 positive samples for Chagas Disease, 10 positive samples for Zika, 10 positive samples for CMV, 10 positive samples for Rubella and 10 positive samples for HIV. No cross-reactivity was observed with samples positive for HTLV, Syphilis, HBsAg, HCV, Chagas Disease, Zika, CMV, Rubella and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE
PRECISION
Repeatability
Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	1.0280	1.8303	0.0210
Standard deviation	0.0061	0.0139	0.0011
Coefficient of variation (%)	0.5890	0.7571	5.0195

Reproducibility
Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	1.0297	1.8051	0.0245
Standard deviation	0.0068	0.0194	0.0020
Coefficient of variation (%)	0.6641	1.0698	8.0353

Considering the variation of all samples tested in the precision study, it was observed that the kit presented a coefficient of variation of less than 20%.

ANALYTICAL SENSITIVITY
The analytical sensitivity of the BIOLISA TOXOPLASMOSIS IgG DBS kit is 13.5 IU/mL.

SENSITIVITY AND CLINICAL SPECIFICITY
The BIOLISA TOXOPLASMOSIS IgG DBS kit was used for the analysis of clinical samples previously characterized and confirmed using another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA TOXOPLASMOSIS IgG DBS kit is 98.38% and the clinical specificity is 98.70%.

BIOLISA TOXOPLASMOSIS IgG DBS	Reference Result		
	Positive	Negative	Total
Positive	61	1	62
Negative	1	76	77
Total	62	77	139

Clinical Sensitivity: 98.38% (61/62) CI 95% = 91.34 to 99.96%
Clinical Specificity: 98.70% (76/77) CI 95% = 92.98 to 99.97%

DIAGNOSTIC MEANING
Toxoplasma gondii is the causative agent of toxoplasmosis. It is an obligate intracellular protozoan that has been found in many species of birds, reptiles and mammals. The agent can be transmitted through organ transplantation, blood and leukocyte transfusion, contact with contaminated cat feces and ingestion of contaminated raw meat. In adults, the infection is usually benign or asymptomatic. However, symptomatic cases, including fatal cases, occur in immunosuppressed patients who have clinical or laboratory evidence of damage to the central nervous system. In children, the risk of fetal infection varies according to the length of pregnancy when the mother becomes infected. In maternal infections that occur during the first trimester, the likelihood of the infection passing to the fetus is lower. However, if transmission occurs, serious outcomes such as miscarriage and hydrocephalus are more likely. Infections acquired later in pregnancy, where fetal transmissions occur more frequently, tend to be less serious, but can still generate congenital manifestations, including brain calcifications and impaired learning. After infection, IgM antibodies appear within 5 days and show reduced levels within a few weeks or months. IgG antibodies generally appear 1 - 2 weeks after infection, reaching peak levels in 6 - 10 weeks and persisting for life.

BIBLIOGRAPHIC REFERENCES
1. Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
2. Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
3. McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
4. Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
5. Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
6. Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
7. Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
8. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
9. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
10. Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)
11. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

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