Bioclin

BIOLISA TOXOPLASMOSE IgG

REF **K127**

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FINALIDADE

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REFERÊNCIAS BIBLIOGRÁFICAS
1. Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
2. Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
3. McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
4. Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
5. Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
6. Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.

7. Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43..
8. TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Feldman, HA. *Toxoplasmosis: An Overview*. Acad. Med. (1974) 50:110-127.
- Remington, JS. *Toxoplasmosis in the Adult*. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. *Toxoplasma in and Around Us*. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.

6. Krick, JA, and Remington, JS. *Toxoplasmosis in the Adult – An Overview*. N. Engl. J. Med. (1978) 298(10):550-553.

7. Bryan, RT, and Wilson, M. *Toxoplasmosis*. Lab Management (1988) 26:40-43..

8. TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.

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
USAGE INSTRUCTIONS

FUNCTION

Test for quantitative and qualitative of IgG antibodies to *Toxoplasma gondii* in serum, plasma or whole blood on filter paper by enzyme immunoassay in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
BIOLISA Toxoplasmosis IgG kit is a solid phase enzyme immunoassay based on the principle of indirect qualitative and quantitative detection of IgG Antibodies to *Toxoplasma gondii* in human serum, plasma and whole blood samples on filter paper. IgG Antibodies to *Toxoplasma gondii*, in the sample bind to Antigen coated on the microplate forming Complexes Antigen-Antibodies IgG. After the initial incubation, the microplate is washed to remove unbound materials. Anti-human IgG Antibodies conjugated to Peroxidase are added to the microplate, which is then incubated. The Anti-human IgG Antibodies conjugated enzyme bind to IgG Antibodies present, linked to antigens. Further washing is performed to remove excess. After this step, the Substrate is added and incubated producing a blue color, which indicates the amount of IgG Antibodies to *Toxoplasma gondii* in the samples. Stop Solution is added to stop the reaction having a color change from blue to yellow, measured in a microplate reader.

REAGENTS

1- Reference Standards (A - E) - Store between 2 and 8°C. Contains: Five (5) flasks (A - E) of Reference Standards containing IgG Antibodies to Anti-*Toxoplasma gondii* in different concentrations in Buffer Solution, surfactant, stabilizers, dye and preservative. **Potentially infective.**

The concentrations of the Reference Standards (A - E) vary with each batch. See vial label.

2- Conjugate - Store between 2 and 8°C. Contains: Buffer Solution, Human Anti-IgG Antibody linked to Peroxidase, surfactant, stabilizers, dye and preservative.

3- Sensitized Plate - Store at 2 to 8°C.

4- Concentrated Washing - Store at 2 to 8°C. Contains: Buffer Solution, surfactant and preservative.

5- Sample Diluent - Store at 2 to 8°C. Contains: Buffer Solution, stabilizers, surfactant, chelator and preservative.

6- Substrate - Store at 2 to 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.

7- Negative Control - Store at 2 to 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**

8- Positive Control - Store between 2 and 8°C. Contains: Antibodies IgG Anti-*Toxoplasma gondii*, Buffer Solution, dye, stabilizers, surfactant and preservative. **Potentially infective.**

9- Stop Solution - Store at 2 to 8°C. Contains: Hydrochloric Acid 1 M.

10- Plate Sealers

11- Filter Paper Elution Buffer - Store at 2 to 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.

12- Negative Extraction Control - Store at 2 to 8°C. Contains: Non-reactive sample for anti-*Toxoplasma gondii* IgG antibodies impregnated in filter paper. **Potentially infective.**

13- Positive Extraction Control - Store at 2 to 8°C. Contains: Reactive sample for anti-*Toxoplasma gondii* IgG antibodies impregnated in filter paper. **Potentially infective.**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 – Reference Standards (A – E)	1 Flask (A – E) x 300 µL.	2 Flasks (A – E) x 300 µL.	5 Flasks (A – E) x 300 µL.
2 - Conjugate	1 Flask x 12mL.	2 Flasks x 12mL.	5 Flasks x 12mL.
3 – Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
4 –Concentrate Washing	1 Flask x 50mL.	2 Flasks x 50mL.	5 Flasks x 50mL.
5 – Sample Diluent	1 Flask x 42mL.	2 Flasks x 42mL.	5 Flasks x 42mL.
6 - Substrate	1 Flask x 12mL.	2 Flasks x 12mL.	5 Flasks x 12mL.
7 –Negative Control	1 Flask x 300 µL.	2 Flasks x 300 µL.	5 Flasks x 300 µL.
8 –Positive Control	1 Flask x 300 µL.	2 Flasks x 300 µL.	5 Flasks x 300 µL.
9 – Stop Solution	1 Flask x 12mL.	2 Flasks x 12mL.	5 Flasks x 12mL.
10 – Plate Sealers	3 Units	6 Units	15 Units
11 – Filter Paper Elution Buffer	1 Flask x 20 mL.	2 Flasks x 20 mL.	5 Flasks x 20 mL.
12 – Negative Extraction Control	1 unit	2 units	5 units
13 – Positive Extraction Control	1 unit	2 units	5 units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

- Operating instructions (manual)

Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes of 5 to 300 µL with lower coefficient of variation than 1,5%.
- 2- Re-pipettor for repetitive pipetting volumes of 300 µL, with lower coefficient of variation than 1,5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Stopwatch or watch.
- 7- Flask to store the washing solution after diluted.
- 8- Distilled water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.
- 11- Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.

4- The water used in material cleaning must to be recent and free of contaminants.

5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.

6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.

7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.

8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.

9- As a safety measure, you should cover the plate during the reaction.

10- You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test;

11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps;

12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

14- Do not use the product in case of damaged packaging.

15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

Whole Blood on Filter Paper

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Samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8°C. Storage for more than 2 years must be carried out the temperature of -20°C.

PROCESS DESCRIPTION

The results of the stability test show that the kit Biolisa Toxplasmose IgG is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Reagent N°4 (Concentrated Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

Samples of Serum or Plasma

Before starting the assay, place all reagents, Reference Standards (A - E), Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

1- Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used Reference Standards (A – E), for the original sealed packaging.

2- Select the first cavity for Blank (OPTIONAL).

3- Pipette 100 µL of Sample Diluent, included in the Blank cavity.

4- Pipette 5 µL of Reference Standards (A – E), Sample and Controls into previously determined wells. Observe the color change of the diluent at the time of sample addition. The color change indicates that the sample was added to the well correctly.

5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

7- Remove the sealing of wells.

8- Discard the contents of the wells by aspiration (Washer) or by decanting (manual).

Use approximately 300 µl of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.

Note: Poor washing and drying can cause inadequate results.

9- Pipette 100 µl of Conjugate into each well except in the Blank cavity.

10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

12- Remove the sealer from plate cavities.

13- Repeat item 8.

14- Pipette 100 µl of Substrate into all wells, even in the cavity for White included in the Blank cavity.

15- Shake gently for ± 30 seconds. Cover wells with plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

17- Remove the sealer from plate cavities.

18- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.

19- Mix gently for ± 30 seconds.

20- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

Whole Blood Sample on Filter Paper

Before starting the assay, place all reagents, Reference Standards (A - E), Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

1- Select the wells to be used considering: Reference Standards (A – E), Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.

2- Select the first cavity for Blank (OPTIONAL).

3- Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Sample on Filter Paper, previously perforated, in the determined wells.

4- Pipette 100 µL of Elution Buffer in Filter Paper, included in the Blank cavity.

5- Pipette 5 µL of Reference Standards (A – E) and Controls into previously determined wells.

6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

7- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

8- Remove the sealing of wells.

9- Discard the contents of the wells by decanting (manual). Remove the paper discs, if necessary, with the aid of a needle or tip.
Use approximately 300 µl of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.

Note: Poor washing and drying can cause inadequate results.

10- Pipette 100 µl of Conjugate into each well, included in the Blank cavity.

11- Mix gently for ± 30 seconds. Cover cavities with plate sealer.

12- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

13- Remove the sealer from plate cavities.

14- Repeat item 9.

15- Pipette 100 µl of Substrate into all wells, included in the Blank cavity

16- Shake gently for ± 30 seconds. Cover wells with plate sealer.

17- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

18- Remove the sealer from plate cavities.

19- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity

20- Mix gently for ± 30 seconds.

21- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Samples of Serum, Plasma and Filter Paper Whole Blood

Verify if the results obtained by the reading of the Controls are compatible to the values below:

ITEM	ABSORBANCE
Blanco	< 0,100
Reference Standard A	< 0,100
Reference Standard B	> 0,2 e < 0,6
Reference Standard C	> B e < D
Reference Standard D	> C e < E
Reference Standard E	> 1,500
Negative Control	< 0,100
Positive Control	> 1,000
Negative Extraction Control	< 0,300
Positive Extraction Control	> 1,000

If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Samples of Serum or Plasma

Cut Off regarded as the mean absorbance obtained with the Reference Standard B.

Example:

ITEM	ABSORBANCE
Reference Standard B	A1 = 0,356
	A2 = 0,352
Cut Off = Average Absorbance from Reference Standard B	Cut Off = (0,356 + 0,352) / 2 Cut Off = 0,354

Calculate the index by dividing the Sample absorbance by the value of the Cut Off.

Example:

ITEM	ABSORBANCE
Sample	0,933
Cut Off Value	0,354
Index = Sample / Cut Off Value	Index = 0,933 / 0,354 Index = 2,635

Note: The data presented in the examples are for illustration only and can not be used to calculate the results.

QUANTITATIVE

A calibration curve is used to determine the concentration of anti-*Toxoplasma gondii* antibodies in unknown samples.

Preparation of Calibration Curve

Serum or Plasma Technique

Record the absorbance obtained in the microplate reader as outlined in example 1. Calculate the averages of duplicates (if they are carried duplicates). Plot the absorbances of each Reference Standard *versus* the corresponding concentration in UI/mL on graph paper (before you plot them on the graph.) Plot the curve.

Note: The concentration of the Reference Standards may vary from batch to batch.

Whole Blood Technique on Filter Paper

Record the absorbances obtained in the microplate reader, as shown in example 1. Calculate the averages of duplicates (if duplicates are performed).

Define the new concentration values for the filter paper technique by multiplying the concentration of each Reference Standard (described on the label) by a factor of 2.5. *

* This value may vary according to the lot.

Plot the average absorbances of each Reference Standard versus the corresponding concentration in IU / mL on graph paper (before plotting them on the graph) plot the curve.

Note: The concentration of the Reference Standards may vary from batch to batch.

Example 1

STANDARDS	ABSORBANCE	AVERAGE ABSORBANCE	CONCENTRACIÓN FOR SERUM AND PLASMA TECHNIQUE	CONCENTRACIÓN FOR TOTAL BLOOD TECHNIQUE ONFILTER PAPER
A	0,006	0,006	0,0	0 x 2,5 = 0,0
	0,007			
B	0,356	0,354	32,0	32 x 2,5 = 80,0
	0,352			
C	0,872	0,872	75,0	75 x 2,5 = 187,5
	0,873			
D	1,521	1,527	150,0	150 x 2,5 = 375,0
	1,534			
E	1,851	1,837	300,0	300 x 2,5 = 750,0
	1,824			

Samples that have absorbance above the Reference Standard E, must be pre-diluted using Sample Diluent and should be retested. The concentration must be multiplied by the dilution factor. Automatic reading and calculation can be performed through the point to point function in appropriate programs of computer.

Note: The data presented in the examples are for illustration only and may not be used in lieu the calibration curve, which must be constructed in the laboratory.

RESULTS

Serum, Plasma and Whole Blood Samples on Filter Paper

RESULTS	QUALITATIVE RESULTS FOR SERUM, PLASMA AND WHOLE BLOOD SAMPLES ON FILTER PAPER	QUANTITATIVE FOR SAMPLES OF SERUM OR PLASMA	QUANTITATIVE FOR SAMPLES OF WHOLE BLOOD ON FILTER PAPER
	INDEX	CONCENTRATION	CONCENTRATION
Negative	≤ 0,90	≤ 28	≤ 72
Positive	≥ 1,0	≥ 32	≥ 80
Indeterminate	0,91 – 0,99	29 - 31	72,1 – 79,9

Note: In case of indeterminate results, the sample must be retested. Samples that results have been obtained repeatedly indeterminate should be retested using an alternative method. If the results remain uncertain, one must collect a new sample in two weeks. Should prevail the result of the last sample collected.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test there should not be based on a single run. This should include confirmation of other tests before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPETIBILITY	SAMPLE		
	1	2	3
Average	1,560	1,217	0,051
Standard Deviation	0,031	0,032	0,003
Coefficient of Variation (%)	2,026	2,694	6,158

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	2,226	2,254	0,045
Standard Deviation	0,106	0,131	0,002
Coefficient of Variation (%)	4,751	5,811	4,660

Clinical sensitivity and specificity

Serum and Plasma Sample

BIOLISA Toxoplasmosis IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA Toxoplasmosis IgG kit is 99,46% and clinical specificity is 99,16%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	BIOLISA Toxoplasmosis IgG
Positive Sample	188	187
Negative Sample	120	119
Total Tested Sample	308	

Clinical Sensitivity: 99,46% (187/188)

Clinical Specificity: 99,16% (119/120)

Whole Blood Samples on Filter Paper

BIOLISA Toxoplasmosis IgG Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA Toxoplasmosis IgG kit is 98,38% and clinical specificity is 98,70%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	BIOLISA Toxoplasmosis IgG
Positive Sample	62	61
Negative Sample	77	76
Total Tested Sample	139	

Clinical Sensitivity: 98,38% (61/62)

Clinical Specificity: 98,70% (76/77)

Linearity

The reaction is able to detect concentrations up to the highest point on the calibration curve. For samples with higher values, dilute the same with Sample Diluent, repeat the dosage and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Toxoplasma gondii is the causative agent of toxoplasmosis. It is an obligate intracellular protozoan that has been found in many species of birds, reptiles and mammals. The agent can be transmitted through organ transplants, blood transfusions and leukocyte, contact with cat feces contaminated and ingestion of contaminated raw meats. In adults, the infection is usually benign or asymptomatic. However, symptomatic cases, including fatal cases occur in patients immunosuppressed patients who have clinical or laboratory evidence of central nervous system damage. In children, the risk of fetal infection varies according to the time of pregnancy when the mother is infected. In maternal infections that occur during the first quarter, the probability of infection through to the fetus is smaller. However, if transmission occurs, severe outcomes such as miscarriage and hydrocephalus are more likely. Infections acquired later in pregnancy, where the transmissions fetal occur most frequently tend to be less severe, but even they can generate congenital manifestations, including brain calcification and impaired learning. After infection, IgM antibodies appear 5 days and had reduced levels within a few weeks or months. IgG antibodies usually appear 1 - 2 weeks after infection, reaching levels peak at 6 - 10 weeks persisting for life.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. *Toxoplasma in and Around Us*. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.

- Krick, JA, and Remington, JS. *Toxoplasmosis in the Adult – An Overview*. N. Engl. J. Med. (1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. *Toxoplasmosis*. Lab Management (1988) 26:40-43.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for Biolisa Toxoplasmosis IgG kit: 10269360196

Review: May/2020

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED