BIOLISA TOXOPLASMOSE IgM DBS

REF **K270**

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FINALIDADE

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- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
- TELELAB, Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med.(1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43..
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Immunoenzymatic test capture (ELISA) for determination of IgM class Antibodies Anti-*T. gondii* in dried blood spot collected on filter paper (DBS). Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA TOXOPLASMOSIS IgM DBS kit is a solid-phase immunoenzymatic assay based on the principle of qualitative detection by capturing Antibodies of the IgM Anti-*Toxoplasma gondii* class in human samples of dried blood collected on filter paper. IgM antibodies present in the sample bind to Anti-IgM Antibodies coated on the microplate forming immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated *Toxoplasma gondii* antigens are added to the microplate which is then incubated. The enzyme-conjugated Antigens bind to the Anti-*Toxoplasma gondii* IgM Antibodies present, bound to the plate coated with Anti-IgM Antibodies. A new wash is carried out to remove the surplus. After this step, the Substrate is added, producing a blue color that indicates the amount of Anti-*Toxoplasma gondii* IgM Antibodies present in the sample. Stopping Solution is added to stop the reaction with a color change from blue to yellow. Color intensity is measured at 450 and 620 nm.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8°C. Contains: Plate sensitized with Anti-IgM Antibodies.
- 2- Conjugate** – Store between 2 and 8°C. Contains: Buffer Solution, *Toxoplasma gondii* Antigens conjugated to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Wash** – Store between 2 and 8°C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Substrate** – Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 5- Stop Solution** – Store between 2 and 8°C. Contains: Hydrochloric Acid 1 M.
- 6- Elution Buffer** – Store between 2 and 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 7- DBS Negative Control** – Store between 2 and 8°C. Contains: Non-reactive sample for Anti-*Toxoplasma gondii* IgM Antibodies impregnated on filter paper. **Potentially infectious.**
- 8- DBS Positive Control** – Store between 2 and 8°C. Contains: Reactive sample for Anti-*Toxoplasma gondii* IgM Antibodies impregnated on filter paper. **Potentially infectious.**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 – Sensitized Plate	1 Unit	2 Units	5 Units
2 - Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3 –Concentrated Wash	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4 – Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
5 - Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6 – Elution Buffer	1 Vial x 20 mL	1 Vial x 40 mL	1 Vial x 100 mL
7 – DBS Negative Control	1 Unit	2 Units	3 Units
8 – DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

Required materials not contained in the Kit:

- 1- Pipettes capable of dispensing volumes from 100 to 500 µL with a coefficient of variation less than 1.5%.
- 2- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (OPTIONAL).
- 3- Microplate washer for technique on filter paper.
- 4- ELISA reader with absorbance capacity at 450 and 630 nm of wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 ± 2 °C.
- 11- Paper shredder (diameter 3 ± 0.2 mm).

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store at 2 to 8°C.
- 4- The water used in material cleaning must be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Hydrochloric Acid, which is a strong acid. Handle it with care.
- 7- All the raw material of product is tested and should be nonreactive for HBsAg, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.
- 8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9- As a safety measure, you should cover the plate during the reaction.
- 10- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.

- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Dried Blood collected on filter paper (DBS) - EDTA or Puncture

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use. Blood samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8 °C. Storage for a period longer than 2 years must be carried out at a temperature of -20°C.

DESCRIPTION OF PROCESS

Stability After Opening

The results of the stability test prove that the BIOLISA TOXOPLASMOSIS IgM DBS kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial N° 3 (Concentrated Wash) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Wash solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, contact SAC (Customer Advice Service).

- Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Separate the wells to be used considering: DBS Positive Control, DBS Negative Control and dried blood spot collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.
- 2- Separate the first cavity for the Blank (OPTIONAL).
- 3- Add a 3mm disk of DBS Positive Control, DBS Negative Control and dried blood spot collected on filter paper, previously perforated, in the determined wells.
- 4- Pipette 100 µL of Elution Buffer into all wells. **Including the Blank well.**
- 5- Gently homogenize for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- After incubation, discard the contents of the wells by aspiration (Washer for filter paper technique). Use approximately 300 µL of **previously prepared** Wash Solution and perform a total of five (5) 5 second shake wash cycles. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause poor results.
- 8- Pipette 100 µL of Conjugate into all wells. **Including the Blank well.**
- 9- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 10- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 11- Remove the plate sealer from the wells.
- 12- Repeat item 7.

- 13- Pipette 100 µL of Substrate into all wells. **Including in the cavity of Blank.**
- 14- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 15- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 16- Remove the plate sealer from the wells.
- 17- Pipette 100 µL of Stop Solution into all wells.
- 18- Gently homogenize for ± 30 seconds.
- 19- Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.150
DBS Negative Control	< 0.150
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the formula below:

Cut Off = DBS Negative Control Mean Absorbance + 0.100.

Example:

ITEM	ABSORBANCE
DBS Negative Control	0.020
	0.022
Cut Off = DBS Negative Control Mean Absorbance + 0.100	((0.020 + 0.022) / 2) + 0.100 = 0.120

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.352
Cut Off Value	0.120
Index = Sample/ Cut Off Value	1.352 / 0.120 = 11.27

Note: Data shown in examples is for illustration only and cannot be used for calculation of results. Each laboratory must validate the Cut Off according to the instrumentation used and the researched population.

INTERPRETATION OF RESULTS

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (No Reagent)	≤ 0.8
Undetermined	Between 0.8 and 1.2
Positive (Reagent)	≥ 1.2

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative, and samples with an index between 0.8 and 1.2 are considered indeterminate. Thus, the sample mentioned in the example, whose absorbance was 1.352 and index 11.27, presents a positive result (index ≥ 1.2).

Notes: In case of indeterminate results, the sample must be reanalyzed in duplicate. Samples that give repeatedly indeterminate results should be retested using an alternate method. If results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected shall prevail. The interpretation of a diagnostic test should not be based on a single test. Additional confirmatory tests must be included before a specimen is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, prior to making a descriptive diagnosis of the disease. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was performed with 48 dried blood samples collected on filter paper negative for Toxoplasmosis IgM, but positive for other infections, in order to assess the possibility of cross-reactivity of these interferents in the result of BIOLISA TOXOPLASMOSIS IgM DBS. Among them 5 samples positive for HTLV, 6 samples positive for Syphilis, 6 samples positive for HBsAg, 5 samples positive for HCV, 7 samples positive for Chagas Disease, 6 samples positive for Zika, 4 samples positive for CMV, 4 samples positive for Rubella and 5 HIV positive samples. No cross-reactivity was observed with positive samples for HTLV, Syphilis, HBsAg, HCV, Chagas Disease, Zika, CMV, Rubella and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

ACCURACY

Repeatability

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	1.3081	0.6429	0.0425
Standard Deviation	0.0480	0.0300	0.0027
Coefficient of Variation (%)	3.6667	4.6622	6.2990

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	1.2985	0.6553	0.0433
Standard Deviation	0.0597	0.0293	0.0030
Coefficient of Variation (%)	4.5930	4.4741	6.9427

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA TOXOPLASMOSIS IgM DBS kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA TOXOPLASMOSIS IgM DBS kit is 98.82% and the clinical specificity is 98.98%.

		Reference Result		
		Positive	Negative	Total
BIOLISA TOXOPLASMOSIS IgM DBS	Positive	84	1	85
	Negative	1	97	98
	Total	85	98	183

Clinical Sensitivity: 98.82% (84/85) CI 95% = 93.60 to 100%

Clinical Specificity: 98.98% (97/98) CI 95% = 94.40 to 100%

DIAGNOSTIC SIGNIFICANCE

Toxoplasma gondii is the causative agent of toxoplasmosis. It is an obligate intracellular protozoan that has been found in many species of birds, reptiles and mammals. The agent can be transmitted through organ transplantation, blood and leukocyte transfusions, contact with contaminated cat feces and ingestion of contaminated raw meat. In adults, the infection is usually benign or asymptomatic. However, symptomatic cases, including fatal cases, occur in immunosuppressed patients who have clinical or laboratory evidence of damage to the central nervous system. After infection, IgM antibodies appear within 5 days and have reduced levels within a few weeks or months. IgG antibodies usually appear 1 - 2 weeks after infection, reaching peak levels at 6 - 10 weeks and persisting for life. In children, the risk of fetal infection varies according to the time of pregnancy, when the mother is infected. In maternal infections that occur during the first trimester, the likelihood of the infection passing to the fetus is lower. However, if transmission occurs, serious outcomes such as miscarriage and hydrocephalus are more likely. Infections acquired later in pregnancy, where fetal transmissions occur more frequently, tend to be less severe, but they can still generate congenital manifestations, including brain calcifications and impaired learning.

BIBLIOGRAPHIC REFERENCES

- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA TOXOPLASMOSIS IgM DBS kit:
10269360387

Review: April/2023

UNIVERSAL SYMBOLOGY

 REF	CATALOG NUMBER	 MADE BY
 LOT	LOT NUMBER	 CONTROL
 MANUFACTURING DATE		 POSITIVE CONTROL
 VALIDITY DATE (last day of the month)		 NEGATIVE CONTROL
 TEMPERATURE LIMIT (store)		 BIOLOGICAL RISK
 CONTENT IS SUFFICIENT FOR <N> TEST		 FLAMMABLE
 SEE INSTRUCTIONS FOR USE		 CORROSIVE
 IN VITRO DIAGNOSTIC PRODUCT		 TOXIC
 KEEP AWAY FROM SUNLIGHT		 DO NOT USE IF PACKAGE IS DAMAGED
 DO NOT REUSE		 PRODUCT STERILIZED
 CAUTION		 DANGER