Bioclin

BIOLISA TOXOPLASMOSE IgM

REF **K126**

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FINALIDADE

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NÚMERO DE TESTES

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- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
- Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
- McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
- Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
- Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
- Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med.(1978) 298(10):550-553.
- Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43..
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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BIOLISA TOXOPLASMOSIS IgM
REF K126
USAGE INSTRUCTIONS

FUNCTION

Test for qualitative determination of IgM antibodies to *Toxoplasma gondii* in serum, plasma or whole blood on filter paper by enzyme immunoassay in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA TOXOPLASMOSIS IgM kit is a solid phase immunoenzymatic assay based on the principle of qualitative detection by capture of IgM Antibodies for *Toxoplasma gondii* in human serum, plasma and whole blood samples on filter paper. IgM antibodies present in the sample bind to the anti-IgM Antibody coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. *Toxoplasma gondii* antigens conjugated to Peroxidase are added to the microplate which is then incubated. Enzyme-conjugated antigens bind to the Anti-*Toxoplasma gondii* IgM Antibodies present, attached to the plate coated with anti-IgM antibodies. New washing is performed to remove surplus. After this step, the Substrate is added and incubated, producing a blue color indicating the amount of Anti-*Toxoplasma gondii* IgM Antibodies present in the sample. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store at 2 to 8°C.
- 2- Conjugate** - Store between 2 and 8°C. Contains: Buffer Solution, Antigens of Toxoplasma gondii conjugated to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store at 2 to 8°C. Contains: Buffer Solution, surfactant and preservative.
- 4- Sample Diluent** - Store at 2 to 8°C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
- 5- Substrate** - Store at 2 to 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store at 2 to 8°C. Contains: Hydrochloric Acid 1 M.
- 7- Negative Control** - Store at 2 to 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: IgM Antibodies *Toxoplasma gondii*, Buffer Solution, dye, stabilizers, surfactant and preservative. **Potentially infective.**
- 9- Plate Sealers.**
- 10- Elution Buffer on Filter Paper** - Store at 2 to 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 11- Negative Extraction Control** - Store at 2 to 8°C. Contains: Non-reactive sample for impregnated IgM anti-*Toxoplasma gondii* antibodies on filter paper. **Potentially infective.**
- 12- Positive Extraction Control** - Store at 2 to 8°C. Contains: Reactive sample for IgM anti-*Toxoplasma gondii* antibodies impregnated in filter paper. **Potentially infective**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 – Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2 - Conjugate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
3 –Concentrate Washing	1 Flask x 50mL	2 Flasks x 50mL	5 Flasks x 50mL
4 – Sample Diluent	1 Flask x 42mL	2 Flasks x 42mL	5 Flasks x 42mL
5 - Substrate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
6 – Stop Solution	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
7 –Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
8 –Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9 – Plate Sealers	3 Units	6 Units	15 Units
10 – Filter Paper Elution Buffer	1 Flask x 20 mL	2 Flasks x 20 mL	5 Flasks x 20 mL
11 – Negative Extraction Control	1 unit	2 units	5 units
12 – Positive Extraction Control	1 unit	2 units	5 units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

Required materials not contained in the Kit:

- 1- Pipette capable of dispensing volumes from 5 to 300 µL with lower coefficient of variation than 1,5%.
- 2- Re-pipettor for repetitive pipetting volumes of 300 µL with lower coefficient of variation than 1,5% or multicontrol pipette (Optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Stopwatch or watch.
- 7- Flask to store the Washing Solution after diluted.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.
- 11- Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store at 2 to 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Chloridric Acid, which is a strong acid. Handle it with care.

7- All the raw material of product is tested and should be nonreactive for HBsAg, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.

8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.

9- As a safety measure, you should cover the plate during the reaction.

10- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.

11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.

12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

14- Do not use the product in case of damaged packaging.

15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

Whole Blood on Filter Paper

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Samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8°C. Storage for more than 2 years must be carried out the temperature of -20°C.

DESCRIPTION OF PROCESS

The results of the stability test show that the kit Biolisa Toxoplasmos IgM is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Reagent N°3 (Concentrated Washing Solution) in 1000 mL of distilled water or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

Samples of Serum or Plasma

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

- 1- Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used, for the original sealed packaging.
- 2- Select the first cavity for Blank (OPTIONAL).
- 3- Pipette 100 µL of Sample Diluent, included in the Blank cavity.
- 4- Pipette 5 µL of Controls and Samples into previously determined wells. Observe the color change of the diluent at the time of sample addition. The color change indicates that the sample was added to the well correctly.

- 5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 7- Remove the sealing of wells.
- 8- Discard the contents of the wells by aspiration (Washer) or by decanting (manual).
- Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
- Note:** Poor washing and drying can cause inadequate results.
- 9- Pipette 100 µL of Conjugate into each well except in the Blank cavity (if you made this option).
- 10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 12- Remove the sealer from plate cavities.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells, included in the Blank cavity.
- 15- Shake gently for ± 30 seconds. Cover wells with plate sealer.
- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 17- Remove the sealer from plate cavities.
- 18- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.
- 19- Mix gently for ± 30 seconds.
- 20- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

Whole Blood Sample on Filter Paper

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

- 1- Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
- 2- Select the first cavity for Blank (OPTIONAL).
- 3- Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Sample on Filter Paper, previously perforated, in the determined wells.
- 4- Pipette 100 µL of Elution Buffer in Filter Paper, included in the Blank cavity.
- 5- Pipette 5 µL of Controls into previously determined wells.
- 6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 7- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 8- Remove the sealing of wells.
- 9- Discard the contents of the wells by decanting (manual). Remove the paper discs, if necessary, with the aid of a needle or tip.
- Use approximately 300 µL of Washing Solution, **previously prepared** to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
- Note:** Poor washing and drying can cause inadequate results.
- 10- Pipette 100 µL of Conjugate into each well, included in the Blank cavity.
- 11- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
- 12- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 13- Remove the sealer from plate cavities.
- 14- Repeat item 9.
- 15- Pipette 100 µL of Substrate into all wells, included in the Blank cavity.
- 16- Shake gently for ± 30 seconds. Cover wells with plate sealer.
- 17- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 18- Remove the sealer from plate cavities.
- 19- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.
- 20- Mix gently for ± 30 seconds.
- 21- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Samples of Serum, Plasma and Whole Blood Sample on Filter Paper
Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0,150
Negative Control	< 0,150
Positive Control	> 1,000
Negative Extraction Control	< 0,300
Positive Extraction Control	> 1,000

If the values are out of the expected values, you must repeat the technique.

CALCULATIONS

QUALITATIVE

Sample of Serum or Plasma

Calculate Cut Off according to the formula below:

Cut Off = (Positive Control Mean Abs. * 0,1) + 0,100*

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N°8)	A1 = 1,780
	A2= 1,700
Cut Off = (Positive Control Mean Abs. x 0,1) + 0,100*	Cut Off = ((1,780 + 1,700/2) x 0,1) + 0,100* Cut Off = 0,274

* This factor may vary with each batch.

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1,245
Cut Off Value	0,274
Index = Sample/ Cut Off Value	Index = 1,245 / 0,274 Index = 4,544

Whole Blood Sample on Filter Paper

Calculate Cut Off according to the formula below:

Cut Off = (Positive Control Mean Absorbance x 0,1) + 0,200*

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N°8)	A1 = 1,542
	A2 = 1,529
Cut Off = (Positive Control Mean Abs. x 0,1) + 0,200*	Cut Off = ((1,542 + 1,529/2) x 0,1 + 0,200* Cut Off = 0,353

* This factor may vary with each batch.

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1,457
Cut Off Value	0,353
Index = Sample/ Cut Off Value	Index = 1,457 / 0,353 Index = 4,127

INTERPRETATION OF RESULTS

Samples of Serum, Plasma and Filter Paper Whole Blood

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (No Reagent)	≤ 0,8
Undetermined	Between 0,8 and 1,2
Positive (Reagent)	≥ 1,2

The results provided by this kit should be interpreted by the responsible medical professional and not the only criterion for determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

IgM antibodies are present in the serum, plasma and whole blood for a short period of time, disappearing up to six months after infection, while IgG antibodies remain present for a long period of time. However, the high sensitivity of the new diagnostic methodologies, in some cases, allows to find very low levels of IgM antibodies, called residual IgM, for a longer period of time. Detection of these residual IgM antibodies hinders clinical interpretation over the period of infection.

In these cases, to confirm the result, an IgG avidity test is recommended.

The interpretation of a diagnostic test, should not be based on a single run. This should include confirmation of other tests before that a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

QUALITY CONTROL

Accuracy

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	1,893	1,658	0,076
Standard Deviation	0,096	0,027	0,008
Coefficient of Variation (%)	5,105	1,642	11,058

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	1,827	1,663	0,079
Standard Deviation	0,016	0,020	0,011
Coefficient of Variation (%)	0,893	1,236	13,947

Clinical Sensitivity and Specificity

Sample of Serum or Plasma

BIOLISA Toxoplasmosis IgM Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA Toxoplasmosis IgM kit is 98,27% and clinical specificity is > 99,9%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	BIOLISA Toxoplasmosis IgM
Positive Sample	58	57
Negative Sample	57	57
Total Tested Sample	115	

Clinical Sensitivity: 98,27% (57/58)

Clinical Specificity: > 99,9% (57/57)

Whole Blood Samples on Filter Paper

BIOLISA Toxoplasmosis IgM Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA Toxoplasmosis IgM kit is 98,11% and clinical specificity is 96,42%.

CLINICAL SENSITIVITY AND SPECIFICITY	Expected Result	BIOLISA Toxoplasmosis IgM
Positive Sample	53	52
Negative Sample	56	54
Total Tested Sample	109	

Clinical Sensitivity: 98,11% (52/53)

Clinical Specificity: 96,42% (54/56)

DIAGNOSTIC SIGNIFICANCE

Toxoplasma gondii is the causative agent of toxoplasmosis. It is an obligate intracellular protozoan that has been found in many species of birds, reptiles and mammals. The agent can be transmitted through organ transplantation, transfusion of blood and leukocytes, contact with feces of contaminated cats and ingestion of contaminated raw meats. In adults, the infection is usually benign or asymptomatic. However, symptomatic cases, including fatal cases, occur in immunosuppressed patients who have clinical or laboratory evidence of damage to the central nervous system. Following infection, IgM antibodies appear within 5 days and have reduced levels within a few weeks or months. IgG antibodies usually appear 1 - 2 weeks post infection, reaching peak levels in 6-10 weeks persisting for life. In children, the risk of fetal infection varies according to the time of pregnancy when the mother is infected. In maternal infections that occur during the first trimester, the probability of infection passing to the fetus is lower. However, if transmission occurs, serious outcomes such as miscarriage and hydrocephalus are more likely. Infections acquired later in pregnancy, where fetal transmissions occur more frequently, tend to be less severe but can still generate congenital manifestations, including cerebral calcifications and poor learning.

NUMBER OF TESTS

Presentation 1 - 96 tests

Presentation 2 - 192 tests

Presentation 3 - 480 tests

BIBLIOGRAPHIC REFERENCES

1. Feldman, HA. Toxoplasmosis: An Overview. Acad. Med. (1974) 50:110-127.
2. Remington, JS. Toxoplasmosis in the Adult. Acad. Med. (1974) 50:211-227.
3. McLeod, R, and Remington JS. In Harrison.s Principles of Internal Medicine. (1980) 879-885.
4. Frenkel, JK. Toxoplasma in and Around Us. Bioscience. (1973) 23:343-352.
5. Feldman, HA. Epidemiology of Toxoplasma Infections, Epid. Rev. (1982) 4:204-213.
6. Krick, JA, and Remington, JS. Toxoplasmosis in the Adult – An Overview. N. Engl. J. Med. (1978) 298(10):550-553.
7. Bryan, RT, and Wilson, M. Toxoplasmosis. Lab Management (1988) 26:40-43.
8. TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454

E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA Toxoplasmosis IgM kit: 10269360194

Review: May/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED