Bioclin

BIOLISA TSH DBS

REF **K278**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of Thyroid Stimulating Hormone or Thyrotropin (TSH) in biological samples of dried blood collected on filter paper (DBS), for the screening of hypothyroidism in adults and neonates, using an enzyme immunoassay test. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or enzyme immunoassay
The BIOLISA TSH DBS kit is a solid phase enzyme immunoassay, based on the “sandwich” principle for the quantitative detection of Thyroid Stimulating Hormone (TSH) in human dried blood samples collected on filter paper (DBS). To achieve this, the Sensitized Plate is coated with specific antibodies against an antigenic site present exclusively on TSH molecules. Firstly, the dried blood sample is added to the microplate, eluted by the Elution Buffer and incubated under agitation. The TSH molecules present in the sample bind to the anti-TSH antibodies adhered to the microplate which, after incubation, is washed to remove unbound materials. The Conjugate, made up of specific anti-TSH antibodies conjugated to peroxidase, is added and incubated, binding to the TSH molecules immobilized by the anti-TSH antibodies on the microplate. A new wash is performed to remove unbound materials. Then, the Substrate is added and incubated, producing a blue color that proportionally indicates the amount of TSH detected in the sample. The Stop Solution is added to stop the reaction, resulting in a change in color from blue to yellow, measured in a microplate reader.

REAGENTS

- 1 - Reference Standards (A – F)** - Store between 2 and 8 °C. Contains: One (1) filter paper impregnated with six (6) units of Reference Standards (A - F), prepared in blood, containing Thyroid Stimulating Hormone (TSH) in different concentrations. Concentrations vary with each batch (from 0 to 200 µUI/mL), see label. **Potentially infectious.**
- 2- Conjugate** - Store between 2 and 8 °C. Contains: Buffer solution containing specific anti-TSH antibodies linked to peroxidase, stabilizers, preservatives, surfactant and dye.
- 3- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with specific anti-TSH antibodies.
- 4- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer solution (phosphate < 0.5 mol/L, potassium chloride < 100 mmol/L, sodium chloride < 5 mol/L), surfactant and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer solution containing urea peroxide, tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1M Hydrochloric Acid.
- 7- Elution Buffer** - Store between 2 and 8 °C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 8- Controls (1-2)** - Store between 2 and 8 °C. Contains: One (1) filter paper impregnated with two (2) units of Controls, prepared in blood, containing Thyroid Stimulating Hormone (TSH) in different concentrations. Concentrations vary with each batch, see label. **Potentially infectious.**

PRESENTATIONS

REAGENTS	1	2	3	4
	96 cavities	192 cavities	480 cavities	960 cavities
1- Reference Standards (A-F)	1 unit	2 units	3 units	5 units
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 100 mL
3- Sensitized Plate	1 unit	2 units	5 units	10 units
4- Concentrated washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL	1 Vial x 500 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 100 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 100 mL
7- Elution Buffer	1 Vial x 20 mL	1 Vial x 40 mL	1 Vial x 100 mL	1 Vial x 200 mL
8- Controls (1-2)	1 unit	2 units	3 units	5 units

EQUIPMENT AND OPERATING SUPPLIES

- Reagents described in the previous table

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 100 to 300 µL with a coefficient of variation of less than 1.5%.
- 2- Repipettor for repetitive pipetting of volumes from 100 to 300 µL with a coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer for filter paper technique.
- 4- ELISA reader with absorbance capacity at 450 and 620 nm wavelength.
- 5- Room temperature incubator with shaking capacity at 650 rpm.
- 6- Paper shredder (diameter 3 ± 0.2 mm).
- 7- Absorbent paper to dry the microcavities.
- 8- Stopwatch or clock.
- 9- Bottle to store the Washing Solution after dilution.
- 10- Distilled or deionized water.
- 11- Quality Control Tools.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8 °C. Transport at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate must only be opened after reaching room temperature. Return the unused microcavity strips to the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free from contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw materials for the product are tested and are non-reactive for HBsAg, Anti-HIV 1&2 and Anti HCV. However, these tests do not offer complete assurance of the absence of infectious agents. The handling of any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take appropriate biosafety precautions when handling these products.

- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 9- As a protective measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
- 11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation stages.
- 12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through a request through SAC (Customer Advisory Service) from Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Dried blood collected on filter paper (DBS) – EDTA or Puncture

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use.

Blood samples dried on filter paper can be stored for up to 35 days at room temperature, as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, samples must remain refrigerated, between 2 and 8 °C. Storage for more than 2 years must be carried out at a temperature of -20 °C.²

PROCESS DESCRIPTION

Stability After Opening

The stability test results prove that the BIOLISA TSH DBS kit is stable after opening for up to 30 days. This stability may vary depending on test and environmental conditions. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of bottle No. 4 (Concentrated Washing) in a 1:20 ratio in distilled or deionized water. Therefore, dilute 50 mL of the Concentrated Washing solution in distilled or deionized water, for a final volume of 1000 mL. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, consult SAC (Customer Advisory Service).

Before starting the assay, place all reagents, samples, Reference Standards and Controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

- 1- Separate the cavities to be used considering: Reference Standards (from A to F), Controls (1 and 2) and dried blood samples on filter paper. It is recommended to test in duplicate. Return unused microplate strips to the original sealed packaging.
- 2- Separate the first cavity for the Blank (optional).
- 3- Add a 3 mm disc of each of the Reference Standards, Controls and dried blood samples, previously cut, into the determined cavities.
- 4- Pipette 100 µL of Elution Buffer into all wells, including the Blank well. Cover the cavities with plate sealer.
- 5- Incubate for 120 minutes at room temperature (approximately 25 °C) and under constant agitation at 650 rpm, which can vary between 550 and 750 rpm.
- Note:** Agitation within the indicated speed is a crucial step for the efficient elution of dried blood samples and necessary for the correct performance of the assay.

- 6- Remove the sealer from the cavities.
- 7- After incubation, discard the contents of the cavities by aspiration in a washing machine suitable for filter paper technique. Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) washing cycles, shaking for three (3) seconds at the end of each cycle. To ensure the plate dries, at the end of washing, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause inadequate results.
- 8- Pipette 100 µL of Conjugate into all wells, including the Blank well.
- 9- Gently homogenize for ± 10 seconds. Cover the cavities with plate sealer.
- 10- Incubate for 60 minutes at room temperature (approximately 25 °C).
- 11- Remove the plate sealer from the cavities.
- 12- Repeat item 7.
- 13- Pipette 100 µL of Substrate into all wells, including the Blank well.
- 14- Gently homogenize for ± 10 seconds. Cover the cavities with plate sealer.
- 15- Incubate for 30 minutes at room temperature (approximately 25 °C), protected from light.
- 16- Remove the sealer from the cavities.
- 17- Pipette 100 µL of Stop Solution into all wells, including the Blank well.
- 18- Gently homogenize for ± 10 seconds.
- 19- Read using a double filter: 450 nm / 620 nm in up to 10 minutes (maximum).

TECHNIQUE CHECK

Check whether the results obtained for reading the Blank and Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.100
Reference Standard A	< 0.100
Reference Standard B	> A and < C
Reference Standard C	> B and < D
Reference Standard D	> C and < E
Reference Standard E	> D and < F
Reference Standard F	> E

After validating the Reference Standards, the Calibration Curve and calculation of the Control concentration must be carried out, as described in the CALCULATIONS item. Check that the results obtained are within the range referenced on the Controls label.

If the absorbance values of the Reference Standards and/or the concentration of the Controls are outside the expected values, the technique must be repeated.

CALCULATIONS

A calibration curve is used to determine the TSH concentration in unknown samples.

Preparation of the Calibration Curve

If they were performed in duplicate, calculate the average absorbances obtained in the microplate reader for each of the Reference Standards (A-F), as in the example:

Standard Reference	Concentration (See Label)	Absorbance	Medium Absorbance
A	0.0	0.039	0.042
		0.045	
B	10.0	0.305	0.299
		0.293	
C	20.0	0.507	0.503
		0.498	
D	50.0	0.912	0.909
		0.905	
E	100.0	1.415	1.407
		1.399	
F	200.0	2.152	2.075
		1.997	

Plot the average absorbances of each Reference Standard (A-F) versus the corresponding concentration in µIU/mL on graph paper or using computer programs. Plot the curve point by point.

Calculation of TSH Concentration of Samples

Calculate the concentration of each sample analyzed from the Calibration Curve made on graph paper, interpolating the absorbance value obtained when reading the sample with the corresponding concentration value (µIU/mL). The calculation of the concentration of the tested samples can also be carried out using suitable computer programs, through linear regression.

The kit's measurement range is between 0.45 and 200.00 µIU/mL (analytical sensitivity and maximum point of the curve, respectively). Results below the lower limit of detection should be expressed as <0.45 µIU/mL. Results above the upper limit of detection must be expressed as >200.00 µIU/mL.

Note: The data presented is for illustration only and cannot be used to replace the calibration curve, which must be constructed in the laboratory. The Reference Standards and Controls provided in the kit must always be tested with each new assay.

REFERENCE VALUES

Adults	Normal	< 4.0 µIU/mL
	Undetermined / Limit	4.0 – 5.5 µIU/mL
	Changed	> 5.5 µIU/mL

Neonates	Normal	< 8.0 µIU/mL
	Undetermined / Limit	8.0 – 12.0 µIU/mL
	Changed	> 12.0 µIU/mL

These values should be used as a guideline, and each laboratory should create its range of reference values according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PROCESS LIMITATIONS

- 1- The BIOLISA TSH DBS product is intended for screening hypothyroidism in adults and neonates. Patients with borderline and altered values in filter paper samples must be confirmed with serum dosage for exact values.
- 2- The use of dried blood samples on filter paper that were collected, dried and stored unsatisfactorily may generate inadequate results.
- 3- It is recommended that the guidelines established for neonatal screening and diagnosis of congenital hypothyroidism in each region be followed.

4- TSH concentrations may vary between different immunoassay manufacturers, due to the diversity of available antigens and antibodies used in the manufacture of each test.

5- The interpretation of a diagnostic test should not be established based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis.

INTERFERENCES

No interference was observed by Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Acid Oxalic 60 mg/dL, C-Reactive Protein 41.2 mg/dL, anti-Streptolysin O 1023 IU/mL and Rheumatoid Factor 1080 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 40 dried blood samples on filter paper with low concentrations of TSH, but with known high concentrations of other analytes. Among them, 10 samples with Human Chorionic Gonadotropin (hCG) up to 50,000 mIU/mL, 10 samples with Follicle Stimulating Hormone (FSH) up to 1,000 mIU/mL, 10 samples with Luteinizing Hormone (LH) up to 1,000 mIU/mL and 10 samples with Prolactin up to 1,000 ng/mL. No false positive results were observed for the above analytes, up to the maximum concentration tested. Despite the results found, the possibility of cross-reactivities with other analytes or with higher concentrations cannot be completely ruled out. The final diagnosis must consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	2.610	0.471	0.083
Standard Deviation	0.229	0.049	0.010
Coefficient of Variation (%)	8.77	10.42	12.49

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	2.603	0.471	0.094
Standard Deviation	0.244	0.046	0.014
Coefficient of Variation (%)	9.36	9.82	14.85

ANALYTICAL SENSITIVITY

The analytical sensitivity of the BIOLISA TSH DBS kit is 0.45 µIU/mL.

LINEARITY

The reaction is linear up to the concentration at the highest point of the calibration curve.

HIGH DOSE PRO-ZONE EFFECT

No high-dose prozone effect was observed at up to 1,000 µIU/mL TSH.

DIAGNOSTIC MEANING

Thyroid Stimulating Hormone (TSH) is a glycoprotein hormone secreted by the anterior pituitary gland. TSH is made up of two subunits, α and β. The α subunit of TSH is similar to the α subunit found in other glycoprotein hormones such as LH, FSH and hCG. However, the β subunit is specific and differs from hormone to hormone. Thyroid hormones (T4 and T3) are secreted and produced by the thyroid gland, under the regulation of TSH. Furthermore, TSH controls iodide transport and stimulates the secretion of thyroglobulin, the main protein component of the thyroid gland, into the bloodstream. Therefore, TSH plays an important role in the proper functioning and development of the thyroid gland. Thyroid hormone concentrations control TSH secretion through negative feedback. For example, in primary hypothyroidism, the serum T4 and T3 level are low while the TSH level is high. Likewise, in primary hyperthyroidism, the serum level of T4 and T3 are high while the level of TSH is low. When the cause is secondary, that is, it is a pituitary dysfunction, T4, T3 and TSH are equally high or low. Therefore, for diagnosis and monitoring of thyroid dysfunctions, it is recommended to test both glycoprotein hormone and T4/T3 hormones. Furthermore, TSH determinations are also useful for monitoring patients receiving thyroxine replacement therapy. Hypothyroidism is a disorder characterized by decreased production of the thyroid hormones T3 and T4. In adults, the main cause is Hashimoto's Disease, a consequence of an autoimmune response against the thyroid. Depression, brachycardia, constipation, irregular menstruation, memory failures, excessive tiredness, muscle pain, dry skin, hair loss, weight gain and increased cholesterol are among the symptoms of hypothyroidism. Congenital hypothyroidism (CH) is a condition where the newborn produces low levels of thyroid hormones due to abnormal development of the thyroid gland. It is estimated that 1 in every 2,000 births have CH. Although there are no symptoms at birth in the vast majority of cases, congenital hypothyroidism can lead to serious problems with the physical and mental development of affected children if not properly diagnosed and treated with hormone replacement. Therefore, investigating CH using the heel prick test is widely recommended to prevent possible consequences of the disease. The quantification of TSH in dried blood samples on filter paper is effective for screening individuals with pathologies where high TSH levels are found, in addition to monitoring the treatment of these conditions.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA TSH DBS kit: 10269360449

Review: March/2024

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER