BIOSTRIP

REF
**K293**

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FINALIDADE

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**REFERÊNCIAS BIBLIOGRÁFICAS**  
**1-** BERG, B. et. al. Guidelines for evaluation of reagent strips. Exemplified by analysis of urine albumin and glucose concentration using visually read reagent strips. Scand J Clin Lab Invest., v.49, n.7, p.689-699, 1989.  
**2-** FREE, A.H.; FREE, H.M. Urinalysis: its proper role in the physician's office. Clin Lab Med., v.6, n.2, p.253-266, 1986.  
**3-** GYURE, W.L. Comparison of several methods for semiquantitative determination of urinary protein. Clin Chem., v.23, n.5, p.876-879, 1977.  
**4-** KOSS, S. et. al. Proteinuria and renal disease: prognostic value of urine dipstick testing for leukocytes. Pediatr Nephrol., v.21, n.4, p.584-587, 2006.  
**5-** WILSON, L.A. Urinalysis. Nurs Stand., v.19, n.35, p.51-54, 2005.  
**6-** STRASINGER,K.S.; LORENZO. M.S. Urinálise e Fluidos Corporais. Livraria Médica Paulista., v.9, n.5, p.57-82, 2009.  
**7-** QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Method for the semi-quantitative determination of Blood (Hb), Bilirubin (BIL), Urobilinogen (URO), Ketone Bodies (CET), Protein (PRO), Nitrite (NIT), Glucose (GLI), pH, Density (DEN) and Leukocytes (LEU) in urine samples, with a high degree of accuracy. Not automated or semi-automated. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

**Methodology:** Colorimetric.

The semi-quantitative determination of each parameter is made according to the following reactions:

**Urobilinogen:** Based on the diazotization reaction of 4-Methoxybenzene diazonium salt with Urobilinogen, in the presence of strong acid, the color scale of the parameter varies from light pink to dark pink, according to the concentration in the urine sample.

**Glucose:** This is a double sequential enzymatic reaction. In the first stage, glucose oxidase catalyzes the reaction between the glucose present in the urine sample and oxygen from the ambient air, producing gluconic acid and hydrogen peroxide. In the second stage, peroxidase catalyzes the reaction between the hydrogen peroxide, formed previously, and the chromogen potassium iodide, giving rise to a compound colored blue for negative samples and varies from green to brown for positive samples.

**Ketone Bodies:** Acetoacetic acid, if present in the urine sample, reacts with sodium nitroprusside, in an alkaline medium, giving rise to a violet color proportional to the concentration of the acid, in positive samples. In its absence, the parameter takes on a beige color, which characterizes negative samples.

**Bilirubin:** Bilirubin, if present in the urine sample, reacts in an acidic environment with an organic diazonium salt, forming a color that varies from yellowish brown in negative samples to violet, which indicates the highest concentration in positive samples.

**Protein:** Protein, if present in the urine sample, promotes an ion exchange with the tetrabromophenol blue indicator, without causing major changes in the pH of the medium. The color scale varies from light yellow, in negative samples, to intense green, in positive samples.

**Nitrite:** Nitrite is a product of the reduction of nitrate, a normal constituent of urine, carried out by certain bacteria. In the samples, nitrite reacts with an organic acid, forming a diazonium salt, which in turn reacts with a chelator, producing a pink compound. Any intensity of pink color is considered positive.

**pH:** Uses a dual indicator system. The indicator methyl red and bromothymol blue are used to give a wide range of colors, covering pH levels between 5.0 and 9.0. The color scale varies from orange to greenish yellow at pH between 5.0 and 6.5, due to the color change of the methyl red indicator, and from green to blue at pH between 6.5 and 9.0, due to the change color caused by the bromothymol blue indicator.

**Blood/Hemoglobin:** The reaction is based on the activity of hemoglobin pseudo-peroxidase which, if present in the urine sample, catalyzes the reaction between the chromogen and organic peroxide. The color scale varies from yellow, for negative samples, to bluish green, for positive samples of higher concentration.

**Density:** This test is based on the change in pKa of the polyelectrolyte, which ionizes releasing hydrogen ions. Low concentrations of these ions indicate low density samples, and the color of the parameter varies from blue to greenish tones. On the other hand, high concentrations of hydrogen ions indicate high density samples and the color of the parameter takes on yellowish tones.

**Leukocytes:** The reaction detects the presence of leukocyte esterases, enzymes found in leukocytes. Leukocytes, if present in the urine sample, are lysed and release esterases that break down an ester to an aromatic and acidic compound. The product of this decomposition reacts with diazonic salt, which in an acidic medium produces a violet color proportional to the amount of leukocytes in the urine. The color scale varies from white, when negative, to violet, when positive.

REAGENTS

**Number 1 - Urine reagent strip** - Store between 15 and 30°C.

PRESENTATION

| Presentation | Reagent N° 1 |
|--------------|--------------|
| 1            | 25           |
| 2            | 50           |
| 3            | 100          |
| 4            | 150          |
| 5            | 200          |

EQUIPMENT AND OPERATIONAL INPUTS

Urine collection bottle, absorbent paper, and timer. They are found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 15 to 30°C. Temperature variations during 3 days during the transport of this product do not affect the quality. After this period, the product must be packaged according to the recommended temperature. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- Carefully read the instructions for use before carrying out the test.
- The expiration date corresponds to the last day of the month marked on the kit box label. Do not use after the expiration date.
- You should avoid exposing the kit to high temperatures, as well as directly in the sun.
- Do not freeze the kit components, as this will cause irreversible deterioration.
- Do not use the kit when any component presents a visual characteristic that does not comply with what is specified in the product's SDS.
- Maintain the usual safety precautions when handling the components. All samples must be handled as potentially infectious materials.
- We recommend applying local, state and federal environmental protection standards so that materials are disposed of in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request through Quibasa's SAC (Customer Advisory Service).

**10-** Use Good Laboratory Practices (GLP) in the conservation, handling and disposal of materials.

SAMPLES

Preferably use freshly collected urine (up to 4 hours old) in a clean, dry bottle, free from soap residue and acids. If the test is carried out 4 hours after collection, store the urine in a refrigerator or refrigerated place with a temperature between 2 and 8°C. To carry out the tests, remove the samples from the refrigerated place and let them reach room temperature. Homogenize the samples manually. Do not centrifuge them.

STABILITY AND HANDLING

Store test strips between 15 and 30°C. Do not reuse them. Discard them after the expiration date. Do not expose them to sunlight and do not remove the desiccant from the packaging. During use, remove the required number of test strips and immediately close the bottle. Darkening or color change of parameters indicates deterioration. If this is evident or if the test results are questionable compared to what was expected, ensure that the test strip stability conditions are adequate.

**Note:** The kit maintains its performance after the first use and is stable, as long as the test strip stability conditions mentioned above are followed.

TECHNIQUE

This procedure must be followed correctly to obtain reliable results.

- Confirm that the product is within the expiration date printed on the label.
- Remove the test strip from the tube and close it immediately.
- Inspect the test strip. Discoloration and darkening in reagent areas may indicate spoilage. In this case, do not use the product.
- Immerse the test strip completely in the well-homogenized fresh urine sample for about 1 second. Make sure all parameters are moistened.
- Remove the test strip from the sample and remove excess urine by touching the test strip to the side of an absorbent paper.
- Compare the results carefully with the color chart on the tube label in a well-lit environment. For the leukocyte parameter, the reading must be taken between 60 and 120 seconds. For other parameters, the reading time must be between 30 and 60 seconds.

**Note:** When reading, keep the strip in a horizontal position to avoid chemical interactions due to possible excess urine. Changes in color along the edges of test areas or after the recommended reading time has elapsed are not of diagnostic significance.

INTERPRETATION OF RESULTS

Results are obtained by directly comparing the test strip to the color chart printed on the tube label. Calculations and laboratory equipment are not necessary, but can be used if semi-automation is chosen.

1- Urobilinogen

**Expected Values:** The normal range of urobilinogen is 0.1 to 1.0 mg/dL. If results exceed a concentration of 2.0 mg/dL, the patient and urine sample should be further evaluated.

**Detection Limits:** The test will detect urobilinogen at such low concentrations, close to 0.1 mg/dL, that in most normal samples the reactive area may appear light pink. In patients with high urobilinogen excretion, the results must be correlated with the Watson-Schwartz procedure using spectrophotometry.

**Test Limitation:** The test area may, in certain concentrations, react with known interfering substances, such as p-aminosalicylic acid. The test is not a reliable method for detecting porphobilinogen.

2- Glucose

**Expected Values:** The normal kidney excretes a small concentration of glucose in the urine. Concentrations of 100 mg/dL may be considered abnormal if found consistently.

**Detection Limits:** Approximately 50 mg/dL of glucose is already detectable. The test is highly specific for glucose. This way, the reagent area will not present a positive result for any other excreted substance, such as lactose, galactose, fructose or salicylate-reducing metabolites and nalidixic acid. False-negative results can be obtained with the presence of levodopa, ascorbic acid, glutathione and dipyrone.

**Test limitation:** A high density (>1.020) with a high pH and the presence of ascorbic acid (above 50 mg/dL) can cause false-negative results in samples that have a low glucose concentration. Moderately high levels of ketone bodies (>40 mg/dL) may present a false negative for a sample containing a low concentration of glucose (100 mg/dL). Test reactivity can be influenced by density and temperature.

3- Ketone Bodies

**Expected Values:** Ketone Bodies should not be detected with this reagent in normal urine samples.

**Detection Limits:** The test has a sensitivity of 0.5 mg/dL of ketone bodies (acetoacetic acid).

**Test Limits:** Positive results (trace or low) may occur with highly pigmented urine samples or those urines that contain large amounts of levodopa metabolites. High density, low pH and phenolphthalein sulfate can lead to false-positive results.

4- Bilirubin

**Expected Values:** Normally, bilirubin is not detectable in urine even with the most sensitive methods. Even low bilirubin concentrations are abnormal, so further investigations should be required.

**Limits of detection:** The test has a sensitivity of 0.5 mg/dL bilirubin.

**Test Limitation:** Urine samples should not be exposed to light for long periods. Drug metabolites, such as Phenazopyridine Hydrochloride and Serenium, which give color at low pH, can cause false-positive results. Indoxyl sulfate may produce a response that varies from orange-yellow to red, and which may interfere with the interpretation of positive or negative bilirubin readings. Ascorbic acid (> 25 mg/dL) may cause a false-negative result.

5- Protein

**Expected Values:** Normal urine samples contain low concentrations of protein (<20 mg/dL), therefore, persistent high levels of protein in the urine may indicate kidney disease or urinary tract disease. Persistent trace or elevated level results indicate significant proteinuria and thus further clinical testing is needed to assess the significance of the results.

**Detection Limits:** This test has a detection limit of 10~15 mg/dL of protein, being more sensitive to albumin

**Test Limitation:** False-positive results can be found in urine with a basic pH (pH 9). Result interpretation is impaired in urine samples that present turbidity.

**6- Nitrite**  
**Expected Values:** Normally nitrite is not detectable in urine.  
**Detection Limits:** Comparing the reagent area against a white background can assist in detecting low levels. The test is specific for nitrite and will not react with any other substance normally excreted in the urine.  
**Test Limitation:** Ascorbic acid (> 25 mg/dL) may lead to a false-negative result in samples containing low nitrite concentration in urine (<0.03 mg/dL). A negative result does not mean that the patient is free from bacteriuria. Pink dots or pink areas should be interpreted as positive results. Negative results may occur when urinary tract infections are not caused by organisms that do not contain nitrate reductase, when urine has not been held in the bladder long enough (four hours or more) for the reduction of nitrate to nitrite to occur, or when nitrate from diet is absent.

**7- pH**  
**Expected Values:** Urine pH values generally range from 5.0 to 9.0.  
**Detection Limits:** The test measures pH values generally within 1 unit in the range of 5.0 to 9.0.  
**Test Limitation:** Excess urine on the test strip can carry acid from the buffer in the protein reagent area to the pH parameter, causing reading errors. This phenomenon is known as "run-over".

**8- Blood/Hemoglobin**  
**Expected Values:** Normally, hemoglobin is not detected in urine (10 RBC/ $\mu$ L; 0.03 mg/dL). When hemoglobin is detectable in urine it may indicate kidney disease or a urinary tract disorder. Blood can often be found in the urine of menstruating women.  
**Limits of Detection:** The test is generally capable of detecting 10 Cells/ $\mu$ L.  
**Limitation of the Test:** The test is slightly more sensitive to free hemoglobin and myoglobin than intact erythrocytes. High density or high protein can reduce the reactivity of the blood test. Microbial peroxidase associated with urinary tract infection may present false-positive results. Ascorbic acid concentrations (>40 mg/dL) may cause false-negative results at low urine levels.

**9- Density**  
**Expected Values:** Normal density is between 1.001 and 1.035.  
**Detection Limits:** The test allows the determination of density in urine between 1.000, 1.005, 1.010, 1.015, 1.020, 1.025, 1.030.  
**Test Limitation:** Highly alkaline samples can lead to a decrease in the result, while highly acidic samples can significantly increase the result.

**10- Leukocyte**  
**Expected Values:** Normally leukocytes are not detectable in urine.  
**Detection Limits:** The test is generally capable of detecting 20~25 WBC/ $\mu$ L with a trace.  
**Limitation of the Test:** The test result may not always be consistent with the number of cells countable by microscopic examination. High glucose concentration, high density, high albumin level, high formaldehyde concentration, or presence of blood may cause decreased test results. High concentration of oxalic acid from oxidizing agents may cause false-positive results.

**QUALITY CONTROL**  
Test strips must be correctly stored and handled before and during testing. The reaction of the reagent strips must be confirmed by testing known negative and positive samples or by multiple analytical controls that contain normal and abnormal amounts of each analytical being tested.

**SENSITIVITY AND SPECIFICITY**  
**Urobilinogen:** 0.1 mg/dL.  
**Glucose:** 50 - 100 mg/dL.  
**Ketone Bodies:** 5 mg/dL (acetoacetic acid).  
**Bilirubin:** 0.5 mg/dL (bilirubin).  
**Protein:** 30 mg/dL (albumin).  
**Nitrite:** 0.05 mg/dL (nitrite ion).  
**pH:** 5.0 to 9.0.  
**Blood:** 10 RBC/ $\mu$ L (0.03 mg/dL hemoglobin and intact erythrocyte).  
**Density:** 1,000 to 1,030.  
**Leukocyte:** 20-25 WBC/ $\mu$ L (intact and lysed).

**LIMITATIONS OF USE**  
The urine strip is a test to semiquantitatively determine the presence of various parameters in urine. Therefore, as with all laboratory tests, definitive diagnosis or therapeutic decisions should not be based on a single result or method. The definitive clinical diagnosis must be made by the doctor after analyzing clinical and laboratory data.

The effects of drugs or other metabolites on individual test strip tests are not known in all cases. Therefore, it is recommended that in case of doubt, the test is repeated after removing the potential interfering agent, such as medication or vitamin supplement, etc.

The correct reading time shown on the tube label is important for the accuracy of the results, any reading outside this will invalidate the test.

Changes in color that appear along the edge of the test area should be ignored; careful removal of excess urine should eliminate this phenomenon.

| Parameter    | Interferers                                           | Interferences          |
|--------------|-------------------------------------------------------|------------------------|
| Urobilinogen | Formalin                                              | False-negative results |
| Glucose      | Increased density associated with elevated pH         | False-negative results |
|              | Ascorbic acid in concentrations greater than 50 mg/dL |                        |
|              | Ketone bodies in concentrations greater than 40 mg/dL |                        |
| Ketones      | Increased density associated with decreased pH        | False-positive results |
| Bilirubin    | Ascorbic acid in concentrations greater than 50 mg/dL | False-negative results |
| Protein      | High pH                                               | False-positive results |
| Nitrite      | Ascorbic acid in concentrations greater than 25 mg/dL | False-negative results |

|            |                                                       |                        |
|------------|-------------------------------------------------------|------------------------|
| pH         | No interferers were identified                        | -                      |
| Blood      | Increased density                                     | False-negative results |
|            | Ascorbic acid in concentrations greater than 40 mg/dL |                        |
|            | Increased protein concentration                       |                        |
| Leukocytes | Increased hemoglobin concentrations                   | False-negative results |
|            | Increased glucose concentrations                      |                        |
|            | Increased density                                     |                        |
|            | Increased albumin concentrations                      |                        |

**RISK MANAGEMENT**  
**Bioclin**, after a detailed review and critical analysis of all known and/or anticipated hazards, concludes that all risks associated with the Biostrip Product have been assessed, that risk reduction measures have been implemented and that the product does not present any major risks. that the benefits obtained from its use; and that, if used by qualified and trained professionals, aware of the precautions described in the products, it will perform its functions with quality, safety and effectiveness.

**CLINICAL IMPORTANCE**  
Urine testing provides the clinician with information on possible pathologies of the urinary tract, as well as about some extra-reputable diseases. Due to its simplicity, low cost and the ease of obtaining the sample for analysis, it is a routine examination.

**BIBLIOGRAPHIC REFERENCES**  
**1-** BERG, B. et. al. Guidelines for evaluation of reagent strips. Exemplified by analysis of urine albumin and glucose concentration using visually read reagent strips. Scand J Clin Lab Invest., v.49, n.7, p.689-699, 1989.  
**2-** FREE, A.H.; FREE, H.M. Urinalysis: its proper role in the physician's office. Clin Lab Med., v.6, n.2, p.253-266, 1986.  
**3-** GYURE, W.L. Comparison of several methods for semiquantitative determination of urinary protein. Clin Chem., v.23, n.5, p.876-879, 1977.  
**4-** KOSS, S. et. al. Proteinuria and renal disease: prognostic value of urine dipstick testing for leukocytes. Pediatr Nephrol., v.21, n.4, p.584-587, 2006.  
**5-** WILSON, L.A. Urinalysis. Nurs Stand., v.19, n.35, p.51-54, 2005.  
**6-** STRASINGER,K.S.; LORENZO. M.S. Urinálise e Fluidos Corporais. Livraria Médica Paulista., v.9, n.5, p.57-82, 2009.  
**7-** QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**  
Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

■ **QUIBASA QUÍMICA BÁSICA Ltda**  
Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**CUSTOMER SERVICE**  
Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

**Registration number of the BIOSTRIP kit at ANVISA:** 10269360444  
**Review:** February/2024

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 LOT NUMBER

 MANUFACTURING DATE

 VALIDITY DATE (last day of the month)

 TEMPERATURE LIMIT (store)

 CONTENT IS SUFFICIENT FOR <N> TEST

 SEE INSTRUCTIONS FOR USE

 IN VITRO DIAGNOSTIC PRODUCT

 KEEP AWAY FROM SUNLIGHT

 DO NOT REUSE

 CAUTION

 MADE BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 FLAMMABLE

 CORROSIVE

 TOXIC

 DO NOT USE IF PACKAGE IS DAMAGED

 PRODUCT STERILIZED

 DANGER