BIOSTRIP VET

REF VET084

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FINALIDADE

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REFERÊNCIAS BIBLIOGRÁFICAS
1- BERG, B. et. al. Guidelines for evaluation of reagent strips. Exemplified by analysis of urine albumin and glucose concentration using visually read reagent strips. Scand J Clin Lab Invest., v.49, n.7, p.689-699, 1989.
2- FREE, A.H.; FREE, H.M. Urinalysis: its proper role in the physician's office. Clin Lab Med., v.6, n.2, p.253-266, 1986.
3- GYURE, W.L. Comparison of several methods for semiquantitative determination of urinary protein. Clin Chem., v.23, n.5, p.876-879, 1977.
4- KOSS, S. et. al. Proteinuria and renal disease: prognostic value of urine dipstick testing for leukocytes. Pediatr Nephrol., v.21, n.4, p.584-587, 2006.
5- WILSON, L.A. Urinalysis. Nurs Stand., v.19, n.35, p.51-54, 2005.
6- STRASINGER,K.S.; LORENZO. M.S. Urinálise e Fluidos Corporais. Livraria Médica Paulista., v.9, n.5, p.57-82, 2009.
7- eCLINPATH- Cornell University College of Veterinary Medicine, 2024. Disponível em: https://edlinpath.com/.
8- Orcurr, C.J. Ferret urogenital diseases. Veterinary Clinics: Exotic Animal Practice, 6(1), 113-138.

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INSTRUCTIONS FOR USE

FUNCTION

Method for the semi-quantitative determination of Blood (Hb), Bilirubin (BIL), Urobilinogen (URO), Ketone Bodies (KET), Protein (PRO), Nitrite (NIT), Glucose (GLU), pH, Density (DEN) and Leukocytes (LEU) in urine samples from animals. Not automated or semi-automated. Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric.

The semi-quantitative determination of each parameter is made according to the following reactions:

Blood/Hemoglobin: The reaction is based on the activity of hemoglobin pseudo-peroxidase which, if present in the urine sample, catalyzes the reaction between the chromogen and organic peroxide. The color scale varies from yellow, for negative samples, to bluish green, for positive samples of higher concentration.

Bilirubin: Bilirubin, if present in the urine sample, reacts in an acidic environment with an organic diazonium salt, forming a color that varies from yellowish brown in negative samples to violet, which indicates the highest concentration in positive samples.

Urobilinogen: Based on the diazotization reaction of 4-Methoxybenzene diazonium salt with Urobilinogen, in the presence of strong acid, the color scale of the parameter varies from light orange to dark orange, according to the concentration in the urine sample.

Ketone Bodies: Acetoacetic acid, if present in the urine sample, reacts with sodium nitroprusside, in an alkaline medium, giving rise to a violet color proportional to the concentration of the acid, in positive samples. In its absence, the parameter takes on a beige color, which characterizes negative samples.

Protein: Protein, if present in the urine sample, promotes an ion exchange with the tetrabromophenol blue indicator, without causing major changes in the pH of the medium. The color scale varies from light yellow, in negative samples, to intense green, in positive samples.

Nitrite: Nitrite is a product of the reduction of nitrate, a normal constituent of urine, carried out by certain bacteria. In the samples, nitrite reacts with an organic acid, forming a diazonium salt, which in turn reacts with a chelator, producing a pink compound. Any intensity of pink color is considered positive.

Glucose: This is a double sequential enzymatic reaction. In the first stage, glucose oxidase catalyzes the reaction between the glucose present in the urine sample and oxygen from the ambient air, producing gluconic acid and hydrogen peroxide. In the second stage, peroxidase catalyzes the reaction between the hydrogen peroxide, formed previously, and the chromogen potassium iodide, giving rise to a compound colored blue for negative samples and varies from green to brown for positive samples.

pH: Uses a dual indicator system. The indicator methyl red and bromothymol blue are used to give a wide range of colors, covering pH levels between 5.0 and 9.0. The color scale varies from orange to greenish yellow at pH between 5.0 and 6.5, due to the color change of the methyl red indicator, and from green to blue at pH between 6.5 and 9.0, due to the change color caused by the bromothymol blue indicator.

Density: This test is based on the change in pKa of the polyelectrolyte, which ionizes releasing hydrogen ions. Low concentrations of these ions indicate low density samples, and the color of the parameter varies from blue to greenish tones. On the other hand, high concentrations of hydrogen ions indicate high density samples and the color of the parameter takes on yellowish tones.

Leukocytes: The reaction detects the presence of leukocyte esterases, enzymes found in leukocytes. Leukocytes, if present in the urine sample, are lysed and release esterases that break down an ester to an aromatic and acidic compound. The product of this decomposition reacts with diazonic salt, which in an acidic medium produces a violet color proportional to the amount of leukocytes in the urine. The color scale varies from white, when negative, to violet, when positive.

REAGENTS

Number 1 - Urine reagent strip - Store between 15 and 30°C.

PRESENTATION

Presentation	Reagent N° 1
1	25
2	50
3	100
4	150
5	200

EQUIPMENT AND OPERATIONAL INPUTS

Urine collection bottle, absorbent paper, and timer. They are found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 15 to 30°C. Temperature variations during 3 days during the transport of this product do not affect the quality. After this period, the product must be packaged according to the recommended temperature. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- Carefully read the instructions for use before carrying out the test.
- The expiration date corresponds to the last day of the month marked on the kit box label. Do not use after the expiration date.
- You should avoid exposing the kit to high temperatures, as well as directly in the sun.
- Do not freeze the kit components, as this will cause irreversible deterioration.
- Do not use the kit when any component presents a visual characteristic that does not comply with what is specified in the product's SDS.
- Maintain the usual safety precautions when handling the components. All samples must be handled as potentially infectious materials.
- We recommend applying local, state and federal environmental protection standards so that materials are disposed of in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request through Quibasa's SAC (Customer Advisory Service).
- Use Good Laboratory Practices (GLP) in the conservation, handling and disposal of materials.

SAMPLES

Preferably use freshly collected urine (up to 4 hours old) in a clean, dry bottle, free from soap residue and acids. If the test is carried out 4 hours after collection, store the urine in a refrigerator or refrigerated place with a temperature between 2 and 8°C. To carry out the tests, remove the samples from the refrigerated place and let them reach room temperature. Homogenize the samples manually. Do not centrifuge them.

STABILITY AND HANDLING

Store test strips between 15 and 30°C. Do not reuse them. Discard them after the expiration date. Do not expose them to sunlight and do not remove the desiccant from the packaging. During use, remove the required number of test strips and immediately close the bottle. Darkening or color change of parameters indicates deterioration. If this is evident or if the test results are questionable compared to what was expected, ensure that the test strip stability conditions are adequate.

Note: The kit maintains its performance after the first use and is stable, as long as the test strip stability conditions mentioned above are followed.

TECHNIQUE

This procedure must be followed correctly to obtain reliable results.

- Confirm that the product is within the expiration date printed on the label.
- Remove the test strip from the tube and close it immediately.
- Inspect the test strip. Discoloration and darkening in reagent areas may indicate spoilage. In this case, do not use the product.
- Immerse the test strip completely in the well-homogenized fresh urine sample for about 1 second. Make sure all parameters are moistened.
- Remove the test strip from the sample and remove excess urine by touching the test strip to the side of an absorbent paper.
- Compare the results carefully with the color chart on the tube label in a well-lit environment. For the leukocyte parameter, the reading must be taken between 60 and 120 seconds. For other parameters, the reading time must be between 30 and 60 seconds.

Note: When reading, keep the strip in a horizontal position to avoid chemical interactions due to possible excess urine. Changes in color along the edges of test areas or after the recommended reading time has elapsed are not of diagnostic significance.

INTERPRETATION OF RESULTS

Results are obtained by directly comparing the test strip to the color chart printed on the tube label. Laboratory calculations and equipment are not required.

1- Blood/Hemoglobin

Expected Values: Normally, hemoglobin is not detected in urine (10 RBC/ μ L; 0.03 mg/dL). When hemoglobin is detectable in urine it may indicate kidney disease or a urinary tract disorder. **Limits of Detection:** The test is generally capable of detecting 10 Cells/ μ L.

Limitation of the Test: The test is slightly more sensitive to free hemoglobin and myoglobin than intact erythrocytes. High density or high protein can reduce the reactivity of the blood test. Microbial peroxidase associated with urinary tract infection may present false-positive results. Ascorbic acid concentrations (>40 mg/dL) may cause false-negative results at low urine levels.

2- Bilirubin

Expected Values: Dogs have a low renal threshold for bilirubin, which is why weak positive reactions (+) can be produced in healthy animals or in more concentrated samples. The appearance of traces of bilirubin (+) can also be observed in the orin of healthy humans. In other animal species, the presence of bilirubin in the urine requires investigating the appearance of liver disease and cholestasis.

Limits of detection: The test has a sensitivity of 0.5 mg/dL bilirubin.

Test Limitation: Urine samples should not be exposed to light for long periods. Drug metabolites, that promotes changes in the color of the urine at low pH, can cause false-positive results. Ascorbic acid (> 25 mg/dL) may cause a false-negative result.

3- Urobilinogen

Expected Values: The normal range of urobilinogen is 0.1 to 1.0 mg/dL. Elevated urobilinogen may occur in cases of liver disease or hemolytic anemia.

Detection Limits: The test will detect urobilinogen at such low concentrations, close to 0.1 mg/dL, that in most normal samples the reactive area may appear light pink.

Test Limitation: False-negative results may occur in collection bottles with formalin residue or in samples collected a long time ago (due to instability of urobilinogen to light or air).

4- Ketone Bodies

Expected Values: In healthy animals, ketone bodies are not detected in urine. Ketonuria is an indicator of an animal's energy imbalance. In different species, ketonuria can indicate different pathological conditions, such as: Diabetes mellitus in dogs and cats, pregnancy toxemia in sheep and camelids, ketosis in cattle, hunger and malnutrition in young animals (ECLINPATH, 2014).

Detection Limits: The test has a sensitivity of 0.5 mg/dL of ketone bodies (acetoacetic acid).

Test Limits: Positive results (trace or low) may occur with highly pigmented urine samples. High density, low pH and phenolphthalein sulfate can lead to false-positive results.

5- Protein

Expected Values: Normal urine samples contain low concentrations of protein (<20 mg/dL), therefore, persistent high levels of protein in the urine may indicate kidney disease or urinary tract disease. Positive reactions can be found in dogs with high density (>1.030), therefore it is essential to simultaneously evaluate these parameters for the correct interpretation of this analyte in urine.

Detection Limits: This test has a detection limit of 10~15 mg/dL of protein, being more sensitive to albumin

Test Limitation: False-positive results can be found in urine with a basic pH (pH 9). Result interpretation is impaired in urine samples that present turbidity.

6- Nitrite

Expected Values: Normally nitrite is not detectable in urine. The detection of nitrite in urine is an indicator of the presence of nitrifying bacteria in the urinary tract.

Detection Limits: Comparing the reagent area against a white background can assist in detecting low levels. The test is specific for nitrite and will not react with any other substance normally excreted in the urine.

Test Limitation: Ascorbic acid (> 25 mg/dL) may lead to a false-negative result in samples containing low nitrite concentration in urine (<0.03 mg/dL). A negative result does not mean that the animal is free from bacteriuria. Pink dots or pink areas should be interpreted as positive results. Negative results may occur when urinary tract infections are not caused by organisms that do not contain nitrate reductase, when urine has not been held in the bladder long enough (four hours or more) for the reduction of nitrate to nitrite to occur, or when nitrate from diet is absent.

7- Glucose

Expected Values: The elimination of glucose in urine is an abnormal finding in animals and requires investigation for the occurrence of kidney damage when present.

Detection Limits: Approximately 50 mg/dL of glucose is already detectable. The test is highly specific for glucose. This way, the reagent area will not present a positive result for any other excreted substance.

Test limitation: False-negative reactions can be observed in conditions such as: increased density (>1.020) and pH of urine and the presence of ascorbic acid above 50 mg/dL when the glucose concentration is low. Other conditions such as ketone body levels > 40 mg/dL and the use of medications such as salicylates and tetracyclines can also lead to false negative results. Test reactivity can be influenced by density and temperature.

8- pH

Expected Values: Animal diet is a determining factor for urine pH. Carnivorous species have a pH between 5.0 and 6.0, while in herbivorous species the urinary pH is alkaline, between 7.7 and 9 (except young lactating animals that have slightly acidic urine).

Detection Limits: The test measures pH values generally within 1 unit in the range of 5.0 to 9.0.

Test Limitation: The time between collection and analysis can promote variations in urinary pH, leading to alkalization. The presence of bacterial contaminants such as *Streptococcus*, *Ureaplasma* or *Proteus spp.* can promote acidification of urine.

9- Density

Expected Values: Density represents the relationship between weight and volume of urine. In healthy dogs, urinary density can vary between 1.015 and 1.045. In cats, density can vary between 1.035 and 1.060. For large animals (cattle and horses), urinary density in healthy individuals is between 1.015 and 1.030. It is recommended that density assessment be associated with conditions such as: degree of hydration and fluid intake, diet, weight, exercise, age, metabolism and climatic conditions in which the animal is located.

Detection Limits: The test allows the determination of density in urine between 1.000, 1.005, 1.010, 1.015, 1.020, 1.025, 1.030.

Test Limitation: Highly alkaline samples can lead to a decrease in the result, while highly acidic samples can significantly increase the result.

10- Leukocyte

Expected Values: Normally leukocytes are not detectable in urine.

Detection Limits: The test is generally capable of detecting 20–25 WBC/µL with a trace.

Limitation of the Test: The test result may not always be consistent with the number of cells countable by microscopic examination. High glucose concentration, high density, high albumin level, high formaldehyde concentration, or presence of blood may cause decreased test results. High concentration of oxalic acid from oxidizing agents may cause false-positive results.

QUALITY CONTROL

Test strips must be correctly stored and handled before and during testing. The reaction of the reagent strips must be confirmed by testing known negative and positive samples or by multiple analytical controls that contain normal and abnormal amounts of each analytical being tested.

SENSITIVITY AND SPECIFICITY

Blood: 10 RBC/µL (0.03 mg/dL Hb and intact erythrocyte).

Bilirubin: 0.5 mg/dL (bilirubin).

Urobilinogen: 0.1 mg/dL.

Ketone Bodies: 5 mg/dL (acetoacetic acid).

Protein: 30 mg/dL (albumin).

Nitrite: 0.05 mg/dL (nitrite ion).

Glucose: 50 - 100 mg/dL.

pH: 5.0 to 9.0.

Density: 1,000 to 1,030.

Leukocyte: 20-25 WBC/µL (intact and lysed).

LIMITATIONS OF USE

The urine strip is a test to semiquantitatively determine the presence of various parameters in urine. Therefore, as with all laboratory tests, definitive diagnosis or therapeutic decisions should not be based on a single result or method. The definitive clinical diagnosis must be made by the veterinary after analyzing clinical and laboratory data.

The effects of drugs or other metabolites on individual test strip tests are not known in all cases. Therefore, it is recommended that in case of doubt, the test is repeated after removing the potential interfering agent, such as medication or vitamin supplement, etc.

The correct reading time shown on the tube label is important for the accuracy of the results, any reading outside this will invalidate the test. Changes in color that appear along the edge of the test area should be ignored; careful removal of excess urine should eliminate this phenomenon.

INTERFERENCES

Parameter	Interferers	Interferences
Blood	Increased density	False-negative results
	Ascorbic acid in concentrations > 40 mg/dL	
	Increased protein concentration	
Bilirubin	Ascorbic acid in concentrations > 50 mg/dL	False-negative results
	Exposure to direct light	
Urobilinogen	Formalin	False-negative results
	Collection time	
Ketones	Increased density associated with decreased pH	False-positive results
	Samples containing pigment	
	Phenophthalein Sulfate	
Protein	High pH	False-positive results
	High density	

Nitrite	Ascorbic acid in concentrations > 25 mg/dL	False-negative results
	Contamination of the sample with content of fecal origin	False-positive results
Glucose	Increased density associated with elevated pH	False-negative results
	Ascorbic acid in concentrations > 50 mg/dL	
	Ketone bodies in concentrations > 40 mg/dL	
	Salicylates and tetracyclines	
pH	Presence of bacteria that produce changes in the pH of the environment	False-positive results
Leukocytes	Increased hemoglobin concentrations	False-negative results
	Increased glucose concentrations	
	Increased density	
	Increased albumin concentrations	

RISK MANAGEMENT

Bioclin, after a detailed review and critical analysis of all known and/or anticipated hazards, concludes that all risks associated with the Biostrip Vet Product have been assessed, that risk reduction measures have been implemented and that the product does not present any major risks. that the benefits obtained from its use; and that, if used by qualified and trained veterinarian, aware of the precautions described in the products, it will perform its functions with quality, safety and effectiveness.

DIAGNOSTIC SIGNIFICANCE

Urine testing provides the clinician with information on possible pathologies of the urinary tract, as well as about some extra-reputable diseases. Due to its simplicity, low cost and the ease of obtaining the sample for analysis, it is a routine examination. Nutritional aspects, the amount of liquid ingested, physical activities before collection, the use of medications and the methodology for obtaining the sample must be considered when interpreting the results by the veterinarian.

BIBLIOGRAPHIC REFERENCES

- 1- BERG, B. et al. Guidelines for evaluation of reagent strips. Exemplified by analysis of urine albumin and glucose concentration using visually read reagent strips. Scand J Clin Lab Invest., v.49, n.7, p.689-699, 1989.
- 2- FREE, A.H.; FREE, H.M. Urinalysis: its proper role in the physician's office. Clin Lab Med., v.6, n.2, p.253-266, 1986.
- 3- GYURE, W.L. Comparison of several methods for semiquantitative determination of urinary protein. Clin Chem., v.23, n.5, p.876-879, 1977.
- 4- KOSS, S. et al. Proteinuria and renal disease: prognostic value of urine dipstick testing for leukocytes. Pediatr Nephrol., v.21, n.4, p.584-587, 2006.
- 5- WILSON, L.A. Urinalysis. Nurs Stand., v.19, n.35, p.51-54, 2005.
- 6- STRASINGER,K.S.; LORENZO. M.S. Urinálise e Fluidos Corporais. Livraria Médica Paulista., v.9, n.5, p.57-82, 2009.
- 7- eCLINPATH- Cornell University College of Veterinary Medicine, 2024. Disponível em: <https://eclinpath.com/>.
- 8- Orcurr, C.J. Ferret urogenital diseases. Veterinary Clinics: Exotic Animal Practice, 6(1), 113-138.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

■ QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Responsible Veterinary technician: Camila Eckstein CRMV/MG: 20611

Registration Number: Product exempt from registration at the Ministry of Agriculture, Livestock and Supply.

Review: April/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC
PRODUCT

TOXIC

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED

DO NOT REUSE

PRODUCT
STERILIZED

CAUTION

DANGER