COVID-19 PCR

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.

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Review: July/2022 (1)

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	4
Equipaments and Operational Inputs	4
Transportation and storage Conditions	4
Special Care	5
Samples	5
Process Description	6
A . RNAExtraction	6
B . Preparation of the reagents	6
C . Dilution of Quantitative Standard	6
D . PCR Preparation	7
E . Thermocycler settings for Real-Time PCR	8
F . Result Validation	10
G. Result Interpretation	11
Process Limitations	12
Product Performance / Quality Control	12
Comparison of Methods and Methodology Specificity	12
Repeatability	13
Reproducibility	13
Clinical Sensitivity	13
Analytical Sensitivity	13
Linearity	14
Diagnostic Significance	14
Bibliographic References	14
Customer Service	14
Universal Symbolology	15

FUNCTION

Quantitative or qualitative detection of SARS-CoV-2 (COVID-19) virus RNA through reverse transcription followed by real-time polymerase chain reaction (RT-PCR). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene COVID-19 PCR** test is a *in vitro* assay based on the quantitative or qualitative detection of SARS-CoV-2 virus RNA by real-time RT-PCR.

The Real-Time RT-PCR method is used to amplify the pathogen RNA.

A thermocycler for Real-Time PCR is used to amplify and detect the fluorescent probe. Software calculates the SARS-CoV-2 virus RNA concentration, expressed in copies/μL, using the standard curve generated from the quantitative standard, contained in the kit.

Reagent	Presentation
	100 Tests
R1	1 x 110 μL*
R2	1 x 110 μL*
R3	2 x 1.1 mL*
R4	2 x 1.5 mL
R5	1 x 500 μL*
R6	1 x 1.5 mL
R7	1 x 2 mL
R8	1 x 150 μL
R9	1 x 110 μL*

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

R1. E PCR Solution: Primer, probe, TRIS-HCl.

R2. RdRp PCR Solution: Primer, probe, TRIS-HCl.

R3. Mix Taq: Polymerases, dNTPs, $MgCl_2$, stabilizers.

R4. Mix Buffer: TRIS-HCl.

R5. Standard A (2×10^5 copies/ μL): Plasmid, TRIS-HCl, EDTA.

R6. Diluent: TRIS-HCl, EDTA

R7. Water: DNase/RNase free water.

R8. Negative Control: TRIS-HCl.

R9. Endogenous PCR Solution: Primer, probe, TRIS-HCl

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item.
- Usage instructions (manual).

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- 1.5 ml microcentrifuge tubes.
- 4- PCR reaction tubes or plate.
- 5- Powder free latex gloves or similar material.
- 6- Microcentrifuge (3,000 - 12,000 rpm).
- 7- Vortex.
- 8- Micropipettes and sterile filter pipette tips (0.5-10 μL , 10-100 μL , 100-1000 μL).
- 9- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C). **After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing.

The kit may be transported between 2 and 30°C for 7 days. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain accurate results.

3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.

4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.

5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.

7- Avoid repeated thawing and freezing of the reagents.

8- Store the RNA samples at -20°C in case they will not be used immediately.

9- Do not use reagents after expiration date.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with RNA samples extracted from sputum, nasopharyngeal or oropharyngeal swab. Samples should be collected according to laboratory recommendations for molecular testing. They must be transported and stored as between 2 to 8°C for up to 3 days¹.

RNA (extracted): Use RNA samples suitable for RT-PCR amplification with adequate purity and concentration. For long-term storage, it is recommended to freeze samples at -70°C. Avoid repeated freezing and thawing of the samples.

PROCESS DESCRIPTION

A. RNA Extraction

The nucleic acids (RNA) must be extracted in accord to the procedure of the chosen extraction kit.

B. Preparation of the Reagents

Note: The **R5** reagent contain RNA/DNA template. It need to be handle (resuspended) in a separated area.

1- Before open, spin (*spin* pulse) each of the following reagents **E PCR Solution (R1)**, **RdRp PCR Solution (R2)**, **Standard A (R5)** and **Endogenous PCR Solution (R9)**

2- Resuspend each vial of **Mix Taq (R3)*** reagent with 1.1mL of **Mix Buffer (R4)** reagent.

*Resuspend the **Mix Taq (R3)** vials as need of use.

3- Resuspend the reagents, **E PCR Solution (R1)**, **RdRp PCR Solution (R2)** and **Endogenous PCR Solution (R9)** with **Water (R7)** reagent, in according to the following table:

Reagent	Presentation
	100 Tests
E PCR Solution (R1)	110 µL
RdRp PCR Solution (R2)	110 µL
Endogenous PCR Solution (R9)	110 µL

C. Dilution of Quantitative Standard*

1- Ressuspend the **Standard A (R5)**** with 500µL do **Diluent (R6)**.

2- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of resuspended **Standard A (R5)**.

3- Pipette 90µL of the **Diluent (R6)** in each microtube and name them as B, C and D respectively.

4- Pipette 10µL of **Standard A (R5)** in microtube B and mix.

5- Replace the tip and pipette 10µL of microtube B in microtube C and mix.

6- Replace the tip and pipette 10µL of microtube C in microtube D and mix.

7- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/ μ L

Standard B - 2×10^4 copies/ μ L

Standard C - 2×10^3 copies/ μ L

Standard D - 2×10^2 copies/ μ L

*Standard curve dilution must be performed for quantitative testing.

**The standard curve must be prepared for the RdRp gene only.

D. PCR Preparation

ATTENTION: The **E PCR Solution (R1)** and **RdRp PCR Solution (R2)** reagents have the same probe (FAM). Therefore, for each sample, must be prepared 2 individual reactions, one reaction for E test and another reaction for RdRp test.

The **Endogenous PCR Solution (R9)** has VIC probe and sufficient volume to realize 100 reactions, so it must only be used in the reaction for **E**.

1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples, Controls and Quantitative Standards to be used.

D1. Preparation of PCR for the E gene and the Endogenous Control

1- Prepare the final PCR solution in accord with reactions number required.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R3)	10 μ L	250 μ L	500 μ L	1 mL
E PCR Solution (R1)	1 μ L	25 μ L	50 μ L	100 μ L
Endogenous PCR Solution (R9)	1 μ L	25 μ L	50 μ L	100 μ L
Water (R7)	3 μ L	75 μ L	150 μ L	300 μ L

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

D2. Preparation of PCR for the RdRp gene without Endogenous Control

1- Prepare the final PCR solution in accord with reactions number required.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R3)	10µL	250µL	500µL	1mL
RdRp PCR Solution (R2)	1µL	25µL	50µL	100µL
Water (R7)	4µL	100µL	200µL	400µL

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

2- Pipette 15µL of the final PCR solution in microtubes or wells previously determined.

3- Add 5µL of extracted RNA sample or 5µL of Quantitative Standards or 5µL of controls in previously determined tubes.

4- Mix well.

5- The total volume of the reaction is 20µL. Each PCR run should include the relevant controls: Negative/Positive Controls and Quantitative Standards.

6- Transport the tubes/plate to the thermocycler and follow the instructions in the section E (Thermocycler Settings for Real-Time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve or Qualitative Test.

Obs: In the case of qualitative testing, Standard A (R5) must be used as Positive Amplification Control.

2- Set the fluorescent probes:

Passive Reference: Equipments that use ROX as a passive reference must be programmed with the option “**NONE**” for passive reference.

Target	Detector	Quencher
E	FAM	NFQ-MGB
RdRp	FAM	
Endogenous Control	VIC	

Note: • Quantitative Standards do not present Endogenous Control (VIC), because it is used for extraction and amplification control of the samples.

• The extracted samples must be marked with the FAM and VIC detectors.

3- Set the quantitative standards*:

Standard A - 2×10^5 copies/ μ L

Standard B - 2×10^4 copies/ μ L

Standard C - 2×10^3 copies/ μ L

Standard D - 2×10^2 copies/ μ L

**Program used in the quantitative test.*

4- Set the cycling conditions:

Stage	Temperature	Time	Cycles
1	55°C	10 minutes	1
2	95°C	3 minutes	1
3	95°C	15 seconds	50
	60°C	60 seconds	

Set “Data Collection” as “stage 3, step 2 (60@0:60)”.

F. Result Validation

1- Standard Curve

Standard Curve	Permitted Range	Amplification / Detection
Correlation Coefficient (R ²)	$0,95 \leq R^2 \leq 1,00$	Valid

Note: The Threshold value can be used in automatic mode, however, it must be checked for the presence of noise that can result in Threshold values outside the exponential phase of the amplification curves, which in this case, must be adjusted manually. In case of doubt in the analysis, contact the SAC (Customer Advisory Service.).

2- Controls

Con- trols	CT		Result	Amplification/ Detection
	FAM	VIC		
Negative	Undetermined concentration or absence of amplification	Undetermined concentration or absence of amplification	Negative	Valid
Positive*	Determined concentration or presence of amplification	Undetermined concentration or absence of amplification	Positive	Valid

** Positive Control (**R5.Standard A**) used in the qualitative test..

3- Samples

SARS-CoV-2 (COVID-19)		Result	Detection
FAM ¹ (E or RdRp)	VIC ² Endogenous Control		
Determined concentration or Presence of amplification	CT ≤ 35	Positive	Valid
	CT > 35	Positive	Invalid*
Undetermined concentration or absence of amplification	CT ≤ 35	Negative	Valid
	CT > 35	Negative	Invalid*

¹COVID-19

For the SARS-CoV-2 (COVID-19) positive diagnosis, it is necessary to detect **2 target genes** (E and RdRp).

²Endogenous Control

The CT value of the Endogenous Control may vary according to the amount of biological material present in the extracted sample.

- The low CT indicates high concentration of biological material in the sample and that the extraction process was efficient.
- The high CT indicates low concentration of biological material in the sample and / or inefficient extraction.

***Note:** The CT values of Endogenous Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Endogenous Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of RNA per reaction. The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005, this means that the sample concentration is 2.0×10^5 copies/ μ L.

Sample Result in Copies/ μ L (FAM)	Copies per Reaction
$\geq 1 \times 10^6$	> 1.000.000
$2 \leq \text{Quantity} \leq 9 \times 10^5$	Quantity obtained
< 2	< 10

The non-detection of the pathogen RNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient. The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

1- Conversion calculation for copies/mL

$$\text{Result in copies/mL} = \frac{\text{Copies}/\mu\text{L} \times \text{Elution volume } (\mu\text{L})^{**}}{\text{Extracted sample volume (mL)}^{**}}$$

*Result of quantification of the sample provided by the equipment in copies/ μL .

**According to the protocol of the extraction kit used.

Example: Elution volume: 30 μL

Sample volume: 200 μL (0.2 mL)

Quantification of the sample: 5000 copies/ μL

$$\text{Result in copies/mL} = \frac{(5000 \text{ copies}/\mu\text{L} \times 30\mu\text{L})}{0.2\text{mL}}$$

$$\text{Result in copies/mL} = 750.000 \text{ copies/mL}$$

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

The heparin use in the samples collection for molecular tests is not recommended, due to the PCR inhibition which can generate false negative results. The samples with presence of PCR inhibitors should be diluted to reduce their effect on the reaction.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene COVID-19 PCR** test was compared with another protocol for the quantification of the COVID-19 virus RNA. 50 analyzes were performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This show that the test has good methodological specificity.

Precision REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Average Concentration (Log)	2.71	3.63
Standard Deviation (Log)	0.07	0.10
Coefficient of variation (%)	2.46	2.77

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

	Day 1	Day 2	Day 3
Average Concentration (Log)	2.71	2.69	2.75
Standard Deviation (Log)	0.07	0.10	0.70
Coefficient of variation (%)	2.46	3.83	2.49

Clinical Sensitivity

The **Bio Gene COVID-19 PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

Method		Clinical Data		Total
Bio Gene COVID-19 PCR	Result	Positive	Negative	
	Positive	24	0	24
	Negative	0	26	26
TOTAL		24	26	50

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 0.61 copies/ μ L (- 1.21 log). This means the assay is capable to detect approximately 0.305 target molecules in 5 μ L RNA extraction product, used on the amplification reaction. The sensitivity was 18.3 copies/mL.*

* The sensitivity may vary according to the used extraction kit.

Linearity

The linearity found in the performance studies was 1.00×10^5 copies/ μL (5.0 log).

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the family Coronaviridae, and they are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnea, myalgia or fatigue. It can cause diseases ranging from cold common to more serious diseases, such as Respiratory Syndrome Middle East (MERS-CoV) and Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with symptoms very similar to viral pneumonia, after sequencing of the respiratory samples was found to be infection by Coronavirus, which was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents, cases of infection have been reported. Therefore, the detection of COVID-19 by real time PCR is very important due to its accuracy and agility.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene COVID-19 PCR: 10269360323

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER

Bioclin · QUIBASA

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130

Tel +55 31 3439 5454 . www.bioclin.com.br

FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422

C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira