BIO GENE

COVID-19 PCR

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.

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Review: April/2020

INDEX

| | |
|---|----|
| Function | 3 |
| Principle of Action | 3 |
| Presentation | 3 |
| Reagents | 4 |
| Equipments and Operational Inputs | 4 |
| Transportation and storage Conditions | 4 |
| Special Care | 5 |
| Samples | 5 |
| Process Description | 6 |
| A . RNAExtraction | 6 |
| B . Preparation of the reagents | 6 |
| C . Dilution of Quantitative Standard | 7 |
| D . PCR Preparation | 7 |
| E . Thermocycler settings for Real-Time PCR | 9 |
| F . Result Validation | 10 |
| G. Result Interpretation | 11 |
| Process Limitations | 12 |
| Product Performance / Quality Control | 12 |
| Comparison of Methods and Methodology Specificity | 12 |
| Repeatability | 13 |
| Reproducibility | 13 |
| Clinical Sensitivity | 13 |
| Analytical Sensitivity | 14 |
| Linearity | 14 |
| Diagnostic Significance | 14 |
| Bibliographic References | 14 |
| Customer Service | 14 |
| Universal Symbology | 15 |

FUNCTION

Quantitative detection of SARS-CoV-2 (COVID-19) virus RNA through reverse transcription followed by real-time polymerase chain reaction (RT-PCR). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene COVID-19 PCR** test is a *in vitro* assay based on the quantitative detection of COVID-19 virus RNA by real-time RT-PCR.

The Real-Time RT-PCR method is used to amplify the pathogen RNA.

A thermocycler for Real-Time PCR is used to amplify and detect the fluorescent probe. Software calculates the COVID-19 virus RNA concentration, expressed in copies/ μL , using the standard curve generated from the quantitative standard, contained in the kit.

| Reagent | Presentation |
|-----------|-------------------------|
| | 100 Tests |
| R1 | 1 x 110 μL * |
| R2 | 1 x 110 μL * |
| R3 | 4 x 600 μL * |
| R4 | 2 x 1.5 mL |
| R5 | 1 x 500 μL * |
| R6 | 1 x 1.5 mL |
| R7 | 1 x 2 mL |
| R8 | 1 x 150 μL |
| R9 | 1 x 110 μL * |

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

R1. E PCR Solution: Primer, probe, TRIS-HCl.

R2. RdRp PCR Solution: Primer, probe, TRIS-HCl.

R3. Mix Taq: Polymerases, dNTPs, MgCl₂, stabilizers.

R4. Mix Buffer: TRIS-HCl.

R5. Standard A (2 x 10⁵ copies/μL): Plasmid, TRIS-HCl, EDTA.

R6. Diluent: TRIS-HCl, EDTA

R7. Water: DNase/RNase free water.

R8. Negative Control: TRIS-HCl.

R9. Endogenous PCR Solution: Primer, probe, TRIS-HCl

EQUIPMENTS AND OPERATIONAL INPUTS**Materials in the kit:**

- Reagents described in the previous item.
- Usage instructions (manual).

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- 1.5 ml microcentrifuge tubes.
- 4- PCR reaction tubes or plate.
- 5- Powder free latex gloves or similar material.
- 6- Microcentrifuge (3,000 - 12,000 rpm).
- 7- Vortex.
- 8- Micropipettes and sterile filter pipette tips (0.5-10μL, 10-100μL, 100-1000μL).
- 9- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C). **After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing.

The kit may be transported between 2 and 30°C for 7 days. Protect from light and avoid moisture.

SPECIAL CARE**1- For professional *in vitro* diagnostic use only.**

- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.
- 4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.
- 5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.
- 6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.
- 7- Avoid repeated thawing and freezing of the reagents.
- 8- Store the RNA samples at -20°C in case they will not be used immediately.
- 9- Do not use reagents after expiration date.
- 10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 12- Do not use the product in case of damaged packaging.
- 13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with RNA samples extracted from sputum, nasopharyngeal or oropharyngeal swab. Samples should be collected according to laboratory recommendations for molecular testing. They must be transported and stored as between 2 to 8°C for up to 3 days¹.

RNA (extracted): Use RNA samples suitable for RT-PCR amplification with adequate purity and concentration.

For long-term storage, it is recommended to freeze samples at -70°C. Avoid repeated freezing and thawing of the samples.

PROCESS DESCRIPTION

A. RNA Extraction

The nucleic acids (RNA) must be extracted in accord to the procedure of the chosen extraction kit.

B. Preparation of the Reagents

Note: The **R5** reagent contain RNA/DNA template. It need to be handle (resuspended) in a separated area.

1- Before open, spin (*spin* pulse) each of the following reagents: **E PCR Solution (R1)**, **RdRp PCR Solution (R2)**, **Standard A (R5)** and **Endogenous PCR Solution (R9)**

2- Resuspend each vial of **Mix Taq (R3)*** reagent with 600µL of **Mix Buffer (R4)** reagent.

*Resuspend the **Mix Taq (R3)** vials as need of use.

3- Resuspend the reagents, **E PCR Solution (R1)**, **RdRp PCR Solution (R2)** and **Endogenous PCR Solution (R9)** with **Water (R7)** reagent, in according to the following table:

| Reagent | Presentation |
|------------------------------|--------------|
| | 100 Tests |
| E PCR Solution (R1) | 110 µL |
| RdRp PCR Solution (R2) | 110 µL |
| Endogenous PCR Solution (R9) | 110 µL |

C. Dilution of Quantitative Standard

- 1- Resuspend the **Standard A (R5)** with 500 μ L of **Diluent (R6)**.
- 2- Select 3 microtubes (not supplied in the kit) suitable for serial dilution of resuspended **Standard A (R5)**.
- 3- Pipette 90 μ L of the **Diluent (R6)** in each microtube and name them as B, C and D respectively.
- 4- Pipette 10 μ L of **Standard A (R5)** in microtube B and mix.
- 5- Replace the tip and pipette 10 μ L of microtube B in microtube C and mix.
- 6- Replace the tip and pipette 10 μ L of microtube C in microtube D and mix.
- 7- At the end of dilution procedure, it will be available 4 standards with the following concentrations:

Standard A - 2×10^5 copies/ μ L

Standard B - 2×10^4 copies/ μ L

Standard C - 2×10^3 copies/ μ L

Standard D - 2×10^2 copies/ μ L

D. PCR Preparation

ATTENTION: The **E PCR Solution (R1)** and **RdRp PCR Solution (R2)** reagents have the same probe (FAM). Therefore, for each sample, must be prepared 2 individual reactions, one reaction for E test and another reaction for RdRp test.

The **Endogenous PCR Solution (R9)** has VIC probe and sufficient volume to realize 100 reactions, so it must only be used in the reaction for **E**.

- 1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples, Controls and Quantitative Standards to be used.

D1. Preparation of PCR for the E gene and the Endogenous Control

1- Prepare the final PCR solution in accord with reactions number required.

| Reagents | 1 Reaction | 25 Reactions | 50 Reactions | 100 Reactions |
|------------------------------|------------|--------------|--------------|---------------|
| Mix Taq (R3) | 10 µL | 250 µL | 500 µL | 1 mL |
| E PCR Solution (R1) | 1 µL | 25 µL | 50 µL | 100 µL |
| Endogenous PCR Solution (R9) | 1 µL | 25 µL | 50 µL | 100 µL |
| Water (R7) | 3 µL | 75 µL | 150 µL | 300 µL |

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

D2. Preparation of PCR for the RdRp gene without Endogenous Control

1- Prepare the final PCR solution in accord with reactions number required.

| Reagents | 1 Reaction | 25 Reactions | 50 Reactions | 100 Reactions |
|------------------------|------------|--------------|--------------|---------------|
| Mix Taq (R3) | 10µL | 250µL | 500µL | 1mL |
| RdRp PCR Solution (R2) | 1µL | 25µL | 50µL | 100µL |
| Water (R7) | 4µL | 100µL | 200µL | 400µL |

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

2- Pipette 15µL of the final PCR solution in microtubes or wells previously determined.

3- Add 5 μ L of extracted RNA sample or 5 μ L of Quantitative Standards or 5 μ L of **Negative Control (R8)** in previously determined tubes.

4- Mix well.

5- The total volume of the reaction is 20 μ L. Each PCR run should include the relevant controls: Negative Control and Quantitative Standards.

6- Transport the tubes/plate to the thermocycler and follow the instructions in the section E (Thermocycler Settings for Real-Time PCR).

E. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

1- Set the type of experiment:

Quantitative Test with Standard Curve.

2- Set the fluorescent probes:

Passive Reference: Equipments that use ROX as a passive reference must be programmed with the option “**NONE**” for passive reference.

| Target | Detector | Quencher |
|--------------------|----------|----------|
| E | FAM | NFQ-MGB |
| RdRp | FAM | |
| Endogenous Control | VIC | |

Note: • Quantitative Standards do not present Endogenous Control (VIC), because it is used for extraction and amplification control of the samples.

• The extracted samples must be marked with the FAM and VIC detectors.

3- Set the quantitative standards:

Standard A - 2 x 10⁵ copies/ μ L

Standard B - 2 x 10⁴ copies/ μ L

Standard C - 2 x 10³ copies/ μ L

Standard D - 2 x 10² copies/ μ L

4- Set the cycling conditions:

| Stage | Temperature | Time | Cycles |
|-------|-------------|------------|--------|
| 1 | 55°C | 10 minutes | 1 |
| 2 | 95°C | 3 minutes | 1 |
| 3 | 95°C | 15 seconds | 50 |
| | 60°C | 60 seconds | |

Set "Data Collection" as "stage 3, step 2 (60@0:60)".

F. Result Validation**1- Standard Curve**

| Standard Curve | Permitted Range | Amplification / Detection |
|-----------------------------------|---------------------------|---------------------------|
| Correlation Coefficient (R^2) | $0,99 \leq R^2 \leq 1,00$ | Valid |

If the R^2 value does not fall within the limits of range, the result must be considered invalid and the test must be repeated.

2- Negative Control

| CT Negative Control | | Result | Amplification/ Detection |
|---------------------|--------------|----------|-----------------------------|
| FAM | VIC | | |
| Undetermined | Undetermined | Negative | Valid |

3- Samples

| COVID-19 | | Result | Detection |
|---------------------------------|--|----------|-----------|
| FAM ¹ (E or RdRp) | VIC ² Endogenous Control | | |
| Determined concentration | CT \leq 35 | Positive | Valid |
| | CT $>$ 35 | Positive | Invalid* |
| Undetermined concentration | CT \leq 35 | Negative | Valid |
| | CT $>$ 35 | Negative | Invalid* |

¹COVID-19

For the COVID-19 positive diagnosis, it is necessary to detect **2 target genes** (E and RdRp).

²Endogenous Control

The CT value of the Endogenous Control may vary according to the amount of biological material present in the extracted sample.

- The low CT indicates high concentration of biological material in the sample and that the extraction process was efficient.
- The high CT indicates low concentration of biological material in the sample and / or inefficient extraction.

***OBS:** The CT values of Endogenous Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Endogenous Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

G. Result Interpretation

This kit is able to detect from 10 to 1.000.000 copies of RNA per reaction. The thermocycler software calculates automatically the samples concentration.

Example: If the program shows a concentration of 2.00E+005, this means that the sample concentration is 2.0×10^5 copies/ μ L.

| Sample Result in Copies/ μ L (FAM) | Copies per Reaction |
|---|---------------------|
| $\geq 1 \times 10^6$ | > 1.000.000 |
| $2 \leq \text{Quantity} \leq 9 \times 10^5$ | Quantity obtained |
| < 2 | < 10 |

The non-detection of the pathogen RNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit. The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient. The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

1- Conversion calculation for copies/mL

$$\text{Result in copies/mL} = \frac{\text{Copies}/\mu\text{L} * \text{Elution volume } (\mu\text{L})^{**}}{\text{Extracted sample volume (mL)}^{**}}$$

*Result of quantification of the sample provided by the equipment in copies/ μL .

**According to the protocol of the extraction kit used.

Example: Elution volume: 30 μL

Sample volume: 200 μL (0.2 mL)

Quantification of the sample: 5000 copies/ μL

$$\text{Result in copies/mL} = \frac{(5000 \text{ copies}/\mu\text{L} \times 30\mu\text{L})}{0.2\text{mL}}$$

$$\text{Result in copies/mL} = 750.000 \text{ copies/mL}$$

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results.

The heparin use in the samples collection for molecular tests is not recommended, due to the PCR inhibition which can generate false negative results.

The samples with presence of PCR inhibitors should be diluted to reduce their effect on the reaction.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene COVID-19 PCR** test was compared with another protocol for the quantification of the COVID-19 virus RNA. 50 analyzes were

performed and the results were evaluated. Comparing the two tests, the difference between the results obtained was ≤ 0.5 log. This shows that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples with different concentrations, obtaining the following results:

| | Sample 1 | Sample 2 |
|-------------------------------------|----------|----------|
| Average Concentration (Log) | 2.71 | 3.63 |
| Standard Deviation (Log) | 0.07 | 0.10 |
| Coefficient of variation (%) | 2.46 | 2.77 |

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample with different concentrations, obtaining the following results:

| | Day 1 | Day 2 | Day 3 |
|-------------------------------------|-------|-------|-------|
| Average Concentration (Log) | 2.71 | 2.69 | 2.75 |
| Standard Deviation (Log) | 0.07 | 0.10 | 0.70 |
| Coefficient of variation (%) | 2.46 | 3.83 | 2.49 |

Clinical Sensitivity

The **Bio Gene COVID-19 PCR** showed clinical sensitivity of 99.9% and a clinical specificity of 99.9%.

| Method | | Clinical Data | | Total |
|--------------------------|----------|---------------|----------|-------|
| Bio Gene COVID-19 PCR | Result | Positive | Negative | |
| | Positive | 24 | 0 | 24 |
| | Negative | 0 | 26 | 26 |
| TOTAL | | 24 | 26 | 50 |

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 0.61 copies/ μ L (- 1.21 log). This means the assay is capable to detect approximately 0.305 target molecules in 5 μ L RNA extraction product, used on the amplification reaction. The sensitivity was 18.3 copies/mL.*

* The sensitivity may vary according to the used extraction kit.

Linearity

The linearity found in the performance studies was 1.00×10^5 copies/ μ L (5.0 log).

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the family Coronaviridae, and they are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnea, myalgia or fatigue. It can cause diseases ranging from cold common to more serious diseases, such as Respiratory Syndrome Middle East (MERS-CoV) and Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with symptoms very similar to viral pneumonia, after sequencing of the respiratory samples was found to be infection by Coronavirus, which was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents, cases of infection have been reported. Therefore, the detection of COVID-19 by real time PCR is very important due to its accuracy and agility.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 031 5454
E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene COVID-19 PCR: 10269360323

UNIVERSAL SYMBOLOGY

| | | | |
|---|---------------------------------------|---|--|
|  | CATALOG NUMBER |  | MANUFACTURED BY |
|  | BATCH CODE |  | CONTROL |
|  | DATE OF MANUFACTURE |  | POSITIVE CONTROL |
|  | USED BY (last day of month) |  | NEGATIVE CONTROL |
|  | TEMPERATURE LIMITATION (store at) |  | BIOLOGICAL RISK |
|  | CONTAINS SUFFICIENT FOR <N> TESTS |  | INFLAMMABLE |
|  | CONSULT INSTRUCTIONS FOR USE |  | CORROSIVE |
|  | IN VITRO DIAGNOSTIC DEVICE |  | POISON |
|  | EUROPEAN AUTHORIZED REPRESENTATIVE |  | CE MARK |
|  | KEEP AWAY FROM SUNLIGHT |  | DO NOT USE IF PACKAGE IS DAMAGED |

BIO GENE

Bioclin · QUIBASA

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130

Tel +55 31 3439 5454 . www.bioclin.com.br

FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422

C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira