COVID-19 Multiplex PCR

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.

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Review: June/2022

INDEX

Function	3
Principle of Action	3
Presentation	3
Reagents	3
Equipments and Operational Inputs	4
Transportation and storage Conditions	4
Special Care	4
Samples	5
Process Description	5
A . RNAExtraction	5
B . Preparation of the reagents	5
C . PCR Preparation	6
D . Thermocycler settings for Real-Time PCR	6
E . Result Validation	7
F . Result Interpretation	9
Process Limitations	10
Cross Reactivity	10
Product Performance / Quality Control	10
Comparison of Methods and Methodology Specificity	10
Repeatability	10
Reproducibility	10
Clinical Sensitivity	11
Analytical Sensitivity	11
Linearity	11
Diagnostic Significance	11
Bibliographic References	12
Customer Service	12
Universal Symbolology	13

FUNCTION

Qualitative detection of SARS-CoV-2 (COVID-19) virus RNA through reverse transcription followed by real-time polymerase chain reaction (RT-PCR). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

The **Bio Gene COVID-19 Multiplex PCR** test is a *in vitro* assay based on the qualitative detection of SARS-CoV-2 virus RNA by real-time RT-PCR.

The Real-Time RT-PCR method is used to amplify the pathogen RNA.

A thermocycler for Real Time PCR is used to amplify and detect the fluorescent probes for each specific target. The software detects the presence of the pathogens.

Reagent	Presentation
	100 Tests
R1	1 x 110 µL*
R2	1 x 1.1 mL*
R3	1 x 1.5 mL
R4	1 x 500 µL*
R5	1 x 1.5 mL
R6	1 x 1.5 mL
R7	1 x 200 µL

*Lyophilized reagents. The volumes described above refer to the final volume after the reagents resuspension, as described in the item PROCESS DESCRIPTION, subitem B (Preparation of the Reagents).

REAGENTS

R1. PCR Solution: Primer, probe, TRIS-HCl.

R2. Mix Taq: Polymerases, dNTPs, MgCl₂, stabilizers.

R3. Mix Buffer: TRIS-HCl.

R4. Positive Control: Plasmid, TRIS-HCl, EDTA.

R5. Diluent: TRIS-HCl, EDTA

R6. Water: DNase/RNase free water.

R7. Negative Control: TRIS-HCl, EDTA.

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the previous item.
- Usage instructions (manual).

Materials needed but not contained in the kit:

- 1- Real-time PCR thermocycler.
- 2- Laminar flow hood.
- 3- 1.5 ml microcentrifuge tubes.
- 4- PCR reaction tubes or plate.
- 5- Powder free latex gloves or similar material.
- 6- Microcentrifuge (3,000 - 12,000 rpm).
- 7- Vortex.
- 8- Micropipettes and sterile filter pipette tips (0.5-10 μ L, 10-100 μ L, 100-1000 μ L).
- 9- Nucleic acid extraction kit.

TRANSPORTATION AND STORAGE CONDITIONS

The required storage temperature is -20°C (-10 to -30°C). **After resuspension of lyophilized reagents, the product is stable for 6 months from the date of resuspension.** Avoid freezing and thawing.

The kit may be transported between 2 and 30°C for 7 days. Protect from light and avoid moisture.

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- Handle and dispose of all biological samples, reagents and materials as it contain infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills and aerosol. Handle and dispose of wastes in accordance with appropriate security procedures.
- 4- It is require skilled personnel for the molecular biology procedures in order to minimize the risk of erroneous results, degradation of nucleic acids contained in the samples or even sample contamination by amplicons.
- 5- It is necessary to have separate areas for the extraction / preparation of amplification reactions and for the amplification / detection of amplified products. Never introduce an amplification product in the area intended for the extraction or preparation of amplification products.

6- All samples and reagents should be handled inside a laminar flow hood. It is necessary to use sterile and DNase/RNase free filter pipettes tips.

7- Avoid repeated thawing and freezing of the reagents.

8- Store the RNA samples at -20°C in case they will not be used immediately.

9- Do not use reagents after expiration date.

10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

11- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

12- Do not use the product in case of damaged packaging.

13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

This kit should be used with RNA samples extracted from sputum, nasopharyngeal or oropharyngeal swab. Samples should be collected according to laboratory recommendations for molecular testing. They must be transported and stored as between 2 to 8°C for up to 3 days¹.

RNA (extracted): Use RNA samples suitable for RT-PCR amplification with adequate purity and concentration. For long-term storage, it is recommended to freeze samples at -70°C. Avoid repeated freezing and thawing of the samples.

PROCESS DESCRIPTION

A. RNA Extraction

The nucleic acids (RNA) must be extracted in accord to the procedure of the chosen extraction kit.

B. Preparation of the Reagents

Note: The **R4** reagent contain RNA/DNA template. It need to be handle (resuspended) in a separated area.

1- Before open, spin (*spin* pulse) each of the following reagents **PCR Solution (R1)** and **Positive Control (R4)**

- 2- Resuspend the vial of **Mix Taq (R2)** reagent with 1.1 mL of **Mix Buffer (R3)** reagent.
- 3- Resuspend the reagent **PCR Solution (R1)** with 110 μL of **Water (R6)** reagent.
- 4- Resuspend the reagent **Positive Control (R4)** with 500 μL of **Diluent (R5)** reagent.

C. PCR Preparation

- 1- Prior beginning the assay, select sufficient number of microtubes/wells to be used on the reaction, in accord to the number of samples and Controls to be used.
- 2- Prepare the final PCR solution in accord with reactions number required.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Mix Taq (R2)	10 μL	250 μL	500 μL	1 mL
PCR Solution (R1)	1 μL	25 μL	50 μL	100 μL

To prepare different reaction number must be multiplied the volume of reagents for a reaction by the number of required reactions.

- 3- Pipette 11 μL of the final PCR solution in microtubes or wells previously determined.
- 4- Add 9 μL of extracted RNA sample or 9 μL of controls in previously determined tubes.
- 5- Mix well.
- 6- The total volume of the reaction is 20 μL . Each PCR run should include the relevant controls: Negative/Positive Controls.
- 7- Transport the tubes/plate to the thermocycler and follow the instructions in the section D (Thermocycler Settings for Real-Time PCR).

D. Thermocycler Settings for Real-Time PCR

Check the operation manual of the Real-time PCR thermocycler to program the experiment.

- 1- **Set the type of experiment:**
Qualitative Test.

2- Set the fluorescent probes:

Passive Reference: Equipments that use ROX as a passive reference must be programmed with the option “**NONE**” for passive reference.

Target	Detector	Quencher
E	FAM	NFQ-MGB
RdRp	ROX	
Endogenous Control	VIC	

Note:

- Controls do not present Endogenous Control (VIC), because it is used for extraction and amplification control of the samples.
- The extracted samples must be marked with the FAM, ROX and VIC detectors.

3- Set the cycling conditions:

Stage	Temperature	Time	Cycles
1	55°C	10 minutes	1
2	95°C	3 minutes	1
3	95°C	15 seconds	50
	60°C	60 seconds	

Set “Data Collection” as “stage 3, step 2 (60@0:60)”.

E. Result Validation

1- Controls

Controls	CT Permitted Range	Amplification/Detection
Positives (FAM e ROX)	$16 \leq CT \leq 21$	Valid
Negative	Undetermined	Valid
Endogenous Control (VIC)	$CT \leq 35$	Valid

2- Samples

SARS-CoV-2 (COVID-19)		Result	Detection
FAM/ROX ¹	VIC ² Endogenous Control		
Presence of amplification	CT ≤ 35	Positive	Valid
	CT > 35	Positive	Invalid*
Absence of amplification	CT ≤ 35	Negative	Valid
	CT > 35	Negative	Invalid*

¹COVID-19

For the SARS-CoV-2 (COVID-19) positive diagnosis, it is necessary to detect **2 target genes** (E and RdRp).

²Endogenous Control

The CT value of the Endogenous Control may vary according to the amount of biological material present in the extracted sample.

- The low CT indicates high concentration of biological material in the sample and that the extraction process was efficient.
- The high CT indicates low concentration of biological material in the sample and / or inefficient extraction.

***Note:** The CT values of Endogenous Control vary with the process conditions, as the extraction efficiency of the DNA/RNA, the concentration of the samples and the thermocycler settings. Therefore, these conditions must be evaluated when the CT values are not appropriated and, in case it is relevant, the results can be validated.

Example: In some cases, the samples with high numbers of copies of DNA / RNA can inhibit the amplification of the Endogenous Control, resulting in CT value outside the optimal range, this does not invalidate the test results.

If the above requirements are not achieved, the assay must be considered invalid and must be repeated.

F. Result Interpretation

SARS-CoV-2 (COVID-19)			Results
Gene E (FAM)	Gene RdRp (ROX)	Endogenous Control (VIC)	
CT Determined	CT Determined	CT ≤35	Positive
CT Determined	CT Undetermined		Inconclusive
CT Undetermined	CT Determined		Inconclusive
CT Undetermined	CT Undetermined		Negative

Determined CT: presence of curve amplification / Undetermined CT: absence of curve amplification.

Obs.: The result is valid as long as all the controls present values between the of permitted range limits and that the aspect of the sample curve tested be similar the curves obtained in the Positive Control (R4) with sigmoid form.

This kit is able to detect from 10 to 1.000.000 copies per reaction.

The non-detection of the pathogen RNA does not exclude the presence of infection when the pathogen title is below the detection limit of this kit.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only parameter for the determination of diagnosis and/or treatment of the patient.

The results obtained must be evaluated considering the clinical data and laboratory tests of the patient.

PROCESS LIMITATIONS

Cross contamination that eventually occurs during the sample collection, processing, transportation and storage may give false results. The samples with presence of PCR inhibitors should be diluted to reduce their effect on the reaction.

CROSS REACTIVITY

A cross-reactivity study was carried out evaluating 50 negative samples for SARS-CoV-2, but positive for other infections. Among them 12 positive samples for Influenza virus, 15 positive samples for Rhinovirus and 13 positive samples for Respiratory Syncytial Virus, and 10 positive samples for endemic Coronavirus. It was not observed cross-reactivity with positive samples for Influenza virus, Rhinovirus, Respiratory syncytial virus and endemic Coronavirus.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The **Bio Gene COVID-19 Multiplex PCR** test was compared with another protocol for the detection of the SARS-CoV-2 virus RNA. 200 analyzes were performed and the results were evaluated. The results obtained among the kits presented 99% agreement. This shows that the test has good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 positive samples, obtaining the following results:

Samples	N° of repetitions	COVID-19
Sample 27	10	100% Positive
Sample 28	10	100% Positive

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 1 sample, obtaining the following results:

Samples	Day 1	Day 2	Day 3
Sample 28	100% Positive	100% Positive	100% Positive

Clinical Sensitivity

The **Bio Gene COVID-19 Multiplex PCR** showed clinical sensitivity of 98% (CI 95%: 94,1% – 100%) and a clinical specificity >99,9%.

Method		Clinical Data		Total
Bio Gene COVID-19 Multiplex PCR	Result	Positive	Negative	
	Positive	49	0	49
	Negative	1	50	51
TOTAL		50	50	100

Analytical Sensitivity

It was found in the performance study that the analytical sensitivity of the kit is 1.20 copies/μL (10 copies/reaction) for the E gene, and 6.55 copies/μL (60 copies/reaction) for the RdRp gene. This means the assay is capable to detect approximately 2000 copies/mL.

* The sensitivity may vary according to the used extraction kit.

Linearity

The linearity found in the performance studies was 1×10^5 copies/μL (5.0 log).

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the family Coronaviridae, and they are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnea, myalgia or fatigue. It can cause diseases ranging from cold common to more serious diseases, such as Respiratory Syndrome Middle East (MERS-CoV) and Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with symptoms very similar to viral pneumonia, after sequencing of the respiratory samples was found to be infection by Coronavirus, which

was called the new coronavirus (2019-nCoV). In all continents, cases of infection have been reported. Therefore, the detection of COVID-19 by real time PCR is very important due to its accuracy and agility.

BIBLIOGRAPHIC REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Corman VM, Landt O, *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
3. Huang C, Wang Y, *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 031 5454

E-mail: sac@bioclin.com.br

ANVISA registration for Bio Gene COVID-19 Multiplex PCR: 10269360363

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER

Bioclin · QUIBASA

Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130

Tel +55 31 3439 5454 . www.bioclin.com.br

FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422

C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira