CÁLCIO ARSENAZO III

REF K305

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INSTRUCTIONS FOR USE

FUNCTION

Quantitative determination of calcium in biological samples of serum, plasma and urine by color reaction. Manual and Automatic application.

PRINCIPLE OF ACTION

Methodology: Colorimetric Endpoint - Arsenazo III

The calcium reacts with arsenazo III in acid environment forming the blue-colored complex, whose intensity is proportional to the calcium concentration in the sample. The absorbance of the reaction product is measured in wavelengths between 600 and 680 nm.

REAGENTS

Number 1: Arsenazo - Store between 2 and 8°C. Contains: Arsenazo III 0.2 mmol/L, 8-Hydroxyquinoline < 20 mmol/L, buffers, surfactants and preservative.

Number 2: Standard - Store between 2 and 8°C. Contains: Calcium Carbonate 10 mg/dL, stabilizer and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2*
1	50 mL	3 mL
2	2 x 50 mL	3 mL
3	5 x 20 mL	3 mL
4	5 x 40 mL	3 mL
5	10 x 40 mL	3 mL

* The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Automated instruments and semi-automatic biochemistry, clock or stopwatch, pipettes, test tubes, a water bath at 37°C, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- The validity of the reactants is printed on the labels of same. However, once opened, the reagent can suffer from contamination during handling and have stability decreased. Important to note changes (Color change, precipitate formation) that may indicate contamination.

3- Reagents N° 1 and N°2 contains Sodium Azide which can be irritating to skin and mucous membranes. Handle with care.

4- To follow strictly the methodology proposed for obtaining accurate results.

5- The water used in material cleaning must be recent and free from contaminants.

6- Columns deionizer saturated releases alkaline water, several ions and oxidizing and reducing agents, which can significantly

alter the results.

7- It is extremely important the thorough cleaning of all material (pipettes, tips, test tubes, cuvettes, etc.) to obtain accurate results. Using solution sulphochromic or aqueous solutions of HCl 50% or detergent nonionic. Rinse the material in running water and rinse with distilled/deionized water. Dry the material under heater protected from direct contact with the metal part of (put the material in glass vats).

8- Do not allow contact of the sample anticoagulant such as EDTA that is a chelat calcium, citrate that complex Ca⁺², oxalate and fluoride that precipitate Calcium.

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.

12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (Heparin): Free of hemolysis. Do not use Oxalate or EDTA as an anticoagulant by interfering the determination of calcium. Separate serum or plasma within one hour after collection, due to increased permeability of red blood cells to calcium. The analyte is stable in serum or plasma for 7 days at 2 to 8°C and 2 months at -20°C.

Urine: Make the urine collection within 24 hours in a recipient free of calcium. Prior to the collection add into container 10 mL of Nitric Acid 50% v/v. When the sample of 24 hours is not acidified during collection, add 20 mL of HCl 6M, homogenize, wait 60 minutes and take the aliquot for the test. The analyte is stable in urine for 4 days at 2 to 8°C and 30 days at -20°C¹³.

PROCESS DESCRIPTION

The calibration stability of the Calcium Arsenazo III kit installed on refrigerated equipment is at least 28 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Select the wavelength from 600 nm to 680 nm.

Mark 3 test tubes: B (Blank), P (Standard) and A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL
Reagent N° 2	---	10 µL	---
Sample	---	---	10 µL

Homogenize and incubate for 2 minutes at 37°C.

Read the absorbance of Standard and Sample at 600nm to 680nm, matching the zero with Blank. The color is stable for 1 hour.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

CALCULATIONS

Serum Calcium:

Calcium (mg/dL) = $\frac{\text{Test Absorbance}}{\text{Standard Absorbance}} \times 10$

As the reaction follows the Beer-Lambert law, the Calibration Factor can be used.

Calibration Factor = $\frac{10}{\text{Standard Absorbance}}$

Calcium (mg/dL) = Absorbance Test x Factor

Urinary Calcium:

Urinary Calcium = $\frac{\text{Urinary Calcium (mg/dL)} \times \text{Volume (mL)}}{(\text{mg/24h}) \quad 100}$

Ionized Calcium (Cal):

Cal (mg/dL) = $\frac{6 \times \text{Ca} - [(0.19 \times \text{P}) + \text{A}]}{(0.19 \times \text{P}) + \text{A} + 6}$

Ca = Calcium Serum (mg/dL)

P = Total Proteins (g/dL)

A = Albumin (g/dL)

INTERFERENTS

Do not use hemolyzed samples.

Icteric and lipemic samples may interfere with measurements.

Proceed to the dosage as follows:

1 - Dosing the sample normally following the described technique. Record the readings.

2 - Add 0.02 mL of the anticoagulant EDTA Bioclin to the blank and test it. Proceed with the new reading hitting the zero with the blank.

3 - Express the difference multiplying the difference between the two readings by the Calibration Factor.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 915B.

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of calcium in healthy populations of male and female.

Calcium (Serum or Plasma)	
Umbilical Cord	8.2 - 11.2 mg/dL
Premature	6.2 - 11.0 mg/dL
0 - 10 days	7.6 - 10.4 mg/dL
Infants	9.0 - 11.0 mg/dL
2 - 12 years	8.8 - 10.8 mg/dL
Adults	8.8 - 11.0 mg/dL
Calcium Ion	
1 - 18 years	4.80 - 5.52 mg/dL
Adults	4.6 - 5.4 mg/dL
Calcium Urine	
Adults	Up to 180 mg/24h (With calcium-restricted diet: 500 mg/24h)
	Up to 280 mg/24h (Without calcium-restricted diet)

Critical values

Calcium (serum or plasma): less than 6.0 mg/dL and higher to 14.0 mg/dL.

To convert the values to mg/dL to mmol/L (SI), multiply by 0.25. These values should be used as reference, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Comparison of Methods

The Calcium Arsenazo III kit was compared to other method for calcium dosage available in the market. 42 analysis were performed and the results were analyzed. The obtained linear equation was Y = 1.0126X - 0.5691, with linear correlation coefficient equal to 0.9934. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	6.28	11.89	9.33
Standard Deviation (mg/dL)	0.19	0.10	0.07
Variation Coefficient (%)	3.04	0.88	0.76

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	6.38	11.85	9.21
Standard Deviation (mg/dL)	0.15	0.05	0.10
Variation Coefficient (%)	2.34	0.40	1.11

SENSITIVITY

The sensitivity was calculated from 40 determinations a sample free of Calcium. The average found, which indicates the detection limit of the method, was 0.022 mg/dL with a standard deviation of 0.018 mg/dL. The sensitivity, which indicates the quantification limit of the method, corresponds the average plus 3 times the standard deviation and is equal to 0.074 mg/dL.

LINEARITY

The reaction is linear up to a concentration of 20 mg/dL. For higher values dilute the sample 1:2 or 1:4 with deionized or distilled water and repeat the determination. Multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Calcium plays an important role in the composition of the skeleton and reactions as essential in the clotting blood, in maintaining the integrity and permeability cell membranes, the stimulation of the skeletal and heart muscles, in neuromuscular conduction. The hypercalcemia can be seen in hyperparathyroidism, malignant tumor, Paget's disease (osteitis deformans), hypervitaminosis D, in cases of immobilization for fractures. Lower serum calcium can be observed in hypoparathyroidism, renal failure, deficiency in vitamin D (rickets), acute pancreatitis, massive blood transfusion.

BIBLIOGRAPHIC REFERENCES

1. CONNERTY, H. V.; BRIGS, A. R., Amer J. Clin. Pathol., 1.966, 45 - 290.
2. GITELLMAN, H. J., Anal. Biochem, 1.967, 18 - 521.
3. HENRY, R. J., Clinical Chemistry Principles and Technics.
4. JANSEN JW, Helbing AR Eur J Clin Chem Clin Biochem 1991;29:197- 201.
5. LEARY NO, Pembroke A, Duggan PF. Clin Chem 1992;38: 904 - 8.
6. MARTIN, E. W., Harzards of Medication, 1971.
7. MORGAN BR, Artiss JD, Zak B. Clin Chem 1993;39:1608-12.
8. POTTGEN P, Davis ER. Clin Chem 1976;22: 1752.
9. STEM, J., LEWIS V. H. P., Clin. Chim. Acta , 1957, 2 - 576.
10. TIETZ,N.W. Tietz Textbook of Clinical Chemistry, Burtis CA, Ashwood ER eds, 2a. edição, Philadelphia: W.B. Saunders Co., 1994.
11. TODD; SANFORD; DAVIDSOHN, Diagnósticos Clínicos, 1983.
12. TONKS, D. B., Clin. Chem., 1.963, 9 - 217.
13. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
14. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil


CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Calcium Arsenazo III kit: 10269360377

Review: March/2025

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

PRODUCT
STERILIZED

DANGER