CÁLCIO ARSENAZO III

REF **K305**

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CALCIUM ARSENAZO IIIREF **K305****INSTRUCTIONS FOR USE****VETERINARY USE**

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Quantitative determination of calcium in biological samples of serum, plasma and urine by color reaction. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION**Methodology:** Colorimetric Endpoint - Arsenazo III

The calcium reacts with arsenazo III in acid environment forming the blue-colored complex, whose intensity is proportional to the calcium concentration in the sample. The absorbance of the reaction product is measured in wavelengths between 600 and 680 nm.

REAGENTS

Number 1: Arsenazo - Store between 2 and 8 °C. Contains: Arsenazo III 0.2 mmol/L, 8-Hydroxyquinoline < 20 mmol/L, buffers, surfactants and preservative.

Number 2: Standard - Store between 2 and 8 °C. Contains: Calcium Carbonate 10 mg/dL, stabilizer and preservative.

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EQUIPMENTS AND OPERATIONAL INPUTS

Automated instruments and semi-automatic biochemistry, clock or stopwatch, pipettes, test tubes, a water bath at 37 °C, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8 °C. The transport at temperatures up to 30 °C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- The validity of the reactants is printed on the labels of same. However, once opened, the reagent can suffer from contamination during handling and have stability decreased. Important to note changes (Color change, precipitate formation) that may indicate contamination.

3- Reagents N° 1 and N°2 contains Sodium Azide which can be irritating to skin and mucous membranes. Handle with care.

4- To follow strictly the methodology proposed for obtaining accurate results.

5- The water used in material cleaning must be recent and free from contaminants.

6- Columns deionizer saturated releases alkaline water, several ions and oxidizing and reducing agents, which can significantly alter the results.

7- It is extremely important the thorough cleaning of all material (pipettes, tips, test tubes, cuvettes, etc.) to obtain accurate results. Using solution sulphochromic or aqueous solutions of HCl 50% or detergent nonionic. Rinse the material in running water and rinse with distilled/deionized water. Dry the material under heater protected from direct contact with the metal part of (put the material in glass vats).

8- Do not allow contact of the sample anticoagulant such as EDTA that is a chelat calcium, citrate that complex Ca²⁺, oxalate and fluoride that precipitate Calcium.

9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

10- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

11- Do not use the product in case of damaged packaging.

12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (Heparin): Free of hemolysis. Do not use Oxalate or EDTA as an anticoagulant by interfering the determination of calcium. Separate serum or plasma within one hour after collection, due to increased permeability of red blood cells to calcium. The analyte is stable in serum or plasma for 7 days at 2 to 8 °C and 2 months at -20 °C.

Urine: Make the urine collection within 24 hours in a recipient free of calcium. Prior to the collection add into container 10 mL of Nitric Acid 50% v/v. When the sample of 24 hours is not acidified during collection, add 20 mL of HCl 6M, homogenize, wait 60 minutes and take the aliquot for the test. The analyte is stable in urine for 4 days at 2 to 8 °C and 30 days at -20 °C¹³.

PROCESS DESCRIPTION

The calibration stability of the Calcium Arsenazo III kit installed on refrigerated equipment is at least 28 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Select the wavelength from 600 nm to 680 nm.

Mark 3 test tubes: B (Blank), P (Standard) and A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL
Reagent N° 2	---	10 µL	---
Sample	---	---	10 µL

Homogenize and incubate for 2 minutes at 37 °C.

Read the absorbance of Standard and Sample at 600nm to 680nm, matching the zero with Blank. The color is stable for 1 hour.

CALCULATIONS**Serum Calcium:**

$$\text{Calcium (mg/dL)} = \frac{\text{Test Absorbance}}{\text{Standard Absorbance}} \times 10$$

As the reaction follows the Beer-Lambert law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{10}{\text{Standard Absorbance}}$$

$$\text{Calcium (mg/dL)} = \text{Absorbance Test} \times \text{Factor}$$

Urinary Calcium:

$$\text{Urinary Calcium (mg/24h)} = \frac{\text{Urinary Calcium (mg/dL)} \times \text{Volume (mL)}}{100}$$

Ionized Calcium (Cal):

$$\text{Cal (mg/dL)} = \frac{6 \times \text{Ca} - [(0.19 \times \text{P}) + \text{A}]}{(0.19 \times \text{P}) + \text{A} + 6}$$

Ca = Calcium Serum (mg/dL)

P = Total Proteins (g/dL)

A = Albumin (g/dL)

INTERFERENTS

Do not use hemolyzed samples.

Lipemic and icteric samples may interfere with dosages. Oxalate or EDTA as an anticoagulant, as they interfere with calcium determination.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 915B.

REFERENCE VALUES

The reference values were determined in serum samples of healthy populations of dogs, cats, cattle and horses, males and females, adults without distinction of race.

Species	Reference Value
Dogs	8.23 - 13.40 mg/dL
Cats	8.04 - 12.95 mg/dL
Cattle	7.72 - 10.81 mg/dL
Horses	8.93 - 12.57 mg/dL

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Calcium Arsenazo III kit was compared to other method for calcium dosage available in the market. 42 analysis were performed and the results were analyzed. The obtained linear equation was $Y = 1.0126X - 0.5691$, with linear correlation coefficient equal to 0.9934. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	6.28	11.89	9.33
Standard Deviation (mg/dL)	0.19	0.10	0.07
Variation Coefficient (%)	3.04	0.88	0.76

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	6.38	11.85	9.21
Standard Deviation (mg/dL)	0.15	0.05	0.10
Variation Coefficient (%)	2.34	0.40	1.11

SENSITIVITY

The sensitivity was calculated from 40 determinations a sample free of Calcium. The average found, which indicates the detection limit of the method, was 0.022 mg/dL with a standard deviation of 0.018 mg/dL. The sensitivity, which indicates the quantification limit of the method, corresponds the average plus 3 times the standard deviation and is equal to 0.074 mg/dL.

LINEARITY

The reaction is linear up to a concentration of 20 mg/dL. For higher values dilute the sample 1:2 or 1:4 with deionized or distilled water and repeat the determination. Multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Calcium is the most abundant mineral in the body and plays a role in important bone functions as well as intracellular, muscular and neurological events. Calcium absorption into the body occurs in the intestine (ileum), mediated by vitamin D and indirectly by parathyroid hormone (PTH). In horses, intestinal calcium absorption is independent of vitamin D, since this species does not have the hydroxylase enzyme, essential in this process. Calcium absorption is directly affected by corticosteroids, stomach acidity, and calcium chelating agents (such as oxalates and phosphates). Calcium is filtered through the glomerulus and approximately 60% of the calcium is reabsorbed and therefore lesions in this tissue significantly affect this process. In horses, excess calcium is eliminated in the urine.

Hypercalcemia is not common in animals, but among domestic species, dogs and horses are the most affected. The main mechanisms leading to hypercalcemia include increased intestinal absorption, osteolysis, reduced renal excretion, and increased binding of calcium to proteins.

Hypocalcemia results in increased muscle and renal excitability and can occur due to: reduced calcium binding to proteins, reduced gastrointestinal absorption, excessive loss of calcium by the kidneys or hormonal imbalance (common in cattle in cases of milk fever).

BIBLIOGRAPHIC REFERENCES

- Connerty HV, Brigs AR. Amer J Clin Pathol. 1.966, 45 - 290.
- Gitellman HJ. Anal Biochem, 1.967, 18 - 521.
- Henry R.J., Clinical Chemistry Principles and Technics.
- Jansen JW, Helbing AR Eur J Clin Chem Clin Biochem 1991;29:197-201.
- Leary NO, Pembroke A, Duggan PF. Clin Chem 1992;38: 904 - 8.
- Martin EW. Harzards of Medication, 1971.
- Morgan BR, Artiss JD, Zak B. Clin Chem 1993;39:1608-12.
- Pottgen P, Davis ER. Clin Chem 1976;22:1752.
- Stem J, Lewis VHP. Clin. Chim. Acta, 1957, 2 - 576.
- Tietz NW. Tietz Textbook of Clinical Chemistry, Burtis CA, Ashwood ER eds, 2a. edição, Philadelphia: W.B. Saunders Co., 1994.
- Todd, Sanford, Davisohn. Diagnósticos Clínicos, 1983.
- Tonks DB. Clin Chem. 1.963, 9 - 217.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
- EClinPahth. Cornell University College of veterinary Medicine. 2013.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <n> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER