CINOMOSE Ag VET FAST

REF VET005

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FINALIDADE

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USE INSTRUCTIONS

FUNCTION

Rapid immunochromatographic test for qualitative determination of canine distemper antigen in ocular or nasal mucosa. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.
Distemper Ag VET FAST is an *in vitro* diagnostic test for rapid and qualitative detection of distemper virus antigens in ocular or nasal mucosa. The cassette consists of a nitrocellulose membrane, impregnated in the test line with antibody against canine distemper virus. The conjugate consists of antibody against distemper virus, conjugated with colloidal gold particles. When the sample is added to the cassette, it migrates by capillary action until it reaches the conjugate. In this region, the antigens present in the sample interact with the conjugate, forming the complex: antigen - conjugate. The sample continues to migrate through the immunochromatographic strip until reaching the test line, where the antigens also interact with the antibody against canine distemper virus immobilized in the test line, forming the complex: immobilized antibody - antigen - conjugate. This shows the formation of a pink-red line on the test line indicating positive result. Sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Number 1 (R1): Cassette - Store between 15 and 30°C. **Do not freeze.** It contains:

- A)** 01 plastic (device) cassette with 01 test strip containing:
- Control region (C): Antibody anti-mouse antibody immobilized.
 - Test region (T): Immobilized canine distemper virus monoclonal antibody.
 - Conjugate: Monoclonal mouse antibody against canine distemper virus, colloidal gold conjugate.

B) 01 sachet of silica

Number 2 (R2): Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Number3 (R3): Swab - Sterile flexible cotton swab.

Collection Accessory:
- Disposable plastic pipette.

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET005-1	10	10 x 0.5 mL
	VET005-2	25	25 x 0.5 mL
	VET005-3*	10	10 x 0.5 mL
	VET005-4*	25	25 x 0.5 mL
	VET005-5	5	5 x 0.5mL
	VET005-6	1	1 x 0.5 mL

* Unit presentation

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. They are in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.**

The cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

- Product for *in vitro* diagnosis, for veterinary use only.**
- Strictly follow the methodology proposed for obtaining results.
- Carry out the cleaning of the eyes and muzzle, avoiding the sampling of pre-existing secretions.
- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
- The sample must be stored in the kit diluent.
- Open the aluminum envelop to remove the cassette only at the time of use.
- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.
- The samples should reach room temperature before starting the test.
- Do not use reagents from different lots or kits.
- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
- The cassette (test device), the sterile flexible cotton swab and the plastic pipette are disposable and cannot be reused.
- The kit materials that had contact with the sample of the patient must have a special destination. The disposal of sterile flexible cotton swab, the plastic pipettes, the cassettes (test devices) and diluent vials should be disposed of in organic waste.
- Secondary packaging (box), Instruction for use, aluminized sachet and silica sachet can be disposed of in trash common.
- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES

Perform the test with ocular (conjunctive) or nasal secretion samples. Test the sample preferably immediately after collection. If necessary, samples can be stored between 2 to 8°C for up to 2 days or at -20°C for longer periods.


DESCRIPTION OF PROCESS
STABILITY AFTER OPEN

The stability of Distemper Ag VET FAST cassettes (test devices) after opening the aluminized sachet is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE

- Carry out the cleaning of the eyes and muzzle, avoiding the sampling of pre-existing secretions.
- Collect the eye discharge or nasal discharge samples using the sample collection swab.
Note: Use a gentle twistin motion during sample collection to maximize sampling.
- Insert the swab into the diluent tube, dissolving the sample.
Note: Use rotary movements of the swab inside the tube and rub it against the tube wall to maximize dissolution of the secretion in the diluent.
- The sample should be at room temperature before starting the test.
- Remove the cassette (test device) from the protective packaging and place it on a clean, dry and level surface
- Use the pipette that comes with the kit to collect the dissolved sample.

- Add 4 drops of sample (80 µL), slowly, into the sample application well of the cassette.
- Wait for the formation of the lines. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF RESULTS

Reactive Test: Formation of two lines: one in the control region (C) and one in the test region (T). Do not interpret after 20 minutes.	
Not Reactive Test: Line formation only in control region (C). Do not interpret after 20 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

TEST LIMITATIONS

This product is intended to detect canine distemper virus infections. A not reactive result can be observed if the amount of analyte present in the sample is below the detection limit of the test. This can happen, depending on disease period in which the sample was collected, the type of sample used in the test and its storage conditions. A false negative result may also occur in case of inefficiency in the sample extraction step in the diluent. Therefore, for the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 108 samples characterized by canine distemper were analyzed; 41 samples were diagnosed as positive and 67 were diagnosed as negative. The results show that the clinical sensitivity of the product is 97.56% and the clinical specificity is 97.01%.

METHOD		REFERENCE		Total
Cassette Distemper Ag	Result	Reactive	Not Reactive	
	Reactive	40	2	
	Not Reactive	1	65	
Total result		41	67	108

Clinical Sensitivity: (40/41): 97.56% (IC 95%: 87.14% – 99.94%)
Clinical Specificity: (65/67): 97.01% (IC 95%: 89.63% – 99.64%)
Precision: (105/108): 97.22%

ANALYTICAL SENSITIVITY

The product was tested with samples of tear secretion, which contained a determined virus load. The product was able to identify with positive result all the concentrations tested with a value equal to or greater than 10³ TCID₅₀/mL. Therefore, the analytical sensitivity (or detection limit) of 10³ TCID₅₀/mL was established.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Result	Found Result
Sample 1	10	Reactive	Reactive
Sample 2	10	Reactive	Reactive

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Result	Found Result
Sample 1	03	10	Reactive	Reactive
Sample 2	03	10	Reactive	Reactive

CLINICAL SIGNIFICANCE

Distemper is a highly contagious viral disease that affects dogs. The infection is frequent in puppies (aged between 3 and 6 months) due to the immaturity of the system the animal's immune system and failures in the vaccine protocol, but it can also occur in adult dogs (over 6 years). Infected dogs can transmit the virus from the time of infection and can persist in shedding the virus for up to 90 days, even when asymptomatic. The infection occurs through the oronasal route and through the ocular mucosa, where high viral loads are found. The incubation period for distemper varies from 1 to 4 weeks or more, and the most common clinical signs are: diarrhea and vomiting, thickening (hyperkeratosis) of cushions and nostrils, skin papules, dry keratoconjunctivitis, nasal discharge, coughing and dyspnea it may progress to pneumonia. Clinical neurological signs such as ataxia (loss of balance), incoordination, paralysis, seizures, and myoclonus may appear in the early or late stages of the disease. Antibody production progressively increases two weeks after infection, however, due to the rapid, progressive and often fatal course of the disease, detection of viral antigen is safer and more effective for diagnosing infection in its early stages. Furthermore, the detection of antibodies can generate false positive results, since the presence of antibodies can be the result of vaccination or maternal antibodies (when a pup). Thus, a diagnosis based on the detection of the viral antigen in ocular and nasal secretions, where high viral loads are found, is indicated.

BIBLIOGRAPHIC REFERENCES

- 1- QUIBASA: Data from the Research and Development Department.
- 2- Santos BM. Cinomose Canina – Revisão de literatura. Universidade castelo branco, Goiânia, agosto 2006.
- 3- Moraes FC. et al. Diagnóstico e controle da cinomose canina. PUBVET, Londrina, v. 7, n. 14, Ed. 237, Art. 1566, Julho, 2013.
- 4- Martella V, Elia G, Buonavoglia C. Canine Distemper Virus. Veterinary Clinics Small Animal Practice, v.32, p.787-797, 2008.
- 5- Krakoecka S. et al. Serologic response to canine distemper viral antigens in gnotobiotic dogs infected with canine Distemper virus. Journal of Infectious Diseases, v. 132, n. 4, p. 384-392, 1975.
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 11/11/2019 with the number 10.282/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: September/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER