# Bioclin

## CK MB UV

REF **K069**

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#### FINALIDADE

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### INSTRUCTIONS FOR USE


#### FUNCTION

Method for determination of the MB fraction Creatine Kinase (CK) in biological samples of serum or plasma (EDTA or Heparin). Kinetic test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Kinetic UV

The process involves the measurement of the activity of CK in the presence of an antibody against fraction M. This antibody completely inhibits the activity of CK-MM and M fraction of CK-MB, without, however, affect the activity of the B subunit of CK-MB and CK-BB. Based on the premise that the BB dimer is practically nonexistent in the peripheral blood, the residual enzyme activity found will be corresponding only to the fraction B of CK-MB. As the M fraction and fraction B have similar activities, we have the activity of CK-MB corresponds to the value found multiplied by two.

CK - Creatine Kinase

G-6-PDH - Glucose-6-Phosphate-Dehydrogenase

HK - Hexokinase

The speed of reduction of NADP<sup>+</sup> to NADPH is proportional to the activity of CK in the sample.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Imidazole Acetate (pH 6,7) 100 mmol/L, Glucose 20 mmol/L, EDTA 2 mmol/L, NADP<sup>+</sup>, 2 mmol/L, Hexokinase 3500 U/L, Magnesium Acetate 10 mmol/L, N-Acetylcysteine 20 mmol/L, Polyclonal Anti CK-M sufficient to inhibit up to 2,000 U/L CK-MM.

**Number 2 - Enzyme - Substrate** - Store between 2 and 8°C. Contains: Glucose-6-Phosphate-Dehydrogenase 2000 U/L, Creatine Phosphate 30 mmol/L, ADP 2 mmol/L, AMP 5 mmol/L, Diadenosine Pentaphosphate 10 mmol/L.

#### PRESENTATION

Presentation	Reagent N°1	Reagent N°2
1	1 x 40 mL	1 x 10 mL
2	2 x 40 mL	2 x 10 mL
3	4 x 40 mL	4 x 10 mL
4	2 x 40 mL	1 x 20 mL
5	4 x 40 mL	2 x 20 mL
6	1 x 20 mL	1 x 5 mL
7	2 x 20 mL	2 x 5 mL
8	3 x 20 mL	3 x 5 mL
9	4 x 20 mL	4 x 5 mL
10	5 x 20 mL	5 x 5 mL
11	6 x 20 mL	6 x 5 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatch, test tubes, Biocal Bioclin, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.  
**3-** Water used in material cleaning must be recent and free of contaminants.

**4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.

**5-** It is important for the proper performance of the test, a strict control of time, temperature and pH.

**6-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**7-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**8-** Do not use the product in case of damaged packaging.

**9-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum or plasma collected with Heparin or EDTA. The sample is stable for 2 days between 15 and 30°C, 7 days between 2 and 8°C and 1 year at -20°C <sup>3</sup>. Do not use hemolyzed sample.

#### PROCESS DESCRIPTION

The calibration stability of the CK MB UV kit installed on refrigerated equipment is at least 13 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### WORKING REAGENT PREPARATION

**The reagents are ready to use in bireagent format and on automated equipment that uses dedicated presentations.**

For use in monoreagent format, prepare the working reagent as follows:

Mix 4 parts of reagent N°1 with 1 part of reagent N°2. Stable for 1 day at temperatures between 15 to 30°C and 15 days between 2 to 8°C.

#### REACTION CONDITIONS

It is indispensable the use of a thermostated cuvette at 37°C, optical path of 1cm and reading at 340 nm.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

#### TECHNIQUE

Bioclin recommends using the kit Biocal as calibrator and as control serum, Biocontrol N and P Bioclin kits.

#### Monoreagent:

Working Reagent	1.0 mL
Sample	40 µL

Homogenize and transfer to a thermostated cuvette at 37°C and wait for 5 minutes. Make the initial reading, triggering the stopwatch. Repeat the readings at 1, 2 and 3 minutes.

To calibrate the bireagent kit, Bioclin recommends using the BIOCAL Bioclin calibrator serum. We also recommend using the Bioclin Biocontrol N and Biocontrol P kits as control serums.

**Check the equipment's programming by calling Customer Service at 0800 031 5454.**

#### CALCULATIONS

Calculate the mean difference in absorbance per minute ( $\Delta A/\text{min}$ .) and use this value to calculate the result.

CK-MB (U/L) 340 nm =  $\Delta A/\text{min} \times 8254$

#### TECHNICAL FOR CK-MB GREATER THAN 20% OF TOTAL CK VALUE

If the measurement of activity of CK-MB exceeds 20% of total CK activity, the presence of macro CK should be suspected and the following technique to determine the value of CK-MB real should be performed.

Example:

CK Total = 150 U/L

20% CK Total = 30 U/L

If CK-MB > 30 U/L, performing technique below.

- Make the CK-NAC (total) and CK-MB before heating the serum.
- Heating the serum at 45°C for 20 minutes.
- Test the CK-MB again.

Value of CK-MB (U/L) = initial CK-MB - CK-MB\*

CK-MB\* = value found after heating the serum

Note: See Procedure Limitations

#### PROCEDURE LIMITATIONS

The method quantifies all CK-BB isoenzyme in serum. Normally the activity of this isoenzyme is small. However, if a significant increase of this isoenzyme occurs, the activity of CK-MB is overrated. Moreover, macromolecular two isoenzymes have been considered as false positive for CK-MB: a macro CK type 1 and type 2 macro CK. The type 1 is a complex of CK-BB or CK-MM bound to IgA or IgG. The macro CK type 2 is an oligomeric complex of mitochondrial origin and is associated with neoplasias. If the measurement of activity of CK-MB exceeds 20% of total CK activity, one should suspect the presence of these isoenzymes macromolecules.

The readings should be taken in a spectrophotometer that satisfies the following conditions:

Wavelength 340 nm  
Semi trajectory of the pass band 10 nm  
Spurious light lower than 0,5%  
10 mm Thermostated cuvette

#### INTERFERENT

No interference was observed by Glucose until 7 g/L, Hemoglobin until 6 g/L, Triglycerides until 8 mmol/L and Sulfasalazine until 2 mg/dL.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The calibration kit can be made using the theoretical calibration factor based on the molar absorptivity of the NADP, or through the calibrator BIOCAL. The Bioclin recommends using BIOCAL calibrator which is traceable to the reference material ERM-AD455 / IFCC and the IFCC reference method.

#### REFERENCE VALUES

The reference values for this method were obtained by determining CK-MB in healthy populations of male and female.

	37°C
Men/Women	0 - 25 U/L

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE

##### ACCURACY

###### Comparison of Methods

The CK MB UV Kit was compared with other methods for CK-MB dosage available in the market. 42 analysis were performed and the results were evaluated. The linear equation obtained was  $Y = 0.996X - 0.3738$ , with linear correlation coefficient equal to 0.9849. With these results we can conclude that the kit shows good methodological specificity.

##### PRECISION

###### Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	20.53	45.45	33.20
Standard Deviation (U/L)	0.51	0.50	0.76
Coefficient of Variation (%)	2.46	1.11	2.28

###### Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	20.69	43.33	31.98
Standard Deviation (U/L)	1.71	3.26	2.13
Coefficient of Variation (%)	8.25	7.52	6.66

#### SENSITIVITY

Sensitivity was calculated from 40 determinations of a CK-MB-free sample on automated equipment. The mean value found, that indicates the method detection limit, was 0.025 U/L with a standard deviation of 0.158 U/L. The sensitivity, which indicates the quantification limit of the method, corresponds to the mean plus 3 times the standard deviation, and is equal to 0.499 U/L.

Using semi-automatic equipment with a photometric resolution of 0.001 as the minimum detectable absorbance, the detection limit of the kit goes to 8.25 U/L. In this case, the dosages may present inaccurate results, which may be detrimental mainly to analyzes of low results.

#### LINEARITY

The method is linear up to 600 U/L. If the variation of molar extinction per minute is larger than 0.150, repeat the test using a dilution of 1:10 with saline. The result obtained is multiplied by 10.

#### DIAGNOSTIC SIGNIFICANCE

Creatine kinase (CK) is an enzyme present in several tissues of the human body and in several cell types. It is found in high concentrations in the skeletal muscle, heart muscle, brain and gastrointestinal tract, and is responsible for catalyzing the conversion of creatine to phosphocreatine, an important energy reservoir for cells.

Its molecular structure consists of two subunits, the M subunit (Muscular type) and the B subunit (Cerebral type), which determine the formation of 3 isoenzymes. Thus, CK is found in dimeric form as CK-MM (muscle), CK-BB (cerebral) and CK-MB (myocardial).

Because of their presence in these organs and tissues, any cell injury may promote increased serum CK levels.

Thus, elevated levels of CK are observed mainly after cardiac muscle trauma and in skeletal muscle diseases. The CK dosage is especially used in conjunction with CK-MB for diagnosis and monitoring of myocardial infarction.

##### Creatine Kinase Isoenzymes (CK)

**CK-MM:** Corresponds to more than 95% of the enzyme present in skeletal muscle and about 70 to 75% of the myocardial enzyme. Elevated values of this isoenzyme are found in hypothyroidism, muscular dystrophy, hypokalemia and recent surgery.

**CK-BB:** Predominant in the brain and lung, and practically non-existent in the blood.

**CK-MB:** Found predominantly in cardiac muscle. This enzyme undergoes elevations in myocardial diseases, such as myocardial infarction and myocarditis. Following the infarction (3 to 5 hours), the serum concentration of this enzyme rises, reaches a peak after 12-24 hours and remains elevated for a short period (2 to 3 days). As it is a specific fraction for the cardiac muscle, this isoenzyme is of great importance for the monitoring of myocardial infarction.

#### BIBLIOGRAPHIC REFERENCES

- 1 - SCANDINAVIAN Socit For Clinical Chemistry Scand. J. Clin. Lab. INVEST., 1974, 33, 291.
- 2 - Ann. Biol. Clin., 1982, 40,99.
- 3 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 4 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

#### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for CK MB UV kit: 10269360073

Review: October/2025

#### UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE  
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT  
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT  
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC  
PRODUCT

TOXIC

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED

DO NOT REUSE

PRODUCT  
STERILIZED

CAUTION

DANGER