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CK NAC UVREF **K010****INSTRUCTIONS FOR USE****VETERINARY USE**

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST,
CONTACT CUSTOMER ADVISORY SERVICE:



SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of the creatine kinase (CK) in biological samples of serum or plasma (EDTA or Heparin). Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION**Methodology:** Kinetic UV

Determination of creatine kinase according to the reactions:

CK - Creatine Kinase

G-6-PDH – Glucose-6-Phosphate-Dehydrogenase

HK - Hexokinase

The reduction speed of NADP⁺ to NADPH is proportional to the activity of the CK in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8 °C. Contains: Imidazole Buffer < 300 mmol/L, Glucose < 100 mmol/L, NADP < 10 mmol/L, Hexokinase < 50 KU/L, N-Acetylcysteine < 100 mmol/L, stabilizers, surfactant and preservative.

Number 2 - Enzyme Substrate - Store between 2 and 8 °C. Contains: Buffer < 300 mmol/L, Glucose-6-Phosphate-Dehydrogenase < 50 KU/L, Creatine Phosphate < 300 mmol/L, ADP < 30 mmol/L, stabilizers, surfactant and preservative.

PRESENTATION

Presentation	Reagent N°1	Reagent N°2
K010-1-VET	1 x 40 mL	1 x 10 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, test tubes, watch or stopwatch, Biocal Bioclin, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8 °C. The transport at temperature up to 30 °C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- It is important for the proper performance of the test, a strict control of time, temperature and pH.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma collected with Heparin or EDTA. The sample is stable for 4 hours between 15 and 30 °C³, 7 days between 2 and 8 °C and 30 days at -20 °C³. Do not use hemolyzed sample or lipemic.

PROCESS DESCRIPTION

The calibration stability of the CK NAC UV kit installed on refrigerated equipment is at least 11 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix 4 parts of reagent N° 1 with 1 part of reagent N° 2. Stable for 2 days at temperatures between 15 to 30 °C and 15 days if kept between 2 to 8 °C. Do not use hemolyzed sample.

REACTION CONDITIONS

It is indispensable the use of a thermostated cuvette at 37°C, light path of 1cm and reading at 340 nm.

TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits.

To do this, mix the reagents according to the information in the table below.

Working Reagent	1.0 mL
Sample	20 µL

Mix and transfer to a thermostated cuvette at 37 °C and wait 2 minutes. Make the initial reading, triggering the timer. Repeat the readings at 1, 2 and 3 minutes.

CALCULATIONS

Calculate the mean difference in absorbance per minute ($\Delta A/\text{min.}$) and use this value to calculate the result.

CK (U/L) 340 nm = $\Delta A/\text{min} \times 8095$

PROCEDURE LIMITATIONS

The readings should be taken in a spectrophotometer that satisfies the following conditions:

Wavelength 340 nm

Semi trajectory of the pass band 10 nm

Spurious light lower than 0.5%

10 mm Thermostated cuvette.

INTERFERENT

No interference was observed by Glucose until 7 g/L, Hemoglobin until 5 g/L and Triglycerides until 7 mmol/L.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the NADP, or through the BIOCAL calibrator which is traceable to the reference material ERM-AD455 / IFCC and the IFCC reference method.

REFERENCE VALUES

The reference values were determined in serum samples of healthy dogs, cats, cattle and horses, males and females, adults and with no predilection for breed.

Specie	Reference Value
Dog	49 - 360.3 U/L
cat	113 - 386 U/L
Cattle	28.67 - 153.1 U/L
Horse	124.3 - 555.3 U/L

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The CK NAC UV kit was compared to other method for CK NAC dosage available in the market. 42 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 1.0045x + 2.5661$, with linear correlation coefficient equal to 0.9994. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	72.55	412.28	235.85
Standard Deviation (U/L)	0.88	2.89	1.64
Variation Coefficient (%)	1.21	0.70	0.70

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	72.65	411.85	234.68
Standard Deviation (U/L)	0.98	1.61	5.24
Variation Coefficient (%)	1.35	0.39	2.23

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of the analyte (CK NAC). The mean concentration of CK NAC UV was 0.225 U/L with standard deviation of 0.480 U/L. The sensitivity, which indicates the limit of detection, which corresponds to the average plus three times the standard deviation is equal to 1.664 U/L.

LINEARITY

The method is linear to 2000 U/L. If the variation of molar extinction per minute is larger than 0,250, repeat the test using a dilution of 1:10 with physiologic serum. The result obtained is multiplied by 10.

DIAGNOSTIC SIGNIFICANCE

The enzyme Creatine Phosphokinase (CPK) or Creatine Kinase (CK) is present in high concentrations in the cytoplasm

of skeletal and cardiac muscle cells and in the brain, and its main function is to make ATP available for cellular metabolism. The enzyme reaches the blood stream due to changes in permeability in the cell membrane, which lead to cellular leakage of CK from the cytoplasm to the extracellular medium.

CK is the most used enzyme as an indicator of injury in muscle cells. Some of the main causes of increased CK include muscular dystrophies, rhabdomyolysis, myositis, Vitamin E and Selenium deficiency, and inoculation of proteolytic snake venoms.

The increase in serum CK can be observed in some specific cases:

- Puppies have higher CK activity when compared to adults;
- Hemolysis causes release of intracellular enzymes and membrane components that cause an artefactual increase in CK
- Physical activity in horses can rapidly increase CK, but its reduction to normal levels occurs in a few hours;
- Intramuscular drug applications;
- Post-surgical conditions
- Boarding for transporting large animals.

In dogs, neurological disorders can lead to an increase in serum CK due to the release of this enzyme by neurons in inflammatory and degenerative processes.

Serum CK reduction is uncommon, but may be related to significant loss of muscle mass.

BIBLIOGRAPHIC REFERENCES

- 1 - Scandivavian Society For Clinical Chemistry Scand. J. Clin. Lab. INVEST., 1974, 33, 291.
- 2 - Ann. Biol. Clin., 1982, 40,99.
- 3 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 4 - Gonzalez FHD, Silva SC. Introdução à bioquímica clínica veterinária. 2ª ed. Porto Alegre: Editora da Universidade Federal do Rio Grande do Sul. c. 8, p. 318-337, 2006.
- 5 - Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals. San Diego: Academic Express, 6. ed., 1998. 936 p.
- 6 - EClinPahth. Cornell University College of veterinary Medicine. 2013.
- 7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 8 - Paltrinieri S, Pintore L, Balducci F, Giordano A, Costabole A, Bernardini M. Serum creatine kinase isoenzymes in dogs with different neurologic disease. Veterinary Clinical Pathology, v.46, n.1, p.91-99, 2017.
- 9 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
 Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOLOGY