CLORETOS MONORREAGENTE

REF K221

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FINALIDADE

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
INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Chloride in biological samples serum or plasma (Heparin), urine, sweat or liquor. Colorimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Colorimetric

In the presence of Chloride ions, the Mercury Thiocyanate, in an acid environment, forms the Mercuric Chloride and Thiocyanate ions. The Thiocyanate ions react with the Iron ions forming a Iron Thiocyanate, of yellow-orange color, who is proportional to the Chloride concentration in the sample.

REAGENTS

Number 1 - Color Reagent - Store between 2 and 30°C. Contains: Mercury Thiocyanate 2,00 mmol/L, Mercury Chloride < 2,00 mmol/L, Iron Nitrate > 15mmol/L, Nitric Acid solubilizer > 0.1% and stabilizers. **Toxic reagent.**
Number 2 - Standard - Store between 2 and 30°C Contains: Chloride 100 mEq/L and preservatives.

PRESENTATION

REAGENT 1	REAGENT 2
1 x 20 mL	1 x 3 mL
2 x 20 mL	1 x 3 mL
3 x 20 mL	1 x 3 mL
4 x 20 mL	1 x 3 mL
1 x 40 mL	1 x 3 mL
2 x 40 mL	1 x 3 mL
3 x 40 mL	1 x 3 mL
4 x 40 mL	1 x 3 mL
1 x 50 mL	1 x 3 mL
2 x 50 mL	1 x 3 mL
3 x 50 mL	1 x 3 mL
4 x 50 mL	1 x 3 mL
1 x 100 mL	1 x 3 mL
2 x 100 mL	1 x 3 mL
3 x 100 mL	1 x 3 mL
4 x 100 mL	1 x 3 mL
1 x 20 mL	-
2 x 20 mL	-
3 x 20 mL	-
4 x 20 mL	-
1 x 40 mL	-
2 x 40 mL	-
3 x 40 mL	-
4 x 40 mL	-
1 x 50 mL	-

2 x 50 mL	-
3 x 50 mL	-
4 x 50 mL	-
1 x 100 mL	-
2 x 100 mL	-
3 x 100 mL	-
4 x 100 mL	-

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch and stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. Can be found at market specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature shall be between 2 and 30°C. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- **For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain exact results
- 3- Water used in material cleaning must be recent and free of contaminant
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significant.
- 5- The serum or plasma shall be prepared in maximum one hour after collection. A superior time produce CO₂ loss who interferes with the chloride distribution between cells and plasma.
- 6- Urine and liquor must be centrifuge
- 7- Use heparinized plasma, because the regular anticoagulants interfere in the determination.
- 8- Handle reagent N°1 carefully - toxic reagent.
- 9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 10- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 11- Do not use the product in case of damaged packaging
- 12- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma (heparinized), sweat or liquor samples. The analyte is stable for 7 days at 2 to 8°C and 1 year at -20°C. ⁶ Hemolyzed and lipemic samples should be avoided as they may cause a false reagent result.
Urine: stable for 3 days at 2 to 8°C. ⁶ **Do not freeze.**

PROCESS DESCRIPTION

The calibration stability of the Chloride Monoreactive kit installed on refrigerated equipment is at least 31 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control sera.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Biocontrol Bioclin Kits.

Mark 3 test tubes with letters B (Blank), S (Sample), T (Standard) and proceed as follows:

	Blank	Standard	Sample
Reagent N°1	1 mL	1 mL	1 mL
Distilled water	0.01 mL	----	----
Sample	----	----	0.01 mL
Reagente N°2	----	0.01 mL	----

For urine samples, dilute 1: 2 with distilled or deionized water (0.1mL urine + 0.1mL water). Multiply the result by 2.
Homogenize and incubate at room temperature for 2 minutes. Determine the absorbance of the Sample and Standard at 450 nm (450 - 505 nm) by zeroing with white. The color is stable for 2 hours.

CALCULATIONS

Chloride (mEq/L) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times 100$

As the reaction follows the Lambert Beer Law, it can use the Calibration Factor:

Calibration Factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$ (100 mEq/L)

Chloride (mEq/L) = Sample Absorbance x Calibration Factor.
Urine (mEq/24 hours) = mEq/L x Volume (liters).
The results are expressed in mEq/L or mmol/L, for the samples (serum, plasma, liquor and sweat) and in mEq/24 hours or mmol/24 hours for urine.

INTERFERENT

No interference was observed by Bilirubin until 40 mg/dL, Hemoglobin until 180 mg/dL and Triglycerides until 1000 mg/dL.

PROCEDURE LIMITATIONS

The use of automatic pipette appropriately calibrated for pipetting the Standard and Sample, in order to minimize problems with inaccurate volume.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 919b.

REFERENCE VALUES

The reference values in mEq/L for this method were obtained by determining Chlorides in healthy male and female populations.
Serum or plasma (All ages): 98 to 112 mEq/L
Urine: 110 to 250 mEq/24 hours
CSF: 118 to 132 mEq/L
Sweat: 5 to 35 mEq/L

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Chloride Monoreactive kit was compared with another commercially available Chloride dosage method. 42 samples were performed and the results were evaluated. The linear equation obtained was $Y = 1.0006X + 1.0085$ with a correlation coefficient of 0.9998. With these results it can be concluded that the kit has good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mEq/L)	82.97	96.61	129.41
Standard Deviation (mEq/L)	0.23	0.38	0.43
Coefficient of Variation (%)	0.27	0.39	0.34

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mEq/L)	82.69	96.61	128.63
Standard Deviation (mEq/L)	1.17	1.35	0.99
Coefficient of Variation (%)	1.42	1.40	0.77

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Chlorides. The average found was 1.174 mEq/L, with standard deviation of 0.124 mEq/L. The sensitivity, which indicates the indicative of methods detection limit, corresponds the average plus 3 times the standard deviation and is equal to 1.546 mEq/L.

LINEARITY

The reaction is linear between 140 mEq/L. For larger values, dilute the serum with distilled or deionized water, repeat the determination and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Chloride is the major of the extra cellular liquid (LEC) anion, being important regulating the water maintenance and distribution, in a anionic and cationic balance of LEC and in osmotic pressure. Its concentration in blood is regulated by the kidneys, suprarenal gland, lungs, skin, intestine gastric tract and blood pH. In several clinical situations, can occur reduction (hypochloremia) or increase (hyperchloremia) of the plasmatic concentration. Hypochloremia - In nephritis associated to pyelonephritis, and patients in Addisonian crisis, the values of Chloride are reduced. Extended vomits or persistent gastric secretion can lead to chlorine loss and low plasmatic values. It can also occur hypochloremia in metabolic alkalosis, bromate intoxication, cranial injuries, kidney absorption malfunction. Hyperchloremia - Pathological causes can be quoted: acute kidney insufficiency, metabolic acidosis, nephrotic syndrome. In insipidus diabetes, in adrenocortical hyperfunction and hyperparathyroidism. In cystic fibrosis, systemic diseases that lead to intestinal obstruction, bile cirrhosis and lung abnormalities, high values of Chloride in sweat are important indicatives of diagnosis

BIBLIOGRAPHIC REFERENCES

- 1 - SCHALES O.; SCHALES S., J. Biol. Chem., 1941, 140, 879.
- 2 - ZALL, D. M.; FISHER, D.; GARNER, M. Q., Anal. Chem., 1.956, 28,1665.
- 3 - BENNINGTON, L. James, Dict & Encycl. of Lab. Med. and Tech., 1.984.
- 4 - BURTIS, A. Carl; Ashwood, R. Edward, Clin., Chem., Tietz Text book of; 2nd. ed., 1.986, 1366.
- 5 - BURTIS, Carl A. Carl; BRUNS, David E. (Org.). TIETZ Fundamentos de Química Clínica e Diagnóstico Molecular. 7. ed. Rio de Janeiro: Elsevier Editora Ltda., 2016. 1078 p.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 7- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made In Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for the Chloride Monoreactive kit: 10269360329

Review: October/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER