COLESTEROL HDL ENZIMÁTICO

REF K015

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FINALIDADE

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1 - ALLAIN, C. C., and al. Clin. Chem., 1974, 20-470.

2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.

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6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31

7- Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA 2001: 285:2486-97.

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of HDL Cholesterol in biological samples serum or plasma (Heparin or EDTA). Colorimetric enzymatic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic - Trinder

The VLDL Lipoproteins (Very Low Density Lipoprotein), and LDL (Low Density Lipoprotein) and the Kilomicrons are precipitate with the mixture of Phosphotungsten Acid and Magnesium Chloride. After centrifugation, the Cholesterol binded to the High Density Lipoprotein (HDL) is determined by enzymatic colorimetric method.

REAGENTS

Number 1 - Standard - Store between 2 and 8°C. Contains: Sodium Chloride < 250 mmol/L, Cholesterol 40 mg/dL, thickener, solubilizer and emulsifier.

Number 3 - Precipitate Reagent - Store between 2 and 8°C. Contains: Phosphotungstic Acid < 3 mmol/L and Magnesium Chloride Hexahydrate < 110 mmol/L.

PRESENTATION

Reagent	Volume
Nº 1	2 mL
Nº 2	25 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Moreagent Cholesterol kit, Biocontrol N and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Water level in water bath should be above the level of reagents in test tubes.

6- The use of Enzymatic Reagent of **Bioclin** Cholesterol Monoreagent kit is of fundamental importance to a good performance of the system.

7- Determine the factor periodically in every batch of the product.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free from hemolysis or plasma, collected with Heparin or EDTA. After collection, the sample should be centrifuged in up to 3 hours.

HDL cholesterol is stable in serum for up to 7 days at 2 to 8°C ⁶ and for 30 days at -20°C.

Samples stored at temperatures between 15 and 30°C have stability of at most 14 hours.

PROCESS DESCRIPTION

TECHNIQUE

Bioclin recommends, use the Monoreagent Cholesterol (Enzymatic Reagent*) kit and as control serum, Biocontrol N and Biocontrol P Bioclin Kits.

Precipitation of VLDL, LDL and Kilomicrons

In a mixer tube pipette:

Sample..... 250 µL
Precipitant Reagent N° 2 250 µL
Agitate in Vortex or manually during 1 minute. Centrifuge in 3500 rpm during 15 minutes. This step is of fundamental importance to obtain the correct result. The limpid supernatant contains HDL Cholesterol and shall be pipette immediately after centrifugation, to avoid wrongly elevated results.

Determination of HDL Cholesterol

Mark 3 test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Supernatant	---	---	50 µL
Reagent Nº 1	---	50 µL	---
Enzymatic Reagent*	1.0 mL	1.0 mL	1.0 mL

Homogenize well and put it in water bath at 37 °C for 5 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), matching the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

Since it is require for the technique the diluted serum at 1:2 with the Precipitate Reagent, it is necessary to correct this dilution to obtain the correct results. For that reason, during the calculations, the standard value is corrected to the double of its real concentration, which means, 80 mg/dL.

$$\text{HDL (mg/dL)} = \frac{\text{Sample Absorbance} \times 80}{\text{Standard Absorbance}}$$

Reaction follows the Lambert-Beer Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{(80 \text{ mg/dL})}{\text{Standard Absorbance}}$$

$$\text{HDL (mg/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

The VLDL and LDL Cholesterol concentration can be calculated according to the Friedewald equation. This calculation could be used for samples which Triglycerides concentration does not pass 400 mg/dL.

$$\text{VLDL Cholesterol} = \text{Triglycerides} / 5$$

$$\text{LDL Cholesterol} = \text{Cholesterol total} - (\text{HDL} + \text{VLDL})$$

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

For lipemic sample or cloudy supernatant, dilute the sample 1:2 with Sodium Chloride 0.85% and repeat the precipitation. Multiply the result obtained by 2. If the supernatant remains cloudy, the sample can not be used.

INTERFERENCES

Some substances, such as Ascorbic Acid, Hemoglobin above 150 mg/dL, Bilirubin above 20 mg/dL and Triglycerides above 600 mg/dL, interfere with the reaction.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

REFERENCE VALUES ⁷

The reference values in mg/dL, for this method, were obtained through the determination of HDL Cholesterol in healthy populations of male and female.

Reference values below are for individuals with 12 hours fasting and follow Adult Treatment Panel III (ATP III).

Adults (≥ 20 years)	Desirable	> 60 mg/dL
	Limit	40 - 60 mg/dL
	Low	< 40 mg/dL

Childrens and Adolescents (2 to 19 years)	Desirable	> 45 mg/dL
---	-----------	------------

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Comparison of Methods

The HDL Cholesterol Enzymatic kit was compared with another commercially available kit of measurement of HDL Cholesterol. 07 analyses were conducted and the results were eva-luated, obtaining the linear equation of Y = 1.029X - 2.114 with a linear correlation coefficient of 0.999. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	35.15	50.70	70.75
Standard Deviation (mg/dL)	0.67	0.73	0.72
Coefficient of Variation (%)	1.91	1.45	1.01

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Concentration (mg/dL)	34.90	50.72	70.77
Standard Deviation (mg/dL)	0.22	0.03	0.03
Coefficient of Variation (%)	0.62	0.06	0.04

SENSITIVITY

The sensitivity was calculated from 20 determnations of a sample free of HDL Cholesterol. The average found was 1.3307 mg/dL, with standard deviation of 0.0005 mg/dL. The sensitivity, that indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 1.3321 mg/dL.

DIAGNOSTIC SIGNIFICANCE

The determination of High Density Lipoproteins (HDL) and Low Density Lipoprotein (LDL) play an important role in prevention and therapy of coronary heart disease. The HDL are capable of captating the Cholesterol of the pheriferical tissue and reconducting to the liver in order to degradate and excrement it. This mechanism denominated reverse transport of Cholesterol, constitute the theory of HDL anti-atherogenicity, where low concetration of HDL are directly related with the incidence of cardiovascular diseases.

BIBLIOGRAPHIC REFERENCES

- 1 - ALLAIN, C. C., and al. Clin. Chem., 1974, 20-470.
2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
3 - TRINDER, P., Ann. Clin. Biochem, 1969, 6-24.
4 - CARL, A. B. And EDWARD, R. A.; Tietz Texbook of Clinical Chem. 2 nd. ed.; 1994, 1002-1081. 5- Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
5- Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
7- Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA 2001: 285:2486-97.
8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for HDL Cholesterol Enzymatic kit:
10269360082

Review: October/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER