COLESTEROL MONORREAGENTE

REF
K306

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FINALIDADE

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PELIGRO

CHOLESTEROL MONOREAGENT

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

REF K306

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br


INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Chlolesterol in biological samples of serum or plasma (Heparin). Enzymatic colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric - COD - PAP (Trinder)
The enzymatic determination of Cholesterol is made accordingly with the following reactions:

The intensity of formed cherry color is directly proportional to the concentration of Cholesterol in the sample.

REAGENTS

Number 1 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Buffer < 150 mmol/L, Phenol < 9 mmol/L, 4-Aminoantipyrine < 3 mmol/L, Cholesterol Oxidase < 5 KU, Lipoprotein Lipase < 5 KU, Peroxidase < 5 KU, chelator, surfactants, stabilizers and preservative.

Number 2 - Standard* - Store between 2 and 8°C. Contains: Cholesterol 200 mg/dL, solubilizers, surfactant and stabilizer.

PRESENTATION

PRESENTATION	REAGENT Nº 1	REAGENT Nº 2*
1	100 mL	3 mL
2	2 x 100 mL	3 mL
3	4 x 100 mL	3 mL
4	5 x 20 mL	3 mL
5	5 x 40 mL	3 mL
6	10 x 40 mL	3 mL
7	2 x 60 mL	3 mL
8	4 x 60 mL	3 mL
9	6 x 60 mL	3 mL
10	8 x 60 mL	3 mL
11	10 x 60 mL	3 mL

*The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath (37°C), watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The

transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1-** For *in vitro* diagnostic use only.
2- Strictly follow the methodology proposed to obtain exact results.
3- Water used in material cleaning must to be recent and free of contaminants.
4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
5- The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.
6- Water level in water bath should be above the level of reagents in test tubes.
7- Plasma collected with anticoagulants such as oxalate, EDTA or citrate, generate wrongly lowered results.
8- Handle reagents with care. Reagent Nº 1 contains Sodium Azide, irritant to skin and mucous.
9- Slightly pink color development in the Reagent Nº 1 does not damage their quality, and stability as long as the corresponding Blank and Standard dosages are used.
10- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
11- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
12- Do not use the product in case of damaged packaging.
13- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolyses or plasma collected with Heparin. Cholesterol is stable in plasma or serum for up to 7 days between 2 and 8°C and 3 months at -20°C. ⁷

PROCESS DESCRIPTION

The calibration stability of the Cholesterol Monoreagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark three test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Reagent Nº2	--	10 µL	--
Reagent Nº 1	1.0 mL	1.0 mL	1.0 mL

Homogenize well and put it in water bath at 37°C for 10 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm),

matching the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

Cholesterol (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}}$ x 200

Since reaction follows the Lambert-Beer Law, the Calibration Factor can be used.

Calibration Factor = $\frac{\text{Standard Concentration (200 mg/dL)}}{\text{Standard Absorbance}}$

Cholesterol (mg/dL) = Sample Absorbance x Calibration Factor

The results shall be express as mg/dL.

INTERFERENCES

The equilibrium of the assay is affected by some interference substances such as the Ascorbic Acid (even in small concentrations), Hemoglobin higher than 150 mg/dL and Bilirubin higher than 20 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 911C.

REFERENCE VALUES ⁸

The reference values below are for subjects without fasting and fasting for 12 hours.

	Fasting	Without fasting
Adults (> 20 years)	<190 mg/dL	<190 mg/dL
Childrens and adolescents (2 to 19 years)	<170 mg/dL	<170 mg/dL

Total Cholesterol values ≥ 310 mg/dL (adults) or ≥ 230 mg/dL (children and adolescents) may be indicative of familial hypercholesterolemia if secondary dyslipidemias are ruled out.

To convert the values from mg/dL into mmol/L (SI) multiply by 0.026.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical

professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Recovery

The recovery analysis was mad through 5 determinations of samples. Accuracies were calculated and were found in good concordance to the reference values, obtaining a recuperation between 95% and 105%.

Comparison of Methods

The Cholesterol Monoreagent kit was compared with another method commercially available to measure the dosage of Cholesterol. 42 analyzes were performed and results were evaluated. The linear equation obtained was $Y = 1.0118X - 0.5703$ with correlation coefficient 0.9932. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	102.75	207.40	386.73
Standard Deviation (mg/dL)	2.82	1.08	3.10
Correlation Coefficient (%)	2.74	0.52	0.80

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	100.48	203.86	385.53
Standard Deviation (mg/dL)	2.23	5.60	2.26
Correlation Coefficient (%)	2.22	2.75	0.59

SENSITIVITY

Sensitivity was calculated through 40 determinations of a sample free of cholesterol. The average found, which indicates the limit of detection, was -1.99 mg/dL with a standard deviation of 0.89 mg/dL, the sensitivity, which indicates the limit of quantification from the method, corresponds the average plus 3 times the standard deviation, and is equal to 0.67 mg/dL.

LINEARITY

Reaction is linear up to 500 mg/dL. For samples with higher values than 500 mg/dL, dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The concentration of plasmatic Cholesterol is influenced by hereditary characters, endocrine function, nutrition and integrity of vital organs such as liver and kidneys. Numerous investigations confirm the relationship between Total Cholesterol and the

development of atherosclerotic coronary disease.

The Cholesterol is increased in diabetes, nephrotic syndrome, cirrhosis bile in hypothyroidism and in hyperlipoproteinemia type IIa, IIb and III. Research shows that high levels of LDL Cholesterol (Cholesterol Bound to Lipoproteins (LP) Low-Density) are closely related with Ischemic Coronary Heart Disease (IHD). Rather, the elevation HDL Cholesterol (Cholesterol Bound to High Density Lipoprotein) represents a protective factor against DCI. Hypertension, smoking, obesity are other causes responsible for atherosclerosis and DCI.

Lower values may occur in the presence of diseases that affect the liver parenchyma (viral hepatitis, toxic hepatitis), occasionally acute infections (pneumonia, typhoid fever), hyperthyroidism, anemia, malnutrition.

BIBLIOGRAPHIC REFERENCES

- 1 - ALLAIN, C.C., and al. Clin. Chem., 1974. 20-470.
- 2 - TONKS, D.B., Qualyt Control in Clinical laboratories, 1983.
- 3 - TRINDER, P., Ann. Clin. bioche.,, 1969, 6-24.
- 4 - HENRY, J.B., Diagnosys and Management by laboratory Methods, 19ª ed., 1996.
- 5 - CARL, A.B. and EDWARDR .A., Tiedz Textbook of Clinical Chem. 2nd ed., 1994,1002-1081.
- 6 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 8 - Faludi AA, Izar MCO, Saraiva JFK, Chacra APM, Bianco HT, Afune Neto A et al. Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose – 2017. Arq Bras Cardiol 2017; 109(2Supl.1):1-76
- 9 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Cholesterol Monoreagent kit:
10269360378

Review: April/2025

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER