COMPLEMENTO C3

REF K076

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1. Clinical Guide to Laboratory Tests, Edited by NW Tiszt W B Saunders Co., Philadelphia, 483, 1983.
2. Carrol MC. Annual Review of Immunology 1998;16:545-568.
3. Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1:16-45.
4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1966;14:401-406.
6. Young DS. Effects of Disease on Clinical Laboratory Test, 3th ed. AACC Pres, 1997.
7. Müller-Eberhard, H.H., Ann. Rev. Biochem. 44, 697 (1975).
8. Lachmann, P.J., Hobart, M.J. and Ashton, W.P. (1973) in Handbook of Experimental Immunology, 2nd Ed., 16, Ed. D.M. Weir, Blackwell Scientific Publications.
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Clinical Guide to Laboratory Tests, Edited by NW Tiszt W B Saunders Co., Philadelphia, 483, 1983.
2. Carrol MC. Annual Review of Immunology 1998;16:545-568.
3. Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1:16-45.
4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1966;14:401-406.
6. Young DS. Effects of Disease on Clinical Laboratory Test, 3th ed. AACC Pres, 1997.
7. Müller-Eberhard, H.H., Ann. Rev. Biochem. 44, 697 (1975).
8. Lachmann, P.J., Hobart, M.J. and Ashton, W.P. (1973) in Handbook of Experimental Immunology, 2nd Ed., 16, Ed. D.M. Weir, Blackwell Scientific Publications.
9. QUIBASA: Datos do Departamento de Pesquisa e Desenvolvimento.

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C3 COMPLEMENT

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST,
CONTACT CUSTOMER ADVISORY SERVICE:

REF **K076**

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

INSTRUCTIONS FOR USE

FUNCTION

Method for quantitative determination of C3 Complement in biological serum or plasma (EDTA or Heparin). Turbidimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Turbidimetry

Turbidimetric assay for quantifying determination of C3 Complement in human serum or plasma. The anti-C3 antibodies form insoluble compounds when combined with the C3 sample of the patient, causing a proportional relationship of the absorbance C3 concentration in the sample, and that can be measured by compared with a calibrator of known concentration of C3.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, PBS buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for C3 Complement < 5 g/L, PBS buffer and preservative.

PRESENTATION

Presentations	Reagent	
	N° 1	N° 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
3	3 x 25 mL	3 x 2 mL
4	4 x 25 mL	4 x 2 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- It is important for the good performance of the test a rigorous control of time and temperature.

5- Reagents n° 1 and 2 contain Sodium Azide, irritant to skin and mucous. Handle carefully.

6- Do not freeze reagents. Freezing of antibody or diluents may affect their functionality.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Fresh serum or plasma obtained with heparin or EDTA. The analyte is stable for 2 days at 2 to 8°C. If you need to store for a longer period, the sample should be frozen.

PROCESS DESCRIPTION

The calibration stability of the C3 Complement kit installed on refrigerated equipment is at least 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Reagents are ready to use.

TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

PROCEDURE LIMITATIONS

Interferences: Bilirubin (40 mg/dL), Hemoglobin (19 g/dL) Rheumatoid Factor (600 UI/mL) do not interfere. Lipids above (1000 mg/dL) interfere.

Note:

Clinical diagnosis should not be made using only the results of a single test, and should consider at the same time all the patient's clinical data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

75 - 135 mg/dL (IFCC)

These values should be used as guidance and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 98 and 104%.

Comparison of Methods

The C3 Complement kit was compared with other commercially available methods for determination of C3 Complement available. 40 tests were performed and the results were evaluated. The linear equation obtained was $Y = 1.034X + 2.318$ and correlation coefficient 0.984. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2
Average Concentration (mg/dL)	96.25	55.72
Standard Deviation (mg/dL)	0.58	0.36
Coefficient of Variation (%)	0.60	0.65

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2
Average Concentration (mg/dL)	95.97	54.78
Standard Deviation (mg/dL)	0.36	0.90
Coefficient of Variation (%)	0.37	1.65

QUIBASA QUÍMICA BÁSICA Ltda
 Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
 Phone.: 0800 0315454 | Email: sac@bioclin.com.br

ANVISA registration for C3 Complement kit:10269360131

Review: November/2025

SENSITIVITY

The sensitivity was calculated from 10 determinations of a sample free the C3 Complement. The average found, that indicates the method detection limit, was 0.22 mg/dL and the standard deviation of 0.04 mg/dL. The sensitivity, that indicates the method quantification limit, corresponds the average plus 3 times the standard deviation, and is equal to 0.35 mg/dL.

LINEARITY

The reaction is linear up to the level of the highest calibrator concentration.

DIAGNOSTIC SIGNIFICANCE

The C3 complement, component with the highest concentration of the whole plasma complement system is activated through the classical pathway and the alternative pathway. It is a protein synthesized by the liver, where bacterial endotoxins can induce their synthesis by monocytes and fibroblasts.

The concentration of C3 increases as a result of a response acute phase (trauma, inflammation or tissue necrosis), obstruction bile and focal glomerulosclerosis.

The C3 concentration may be impaired as a result of a genetic defect, the increased risk of infection particularly with encapsulated bacteria or a deficiency that causes acquired vascular problems and severe infections.

BIBLIOGRAPHIC REFERENCES

1. Clinical Guide to Laboratory Tests, Edited by NW Tiszt W B Saunders Co., Philadelphia, 483,1983.
2. Carrol MC. Annual Review of Immunology 1998;16:545-568.
3. Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993; Vol1:16-45.
4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1966;14:401-406.
6. Young DS. Effects of Disease on Clinical Laboratory Test, 3th ed. AACC Pres, 1997.
7. Müller-Eberhard, H.H., Ann. Rev. Biochem. 44, 697 (1975).
8. Lachmann, P.J., Hobart, M.J. and Ashton, W.P. (1973) in Handbook of Experimental Immunology, 2nd Ed., 16, Ed. D.M. Weir, Blackwell Scientific Publications.
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

UNIVERSAL SYMBOLOGY