COMPLEMENTO C4

REF K077

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1. Clinical Guide to Laboratory Tests, Edited by NW Tiszt W B Saunders Co., Philadelphia, 483, 1983.
2. Carrol MC. Annual Review of Immunology 1998; 16:545-568.
3. Lambris J D. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993;

Vol1:16-45.

4. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.

5. Dati F et al. Eur J Clin Chem Clin Biochem 1966; 14:401-406.

6. Young DS. Effects of Disease on Clinical Laboratory Test, 3th ed. AACC Pres, 1997.

7. Müller-Eberhard, H.H., Ann. Rev. Biochem. 44, 697 (1975)

8. Lachmann, P.J., Hobart, M.J. and Ashton, W.P. (1973) in Handbook of Experimental Immunology, 2nd Ed., 16, Ed. D.M. Weir, Blackwell Scientific Publications

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REF K077

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method of quantitative determination of C4 Complement in biological samples serum or plasma (EDTA or Heparin). Turbidimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Turbidimetry

Turbidimetric assay for quantification of C4 Complement human serum or plasma. The anti-C4 form insoluble compounds when combined with the C4 sample of the patient, causing a proportional relationship of the absorbance C4 concentration in the sample, and that can be measured by compared with a calibrator of known concentration of C4.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, Buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for C4 Complement < 5g/L, Buffer and preservative.

PRESENTATION

Presentation	Reagent	
	Nº 1	Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
3	3 x 25 mL	3 x 2 mL
4	4 x 25 mL	4 x 2 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

This kit is specific for use in automated equipments.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Reagents nº 1 and 2 contain Sodium Azide, irritant to skin and mucous. Handle carefully.

5- Do not freeze the reagents. The freezing of Antiserum or Buffer can affect their performance.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma obtained with heparin or EDTA. This sample is stable for 2 days if kept in temperatures between 2 and 8°C or 3 months at -20°C. Do not use samples highly hemolyzed or lipemic.

PROCESS DESCRIPTION

The calibration stability of the C4 Complement kit installed on refrigerated equipment is up to 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Reagents are ready to use.

TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

PROCEDURE LIMITATIONS

Interferences: Bilirubin (40 mg/dL), hemoglobin (19 g/dL) rheumatoid factor (600 UI/mL) do not interfere. Lipids above (1000 mg/dL) interfere.

Note: Clinical diagnosis should not be made using only the results of a single test, and should consider at the same time all the patient's clinical data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

9-36 mg/dL (IFCC)

These values should be used as guidance and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 96 and 102%.

Comparison of Methods

The C4 Complement kit was compared with other commercially available methods for measurement of C4. 40 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1.069X - 0.693$ and the correlation coefficient 0.991. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2
Concentration (mg/dL)	17.44	9.72
Standard Deviation (mg/dL)	0.08	0.09
Coefficient of Variation (%)	0.48	0.95

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2
Concentration (mg/dL)	17.49	9.80
Standard Deviation (mg/dL)	0.10	0.14
Coefficient of Variation (%)	0.60	1.38

SENSITIVITY

The sensitivity was calculated from 10 determinations of a sample free the C4 Complement. The average found, that indicates the method detection limit, was 0.126 mg/dL and the standard deviation of 0.021 mg/dL. The sensitivity, that indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.188 mg/dL.

LINEARITY

The reaction is linear up to the level of the highest calibrator concentration.

DIAGNOSTIC SPECIFICITY

Complement C4, the second reactive component of the classical pathway of complement activation. It is a protein synthesized by the liver, or can be synthesized by monocytes or other tissue. The C4 concentration in plasma increases as a result of an acute phase response (trauma, inflammation or necrosis tissue). A complete genetic deficiency induces the loss of C4 concentration in plasma, associated with high prevalence of autoimmune diseases or collagen vascular diseases, particularly Lupus Erythematosus (SLE). Its concentration can also decrease as a consequence of consumption in the formation of immune complexes.

BIBLIOGRAPHIC REFERENCES

1. Clinical Guide to Laboratory Tests, Edited by NW Tiszt W B Saunders Co., Philadelphia, 483,1983.
2. Carrol MC. Annual Review of Immunology 1998;16:545-568.
3. Lambris JD. Cruse JM Lewis RE Jr (eds): Complement Today. Complement Profiles. Basel, Karger, 1993;

Vol1:16-45.

4. Pesce A J and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, 1987.
5. Dati F et al. Eur J Clin Chem Clin Biochem 1966;14:401-406.
6. Young D S. Effects of Disease on Clinical Laboratory Test, 3th ed. AACC Pres, 1997.
7. Müller-Eberhard, H.H., Ann. Rev. Biochem. 44, 697 (1975)
8. Lachmann, P.J., Hobart, M.J. and Ashton, W.P. (1973) in Handbook of Experimental Immunology, 2nd Ed., 16, Ed. D.M. Weir, Blackwell Scientific Publications
9. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for C4 Complement kit: 10269360132

Review: March/2026

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER