Bioclin

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REF **K222**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Creatinine in biological samples of serum, plasma (Heparin, Oxalate, Fluoride, Citrate or EDTA) or urine. Colorimetric Kinetic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Kinetic

Fixed Time Kinetic: Creatinine reacts with Picrate alkaline in half-tamponed, obtaining a chromogen whose absorbance is proportional to the concentration of Creatinine in the sample. The nonspecific chromogens are eliminated by a pre-reading because they have immediate formation.

REAGENTS

Number 1 - Buffer - Store between 15 and 30°C. Contains: Sodium Hydroxide < 500 mmol/L and Sodium Carbonate < 75 mmol/L.

Number 2 - Picric Acid - Store between 15 and 30°C. Contains: Picric Acid < 60 mmol/L.

Number 3 - Standard - Store between 15 and 30°C. After handling, store in refrigerator between 2 and 8°C to prevent evaporation. Contains: Creatinine 4,0 mg/dL, stabilizer and preservative.

PRESENTATION

Reagent N°1	Reagent N°2	Reagent N°3
1 x 240 mL	1 x 60 mL	1 x 3 mL
4 x 40 mL	2 x 20 mL	1 x 3 mL
4 x 60 mL	3 x 20 mL	1 x 3 mL
4 x 20 mL	4 x 5 mL	1 x 3 mL
4 x 32 mL	2 x 16 mL	1 x 3 mL
1 x 40 mL	1 x 10 mL	1 x 3 mL
2 x 40 mL	2 x 10 mL	1 x 3 mL
4 x 40 mL	4 x 10 mL	1 x 3 mL
4 x 60 mL	4 x 15 mL	1 x 3 mL
1 x 160 mL	1 x 40 mL	1 x 3 mL
1 x 240 mL	1 x 60 mL	-
4 x 40 mL	2 x 20 mL	-
4 x 60 mL	3 x 20 mL	-
4 x 20 mL	4 x 5 mL	-
4 x 32 mL	2 x 16 mL	-
1 x 40 mL	1 x 10 mL	-
2 x 40 mL	2 x 10 mL	-
4 x 40 mL	4 x 10 mL	-
4 x 60 mL	4 x 15 mL	-
1 x 160 mL	1 x 40 mL	-

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, watch or stopwatch, pipettes, test tubes, Biocontrol N Bioclin and Biocontrol P Bioclin. These items are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma collected with Heparin, Oxalate, Fluoride, Citrate or EDTA. The compound is stable for 07 days between 2 to 8°C and 3 months at -20°C.

Isolated urine: The urine should be collected and preserved in refrigerator until the time of dosage.

24 hours Urine: The urine should be collected within a period of 24 hours and preserved in refrigerator until the time of dosage.

PROCESS DESCRIPTION

The calibration stability of the Creatinine Automated kit installed in refrigerated equipment is at least 7 days. This stability may vary according to the test conditions, equipment and environment. Therefore, it is suggested to monitor the product's performance using control serum.

PREPARATION OF WORKING REAGENT

The reagents are ready to use in bireagent format and on automated equipment that uses dedicated presentations.

For use in monoreagent format, prepare the working reagent as follows:

Mix 4 parts of Reagent N°1 with 1 part of Reagent N°2. Store in amber bottle. The reagent is stable for 24 hours at 15 to 30°C and 07 days at 2 to 8°C. Homogenize the reagent of work before starting the technique.

TECHNIQUE

We also recommend using the Biocontrol N and Biocontrol P Bioclin kits as control serums.

Monoreagent:

Add 100 µL of Standard or Sample (serum, plasma or urine) to 1 mL of Working Reagent (previously homogenized), mix and immediately transfer to a thermostated cuvette 37°C. Measure absorbance of Standard and Sample at 510 nm (490 - 520nm) at 30 and 90 seconds

For the calibration of the kit in bireagent format, Bioclin recommends using the BIOCAL Bioclin calibrating serum.

Check the equipment's programming by calling Customer Service at 0800 031 5454.

Urine Dosage

Take note of the volume collected in mL, centrifuge an aliquot of urine sample to be measured and make a 1:25 dilution (0.1 mL of urine and 2.4 mL of distilled water). Then make the estimation as in technique described above. Multiply the result by 25.

CALCULATIONS

ΔA of Standard = Abs. 90 seconds - Abs. 30 seconds
or Sample

$$\text{Uncompensated Creatinine (mg/dL)} = \frac{\Delta A \text{ of Sample} \times \text{Standard Concentration (4mg/dL)}}{\Delta A \text{ of Standard}}$$

Since the reaction follows the Lambert-Beer Law, the Calibration Factor may be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (4mg/dL)}}{\Delta A \text{ of Standard}}$$

$$\text{Uncompensated Creatinine (mg/dL)} = \Delta A \text{ of Sample} \times \text{Calibration Factor}$$

$$\text{Creatinine in Urine (mg/24h)} = \frac{\text{mg/dL} \times \text{Volume of 24h (mL)}}{100}$$

COMPENSATION FACTOR:

The Jaffe method suffers significant interference from other plasma compounds, mainly proteins, which provide a positive bias in serum creatinine concentrations. The use of the Compensation Factor minimize error's of agreement with other more specific techniques and must be applied as described below:

$$\text{Compensated (mg/dL) Creatinine} = \text{Uncompensated (mg/dL) Creatinine} - \text{Compensating Factor (0.30 mg/dL)}$$

Do not apply the Compensation Factor for urine samples.

CLEARANCE

Hydrate the patient, with at least 600 mL of water and then guide you to empty the bladder. Creatinine clearance should be performed on a sample taken at 4, 12 or 24 hours. In any time of the chosen time interval, collecting a sample of blood. Measure the total volume of urine and calculate the volume per minute (VM).

$$\text{VM} = \frac{\text{Urinary volume in determined time (mL)}}{\text{Determined time (in minutes)}}$$

Dose the serum (applying the compensation factor) and urine creatinine following the methodologies proposed.

To calculate creatinine clearance:

U = Urinary Creatinine (mg/dL)

S = Compensated Creatinine (mg/dL)

VM = Volume/minute

A = Value of the patient's body surface area in m²

$$\text{Creatinine Clearance (mL/min)} = \frac{U \times \text{VM}}{S}$$

Note: This value should be corrected for patient's body surface. Using the nomogram of height-weight correlation. Multiply the value of clearance by 1.73 and divide by the surface of patient's body.

$$\text{Creatinine Clearance (mL/min/1.73 m}^2\text{)} = \frac{U \times \text{VM} \times 1.73}{S \times A}$$

Example: Clearance (24-hour urine)

U = Urinary Creatinine = 110 mg/dL

S = Compensated Creatinine = 0.9 mg/dL

24-Hour urine volume = 1480 mL in 24 hours (1440 minutes) de 24 hours

$$\text{Volume / minute} = \text{VM} = \frac{1480}{1440} = 1,0277$$

Patient's weight = 60 kg
Patient height = 165 cm
Body surface area = 1,66 m²

$$\text{Clearance} = \frac{110}{0.9} \times 1,0277 = 125,60 \text{ (mL/min)}$$

$$\text{Clearance (mL/min/1.73 m}^2\text{)} = \frac{110}{0.9} \times 1,0277 \times \frac{1,73}{1,66} = 130,889$$

Glomerular Filtration Rate (GFR) according to CKD-EPI (2021)

The equations below are as recommended by the National Kidney Foundation and are designed to estimate Glomerular Filtration Rate (GFR) in adults using creatinine values traceable to NIST reference materials assigned by the IDMS method, taking into account serum creatinine, sex, and age.

The results are expressed in mL/min/1.73m².

	Serum creatinine (mg/dL)	Equation
Women	≤0.7	GFR = 142 × (Scr/0.7) ^{-0.241} × (0.9938) ^{Age} × 1.012
	>0.7	GFR = 142 × (Scr/0.7) ^{-1.200} × (0.9938) ^{Age} × 1.012
Men	≤0.9	GFR = 142 × (Scr/0.9) ^{-0.302} × (0.9938) ^{Age}
	>0.9	GFR = 142 × (Scr/0.9) ^{-1.200} × (0.9938) ^{Age}

Scr = Compensated Serum Creatinine.

Values above 90 mL/min/1.73m² are released as greater than 90 mL/min/1.73m².

PROCEDURE LIMITATIONS

This process can only be accomplished with a spectrophotometer cuvette at 37°C, due to the small time interval reaction. **The methodology described should not be used for extremely lipemic samples.**

INTERFERENT

Triglyceride concentrations up to 900 mg/dL, Bilirubin up to 5 mg/dL and Hemoglobin up to 180 mg/dL do not interfere significantly with the reaction.

Samples with Bilirubin, Hemoglobin and Triglycerides in concentrations greater than those above mentioned can be diluted in NaCl 150 mmol/L (0.85%) before performing the tests to reduce interferences. The result obtained should be multiplied by the dilution factor.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard and the BIOCAL calibrator are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 914, making it equivalent to the IDMS method

(Isotope Dilution Mass Spectrometry), as recommended by the NKDEP (National Kidney Disease Education Program).

REFERENCE VALUES

The reference values in mg/dL, for this method were obtained through the determination of uncompensated Creatinine in healthy populations of male and female.

Serum/Plasma ⁵	Adult (Male)	0.70 - 1.20 mg/dL
	Adult (Female)	0.50 - 1.00 mg/dL
	Newborn	0.31 - 0.92 mg/dL
	2 weeks - 1 year	0.16 - 0.39 mg/dL
	1 - 3 years	0.17 - 0.35 mg/dL
	3 - 5 years	0.26 - 0.42 mg/dL
	5 - 7 years	0.29 - 0.48 mg/dL
	7 - 9 years	0.34 - 0.55 mg/dL
	9 - 11 years	0.32 - 0.64 mg/dL
	11 - 13 years	0.42 - 0.71 mg/dL
13 - 15 years	0.46 - 0.81 mg/dL	
Urine	Adult (Male)	1500 - 2500 mg/24h or 21 - 26 mg/kg/24h
	Adult (Female)	800 - 1500 mg/24h or 16 - 22 mg/kg/24h
	Children 2 - 3 years	6 - 22 mg/kg/24h
	Children > 3 years	12 - 30 mg/kg/24h
Creatinine clearance	Adult (Male)	97 - 137 mL/min/1.73m ²
	Adult (Female)	88 - 128 mL/min/1.73m ²
	Children	70 - 140 mL/min/1.73m ²

mg/kg weight = mg/24h divided by body weight

To convert the values to mg/dL to mmol/L (SI), multiply by 0.0884. These values should be used as guidelines, each laboratory should create their own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Creatinine Automated kit was compared with other methods commercially available for measurement of Creatinine. 40 analyses were performed and the results were evaluated. The linear equation was obtained $y = 0.9913x - 0.0248$ and correlation coefficient 0.9991. With these results, we can conclude that the kit shows good methodology specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Samples		
	1	2	3
Average Concentration	1.01	1.94	5.55
Standard Deviation	0.01	0.01	0.04
Variation Coefficient (%)	1.20	0.74	0.72

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Samples		
	1	2	3
Average Concentration	1.02	1.95	5.60
Standard Deviation	0.01	0.02	0.05
Variation Coefficient (%)	1.34	0.97	0.87

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of Creatinine. The average found that indicates the method detection limit was -0.009 mg/dL, with a standard deviation 0.011 mg/dL. The sensitivity, that indicates the method quantification limit, corresponds the average plus 3 times the standard deviation, and is equal to 0.034 mg/dL.

LINEARITY

The reaction color is linear to the concentration of 12 mg/dL. For higher values dilute the sample with Sodium Chloride 0.85% and repeat the determination. Multiply the value obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Creatinine is synthesized primarily in the liver and kidneys, is produced and excreted at a constant rate independent of diet, degree of hydration and protein metabolism. With renal failure, the laboratory data provide insurance rates as the ability kidneys to excrete, secrete and reabsorb creatinine.

Creatinine is elevated in acute and chronic renal failure, in urinary tract obstruction, congestive heart failure, in dehydration, shock and in diabetes *mellitus*.

Decreased levels may be observed in muscular dystrophies, malnutrition, reduced muscle mass and severe renal disease.

BIBLIOGRAPHIC REFERENCES

- 1 - SLOT, C.; Sand. J. Clin. Lab. Invest. 17:381, 1965.
- 2 - HENRY, R. J.: Clinical Chemistry - Principles and Technics, 4 Ed., Harper and Row. New York, 1964.
- 3 - TODD, SANFORD, DAVIDSOHN; Diagnósticos Clínicos, 16 Ed., Editora Manole, 1983.
- 4 - CARL, A. B. and EDWARD, R. A.: Tietz Textbook of Clinical Chem 2nd ed., 1994, 1531-1539.
- 5 - CERIOTTI, F. et al. Reference intervals for serum creatinine concentrations: Assessment of available data for global application. Clinical Chemistry, 2008.
- 6 - Inker LA, Eneanya ND, Coresh J, Tighiouart H, Wang D, Sang Y, et al. New creatinine- and cystatin C-based equations to estimate GFR without race. N Engl J Med. 2021;385(19):1737-49
- 7- National Kidney Foundation. CKD-EPI Creatinine Equation (2021). Disponivel em: <https://www.kidney.org/content/ckd-epi-creatinine-equation-2021>.
- 8- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Creatinine Automated kit: 10269360328

Review: October/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER