# Bioclin

## CREATININA CINÉTICA

REF **K067**

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#### FINALIDADE

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### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of Creatinine in biological samples of serum, plasma (Heparin, Oxalate, Fluoride, Citrate or EDTA) and urine. Colorimetric Kinetic test, for *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Colorimetric Kinetic

Fixed Time Kinetic: Creatinine reacts with Picrate alkaline in half-tamponed, obtaining a chromogen whose absorbance is proportional to the concentration of Creatinine in the sample. The nonspecific chromogens are eliminated by a pre-reading because they have immediate formation.

#### REAGENTS

**Number 1 - Buffer** - Store between 15 and 30°C. Contains: Sodium Hydroxide < 500 mmol/L and Sodium Carbonate < 75 mmol/L.

**Number 2 - Picric Acid** - Store between 15 and 30°C. Contains: Picric Acid < 60 mmol/L.

**Number 3 - Standard** - Store between 15 and 30°C. After handling, store in refrigerator between 2 and 8°C to prevent evaporation. Contains: Creatinine 3.0 mg/dL, stabilizer and preservative.

#### PRESENTATION

Presentation	Reagent N°1	Reagent N°2	Reagent N°3
1	1 x 100 mL	1 x 100 mL	1 x 3 mL
2	2 x 100 mL	2 x 100 mL	1 x 3 mL
3	2 x 200 mL	2 x 200 mL	1 x 3 mL
4	1 x 40 mL	1 x 40 mL	1 x 3 mL
5	2 x 40 mL	2 x 40 mL	1 x 3 mL
6	3 x 40 mL	3 x 40 mL	1 x 3 mL
7	4 x 40 mL	4 x 40 mL	1 x 3 mL
8	5 x 40 mL	5 x 40 mL	1 x 3 mL
9	6 x 40 mL	6 x 40 mL	1 x 3 mL
10	1 x 60 mL	1 x 60 mL	1 x 3 mL
11	2 x 60 mL	2 x 60 mL	1 x 3 mL
12	3 x 60 mL	3 x 60 mL	1 x 3 mL
13	4 x 60 mL	4 x 60 mL	1 x 3 mL
24	5 x 20 mL	5 x 20 mL	1 x 3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

#### SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum or plasma collected with Heparin, Oxalate, Fluoride, Citrate or EDTA Bioclin. The compound is stable for 07 days between 2 and 8°C. Isolated urine: The urine should be collected and preserved in refrigerator until the time of dosage.

24 hours Urine: The urine should be collected within a period of 24 hours and preserved in refrigerator until the time of dosage.

#### PROCESS DESCRIPTION

The calibration stability of the Creatinine Kinetic kit installed on refrigerated equipment for up to 8 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### PREPARATION OF WORKING REAGENT

Mix 1 part of Reagent N° 1 with 1 part of Reagent N° 2. Store in amber bottle. Stable 24 hours between 15 and 30°C. Homogenize the reagent of work before starting the technique.

#### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Add 100 µL of Standard or Sample (serum, plasma or urine) to 1 mL of Working Reagent (previously homogenized), mix and immediately transfer to a thermostated cuvette 37°C. Measure absorbance of Standard and Sample at 510 nm (490 - 520nm) at 30 and 90 seconds.

#### Urine Dosage

Take note of the volume collected in mL, centrifuge an aliquot of urine sample to be measured and make a 1:25 dilution (0,1 mL of urine and 2,4 mL of distilled water). Then make the estimation as in technique described above. Multiply the result by 25.

#### CALCULATIONS

$\Delta A$  of Standard = Abs. 90 seconds - Abs. 30 seconds  
or Sample

$$\text{Creatinine (mg/dL)} = \frac{\Delta A \text{ of Sample} \times \text{Standard Concentration (3mg/dL)}}{\Delta A \text{ of Standard}}$$

Since the reaction follows the Lambert-Beer Law, the Calibration Factor may be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (3mg/dL)}}{\Delta A \text{ of Standard}}$$

$$\text{Creatinine (mg/dL)} = \Delta A \text{ of Sample} \times \text{Calibration Factor}$$

$$\text{Creatinine in Urine (mg/24h)} = \frac{\text{mg/dL} \times \text{Volume of 24h (mL)}}{100}$$

#### CLEARANCE

Hydrate the patient, with at least 600 mL of water and then guide you to empty the bladder. Creatinine clearance should be performed on a sample taken at 4, 12 or 24 hours. In any time of the chosen time interval, collecting a sample of blood. Measure the total volume of urine and calculate the volume per minute (VM).

$$\text{VM} = \frac{\text{Urinary volume in determined time (mL)}}{\text{Determined time (in minutes)}}$$

Dose the serum and urine creatinine following the methodologies proposed.

To calculate creatinine clearance:

U = Urinary Creatinine (mg/dL)

S = Serum Creatinine (mg/dL)

VM = Volume/minute

A = Value of the patient's body surface area in m<sup>2</sup>

$$\text{Creatine Clearance (mL/min)} = \frac{U \times \text{VM}}{S}$$

Note: This value should be corrected for patient's body surface. Using the nomogram of height-weight correlation.

Multiply the value of clearance by 1.73 and divide by the surface of patient's body.

$$\text{Creatine Clearance (mL/min/1.73 m}^2\text{)} = \frac{U \times \text{VM} \times 1.73}{S \times A}$$

Example: Clearance (24-hour urine)

U = Urinary Creatinine = 110 mg/dL

S = Serum Creatinine = 0,9 mg/dL

24-Hour urine volume = 1480 mL in 24 hours (1440 minutes)

$$\text{Volume/minute} = \text{VM} = \frac{1480}{1440} = 1.0277$$

Patient's weight = 60 kg

Patient height = 165 cm

Body surface area = 1.66 m<sup>2</sup>

$$\text{Clearance (mL/min)} = \frac{110 \times 1.0277}{0.9} = 125.60$$

$$\text{Clearance (mL/min./1.73 m}^2\text{)} = \frac{110 \times 1.0277 \times 1.73}{0.9 \times 1.66} = 130.889$$

#### PROCEDURE LIMITATIONS

This process can only be accomplished with a spectrophotometer cuvette at 37°C, due to the small time interval reaction. **The methodology described should not be used for extremely lipemic samples.**

**INTERNAL QUALITY CONTROL**

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**TRACEABILITY**

The kit's standard and the BIOCAL calibrator are traceable to the reference material NIST (National Institute of Standards and Technology) SRM 914, making it equivalent to the IDMS method (Isotope Dilution Mass Spectrometry), as recommended by the NKDEP (National Kidney Disease Education Program).

**REFERENCE VALUES**

The reference values in mg/dL, for this method were obtained through the determination of Creatinine in healthy populations of male and female.

<b>Serum/Plasma<sup>5</sup></b>	<b>Adult</b>	0.4 - 1.4 mg/dL
	<b>Newborn</b>	0.31 - 0.92 mg/dL
	<b>2 weeks - 1 year</b>	0.16 - 0.39 mg/dL
	<b>1 - 3 years</b>	0.17 - 0.35 mg/dL
	<b>3 - 5 years</b>	0.26 - 0.42 mg/dL
	<b>5 - 7 years</b>	0.29 - 0.48 mg/dL
	<b>7 - 9 years</b>	0.34 - 0.55 mg/dL
	<b>9 - 11 years</b>	0.32 - 0.64 mg/dL
	<b>11 - 13 years</b>	0.42 - 0.71 mg/dL
	<b>13 - 15 years</b>	0.46 - 0.81 mg/dL
<b>Urine</b>	<b>Adult (Male)</b>	1500 - 2500 mg/24h or 21 - 26 mg/kg/24h
	<b>Adult (Female)</b>	800 - 1500 mg/24h or 16 - 22 mg/kg/24h
	<b>Children 2 - 3 years</b>	6 - 22 mg/kg/24h
	<b>Children &gt; 3 years</b>	12 - 30 mg/kg/24h
<b>Creatinine clearance</b>	<b>Adult (Male)</b>	97 - 137 mL/min/1.73m <sup>2</sup>
	<b>Adult (Female)</b>	88 - 128 mL/min/1.73m <sup>2</sup>
	<b>Children</b>	70 - 140 mL/min/1.73m <sup>2</sup>

mg/kg weight = mg/24h divided by body weight

To convert the values to mg/dL to mmol/L (SI), multiply by 0.0884.

These values should be used as guidelines, each laboratory should create their own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****ACCURACY****Comparison of Methods**

The Creatinine Kinetics kit was compared with other methods commercially available for measurement of Creatinine. 42 analyses were performed and the results were evaluated. The linear equation was obtained  $Y = 1.0432X - 0.0708$  and correlation coefficient 0.9948. With these results, we can conclude that the kit shows good methodology specificity.

**PRECISION****Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Samples		
	1	2	3
<b>Average Concentration</b>	0.77	4.31	2.57
<b>Standard Deviation</b>	0.05	0.11	0.05
<b>Variation Coefficient (%)</b>	6.18	2.63	1.85

**Reproducibility**

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Samples		
	1	2	3
<b>Average Concentration</b>	0.77	4.34	2.57
<b>Standard Deviation</b>	0.01	0.05	0.02
<b>Variation Coefficient (%)</b>	0.65	1.20	0.59

**SENSITIVITY**

Sensitivity was calculated from 40 determinations of a sample free of Creatinine. The average found was 0.105 mg/dL, with a standard deviation 0.010 mg/dL. The sensitivity, which indicates the method detection limit, corresponds to the average plus 3 times the standard deviation, and is equal to 0.135 mg/dL.

**LINEARITY**

The reaction color is linear to the concentration of 10 mg/dL. For higher values dilute the sample with Sodium Chloride 0.85% and repeat the determination. Multiply the value obtained by the dilution factor.

**DIAGNOSTIC SIGNIFICANCE**

Creatinine is synthesized primarily in the liver and kidneys, is produced and excreted at a constant rate independent of diet, degree of hydration and protein metabolism. With renal failure, the laboratory data provide insurance rates as the ability kidneys to excrete, secrete and reabsorb creatinine.

Creatinine is elevated in acute and chronic renal failure, in urinary tract obstruction, congestive heart failure, in dehydration, shock and in diabetes *mellitus*.

Decreased levels may be observed in muscular dystrophies, malnutrition, reduced muscle mass and severe renal disease.

**BIBLIOGRAPHIC REFERENCES**

- 1 - SLOT, C.; Sand. J. Clin. Lab. Invest. 17:381, 1965.
- 2 - HENRY, R. J.: Clinical Chemistry - Principles and Technics, 4 Ed., Harper and Row. New York, 1964.
- 3 - TODD, SANFORD, DAVIDSOHN; Diagnósticos Clínicos, 16 Ed., Editora Manole, 1983.
- 4 - CARL, A. B. and EDWARD, R. A.: Tietz Textbook of Clinical Chem 2nd ed., 1994, 1531-1539.
- 5 - CERIOTTI, F. et al. Reference intervals for serum creatinine concentrations: Assessment of available data for global application. Clinical Chemistry, 2008.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**CUSTOMER SERVICE**

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Creatinine Kinetic kit: 10269360201

Review: August/2022

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <math>n</math>-TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER