DÍMERO D

REF K290/K291/K292

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- RAMZI, D.W., et al. American Family Physician 69(12): 2829-2836, 2004.
- WELLS, P.S., et al. New England Journal of Medicine 349: 1227-1235, 2003.
- WELLS, P.S., et al. Annals of Internal Medicine 135: 98-107, 2001.
- HUMPHREYS, C.W., et al. Thrombosis Research 113(5): 275-82, 2004.
- ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Lower-extremity Deep Vein Thrombosis. Annals of Emergency Medicine 41: 124-135, 2003.
- ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Pulmonary Embolism. Annals of Emergency Medicine 41(2): 257-270, 2003.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Leonard Naymagon, et al. Thromb Res. 2020 Dec; 196: 99-105.
- Linkins LA, Takach Lapner S. Int J Lab Hematol. 2017 May;39 S1:98-103.
- Suzuki K, Wada H, Imai H, et al. J Thromb Haemost. 2018;16: 1442-1444.
- Toh CH, Hoots WK. J Thromb Haemost. 2007;5: 604-606.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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REF **K290/K291/K292**

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Test for the quantitative determination of D-Dimer (D-D) in a biological sample of citrated plasma using an immunoturbidimetric test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry.

The reaction allows quantifying, using a turbidimetric method, the concentration of D-D present in the sample. The suspended latex particles coated with Anti-D-Dimer antibodies are agglutinated by Dimer D in the plasma. The agglutination process causes an increase in optical density. The intensity of the scattered light is proportional to the D-D concentration.

REAGENTS

Reagent Number 1 - Buffer - Store at 2 to 8 ° C. Contains: Buffer <200 mmol / L, stabilizers, surfactant, and preservative.

Reagent Number 2 - Latex D Dimer - Store at 2 to 8 ° C. Contains: Suspension of Latex Particles sensitized with Anti-Dimer D antibodies, stabilizers, and preservatives.

Reagent Number 3 - Calibrator - Store at 2 to 8 ° C. Contains: D-dimer, stabilizers, and preservatives. Potentially infective.

Reagent Number 4 - Control Level 1 - Store at 2 to 8 ° C. Contains: D-dimer, stabilizers, and preservatives. **Potentially infective.**

Reagent Number 5 - Control Level 2 - Store at 2 to 8 ° C. Contains: D-dimer, stabilizers, and preservatives. **Potentially infective.**

Attention: The D-dimer concentration for the calibrator varies by lot - see the vial label.

The controls concentrations vary on each lot. **Refer to table sent along with the product.**

Note: Reagents, controls, and calibrators may be marketed separately, and the reagent number may vary depending on the presentation.

PRESENTATION

Reagent	Presentation			
	K290-1	K291-1	K292-1	K290-2
Buffer	1 x 20 mL	-	-	1 x 20 mL
Latex D Dimer	1 x 10 mL	-	-	1 x 10 mL
Calibrator	1 x 1 mL	1 x 1 mL	-	-
Control Level 1	1 x 1 mL	-	1 x 1 mL	-
Control Level 2	1 x 1 mL	-	1 x 1 mL	-

EQUIPMENT AND OPERATIONAL INPUTS

Automated biochemistry analyzer or automated coagulometer are found in the specialized market of articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30 ° C should not exceed 7 (seven) days. Keep out of the light and avoid moisture. Do not freeze.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The water used to clean the material must be fresh and free from contaminants.
- 4- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents, significantly altering the results.
- 5- The reagents were tested for anti-HIV, anti-HCV, and HBs antigen antibodies using the latest methods and show negative results. The risk of infection cannot be excluded, and the reagent must be handled with the same care observed for the patient's plasma. Potentially infective.
- 6- To avoid possible analytical errors, the distilled water used in the product's reconstitution should be of adequate quality. Thus, it is recommended that the water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- 7- Make a new calibration curve whenever necessary and for each batch of the product.
- 8- We recommend applying local, state, and federal environmental protection standards so that the disposal of reagents and biological material is carried out following current legislation.
- 9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Service of Customer Service) of Quibasa.
- 10- Do not use the product in case of damage to the packaging.
- 11- It is essential that the instruments and equipment used are appropriately calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma citrated

Anticoagulant: Use the Bioclin A001 Citrate kit as an anticoagulant according to the instructions for use, or 3.8% Sodium Citrate (0.130 M) in the proportion of 9 parts of blood to 1 part of anticoagulant (as an example: 4.5 mL of blood + 0.5 mL of anticoagulant).

Collection: Blood must be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, bubble formation, and tissue fluid aspiration.

Preparation: The blood must be homogenized mixed with the anticoagulant immediately after collection. Centrifuge immediately after collection, at high speed (approximately 3000 rpm), for 15 minutes.

Stability

The samples must be tested within 3 hours after collection. The analyte is stable for 8 hours at 20 to 25 ° C, 4 days at 2 to 8 ° C, and 1 month at -20 ° C.

PROCESS DESCRIPTION

The calibration stability of the DIMER D kit installed may vary according to the conditions of the test, the analyzer, and the environment. Therefore, it is suggested to monitor the performance of the product using controls.

TECHNIQUE

The reagents for the "D-DIMER" presentations are ready for use. For D-Dimer Buffer and Latex reagents from the "D-DIMER" presentations, use the D-Dimer Calibrator and the D-Dimer Controls, levels 1 and 2.

The "D-Dimer" presentation kits can be used in biochemical analyzers or automated coagulometers.

Check the application for the specific equipment on the website www.bioclin.com.br or through the Customer Service.

D-DIMER CALIBRATION - To obtain the calibration curve for the "D-DIMER" presentations, a serial dilution of the calibrator must be carried out, as described below or in accordance with the equipment dilution table:

Dilución	1:1	1:2	1:4	1:8	1:16	1:32
D Dimer Calibrator	250 µL	250 µL	250 µL of 1:2 dilution	250 µL of 1:4 dilution	250 µL of 1:8 dilution	250 µL of 1:16 dilution
Distilled Water (µL)	0	250	250	250	250	250
Concentration (ng/mL)	8,000	4,000	2,000	1,000	500	250

The calibration curve will consist of 6 points, with the last point of the curve having a concentration equal to 250 ng/mL.

INTERFERENCES

Samples with a concentration of up to 20 mg/dL of Bilirubin, 2000 mg/dL of Triglycerides, and 500 mg/dL of Hemoglobin do not produce significant interference.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program in which procedures, standards, limits, and tolerance for variations are established. It is important to note that all measurement systems have a characteristic analytical variability, which the laboratories themselves must monitor. It is recommended to use controls, which allow evaluating the dosages' precision and accuracy.

REFERENCE VALUES

The present method's reference values were obtained by determining the D-Dimer in healthy male and female populations.

Plasma: Less than 500 ng / mL (0.5 µg / mL)

Conversion: µg / mL = ng / mL ÷ 1000

These values should be used as a guide, and each laboratory should reestablish its reference values according to the population served.

Note: The medical professional must always interpret the results in conjunction with the patient's medical history, physical examination, and other results.

PRODUCT PERFORMANCE

ACCURACY

Comparison of methods and methodological specificity
The DIMER D kit was compared with another commercially available method for dosing D-Dimer. Twenty three analyzes were performed, and the results were evaluated. The linear equation obtained was $y = 1.6361x - 669.57$, and the correlation coefficient was 0.994. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (ng/mL)	269.98	488.95	2322.8
Standard deviation (ng/mL)	14.678	27.355	31.94
Coefficient of variation (%)	5	5.6	1.4

The product presents results within the expected variation range (coefficient of variation $\leq 10\%$).

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (ng/mL)	263.58	748.37
Standard deviation (ng/mL)	17.452	21.629
Coefficient of variation (%)	7	3

The product presents results within the expected variation range (coefficient of variation $\leq 10\%$).

SENSITIVITY

The analytical sensitivity of the D Dimer kit is 220 ng/mL.

PROZONE EFFECT

No prozone effect was observed up to a concentration of 10,000 ng/mL.

LINEARITY

The assay determines D-Dimer concentrations within the range of 220 ng / mL to the calibration curve's highest point.

DIAGNOSTIC MEANING

D-dimer is generated as a product of the fibrinolysis process from fibrin degradation by plasmin. Elevated D-dimer levels are found in clinical conditions, such as deep vein thrombosis (DVT), pulmonary embolism and disseminated intravascular coagulation (ICD), sepsis, postoperatively. D-dimer levels increase during pregnancy, and elevated levels are associated with complications during pregnancy.

BIBLIOGRAPHIC REFERENCES

1. FEDULLO, P.F., et al. New England Journal of Medicine 349: 1247-1256, 2003.
2. S.Z. GOLDBERGER. New England Journal of Medicine 339: 93-104, 1998.
3. KLINE, J.A., et al. Journal of Thrombosis and Haemostasis 2(8):1247-1255, 2004.

4. RAMZI, D.W., et al. American Family Physician 69(12): 2829-2836, 2004.
5. WELLS, P.S., et al. New England Journal of Medicine 349: 1227-1235, 2003.
6. WELLS, P.S., et al. Annals of Internal Medicine 135: 98-107, 2001.
7. HUMPHREYS, C.W., et al. Thrombosis Research 113(5): 275-82, 2004.
8. ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Lower-extremity Deep Vein Thrombosis. Annals of Emergency Medicine 41: 124-135, 2003.
9. ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Pulmonary Embolism. Annals of Emergency Medicine 41(2): 257-270, 2003.
10. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
11. Leonard Naymagon, et al. Thromb Res. 2020 Dec; 196: 99-105.
12. Linkins LA, Takach Lapner S. Int J Lab Hematol. 2017 May;39 S1:98-103.
13. Suzuki K, Wada H, Imai H, et al. J Thromb Haemost. 2018;16: 1442-1444.
14. Toh CH, Hoots WK. J Thromb Haemost. 2007;5: 604-606.
15. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY WARRANTY

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The reagents' quality is assured until the expiration date stated on the presentation packaging when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Tel.: 0800 031 5454 | E-mail: sac@bioclin.com.br

ANVISA registration for D DIMER kit: 10269360344

Review: October/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER