TÉCNICA

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TECHNIQUE

Buffer and Latex D-Dimer reagents are ready for use. For the use of the D-Dimer buffer and latex reagents, use the 6 levels of D-Dimer calibration as a calibrator and control the 2 levels of D-Dimer controls.

The kit is only intended for use in automatic biochemical analyzers. **Check the application for the specific equipment on the website www.bioclin.com.br or through the Custome Service.**

Calibration

To perform the calibration curve, use the D-Dimer CALIBRATOR reagents from levels 1 to 6.

INTERFERENCES

Samples with a concentration of up to 20 mg/dL of Bilirubin, 2000 g/dL of Triglycerides, and 500 mg/dL of Hemoglobin do not produce significant interference.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program in which procedures, standards, limits, and tolerance for variations are established. It is important to note that all measurement systems have a characteristic analytical variability, which the laboratories themselves must monitor. It is recommended to use controls, which allow evaluating the dosages' precision and accuracy.

REFERENCE VALUES

The present method's reference values were obtained by determining the D-Dimer in healthy male and female populations.

Plasma: Less than 500 ng/mL (0.5 µg/mL)
Conversion: µg/mL = ng/mL + 1000

These values should be used as a guide, and each laboratory should reestablish its reference values according to the population served. **Note:** The medical professional must always interpret the results in conjunction with the patient's medical history, physical examination, and other results.

The responsible medical professional must interpret the results provided by this kit, and it is not the only criterion for determining the diagnosis and treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Comparison of methods and methodological specificity The D-Dimer kit was compared with another commercially available method for dosing D-Dimer. Twenty analyzes were performed, and the results were evaluated. The linear equation obtained was Y = 0.7928X - 135.9, and the correlation coefficient was 0.997. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION
REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (ng/mL)	1190	14360
Standard Deviation (ng/mL)	40	40
Coefficient of Variation (%)	3.36	0.29

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- FEDULLO, P.F., et al. New England Journal of Medicine 349: 1247-1256, 2003.
- S.Z. GOLDBABER. New England Journal of Medicine 339: 93-104, 1998.
- KLIN, J.A., et al. Journal of Thrombosis and Haemostasis 2(8):1247-1255, 2004.
- RAMZI, D.W., et al. American Family Physician 69(12): 2829-2836, 2004.
- WELLS, P.S., et al. New England Journal of Medicine 349: 1227-1235, 2003.
- WELLS, P.S., et al. Annals of Internal Medicine 135: 98-107, 2001.
- HUMPHREYS, C.W., et al. Thrombosis Research 113(5): 275-82, 2004.
- ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Lower-extremity Deep Vein Thrombosis. Annals of Emergency Medicine 41: 124-135, 2003.
- ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Pulmonary Embolism. Annals of Emergency Medicine 41(2): 257-270, 2003.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

REPRODUCTIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (ng/mL)	1190	14380
Standard Deviation (ng/mL)	50	40
Coefficient of Varation (%)	4.22	0.28

SENSITIVITY

The analytical sensitivity of the D-Dimer kit is 170 ng/mL.

Prozone Effect

No prozone effect was observed up to a concentration of 50,000ng/mL.

LINEARITY

The operating range for the D-Dimer reagent was designed to be 170 to 24,891 ng/mL. However, by itself deal with a method with multipoint calibration, it should be considered the concentration of the last calibration point as linearity of the kit. Samples with concentration above the last point of calibration should be diluted.

DIAGNOSTIC SIGNIFICANCE

D-Dimer is generated as a product of the fibrinolysis process from fibrin degradation by plasmin. Elevated D-Dimer levels are found in clinical conditions, such as deep vein thrombosis (DVT), pulmonary embolism and disseminated intravascular coagulation (CiD), sepsis, postoperatively. D-Dimer levels increase during pregnancy, and elevated levels are associated with complications during pregnancy

BIBLIOGRAPHIC REFERENCES

- FEDULLO, P.F., et al. New England Journal of Medicine 349: 1247-1256, 2003.
- S.Z. GOLDBABER. New England Journal of Medicine 339: 93-104, 1998.
- KLIN, J.A., et al. Journal of Thrombosis and Haemostasis 2(8):1247-1255, 2004.
- RAMZI, D.W., et al. American Family Physician 69(12): 2829-2836, 2004.
- WELLS, P.S., et al. New England Journal of Medicine 349: 1227-1235, 2003.
- WELLS, P.S., et al. Annals of Internal Medicine 135: 98-107, 2001.
- HUMPHREYS, C.W., et al. Thrombosis Research 113(5): 275-82, 2004.
- ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Lower-extremity Deep Vein Thrombosis. Annals of Emergency Medicine 41: 124-135, 2003.
- ACEP Clinical Policy; Critical Issues in the Evaluation and Management of Adult Patients Presenting with Suspected Pulmonary Embolism. Annals of Emergency Medicine 41(2): 257-270, 2003.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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TÓXICO

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The reagents' quality is assured until the expiration date stated on the presentation packaging when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for D-Dimer kit: 10269360344

Review: July/2021