DESIDROGENASE LÁTICA LDH UV

REF K014

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FINALIDADE

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REF K014

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Lactate Dehydrogenase in biological samples serum, plasma, urine or cerebrospinal fluid. Kinetic test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: UV Kinetic - DGKC

The Lactate Dehydrogenase (LDH) catalyzes the reduction of Pyruvate with NADH, and NAD<sup>+</sup> to Lactate. The catalytic concentration is determined by the speed decomposition of NADH, measured by the drop in absorptivity at 340 nm.

REAGENTS

**Number 1 – Buffered Substrate** - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, Pyruvate < 6 mmol/L, stabilizer and preservative.

**Number 2 - Coenzyme** - Store between 2 and 8°C. Contains: Buffer < 100 mmol/L, NADH < 5,0 mmol/L and preservative.

PRESENTATION

| Presentation | Reagent Nº 1 | Reagent Nº 2 |
|--------------|--------------|--------------|
| 1            | 54 mL        | 6 mL         |
| 2            | 1 x 40 mL    | 1 x 10 mL    |
| 3            | 2 x 40 mL    | 2 x 10 mL    |
| 4            | 4 x 40 mL    | 4 x 10 mL    |
| 5            | 2 x 40 mL    | 1 x 20 mL    |
| 6            | 4 x 40 mL    | 2 x 20 mL    |
| 7            | 4 x 27 mL    | 1 x 12 mL    |
| 8            | 6 x 36 mL    | 6 x 4 mL     |
| 9            | 6 x 36 mL    | 2 x 12 mL    |
| 10           | 3 x 36 mL    | 3 x 4 mL     |
| 11           | 3 x 36 mL    | 1 x 12 mL    |
| 12           | 6 x 36 mL    | 1 x 24 mL    |
| 13           | 1 x 20 mL    | 1 x 5 mL     |
| 14           | 2 x 20 mL    | 2 x 5 mL     |
| 15           | 3 x 20 mL    | 3 x 5 mL     |
| 16           | 4 x 20 mL    | 4 x 5 mL     |
| 17           | 5 x 20 mL    | 5 x 5 mL     |
| 18           | 6 x 20 mL    | 6 x 5 mL     |

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatch, test tubes, Biocal Bioclin, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** Water used in material cleaning must to be recent and free of contaminants.

**4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

**5-** It is important, for the good development of the test, a rigorous control of time and temperature.

**6-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**7-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**8-** Do not use the product in case of damaged packaging.

**9-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma collected with EDTA or Heparin, cerebrospinal fluid and urine. LDH is stable in serum and plasma samples for 2 days at 2 to 8°C and 20 days at - 20°C <sup>6</sup>. Samples of cerebrospinal fluid and urine must be stored at 2 to 8°C and dosed within 24 hours.

PROCESS DESCRIPTION

The calibration stability of the Lactate Dehydrogenase LDH UV kit installed on refrigerated equipment is up to 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix four (4) parts of Reagent Nº 1 with one (1) part of Reagent Nº 2. Working Reagent is stable during 14 days, if kept in temperatures between 2 and 8°C.

REACTION CONDITIONS

It is indispensable to use thermostated cuvette at 37°C, and 1 cm optical path reading at 340 nm.

TECHNIQUE

Bioclin recommends using the kit Biocal as calibrator and as control serum, Biocontrol N and P Bioclin Kits.

Add 20 µL of Sample to 1,0 mL of Working Reagent, mixing and transferring to a thermostated cuvette at 37°C and wait 1 minute. Make the initial reading simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the mean differences in absorbance per minute (ΔA/min.) and use it to calculate the result.

CALCULATIONS

LDH (U/L) = ΔA/min. x 8016

Results are expressed as U/L.

PROCEDURE LIMITATIONS

The specifications below refers to semi-automated equipments:

The kinetic method is based on the absorptivity and by this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions: Wavelength 340 nm.

Semi trajectory of the pass band 10 nm

Stray light less than 0,5%

1cm thermostated cuvette

INTERFERENCES

Hemolysis or late separation of serum could cause elevated results due to a high concentration of LDH in erythrocytes. Lipemia (Tryglicerides > 1000 mg/dL) and Bilirubin (> 20 mg/dL) could result in falsely elevated results. No interference was observed by Ascorbic Acid until 60 mg/dL and Rheumatoid Factor until 300 UI/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The calibration kit can be made using the theoretical calibration factor based on the molar absorptivity of NADH, or by calibrator BIOCAL.

REFERENCE VALUES

The reference values in U/L, for this method were obtained through the determination of LDH in healthy populations of male and female.

| Children and Adolescents |               |
|--------------------------|---------------|
| 1 - 3 years              | 490 - 730 U/L |
| 4 - 9 years              | 320 - 520 U/L |
| 10 - 13 years            | 250 - 500 U/L |
| Adults                   |               |
| 200 - 480 U/L            |               |

These values should be use as reference, each laboratory must create their own range of values for reference, according to the served population.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The LDH UV Lactic Dehydrogenase kit was compared with another commercially available method for measuring LDH. 42 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 1.0154X + 3.4167$  and the correlation coefficient was 0.9992. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

| Repeatability                | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Average Concentration (U/L)  | 212.65   | 142.73   | 537.45   |
| Standard Deviation (U/L)     | 9.74     | 6.81     | 22.68    |
| Coefficient of Variation (%) | 4.58     | 4.77     | 4.22     |

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

| Reproducibility              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Average Concentration (U/L)  | 212.51   | 139.66   | 552.79   |
| Standard Deviation (U/L)     | 7.56     | 3.73     | 23.46    |
| Coefficient of Variation (%) | 3.56     | 2.67     | 4.24     |

SENSITIVITY

The sensitivity was calculated from 40 determinations from a sample free of LDH. The average found was 3.95 U/L with a standard deviation of 3.99 U/L. The sensitivity, that indicates the method detection limit corresponds the average plus 3 times the standard deviation and is equal to 15.91 U/L.

LINEARITY

The reaction is linear to the concentration of 2000 U/L. For an average change in absorbance at 340 nm greater than 0.12, repeat the determination, diluting the sample with 0.85% NaCl. Multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Lactate Dehydrogenase (LDH) is present in all body cells, and the highest concentrations are in liver, heart, kidney, muscle skeletal and erythrocytes. The concentration of LDH in serum or plasma is increased in patients with liver diseases, renal disorders, myocardial infarction, many malignant diseases, progressive muscular dystrophy and in any case of hemolysis. The diagnosis should not be done taking into account only the result of a single assay of LDH, but it should be integrate clinical data with the laboratory.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until piration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Lactate Dehydrogenase LDH UV kit: 10269360074

Review: August/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE  
(last day of the month)

TEMPERATURE LIMIT  
(store)

CONTENT IS SUFFICIENT  
FOR <N> TEST

SEE INSTRUCTIONS  
FOR USE

IN VITRO DIAGNOSTIC  
PRODUCT

KEEP AWAY  
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF  
PACKAGE IS  
DAMAGED

PRODUCT  
STERILIZED

DANGER