DIROFILARIOSE Ag VET FAST

REF VET059

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Immunochromatographic test for detection of *Dirofilaria immitis* antigens in samples of whole blood, serum or plasma from canines. For *in vitro* diagnostics only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
 The Heartworm Ag VET FAST is an immunochromatographic test for detecting *Dirofilaria immitis* antigens in samples of whole blood, serum or plasma from dogs. After adding the sample to the cassette, the *Dirofilaria immitis* antigens interact with the conjugate, consisting of an anti-*D. immitis* conjugated to colloidal gold, forming the antigen-conjugate complex. The complex migrates by capillarity through the immunochromatographic strip until it reaches the test line region, consisting of anti-*D immitis* antibody. In this region, the *Dirofilaria immitis* antigens, bound to the conjugate, interact with the antibodies immobilized on the membrane, forming the antibody-antigen-conjugate complex. As a result, a pink line forms, indicating a positive result. The complex continues to migrate along the test strip until it reaches the control line, made up of anti-conjugate antibody. In this way, the formation of a second pink line (control line) is observed, indicating that the test presented a valid result.

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 30 °C. **Do not freeze.** Contains:

- A)** 01 plastic cassette with 01 test strip containing:
 - Control Region (C): anti-conjugate antibody.
 - Test Region (T): anti-*Dirofilaria immitis* antibody.
 - Conjugate: anti-*Dirofilaria immitis* antibody conjugated to colloidal gold

B) 01 sachet of silica

Number 2 (R2) Diluent - Store between 2 and 30 °C. Do not freeze. Contains: Buffered solution and preservative.

Accessories: Disposable plastic pipette, Collection tube with EDTA*, Alcoholic sachet* and Lancet*.

* Items forwarded according to presentation.

PRESENTATION

Reagent	Presentation	Number of Cassettes per Pack	Diluent
Cassette + Diluent	VET059-1	10	1 x 2.5 mL
	VET059-5	5	1 x 1.25 mL
	VET059-3	25	1 x 5.0 mL
	VET059-7*	10	10 x 0.5 mL
	VET059-8*	25	25 x 0.5 mL

*Presentations composed of individualized kits.

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table
 - Instructions for use (manual)

Required materials not contained in the kit:

Clock or stopwatch, PPE (gloves, mask, protective glasses) and trash for disposal. All listed inputs can be found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage and transport temperature of the product must be between 2 and 30 °C. Keep away from direct sunlight. **Do not freeze.**

The cassette is sensitive to air humidity. Therefore, it must always be kept in its primary packaging (aluminized sachet) intact, until the moment of use.

SPECIAL CARES

- 1- Product for *in vitro* diagnosis, for veterinary use only.**
2- Strictly follow the proposed methodology to obtain exact results.
3- The cassettes (test devices) do not contain contaminating substances, but as the samples are biological, safety standards and routines must be strictly followed. The use of disposable gloves and other personal protective equipment is essential.
4- Use suitable and suitable containers to store biological samples.
5- Open the aluminum envelope to remove the cassette (test device) only when using it.
6- The kit must not be frozen or exposed to high temperatures, as this will cause deterioration of its components.
7- The reagents, as well as the samples, must reach room temperature before starting the test.
8- Do not use reagents from different batches or kits.
9- Do not use a kit after its expiration date. The expiration date is indicated on the product label.
10- All kit components that come into contact with the sample are disposable and should not be reused.
11- The EDTA collection tube accessory that accompanies some presentations must be used exclusively with this product.
12- The components that came into contact with the patient's sample and the diluent that accompanies the kit must be discarded as biological waste. The kit's lancet, when present, must be discarded in sharps waste.
13- Secondary packaging (box), Instruction for Use, aluminum envelope and silica sachet can be discarded in general trash.
14- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
15- To obtain information related to biosafety or in the event of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through a request through the SAC (Customer Service Advisory) from Quibasa.
16- Do not use the product if the packaging is damaged.

SAMPLES

Whole Blood





Collect the blood sample in a tube containing EDTA, Heparin or Citrate, using standard laboratory procedures. Whole blood can be stored at 2 to 8 °C and tested up to 24 hours after collection. **Do not freeze.**

Serum and Plasma


They must be separated from the blood interface after centrifugation. To obtain serum, wait for the blood to clot before centrifuging the sample. To obtain plasma, collect it using EDTA, Heparin or Citrate as an anticoagulant. Plasma or serum can be stored for up to 3 days at 2 to 8 °C and 6 months at -20 °C.

PROCESS DESCRIPTION

- 1-** The sample must be at a temperature between 15 and 30°C before starting the test.
2- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface, and identify it appropriately.
3- Collect the sample using the disposable plastic pipette provided in the kit.
4- Press the pipette bulb and touch it to the sample. Relieve the pressure on the bulb so that the sample can be aspirated. Aspirate up to the trace to obtain a volume of 20 uL.
5- To dispense the sample, press the bulb again, touching the lower part of the pipette to the test membrane in the sample application hole.
Note: Alternatively, an automatic pipette can be used to collect and add the sample.
6- Hold the Diluent bottle (Reagent No. 2) vertically and apply 2 drops to the sample application hole (70 µL).
7- Wait for the lines to form. Interpret the results within 15 to 20 minutes. **Do not interpret after 20 minutes.**

INTERPRETATION OF THE RESULT

Reagent Test: Formation of a red line in the test region (T) and another line in the control region (C) in the first 15 to 20 minutes. Do not interpret after 20 minutes.	
Non-Reactive Test: Formation of a red line in the control region (C) and complete absence of a red line in the Test region (T). Do not interpret after 20 minutes.	
Invalid Result: Complete absence of line in the control region (C) with or without a red line in the test region (T). Test the sample again.	 

Notes:
1- The line in the control region may appear before 15 minutes of test incubation, and this does not mean that results can be interpreted ahead of time.
2- The results should not be interpreted after 20 minutes to avoid incorrect results.
3- The results provided by this kit must be interpreted by the responsible veterinary professional and are not the only criteria for determining the diagnosis and/or treatment of the patient.

PROCESS LIMITATIONS

1- This product is intended for screening *Dirofilaria immitis* antigens.
2- It is important to use the correct sample volume, as lower or higher volumes may cause erroneous results.
3- A negative result may be observed if the amount of analyte present in the sample is lower than the detection limit of the test. This can happen, depending on the period of the disease in which the sample was collected, the type of sample used in the test and its storage conditions.
4- For the correct diagnosis, the professional must take into account not only the result of this test, but also the clinical signs, health conditions of the animal and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls that allow the precision and accuracy of dosages to be assessed.

PRODUCT PERFORMANCE

PRECISION
Repeatability
Repeatability was calculated from 10 successive determinations, using 3 samples. The following results were obtained:

Repeat-ability	N° of Repeti-tions	Expected outcome	Result Found
Sample 1	10	Positive	Positive
Sample 2	10	Positive	Positive
Sample 3	10	Negative	Negative

Reproducibility
Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples. The following results were obtained:

Reproduc-ibility	N° of Days	N° of Rep-itations	Expected outcome	Result Found
Sample 1	3	10	Positive	Positive
Sample 2	3	10	Positive	Positive
Sample 3	3	10	Negative	Negative

CLINICAL SENSITIVITY AND SPECIFICITY
In HEARTWORM Ag VET FAST, clinical samples were analyzed in comparison with another commercially available immunochromatography method. The results show that the clinical sensitivity of the kit is 98.7% and the clinical specificity is 99.1%.

METHOD	REFERENCE			Total
	Result	Positive	Negative	
	Positive	75	1	76
	Negative	1	105	106
Result Total		76	106	182

Clinical Sensitivity: 98.7% (75/76) - IC 95%: 92.9% - 99.8%
Clinical Specificity: 99.1% (105/106) - IC 95%: 94.9% - 99.8%
Precision: 98.9% [(75+105) / (76+106)]

CLINICAL SIGNIFICANCE

Heartworm disease is a parasitic zoonosis, caused by nematodes of the genus *Dirofilaria*, which affect dogs and cats, humans and several species of wild animals. Transmission occurs through mosquito bites. The first stage larvae are acquired by the mosquito after a blood meal, and development occurs inside them until they reach the infective larval stage (L3), when they are inoculated into the mammalian host. The larvae then reach the bloodstream and reach the heart through the veins and arteries. The main species that affects dogs is *Dirofilaria immitis*, and its infection is associated with eosinophilia and neutrophilia, which accumulate in the blood vessels and kidneys causing inflammation. Its presence in the heart chambers and blood vessels promotes cardiovascular and pulmonary changes. Clinical symptoms of cough, dyspnea, heart and lung sounds, syncope and intense apathy. Diagnosis is generally made based on clinical signs, however the high frequency of asymptomatic animals may compromise the identification of infected dogs. Thus, the detection of circulating *Dirofilaria* antigen in blood, serum or plasma samples can favor early diagnosis, even in animals without apparent signs of the disease.

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QUALITY GUARANTEE
Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE
Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product Licensed at the Ministry of Agriculture, Livestock and Supply since 09/04/2022 under number 10.528/2022

Technical Responsible: Dr. Camila Eckstein (CRMV/ MG 20.611)

Revision: February/2024

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER