Bioclin

FATOR REUMATÓIDE

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method of quantitative determination of Rheumatoid Factor (RF) in biological samples of serum or plasma. Immunoturbidimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

The reaction permits you to quantify, using the turbidimetric method, the concentration of RF present in sample. The suspended latex particles covered with human gamma globulin are agglutinated by Rheumatoid Factors in serum. The agglutination process causes an increase in optical density. The intensity of scattered light is proportional to the concentration of RF.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.

Number 2 - RF Latex - Store between 2 and 8°C. Contains: Suspension of Latex particles sensitized with human Gamma Globulin, stabilizers and preservative.

Number 3 - Calibrator - Store between 2 and 8°C. Contains: Rheumatoid Factor, stabilizers and preservative. **Potentially infectious.**

Warning: The concentration of Rheumatic Factor varies according to the lot - see the bottle label.

PRESENTATION

Presentation	Reagent		
	Nº 1	Nº 2	Nº 3
1	18 mL	2 mL	1.0 mL
2	45 mL	5 mL	1.0 mL
3	2 x 45 mL	2 x 5 mL	2 x 1.0 mL
4	20 mL	5 mL	1.0 mL
5	2 x 20 mL	2 x 5 mL	1.0 mL
6	3 x 20 mL	3 x 5 mL	2 x 1.0 mL
7	4 x 20 mL	4 x 5 mL	2 x 1.0 mL
8	40 mL	10 mL	1.0 mL
9	2 x 40 mL	2 x 10 mL	2 x 1.0 mL
10	3 x 40 mL	3 x 10 mL	3 x 1.0 mL
11	4 x 40 mL	4 x 10 mL	4 x 1.0 mL
12	60 mL	15 mL	2 x 1.0 mL
13	2 x 60 mL	2 x 15 mL	3 x 1.0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagent Nº 3 should be handled cautiously since it is subject to biological contamination.

6- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.

7- The reagent has been tested for HIV and HBsAg using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma obtained free of hemolysis and severe lipemia. The analyte is stable between 2 and 8°C up to 2 days and 3 months at -20 °C¹¹.

Samples with presence of fibrin should be centrifuged before testing.

PROCESS DESCRIPTION

The calibration stability of the Rheumatoid Factor kit installed on refrigerated equipment is 10 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

The reagents are ready to use. Thoroughly mix the Reagent Nº 2 before use.

For calibration of the reaction, utilize preferentially the kit's internal calibrator or the Bioclin Multical kit. For control use the Bioclin Multicontrol Kit.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

PROCEDURE LIMITATIONS

The performance of the product depends directly on the performance of the automated equipment used. Therefore, to ensure best results of the test is very important to check periodically the optical and pipetting system of the equipment.

INTERFERENT

No interference was observed for Bilirubin up to 20 mg/dL, Hemoglobin up to 10 g/L and Lipids up to 10 g/L.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 64/002 (Rheumatoid Arthritis Serum).

REFERENCE VALUES

The reference values for this method were obtained through determination of Rheumatoid Factor in healthy male and females populations. Less than 20 UI/mL.

These values should be used as guidance and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 95 and 107%.

Comparison of Methods

The Rheumatoid Factor t kit was compared with other commercially available methods for determination of Rheumatic Factor available. 29 analyses were performed and the results were evaluated. We obtained the linear equation $Y = 0.941X + 1.360$ and correlation coefficient 0.995. With these results we can conclude that the kit good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	48.57	130.25	75.51
Standard Deviation (UI/mL)	0.17	0.77	0.27
Coefficient of Variation (%)	0.34	0.59	0.36

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	46.00	126.44	71.93
Standard Deviation (UI/mL)	2.39	3.98	3.59
Coefficient of Variation (%)	5.20	3.15	4.99

SENSITIVITY

Sensitivity was calculated from 10 determinations of a sample Rheumatic Factor. The average found was 0.711 UI/mL with a standard deviation of 1.513 UI/mL. The sensitivity, which indicates the method detection limit, corresponds to the average plus 3 times the value of the standard deviation and is equal to 5.250 UI/mL.

LINEARITY

The reaction is linear until the highest point in the calibrator curve. For samples with higher values than it is recommended to dilute the sample with Sodium Chloride 0.85% repeat the dosage and multiply the result by the dilution factor.

HIGH DOSE HOOK EFFECT

Hook effect wasn't observed with a high concentration of Rheumatoid Factor up to 800 UI/mL.

DIAGNOSTIC SIGNIFICANCE

Rheumatoid factor refers to a group of macroglobulins (antiglobulin) that react with the Fc fragment of immunoglobulins IgG. In rheumatoid arthritis the finding of Rheumatoid Factor represents the most important serological data, present in about 75% of patients. Its plasma level can be significantly increased in old age, connective tissue diseases, liver diseases syphilis, tuberculosis, leprosy, bacterial endocarditis, mononucleosis, sarcoidosis, leishmaniasis, rubella, cancer, parasitic infestations, blood transfusions, renal transplant and Sjogren syndrome.

BIBLIOGRAPHIC REFERENCES

- 1 - PESCE, A. J.; KAPLAN, L. A.: Methods in Clinical Chemistry the C. V. Mosby Company, 1987.
- 2 - HELLSING, K., Profides in the Biological Fluids, 1973, 23, 579.
- 3 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976, 22,657.
- 4 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980;26, 1459.
- 5 - GALVIN, J. P. et al; Clin. Lab Assays 73, 4th, 1983.
- 6 - PASSING, H.; BABLOCK, W., J., Clin. Chem. Clin. Biochem 21, 709; 1983.
- 7 - NOKAMA, Y.; NAKAMURA, R. M., J. Clin. Lab. Anal., 1:15; 1987.
- 8 - OKOMURA, J. M.; MIYAGI, J. M.; TERADA, K. J., Clin. Lab. Anal. 4, 231; 1990.
- 9 - MOORE, T. L., Clin. Biochem, 26; 75-84; 1993.
- 10 - SINGER, J. M. Et al Amer J Med 21, 888; 1952.
- 11 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 12 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | Email: sac@bioclin.com.br

Rheumatic Factor kit registration number at ANVISA: 10269360156

Review: August/2022

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER