Fe AUTOMAÇÃO

REF K308

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FINALIDADE

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PRECAUCIÓN

PELIGRO

FUNCTION
Method for determination of Iron Serum in human samples of serum and plasma (Heparin). Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION
Methodology: Modified Goodwin/Ferrozine
Iron is released from Transferrin in an acid enviroment and lowered to its ferrous state by the action of Hydroxylamine. Subsequently, reacts with Ferrozine leading to the formation of a violet colored complex.

REAGENTS
Number 1 - Reduction Buffer - Store between 2 and 8°C. Contains: Sodium Acetate < 900 mmol/L, Hydroxylamine < 800 mmol/L, surfactants and stabilizer.
Number 2 - Color Reagent - Store between 2 and 8°C. Contains: Ferrozine < 25 mmol/L, Hydroxylamine < 20 mmol/L and stabilizers.
Number 3 - Standard* - Store between 2 and 8°C. Contains: Iron 200 µg/dL, stabilizer and preservative.

PRESENTATION

Reagent N° 1	Reagent N° 2	Reagent N° 3*
1 x 40 mL	1 x 10 mL	1 x 3 mL
2 x 40 mL	2 x 10 mL	1 x 3 mL
4 x 40 mL	4 x 10 mL	1 x 3 mL
6 x 40 mL	6 x 10 mL	1 x 3 mL
8 x 40 mL	8 x 10 mL	1 x 3 mL
2 x 40 mL	1 x 20 mL	1 x 3 mL
4 x 40 mL	2 x 20 mL	1 x 3 mL
1 x 20 mL	1 x 5 mL	1 x 3 mL
2 x 20 mL	2 x 5 mL	1 x 3 mL
3 x 20 mL	3 x 5 mL	1 x 3 mL
4 x 20 mL	4 x 5 mL	1 x 3 mL
5 x 20 mL	5 x 5 mL	1 x 3 mL

*The R3 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS
Biochemistry automated equipment or semi-automated, water bath, pipettes, test tubes, Biocontrol N Bioclin and Biocontrol P Bioclin. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS
Store temperatures must be between 2 and 8°C. The transport at temperatures between 15 and 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE
1- For *in vitro* diagnostic use only.
2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.
4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
5- Hemolysis, even if discrete, can interfere with dosage.
6- Do not use plasma.
7- Every material used for analyses must be previously immersed in a recent HCL (7%) solution for 30 minutes. Rinse it with deionized water.
8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
10- Do not use the product in case of damaged packaging.
11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES
Serum and plasma (Heparin) obtained free of hemolysis, stable for 21 days between 2 and 8°C and 1 year at -20°C. ⁷

PROCESS DESCRIPTION
The calibration stability of the Fe Automated kit installed on refrigerated equipment is up to 19 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE
Bioclin recommends, as control serum, Biocontrol N Bioclin and Biocontrol P Bioclin Kit.

Mark 3 test tubes with letters B (Blank), P (Standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Reagent N°1	0.8 mL	0.8 mL	0.8 mL
Distilled Water	0.1 mL	----	----
Sample	----	----	0.1 mL
Reagent N°3	----	0.1 mL	----

Smoothly homogenize and determine the absorbances of the Sample and Standard at 560 nm (540 - 580 nm), hitting the zero with the distilled water. This will be the absorbance A₁.

And then add:

	Blank	Standard	Sample
Reagent N°2	0.2 mL	0.2 mL	0.2 mL

Smoothly homogenize and incubate in water-bath at 37°C, per 5 minutes. Determinate the absorbances of Sample and the Standard in 560 nm (540 - 580 nm), matching the zero with the Blank. This will be the absorbance of the sample A₂.

CALCULATIONS
Due to the dilution, the A₁ reading of the Standard and Sample should be corrected for the final volume, using the factor of 0.82. Thus, A_{1CORRECTED} is obtained.

Standard-A_{1CORRECTED} = Standard A₁ x 0.82
Sample-A_{1CORRECTED} = Sample A₁ X 0.82

Iron Serum (µg/dL) =
$$\frac{\text{Sample (A}_2 - \text{A}_{1\text{CORRECTED}}) \times 200}{\text{Standard (A}_2 - \text{A}_{1\text{CORRECTED}})}$$

The color reaction follows the Lambert Beer Law, it can use the Calibration Factor:

Calibration Factor =
$$\frac{200}{\text{Standard (A}_2 - \text{A}_{1\text{CORRECTED}})}$$

Iron Serum (µg/dL) =
$$\frac{\text{Sample (A}_2 - \text{A}_{1\text{CORRECTED}}) \times \text{Calibration Factor}}{\text{Factor}}$$

The results are expressed in µg/dL.

PROCEDURE LIMITATIONS
The use of water contaminated by Iron increase white color intensity, changing the results.

INTERFERENCE
Triglycerides up to 1000 mg/dL, Bilirubin up to 20 mg/dL and Copper up to 7.5 mg/dL did not interfere with the results. Some drugs can raise levels of Serum Iron such as Chloramphenicol, Estrogens, Ethanol, Methotrexate and Anovulatory. Other drugs may decrease the levels of Serum Iron such as: Acetylsalicylic Acid (high doses), Allopurinol, Cortisone and Corticotrophin. Samples of patients on hemodialysis, using Heparin, or with presence of fibrinogen are not interfered with in the analysis.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY
The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 937.

REFERENCE VALUES
The reference values in µg/dL, for this method, were obtained through the determination of Iron in healthy populations of male and female.

Iron Serum (µg/dL)		
Newborn		100 - 250
Infant		40 - 100
Children		50 - 120
Adult	Male	65 - 170
	Female	50 - 170

To convert the values from µg/dL to mmol/L, multiply by 0.179.
These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE ACCURACY
Comparison of Methods and Methodological Specificity
The Fe Automated kit was compared with other method for determination of Iron dosage commercially available. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0.9744X + 1.6825$ and correlation coefficient 0.999. With these results we can conclude that the kit shows good methodological specificity.

PRECISION
Repeatability
The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration Promedio (µg/dL)	72.90	222.73	27.90
Standard Deviation (µg/dL)	1.43	0.93	0.63
Coefficiente de Variación (%)	1.96	0.42	2.27

Reproducibility
The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Concentration Promedio (µg/dL)	80.48	232.46	30.46
Standard Deviation (µg/dL)	6.57	8.43	2.22
Coefficiente de Variación (%)	8.16	3.63	7.28

SENSITIVITY
The sensitivity was calculated from 40 determinations from a sample free of iron. The average found, which indicates the detection limit, was 1.48 µg/dL with standard deviation of 0.51 µg/dL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 2.99 µg/dL.

LINEARITY
The reaction is linear up to 1500 µg/dL. For larger values than 1500 µg/dL, dilute the serum with distilled water, repeat dosage and multiply the result by a factor of dilution.

DIAGNOSTIC SIGNIFICANCE
Iron is the most important transition metal and abundant in the human body. It is essential to the synthesis of hemoglobin and myoglobin, besides being present in several biochemical reactions of metabolism. The Serum Iron is elevated in sideroblastic anemia, hemolytic anemias, hereditary anemias (Thalassemia Major), acute hepatitis, chronic liver diseases, hemochromatosis and massive blood transfusions.
Low Serum Iron levels are found in chronic infection diseases, neoplasms, parasites, menstrual period, large bleeding, pregnancy and dietary deficiencies.

BIBLIOGRAPHIC REFERENCES
1 - PETERS T., et coll.: J. Lab and Clin Med. 48:280, 1956.
2 - HENRY, R. J., Cannon. D. C. and Winkleman, J. W., eds.: Clinical Chemistry: Principles and Technics, Haper and Row Publishers, 1974.
3 - ZAK, B.: Clin Chim. Acta 3: 328, 1958.
4 - BRITTENHAM, G.: Clin Chim. Acta. 91: 203, 1979.
5 - TONKS, D. B.: Chim. Acta 2: 393, 1983.
6 - BURTIS, C.A., ASHWOOD, E.R., Tietz Textbook of Clinical Chemistry, 2ª Edição, Philadelphia: W.B. Saunders, 1986:2175-2211.
7- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
8- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Fe Automated Family kit: 10269360375

Review: April/2025

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 LOT NUMBER

 MANUFACTURING DATE

 VALIDITY DATE
(last day of the month)

 TEMPERATURE LIMIT
(store)

 CONTENT IS SUFFICIENT
FOR <N> TEST

 SEE INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC
PRODUCT

 KEEP AWAY
FROM SUNLIGHT

 DO NOT REUSE

 CAUTION

 MADE BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 FLAMMABLE

 CORROSIVE

 TOXIC

 DO NOT USE IF
PACKAGE IS
DAMAGED

 PRODUCT
STERILIZED

 DANGER