FIV Ac / FeLV Ag COMBO VET FAST

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FINALIDADE

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REF VET013

INSTRUCTIONS FOR USE

FUNCTION
Fast combined immunochromatographic test for the qualitative determination of antibodies against the virus of feline immunodeficiency (FIV) and feline leukemia virus antigen (FeLV) in feline serum, plasma or whole blood samples. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION
Methodology: Immunochromatography.
The FIV Ac / FeLV Ag Combo VET FAST is an *in vitro* diagnostic test for the rapid and qualitative detection of total antibodies against feline immunodeficiency virus (FIV Ac) and feline leukemia virus antigen (FeLV Ag) in feline blood, serum or plasma.

The test consists of two independent cassettes, linked by a support (combo cassette): one detects anti-FIV antibodies and the other detects the FeLV antigen. Each cassette consists of a nitrocellulose membrane, impregnated in the test line with the feline immunodeficiency virus antigen (FIV cassette) and antibodies against feline leukemia virus (FeLV cassette). The conjugate consists of FIV antigen (FIV cassette) and anti-FeLV antibody (FeLV cassette); both conjugated with colloidal gold particles.

When the sample is added to the cassettes, it migrates by capillary action through the test strip until it reaches the conjugate. The following occurs in this region:

In the FIV cassette: Anti-FIV antibody, if present in the sample, interact with the conjugate, formed of colloidal gold conjugated FIV antigen. Therefore, formation of the complex occurs: anti-FIV antibody – FIV antigen with colloidal gold. The sample continues to migrate through the test strip until it reaches the test line, where the anti-FIV antibody also interacts with the immobilized FIV antigen in the test line, forming the complex: conjugate - anti-FIV antibody - antigen immobilized. This shows the formation of a pink-red line on the test line; indicating positive result.

In the FeLV cassette: the FeLV antigen, possibly present in the sample, interacts with the conjugate, formed of colloidal gold conjugated anti-FeLV antibodies. This is how complex formation occurs: FeLV antigen - colloidal gold conjugated anti-FeLV antibody. The sample continues to migrate through the test strip until it reaches the test line, where the FeLV antigen also interacts with the anti-FeLV antibody immobilized in this area, forming the complex: conjugate - FeLV antigen - immobilized anti-FeLV antibody. This shows the formation of a pink-red line on the test line; indicating positive result.

In both cassettes, sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS
Number 1 (R1): Cassette - Store between 15 and 30°C. **Do not freeze.** It contains:

- A) FIV cassette:**
01 plastic cassette with 01 test strip containing:
- Control Region (C): Anti-conjugate antibody.
- Test Region (T): Immobilized FIV antigen.
- Conjugate: FIV antigen, conjugate to colloidal gold particles.

- B) FeLV cassette:**
01 plastic cassette with 01 test strip containing:
- Control Region (C): Anti-conjugate antibody.
- Test Region (T): Immobilized monoclonal antibody against feline leukemia virus.
- Conjugate: Anti-FeLV antibody, conjugate to colloidal gold particles.

C) 01 Sachet of Silica

Number 2 (R2): Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

- Collection Accessories:**
- Dispisable Plastic Pipette
- Collection tube with EDTA

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET013-1	10	1 x 2,5 mL
	VET013-2	25	1 x 5,0 mL
	VET013-3*	10	10 x 0,5 mL
	VET013-4*	25	25 x 0,5 mL
	VET013-5	5	1X1,25 mL

*Individual presentations


EQUIPMENTS AND OPERATIONAL INPUTS
Clock or stopwatch. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS
The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.**
The combo cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

- SPECIAL CARE**
1- Product for *in vitro* diagnosis, for veterinary use only.
2- Strictly follow the methodology proposed for obtaining results.
3- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
4- Use clean, dry containers to package biological samples.
5- Open the aluminum envelop to remove the combo cassette only at the time of use.
6- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.
7- The samples should reach room temperature before starting the test.
8- Do not use reagents from different lots or kits.
9- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
10- The combo cassette (test devices) and the plastic pipette are disposable and cannot be reused.
11- The EDTA collector tube accessory must be used exclusively with this product.
12- Blood samples with clots, microcoagulants and fibrin mesh formation can cause the filter mesh to clog, making it impossible for the sample to pass through the test membrane correctly.
13- The kit materials that had contact with the sample of the patient must have a special destination. The disposal of plastic pipette, the cassettes (test devices) and diluent vials should be disposed of in organic waste.
14- Secondary packaging (box), Instruction for use, envelope aluminum and silica sachet can be disposed of in trash common.
15- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.
16- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES
Perform the test with serum, plasma or whole blood.
Serum: Allow blood to clot and then centrifuge to separate serum.
Plasma: Collect whole blood in a tube containing Heparin, Citrate, or EDTA. Centrifuge the blood and separate the plasma.
Whole Blood: Collect the blood sample in a tube containing EDTA, Heparin, or Citrate using standard laboratory procedures.
Whole blood samples should be analyzed within 3 days, can be stored at 2 and 8°C and never be frozen.
Serum and plasma samples should be analyzed immediately after collection, but if necessary, they can be stored at 2 and 8°C for up to 7 day sor at -20°C for longer periods.

- DESCRIPTION OF PROCESS**
1- The sample should be at room temperature before starting the test.
2- Remove the combo cassette (test devices) from the protective packaging and place it on a clean, dry and level surface.
3- Collect the sample using the disposable plastic pipette that accompanies the kit.
4- Press the pipette bulb and touch it to the sample. Relieve the pressure on the bulb so that the sample can be aspirated. Aspirate up to the trace to obtain a volume of 20 µL. To dispense the sample, press the bulb again, touching the lower part of the pipette to the test membrane in the sample application hole. Perform this procedure on both test cassettes (FIV Ac and FeLV Ag).
Note: Alternatively, an automatic pipette can be used for sample collection and addition.
5- Hold the diluent bottle (Reagent Nº 2) in a vertical position and apply 2 drops (70 µL) in the sample application hole of each combo device (2 drops in FIV Ac and 2 drops in FeLV Ag).
6- Wait for the lines to form. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

Reagent Test: Formation of two lines: one in the control region (C) and one in the test region (T). Any color intensity (light or dark) means the result is reactive. Do not interpret after 20 minutes.	
No Reactive Test: Line formation only in control region (C). Do not interpret after 20 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

Note: Each cassette must be interpreted individually.

TEST LIMITATIONS
This product is intended for detection of infection by feline immunodeficiency virus (FIV - detection of Antivirus Antibodies) and by feline leukemia virus (FeLV - detection of viral Antigen). A no reactive result can be observed if the amount of analyte present in the sample is below the detection limit of the test. This may happen depending on the period of the disease in which the sample was collected, the type of sample used in the test or the storage conditions. For FIV, negative reactions can occur especially in the first few weeks after infection or in immunosuppressed animals. Cats in the regressive phase of FeLV do not have the virus circulating and present negative results when detecting the viral antigen. In these cases, consult the responsible veterinarian and consider repeating the test after a few weeks. The definitive diagnosis of FIV and FeLV requires a comprehensive assessment of several factors by the veterinarian.

Isolated results should not be used as a criterion for definitive diagnosis, requiring consideration of factors such as the animal's health condition at the time of sample collection (asymptomatic or symptomatic feline) and the use of complementary diagnostic tools that may be available for diagnosis. confirmation of results (AAFP, 2020).

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 102 samples of FIV Ac were analyzed, which 41 samples were diagnosed as reagent and 61 were diagnosed as no reactive. The results show that the clinical sensitivity of the product is 97.56% and the clinical specificity is 96.72%.

METHOD		REFERENCE		Total
FIV Ac Cassette	Result	Reagent	No Reactive	
	Reagent	40	2	
	No Reactive	1	59	
Total Result		41	61	102

Clinical Sensivity:(40/41): 97.56% (IC 95% 92.9 - 100%)
Clinical Specificity: (59/61): 96.72% (IC 95% 92.3 - 100%)
Precision: (99/102): 97.05%

A total of 104 samples of FeLV Ag were analyzed as well, which 46 samples were diagnosed as reactive and 58 were diagnosed as no reactive. The results show that the clinical sensitivity of the product is >99% and the clinical specificity is 98.27%.

METHOD		REFERENCE		Total
FeLV Ag Cassette	Result	Reagent	No Reactive	
	Reagent	46	1	
	No Reactive	0	57	
Total Result		46	58	104

Clinical Sensivity: (46/46): >99.00% (IC 95% 95% - 100%)
Clinical Specificity: (57/58): 98.27% (IC 95% 94.87% - 100%)
Precision: (103/104): 99.04%

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 2 different samples of FIV Ac and 2 different samples of FeLV Ag. The following results was obtained:

Samples	N° of Repetitions	Expected Result	Found Result
FIV Ac Sample 1	10	Reagent	Reagent
FIV Ac Sample 2	10	Reagent	Reagent
FeLV Ag Sample 1	10	Reagent	Reagent
FeLV Ag Sample 2	10	Reagent	Reagent

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples of FIV Ac and 2 different samples of FeLV Ag. The following results were selected:

Samples	N° of Days	N° de Repetitions	Expected Result	Found Result
FIV Ac Sample 1	03	10	Reagent	Reagent
FIV Ac Sample 2	03	10	Reagent	Reagent
FeLV Ag Sample 1	03	10	Reagent	Reagent
FeLV Ag Sample 2	03	10	Reagent	Reagent

CLINICAL SIGNIFICANCE

The Feline Immunodeficiency Virus (FIV) and the Feline Leukemia Virus (FeLV) are the cause of important diseases that reduce the life expectancy of cats. FIV transmission occurs mainly through bites and scratches, while FeLV transmission occurs through contact with secretions such as saliva and sharing objects.

In FIV infection, the main clinical signs observed are fever and generalized lymphadenopathy, chronic opportunistic infections with a high rate of recurrence, neurological alterations and tumors. After infection, there is a significant and persistent increase in antibodies against the virus, which makes antibody detection a good diagnostic tool. In FeLV infection, the appearance of tumors is common and animals often become susceptible to secondary infections with low efficacy in its treatment. The production of antibodies against FeLV is deficient and nonspecific, while the circulation of viral antigen is high and persistent, making detection of viral antigen a safe and effective tool for diagnosis.

It is common for asymptomatic cats to occur FIV and FeLV, so the diagnosis is essential for the application of the appropriate behavior with the animal and for the control of the spread of these diseases.

BIBLIOGRAPHIC REFERENCES

- 1- QUIBASA: Data from the Research and Development Department.
- 2- MARTINS ES, TOGNOLI GK, ILHA PHO, et al. Prevalência de imunodeficiência viral felina e leucemia viral felina no Distrito Federal. Archives of Veterinary Science, v.17, n.S1, p.274-276, 2012.
- 3- HARBOUR DA, CANEY SMA, SPARKES AH. Feline Immunodeficiency Virus Infection. p. 607-622. Cap.24. IN: Chandler EA, Gaskeç CJ, Kaskel RM. Feline Medicine and Therapeutics. Blackwell: Wiley. 2004.
- 4- BIEZUS G, MACHADO G, FERIAN PE, COSTA UM, et al. Prevalence of and factors associated with feline leukemia virus (FeLV) and feline immunodeficiency virus (FIV) in cats of the state of Santa Catarina, Brazil. Comparative Immunology, Microbiology and Infectious Disease. v.63, p.17-21, 2019.
- 5- AAFP (American Association of Feline Practioners). Feline Retrovirus testing and Management and Guidelines. Journal of Feline Medicine and Surgery, v.22, p. 5-30, 2020.

QUALITY GUARANTEE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Advisory Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 11/11/2019 with the number 10.283/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: February/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER