## FOSFATASE ALCALINA IFCC

REF K224

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## FINALIDADE

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## REFERENCIAS BIBLIOGRAFICAS

1- SCHUMANN, Gerhard; KLAUKE, Rainer; CANALIAS, Francesca; BOSSERT-REUTHER, Steffen; FRANCK, Paul F.h.; GELLA, F.-javier; JØRGENSEN, Poul J.; KANG, Dongchon; LESSINGER, Jean-marc; PANTEGHINI, Mauro. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 9: reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clinical Chemistry And Laboratory Medicine, [s.l.], v. 49, n. 9, p. 1439-1446, 1 jan. 2011. Walter de Gruyter GmbH. <http://dx.doi.org/10.1515/cclm.2011.621>.

2- GARDNER, M. D. et al – Frequency distribution and 'reference values' of plasma alkaline phosphatase (EC3.1.3.1) activity in the adult population of Scottish new town – J. Clin. Pathol. 31/12:1202 (1978).

3- Thomas, L; et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in seru. J Lab Med 2005;29(5):301–308

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### INSTRUCTIONS FOR USE

#### FUNCTION

Method for determination of Alkaline Phosphatase. In serum or plasma (Heparin) samples. Kinetic test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Kinetic IFCC.

The Alkaline Phosphatase catalyzes the transfer of Phosphate group from the substrate p-Nitrophenylphosphate (pNPP) to 2-Amino-2-Methyl-1-Propanol (AMP), releasing p-Nitrophenol according to the equation below:

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: 2-Amino 2-Methyl Propanol < 1.0 mol/L, activators, stabilizer and preservative.

**Number 2 - Substrate** - Store between 2 and 8°C. Contains: p-NFF Substrate p-Nitrophenylphosphate.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 40 mL	1 x 10 mL
2	2 x 40 mL	2 x 10 mL
3	3 x 40 mL	3 x 10 mL
4	4 x 40 mL	4 x 10 mL
5	2 x 40 mL	1 x 20 mL
6	4 x 40 mL	2 x 20 mL
7	6 x 36 mL	2 x 12 mL
8	1 x 20 mL	1 x 5 mL
9	2 x 20 mL	2 x 5 mL
10	3 x 20 mL	3 x 5 mL
11	4 x 20 mL	4 x 5 mL
12	1 x 60 mL	1 x 15 mL
13	2 x 60 mL	2 x 15 mL
14	3 x 60 mL	3 x 15 mL
15	4 x 60 mL	4 x 15 mL
16	1 x 120 mL	2 x 15 mL
17	2 x 120 mL	4 X 15 mL
18	1 x 120 mL	1 x 30 mL
19	2 x 120 mL	2 x 30 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, atch or stopwatch, test tubes, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin They can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.
- It is important, for the good development of the test, a rigorous control of time and temperature.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum or plasma (Heparin) free of hemolysis. The analyte is stable for 7 days if kept in temperature between 2 and 8°C and 30 days if kept at 10°C negative.

#### PROCESS DESCRIPTION

The calibration stability of the Alkaline Phosphatase IFCC kit installed on refrigerated equipment is up to 3 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### WORKING REAGENT PREPARATION

**The reagents are ready to use in bireagent format and on automated equipment that uses dedicated presentations.**

For use in monoreagent format, prepare the working reagent as follows:

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

Mix four (4) parts of Reagent N° 1 with one (1) part of Reagent N° 2. Working Reagent is stable during 14 days, if kept between 2 and 8°C. Store away from light.

#### REACTION CONDITIONS

It is indispensable conditions to use a thermostated cuvette at 37°C, 1 cm optical path and reading at 405 nm.

#### TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin

#### Monoreagent:

Kits. Add 20 µL of Sample to 1.0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37°C and wait for 1 minute. Make the initial Reading simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the mean differences in absorbance per minute (Abs/min.) and use it to calculate the result.

To calibrate the bireagent kit, Bioclin recommends using the BIOCAL Bioclin calibrator serum. We also recommend using the Bioclin Biocontrol N and Biocontrol P kits as control serums.

**Check the equipment's programming by calling Customer Service at 0800 031 5454.**

#### CALCULATIONS

ALP (U/L) = Abs/min. x 2757  
Results are expressed as U/L.

#### LIMITATIONS OF THE PROCEDURE

The kinetic method is based on the molar absorptivity. By this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions:  
Wavelength 405 nm  
Semi trajectory of the pass band 10 nm  
Stray light less than 0.5%  
1cm thermostated cuvette.

#### INTERFERENCES

Do not use hemolyzed or lipemic samples, as the results obtained may be altered. Triglycerides above 1800mg / dL lead to falsely high results.  
Bilirubin values up to 35mg / dL do not interfere with the reaction. Hemoglobin values up to 40mg/dL also do not interfere with the reaction.  
Collection tubes containing citrate, oxalate, EDTA and fluoride should be avoided, as these are inhibitors of alkaline phosphatase activity.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to

variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

## TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the p-Nitrophenol, or through the Biocal calibrator. Bioclin recommends the usage of the Biocal calibrator which is also traceable to the molar absorptivity of the p-Nitrophenol.

## REFERENCE VALUES

The reference values in U/L, for this method were obtained through the determination of Alkaline Phosphatase in healthy populations.

Adult (men): 40 - 130 U/L

Adult (women): 35 - 105 U/L

Children: 75 to 390 U/L

(Including teenager up to 15 years old)

These values should be used as reference, being that every laboratory must create their own range of values for reference, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional and are not the sole criterion for determining the patient's diagnosis and/or treatment.

## PRODUCT PERFORMANCE QUALITY CONTROL ACCURACY

### Comparison Of Methods And Specificity Methodology

The Alkaline Phosphatase IFCC was compared with other commercially available methods for measurement of Phosphatase. 42 analyzes were performed and the results were evaluated. The linear equation obtained was  $y = 1.1027x - 0.2651$  and the correlation coefficient 0.9999. With these results, we can conclude that the kit shows good methodological specificity.

## PRECISION

### Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
<b>Average Concentration (U/L)</b>	153.1	145.9	327.9
<b>Standard Deviation (U/L)</b>	1.7	1.8	2.6
<b>Coefficient of Variation (%)</b>	1.14	1.22	0.78

### Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
<b>Average Concentration (U/L)</b>	154.5	147.4	328.9
<b>Standard Deviation (U/L)</b>	1.7	1.4	0.9
<b>Coefficient of Variation (%)</b>	1.13	0.92	0.26

## SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Phosphatase. The average found by the method detection limit was 0.83 U/L. The sensitivity, that indicates the method quantification limit, corresponds to 3 times the standard deviation is equal to 2.48 U/L.

## LINEARITY

The reaction is linear up to the concentration of 1500 U/L. For sample with value above 1500 U/L, we recommend diluting the sample with NaCl 0.85%. Multiply the results obtained by the dilution factor.

## DIAGNOSTIC SIGNIFICANCE

Alkaline Phosphatase is present in high concentrations in bone, liver, intestine and placenta. Constitutes a useful indicator of liver and bone diseases associated with osteoblastic hyperactivity. The conditions under which Alkaline Phosphatase is increased are: the liver and biliary tract, metastasis liver and bone metastasis, in acromegaly, in hyperthyroidism, in rickets, infectious mononucleosis, Paget's disease (osteoblastic activity) and in physiologic bone growth. Decreased values of Alkaline Phosphatase may be observed in hypothyroidism, growth retardation in children, hypophosphatasia (inborn metabolic error) and severe malnutrition/placenta.

## BIBLIOGRAPHIC REFERENCES

- SCHUMANN, Gerhard; KLAUKE, Rainer; CANALIAS, Francesca; BOSSERT-REUTHER, Steffen; FRANCK, Paul F.h.; GELLA, F.-javier; JØRGENSEN, Poul J.; KANG, Dongchon; LESSINGER, Jean-marc; PANTEGHINI, Mauro. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 9: reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clinical Chemistry And Laboratory Medicine, [s.l.], v. 49, n. 9, p. 1439-1446, 1 jan. 2011. Walter de Gruyter GmbH. <http://dx.doi.org/10.1515/cclm.2011.621>.
- GARDNER, M. D. et al – Frequency distribution and 'reference values' of plasma alkaline phosphatase (EC3.1.3.1) activity in the adult population of Scottish new town – J. Clin. Pathol. 31/12:1202 (1978).
- Thomas, L; et al. Consensus of DGKL and VDHG for interim reference intervals on enzymes in serum. J Lab Med 2005;29(5):301–308
- Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

## QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

## QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP: 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 31 3439.5454 | E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)

CNPJ: 19.400.787/0001-07 - Made in Brazil

## CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for ALKALINE PHOSPHATASE

IFCC kit: 10269360333

Review: September/2025

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE (last day of the month)

TEMPERATURE LIMIT (store)

CONTENT IS SUFFICIENT FOR <N> TEST

SEE INSTRUCTIONS FOR USE

IN VITRO DIAGNOSTIC PRODUCT

KEEP AWAY FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF PACKAGE IS DAMAGED

PRODUCT STERILIZED

DANGER