FOSFATASE ALCALINA IFCC

[REF] K224

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FINALIDADE

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1 - SCHUMANN, Gerhard; KLAUKE, Rainer; CANALIAS, Francesca; BOSSERT-REUTHER, Steffen; FRANCK, Paul F.h.; GELLA, F.-javier; JØRGENSEN, Poul J.; KANG, Dongchon; LESSINGER, Jean-marc; PANTEGHINI, Mauro. *IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 9: reference procedure for the measurement of catalytic concentration of alkaline phosphatase*. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clinical Chemistry And Laboratory Medicine, [s.l.], v. 49, n. 9, p. 1439-1446, 1 jan. 2011. Walter de Gruyter GmbH. <http://dx.doi.org/10.1515/cclm.2011.621>.
2 - GARDNER, M. D. *et al – Frequency distribution and ‘reference values’ of plasma alkaline phosphatase* (EC3.1.3.1) activity in the adult population of Scottish new town – J. Clin. Pathol. 31/12:1202 (1978).
3 - Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

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1 - SCHUMANN, Gerhard; KLAUKE, Rainer; CANALIAS, Francesca; BOSSERT-REUTHER, Steffen; FRANCK, Paul F.h.; GELLA, F.-javier; JØRGENSEN, Poul J.; KANG, Dongchon; LESSINGER, Jean-marc; PANTEGHINI, Mauro. *IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 9: reference procedure for the measurement of catalytic concentration of alkaline phosphatase*. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clinical Chemistry And Laboratory Medicine, [s.l.], v. 49, n. 9, p. 1439-1446, 1 jan. 2011. Walter de Gruyter GmbH. <http://dx.doi.org/10.1515/cclm.2011.621>.
2 - GARDNER, M. D. et al – *Frequency distribution and 'reference values' of plasma alkaline phosphatase* (EC3.1.3.1) activity in the adult population of Scottish new town – J. Clin. Pathol. 31/12:1202 (1978).
3 – Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS


FUNCTION

Method for determination of Alkaline Phosphatase. In serum or plasma (Heparin) samples.Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Kinetic IFCC

The Alkaline Phosphatase catalyzes the transfer of Phosphate group from the substrate p-Nitrophenylphosphate (pNPP) to 2-Amino-2-Methyl-1-Propanol (AMP), releasing p-Nitrophenol according to the equation below:

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: 2-Amino 2-Methyl Propanol < 1.0 mol/L, activators, stabilizer and preservative.

Number 2 - Substrate - Store between 2 and 8°C. Contains: p-NFF Substrate p-Nitrophenylphosphate

PRESENTATION

Presentation	Reagent N°1	Reagent N°2
1	R1: 1 x 40 mL	R2: 1 x 10 mL
2	R1: 2 x 40 mL	R2: 2 x 10 mL
3	R1: 3 x 40 mL	R2: 3 x 10 mL
4	R1: 4 x 40 mL	R2: 4 x 10 mL
5	R1: 2 x 40 mL	R2: 1 x 20 mL
6	R1: 4 x 40 mL	R2: 2 x 20 mL
7	R1: 6 x 36 mL	R2: 2 x 12 mL
8	R1: 1 x 20 mL	R2: 1 x 5 mL
9	R1: 2 x 20 mL	R2: 2 x 5 mL
10	R1: 3 x 20 mL	R2: 3 x 5 mL
11	R1: 4 x 20 mL	R2: 4 x 5 mL
12	R1: 1 x 60 mL	R2: 1 x 15 mL
13	R1: 2 x 60 mL	R2: 2 x 15 mL
14	R1: 3 x 60 mL	R2: 3 x 15 mL
15	R1: 4 x 60 mL	R2: 4 x 15 mL
16	R1: 1 x 120 mL	R2: 2 x 15 mL
17	R1: 2 x 120 mL	R2: 4 X 15 mL
18	R1: 1 x 120 mL	R2: 1 x 30 mL
19	R1: 2 x 120 mL	R2: 2 x 30 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, atch or stopwatch, test tubes, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures between 15 and 30°C should not exceed 72 (seventy two) hours. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional in vitro diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain exact results.
- 3-** Water used in material cleaning must to be recent and free of contaminants.
- 4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers, that may alter the results significantly.
- 5-** It is important, for the good development of the test, a rigorous control of time and temperature.
- 6-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8-** Do not use the product in case of damaged packaging.
- 9-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma (Heparin) free of hemolysis. The analyte is stable for 07 days if kept in temperature between 2 and 8°C and 30 days if kept at 10°C negative.

PROCESS DESCRIPTION

The calibration stability of the Alkaline Phosphatase IFCC kit installed on refrigerated equipment is at least 3 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORKING REAGENT

Mix four (4) parts of Reagent N° 1 with one (1) part of Reagent N° 2. Working Reagent is stable during 14 days, if kept between 2 and 8°C. Store away from light.

REACTION CONDITIONS

It is indispensable conditions to use a thermostated cuvette at 37°C, 1 cm optical path and reading at 405 nm.

TECHNIQUE

Bioclin recommends using the kit Biocal Bioclin as calibrator and as control serum, Biocontrol N and P Bioclin Kits. Add 20 µL of Sample to 1.0 mL of Working Reagent, mix and transfer to a thermostated cuvette at 37°C and wait for 1 minute. Make the initial Reading simultaneously starting the stopwatch. Repeat readings after 1, 2 and 3 minutes. Calculate the mean differences in absorbance per minute (Abs/min.) and use it to calculate the result

CALCULATIONS

Bi reagent
ALP (U/L) = Abs/min. x 2757
Monoreagent
ALP (U/L) = Abs/min. x 3296
Results are expressed as U/L

LIMITATIONS OF THE PROCEDURE

The kinetic method is based on the molar absorptivity. By this reason, the readings must be conducted in a spectrophotometer that satisfies the following conditions: Wavelength 405 nm Semi trajectory of the pass band 10 nm Stray light less than 0.5% 1cm thermostated cuvette.

INTERFERENCES

Do not use hemolyzed or lipemic samples, as the results obtained may be altered. Triglycerides above 1800mg / dL lead to falsely high results. Bilirubin values up to 35mg / dL do not interfere with the reaction. Hemoglobin values up to 40mg/dL also do not interfere with the reaction. Collection tubes containing citrate, oxalate, EDTA and fluoride should be avoided, as these are inhibitors of alkaline phosphatase activity.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the p-Nitrophenol, or through the Biocal calibrator. Bioclin recommends the usage of the Biocal calibrator which is also traceable to the molar absorptivity of the p-Nitrophenol.

REFERENCE VALUES

The reference values in U/L, for this method were obtained through the determination of Alkaline Phosphatase in healthy populations.

Adults: 27 to 100 U/L

Children: 75 to 390 U/L

(Including teenager up to 15 years old)

These values should be use as reference, being that every laboratory must create their own range of values for reference, according to the population served.

PRODUCT PERFORMANCE QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND SPECIFICITY METHODOLOGY

The Alkaline Phosphatase IFCC was compared with other commercially available methods for measurement of Phosphatase. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $y = 1.1027x - 0.2651$ and the correlation coefficient 0.9999. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	153.1	145.9	327.9
Standard Deviation (U/L)	1.7	1.8	2.6
Coefficient of Variation (%)	1.14	1.22	0.78

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	154.5	147.4	328.9
Standard Deviation (U/L)	1.7	1.4	0.9
Coefficient of Variation (%)	1.13	0.92	0.26

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free the Phosphatase. The standard deviation was 0.83 U/L. The sensitivity, that indicates the method detection limit, corresponds to 3 times the standard deviation is equal to 2.48 U/L.

Linearity

The reaction is linear up to the concentration of 1500 U/L. For sample with value above 1500 U/L, we recomend diluting the sample with NaCl 0.85%. Multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Alkaline Phosphatase is present in high concentrations bone, liver, intestine and placenta. Constitutes a useful indicator of liver and bone diseases associated with osteoblastic hyperactivity. The conditions under which Alkaline Phosphatase is increased are: the liver and biliary tract, metastasis liver and bone metastasis, in acromegaly, in hyperthyroidism, in rickets, infectious mononudeosis, Paget's disease (osteoblastic activity) and in physiologic bone growth. Decreased values of Alkaline Phosphatase may be observed in hypothyroidism, growth retardation in children, hypophosphatasia (inborn metabolic error) and severe malnutritionplacenta.

BIBLIOGRAPHIC REFERENCES

1 - SCHUMANN, Gerhard; KLAUKE, Rainer; CANALIAS, Francesca; BOSSERT-REUTHER, Steffen; FRANCK, Paul F.h.; GELLA, F.-javier; JØRGENSEN, Poul J.; KANG, Dongchon; LESSINGER, Jean-marc; PANTEGHINI, Mauro. *IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 9: reference procedure for the measurement of catalytic concentration of alkaline phosphatase*. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clinical Chemistry And Laboratory Medicine, [s.l.], v. 49, n. 9, p. 1439-1446, 1 jan. 2011. Walter de Gruyter GmbH. <http://dx.doi.org/10.1515/cclm.2011.621>.

2 - GARDNER, M. D. et al – *Frecuency distribution and 'reference values' of plasma alkaline phosphatase* (EC3.1.3.1) activity in the adult population of Scottish new town – J. Clin. Pathol. 31/12:1202 (1978).

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

MANUFACTURER DATA

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454.
E-mail: sac@bioclin.com.br

ANVISA registration for Alkaline Phosphatase IFCC kit:
10269360333

Review: February/2021

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED