FÓSFORO UV

REF K309

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USAGE INSTRUCTIONS

FUNCTION

Method of determination of Phosphorus in biological samples of serum, plasma (EDTA o Citrate) and urine. UV ending point, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: UV Ending Point

The inorganic Phosphorus in acid reacts with the Molybdate, forming Phosphomolybdate, which intensity of color developed is proportional to the concentration of Phosphorus present in the sample.

REAGENTS

Number 1 - Molybdate - Store between 15 and 30°C. Contains: Ammonium Molybdate < 5,0 mmol/L, surfactant, activators, stabilizer and preservative.

Number 2 - Standard - Store between 15 and 30°C. Contains: Phosphorus 4mg/dL and preservatives.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2*
1	100 mL	3 mL
2	5 x 20 mL	3 mL
3	5 x 40 mL	3 mL
4	10 x 40 mL	3 mL

* The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter with 340 nm reading, watch or stopwatch, pipettes, test tubes, water bath, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperatures should be between up to 30°C, should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagent N°2 contain Sodium Azide which can be irritating to skin and mucous membranes. Handle with care.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma obtained with EDTA or Citrate and 24 hour urine. The sample is stable for 7 days if kept in temperature between 2 and 8°C. Serum or plasma must be separated at least an hour after collected, due to liberation of hematic phosphorus. Acidification with HCl maintains urine in refrigerator for 15 days.

PROCESS DESCRIPTION

The calibration stability of the Phosphorus UV kit installed on refrigerated equipment is 07 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), A (Sample), S (Standard):

	Blank	Standard	Sample
Sample	---	---	10 µL
Reagent Nº 2	---	10 µL	---
Reagent Nº 1	1,0 mL	1,0 mL	1,0 mL

Homogenize and place it in a water bath at 37°C for 5 minutes. Make the readings in 340 nm (334 - 360 nm), adjusting the zero with test tube "B" (Blank). Color reaction is stable for 20 minutes.

For urine dosage: Homogenize, separate 10 mL and adjust the pH between 1 and 3 with concentrated HCl. When performing the technique, dilute the urine 1:10 with

distilled or deionized water. After that perform the dosage according to the technique above described. Multiply the results found by 10.

CALCULATIONS

Phosphorus (mg/dL) = $\frac{\text{Sample Absorbance} \times 4}{\text{Standard Absorbance}}$

Since the reaction follows the Lambert-Beer law, the Calibration Factor can be used.

Calibration Factor = $\frac{\text{Standard Concentration (4 mg/dL)}}{\text{Standard Absorbance}}$

Phosphorus (mg/dL) = Sample Absorbance x Calibration Factor

24 hour urine

Urine (mg/24h) = $\frac{(\text{mg/dL}) \times \text{Urinary Volume (mL)}}{100}$

INTERFERENT

Hemolysis, lipemia or icterus are factors that may lead to falsely elevated results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 186G.

REFERENCE VALUES

The reference values in mg/dL, for this method, were obtained through the determination of phosphorus in healthy populations of male and female.

Serum	Up to 10 days	4,5 - 9,0 mg/dL
	10 days to 2 years	4,5 - 6,7 mg/dL
	2 to 12 years	4,5 - 5,5 mg/dL
	Over 12 years	2,5 - 4,8 mg/dL
	Adults	2,5 - 4,8 mg/dL
Urine	340 - 1000 mg/24h	

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

To convert the values from mg/dL to mmol/L (SI), multiply by 0.322.

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Comparison of Methods

The Phosphorus UV kit was compared with other methods for measurement of phosphorus commercially available. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0.9986X + 0.0074$, with correlation coefficient 0.9989. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION
Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	2.33	6.36	4.38
Standard Deviation (mg/dL)	0.07	0.11	0.12
Coefficient of Variation (%)	3.08	1.78	2.71

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	2.38	6.41	4.34
Standard Deviation (mg/dL)	0.05	0.05	0.04
Coefficient of Variation (%)	1.95	0.80	0.85

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Phosphorus. The average found, which indicates the detection limit, was 0.091 mg/dL, with standard deviation of 0.006 mg/dL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.109 mg/dL.

LINEARITY

The color of the reaction is linear up to the concentration of 15 mg/dL. For higher values, dilute the sample with distilled or deionized water, perform a new determination and multiply the value obtained by the dilution factor.

DIAGNOSTIC SPECIFICITY

Phosphorus homeostasis is maintained primarily by renal function, intestinal absorption, bone metabolism and function parathyroid gland. Phosphorus has a role in bone composition, the acid-base balance, metabolism of fats and carbohydrates, in addition to being presente the formation of compounds such as nucleic acids.

We can observe the increased values of phosphorus in the following conditions: renal failure (renal tuberculosis, chronic glomerulonephritis, hydronephrosis, pyelonephritis) hypoparathyroidism, excessive vitamin D.

The decreased level of phosphorus comes from many causes as: vitamin D deficiency (rickets, osteomalacia) intake of antacids, liver disease, syndrome Fanconi, hyperparathyroidism and use of thiazide diuretics.

BIBLIOGRAPHIC REFERENCES

1 - BAGINSKI, E. S.: Amer J. Med. Tec., 35: 475, 1969.
2 - GODWIN, J. F., Clin. Chem., 16: 776, 1970.
3 - TONKS, D. B., Clin. Chem., 9 – 217, 1963.
4 - CHRISTIAN, D. G.: Amer. J. Clin. Pathol. 54: 118, 1979.
5 - HENRY, R. J.: Clinical Chemistry - Principles and Technics, 2 Ed., Haper-Row, York, 1974.
6 - TODD; SANFORD; DAVIDSOHN, Diagnósticos Clínicos, 1983.
7 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP: 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Phosphorus UV kit: 10269360379

Review: April/2025

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 LOT NUMBER

 MANUFACTURING DATE

 VALIDITY DATE (last day of the month)

 TEMPERATURE LIMIT (store)

 CONTENT IS SUFFICIENT FOR <N> TEST

 SEE INSTRUCTIONS FOR USE

 IN VITRO DIAGNOSTIC PRODUCT

 KEEP AWAY FROM SUNLIGHT

 DO NOT REUSE

 CAUTION

 MADE BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 FLAMMABLE

 CORROSIVE

 TOXIC

 DO NOT USE IF PACKAGE IS DAMAGED

 PRODUCT STERILIZED

 DANGER