Bioclin

FRUTOSAMINA

REF **K135**

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FINALIDADE

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1. Armbuster DA, Fructosamine: Structure, Analysis and Clinical Usefulness. Clin. Chem. 1987; 33 (12): 2153-2163
2. Howey JEA, Browning MCK, Fraser CG. Assay of serum fructosamine that minimizes standardization and matrix problems: Use to acess components of biological variation. Clin. Chem. 1987; 33: 269-272
3. BURTIS, C.A., ASHWOOD, E.R., Tietz Clinical Chemistry, W. B. Saunders Company. 2nd edition, 1994; 986-988.
4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Method for the quantitative determination of Fructosamine in biological samples of serum. Colorimetric, for *In vitro* Diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric

The fructosamine present in the serum is formed from the link from the glucose to the ammine group of the protein forming the Schiff base, that after a molecular reorienting becomes ketoamine stable (fructosamine). Under alkaline conditions the fructosamine is converted to a enol shape, that reduces the nitrotetrazolium blue to a purple complex. The difference between the absorbances read in the 540 nm (530-550 nm) filter, is proportional to the fructosamine concentration.

REAGENTS

Number 1: Buffer - Store between 2 and 8°C. Contains: Carbonate Buffer < 20 mol/L, detergent and preservative.
Number 2: Color Reagent - Store between 2 and 8°C. Contains: Carbonate Buffer < 20 mol/L, Nitrotetrazolium Blue < 10 mmol/L and preservative.
Number 3: Calibrator - Store between 2 and 8°C. Contains: Fructosamine.
Warning: Concentration of Fructosamine may vary according to the lot – see bottle label.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2	Reagent Nº 3
1	15 mL	5 mL	2 mL
2	2 x 15 mL	2 x 5 mL	2 mL
3	3 x 15 mL	3 x 5 mL	2 mL
4	4 x 15 mL	4 x 5 mL	2 mL
5	30 mL	10 mL	2 mL
6	2 x 30 mL	2 x 10 mL	2 mL
7	3 x 30 mL	3 x 10 mL	2 mL
8	4 x 30 mL	4 x 10 mL	2 mL
9	60 mL	20 mL	2 mL
10	2 x 60 mL	2 x 20 mL	2 mL
11	3 x 60 mL	3 x 20 mL	2 mL
12	4 x 60 mL	4 x 20 mL	2 mL
13	45 mL	15 mL	2 mL
14	2 x 45 mL	2 x 15 mL	2 mL
15	3 x 45 mL	3 x 15 mL	2 mL
16	4 x 45 mL	4 x 15 mL	2 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.
Note: After opening the reagents are stable for 1 month.

SPECIAL CARE

- 1 - For *in vitro* diagnostic only.**
2 - Strictly follow the methodology proposed to obtain exact results.
3 - Water used in material cleaning must be recent and free of contaminants.
4 - Although the substrate contains preservative, must avoid bacterial contamination.
5 - Reagents Nº1 and Nº 2 contain sodium azide, irritating to mucous and skin. Handle with care.
6 - The samples must be treated as potentially infectious.
7 - Do not ingest. Avoid contact with eyes and skin.
8 - We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
9 - To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

SAMPLES

Serum. The analyte is stable for 2 weeks between 2 and 8°C⁴ and 30 days at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Fructosamine kit installed on refrigerated equipment is 07 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

The kit is indicated for biochemistry automated analyzers only. **Verify the programming for the instrument in our website www.bioclin.com.br or through our SAC.**

RESULTS

The results are expressed in µmol/L.

PROCEDURE LIMITATIONS

Test periodically and calibrate to detect possible changes in the response of the instrument.

INTERFERENCE

The following interference substances were tested did not suffer interference from the substances with the top concentrations of:
Ascorbic Acid 10 mg/dL
Bilirubin 30 mg/dL
Hemoglobin 2000 mg/dL
Triglycerides 2000 mg/dL

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in µmol/L, for this method, were obtained through the determination of fructosamine in healthy populations of male and female.
Fructosamine..... 205 - 285 µmol/L
These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

PRODUCT PERFORMANCE ACCURACY

Comparison of Methods

The Fructosamine Kit was compared with another method commercially available to measure of fructosamine. 07 analysis were conducted and the results were evaluated. The linear equation was obtained Y = 0.996X + 2.305, with a correlation coefficient equal to 0.999. With these results we can conclude that the kit has good methodological specificity.

PRECISION

Repeatability

Repeatability refers to 20 successive measurements of fructosamine were performed with three samples with different concentrations obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (μmol/L)	125.00	241.80	314.85
Standard Deviation (μmol/L)	0.79	0.77	0.81
Coefficient of Variation (%)	0.64	0.32	0.26

Reproducibility

Reproducibility refers to 20 measurements of fructosamine, performed for 3 different days with 3 samples of different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (μmol/L)	124.62	242.17	314.77
Standard Deviation (μmol/L)	0.54	0.32	0.19
Coefficient of Variation (%)	0.43	0.13	0.06

Sensitivity

The sensitivity of the Fructosamine kit was calculated from 20 determinations of a sample free of fructosamine. The average found was of 9.785 μmol/L with standard deviation of 0.067 μmol/L. The sensitivity, which indicates the Method Detection Limit, corresponds to 3 times the standard deviation which is equal to 9.986 μmol/L.

Linearity

The Fructosamine Kit is linear up to 1000 μmol/L.

Diagnostic Significance

Fructosamine is the generic name given to all glycosylated proteins, from which the bigger parcel is of the glycosylated albumin, Who is the biggest plasmatic protein mass after the hemoglobin. The serum level of the fructosamine represents the average value of the blood glucose in the time of twoor three weeks, inferior to the value of the glycohemoglobin and identical to the glycated proteins. The determination of fructosamine for monitoring the glucose metabolism in patients with diabetes, specially the mellitus type II and also adequate to the monitoring and efficiency of the drugs. The fructosamine is elevated in every case of diabetes under inadequate metabolic control, being related to the hyperglycemia and had also been observed that the values return to the reference values twenty days after the stabilization of the glycemia in adequate levels. The reduction of the fructosamine is observed in patients with great loss of albumin or in diseases that increase the proteic catabolism.

Bibliographic References

- 1. Armbuster DA, Fructosamine: Structure, Analysis and Clinical Usefulness. Clin. Chem. 1987; 33 (12): 2153-2163
- 2. Howey JEA, Browning MCK, Fraser CG. Assay of serum fructosamine that minimizes standardization and matrix problems: Use to acess components of biological variation. Clin. Chem. 1987; 33: 269-272
- 3. BURTIS, C.A., ASHWOOD, E.R., Tietz Clinical Chemistry, W. B. Saunders Company. 2nd edition, 1994; 986-988.
- 4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

Quality Assurance

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

Customer Service

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Fructosamine kit: 10269360232

Review: July/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGEDPRODUCT
STERILIZED

DANGER