Bioclin

GAMA GT CINÉTICO

REF

K080

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
K080

USAGE INSTRUCTIONS

FUNCTION
Method for determination of Gamma GT in biological samples of serum or plasma (EDTA or Heparin). Kinetic test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION
Methodology: Modified Szasz / IFCC Standard
This method is based on IFCC (International Federation of Clinical Chemistry), in which the γ -Glutamyltransferase catalyses the transfer of the substrate's Glutamyl radical to a Glycylglycine forming a γ -Glutamylglycylglycine and 5-Amino-2-Nitrobenzoate. The 5-Amino-2-Nitrobenzoate formation speed in a 405 nm is proportional to the activity of the enzyme.

REAGENTS
Number - Buffer - Store between 2 and 8°C. Contains: Buffer < 300 mmol/L, Glycyl Glycine > 50 mmol/L stabilizer and preservative.
Number 2 - Substrate - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, Gamma Glutamyl Carboxy p-Nitroanilide < 100 mmol/L, stabilizers and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 24 mL	1 x 6 mL
2	2 x 24 mL	2 x 6 mL
3	4 x 24 mL	4 x 6 mL
4	1 x 40 mL	1 x 10 mL
5	2 x 40 mL	2 x 10 mL
6	4 x 40 mL	4 x 10 mL
7	2 x 40 mL	1 x 20 mL
8	4 x 40 mL	2 x 20 mL

EQUIPMENTS AND OPERATIONAL INPUTS
Thermostated spectrophotometer, pipettes and test tubes, watch or stopwatch, Biocontrol N, Biocontrol P Bioclin and Biocal Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS
The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE
1- For *in vitro* diagnostic use only.
2- Strictly follow the methodology proposed to obtain exact results.
3- Water used in material cleaning must to be recent and free of contaminants.
4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
5- The reagents N° 1 and N° 2 contain Sodium Azide, irritant to skin and mucous membranes. Handle them with care.

6- It is important for good test performance, a rigorous controlling time and temperature.
7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
9- Do not use the product in case of damaged packaging.
10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES
Serum or plasma collected with Heparin or EDTA, obtained free of hemolysis. The enzyme is stable for 7 days if kept between 2 and 8°C and 1 year at -20°C.⁷ Anticoagulant containing Citrate, Fluoride or Oxalate may inhibit Gamma GT activity.

PROCESS DESCRIPTION
The stability of the Gamma GT Kinetic product after opening (In Use / On Board – Product installed in equipment) is 30 days. This stability may vary according to test conditions. Therefore, it is suggested to monitor the performance of the product using control serum.

PREPARATION OF WORKING REAGENT
Mix 4 mL of Reagent N° 1 with 1 mL of Reagent N° 2. The Working Reagent is stable for 5 days between 15 and 30°C and 3 weeks between 2 and 8°C.

CONDITIONS OF REACTION
It is an indispensable condition using a thermostated cuvette at 37°C with 1 cm optical path and reading at 405 nm.

TECHNIQUE
Bioclin recommends using the kit as control serum for the Biocontrol N and P Bioclin Kits and as a calibrator for the kit Biocal Bioclin.

Add 50 μ L of sample to 1 mL of Working Reagent, mix and transfer to thermostated cuvette at 37°C and wait 1 minute. Make the initial reading, simultaneously starting the stopwatch. Repeat the readings at 1, 2 and 3 minutes. Calculate the average differences in absorbance per minute ($\Delta A/\text{min.}$) and use this value to calculate the result.

CALCULATIONS
Considering the optical path of 1 cm and reaction temperature 37°C the calculation will be:
Gamma GT (U/L) = $\Delta A/\text{min.}$ x Calibration Factor

Filter (nm)	Calibration Factor
405	2638

PROCEDURE LIMITATIONS
The kinetic method is based on the molar absorptivity of 5-Amine-2-Nitrobenzoate, for this reason the readings should be performed in a spectrophotometer that satisfies the following conditions:
Wavelength: 405 nm
Semi trajectory of the pass band: 10 nm
Stray light less than 0.5%
10 mm thermostated cuvette.

INTERFERENTS
Some drugs like Alcohol, Phenytoin, Phenobarbital and Sulfonamides can raise the values of Gamma GT. Contraceptives can reduce Gamma GT levels. Certain sera controls have stabilizers and/or preservatives that can interfere with the linearity, causing variations in the ΔA . This does not occur with biological samples.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY
The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the 5-Amine-2-Nitrobenzoate, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material ERM-AD452 / IFCC and the IFCC reference method.

REFERENCE VALUES
The following reference values were obtained in healthy populations. The dosage was carried out in the temperature of 37°C.

	Male (U/L)	Female (U/L)
0 - 6 months	12 - 122	15 - 132
6 - 12 months	1 - 39	1 - 39
1 - 12 years	3 - 22	4 - 22
13 - 18 years	2 - 42	4 - 24
Adults	$\leq$ 55	$\leq$ 38

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE ACCURACY
Recovery
The recovery analysis was performed with 05 determinations samples. Accuracies were calculated, and met in good agreement with the reference values, obtaining a recovery between 98 and 104%.

Comparison of Methods
The Gamma GT Kinetic kit was compared with other method for determination of Gamma GT. 42 tests were performed and the results were evaluated. The linear equation obtained was $Y = 1.0248X - 1.9615$ and the correlation coefficient obtained 0.9982. With these results we can conclude that the kit shows good methodological specificity.

PRECISION
Repeatability
The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	25.83	108.48	66.53
Standard Deviation (U/L)	1.74	4.95	2.10
Coefficient of Variation (%)	6.73	4.56	3.16

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (U/L)	25.43	106.96	66.17
Standard Deviation (U/L)	0.38	2.01	1.40
Coefficient of Variation (%)	1.48	1.88	2.11

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Gamma GT. The average was 1.950 U/L with a standard deviation of 0.221 U/L. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 2.612 U/L.

Linearity

The reaction is linear to the concentration of 700 U/L. For a variation in absorbance ≥ 0.200, repeat the determination, diluting the sample 1:10 with a solution of Sodium Chloride 9 g/L. Multiply the results obtained by 10.

Diagnostic Significance

Gamma Glutamyl Transferase is a serum enzyme that originates Hepato Biliary System primarily (SHB). In diseases of SHB, as in biliary obstruction intra and extra-hepatic, levels serum values reach 5 to 30 times above the Limits Reference Method (LRM). Its sensitivity in detecting obstructive jaundice, cholangitis and cholecystitis is higher than that of enzymes such as Alkaline Phosphatase (ALP) Arilamidase (Leucineaminopeptidase - LAP) and 5' Nucleotidase - NTP (5' N). Infectious hepatitis occur in increments of 2 to 5 times above the LRM. In pancreatic diseases associated with obstruction hepato biliary the increases are the order of 5 to 15 times higher than the LRM. Bone diseases in the values of Gamma GT are within the LRM is that the joint analysis of Gamma GT and ALP is a valuable asset to distinguish liver disease of bone disease. In muscle disease and renal failure levels are normal, occurring small increase in nephrosis lipid untreated. In the infarct, the values are within the LRM, except when accompanied by liver involvement. The use of drugs as Dilatim (Phenyl Hydantoin), Phenobarbital, and such as Alcohol, increase the level of the enzyme, due to the effect toxic liver. As for Alcohol, is highly significant Gamma GT determination in monitoring the treatment of alcoholic patients.

Bibliographic References

- 1 - J., Biol. Chem., 1939, 128, 537 - 550.
- 2 - DIMOV, D. M.; KULHANEZ, V., Clin. Chim., Acta., 1.967, 16,271 - 277.
- 3 - NAFTALIN, L.; SEXTOM, M.; WHITAKER, S. F.; TRACEY, D., Clin., Chim., Acta., 1.969, 26, 293 - 296.
- 4 - SZASZ, G., Clin., Chem., 1.969, 15, 124 - 136.
- 5 - IFCC Reference Procedure for the Measurement og Catalytic Concentration of γ-Glutamyltransferase. Clin Chem Lab Med 2002; 40:734-38.
- 6 - SOLDIN, S.J., BRUGNARA C., WONG, E.C. Pediatric Reference Intervals, 5.ed.Washington: AACC Press, 2005.p.98-99.
- 7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

Quality Assurance

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

Customer Service

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Gamma GT Kinetic kit: 10269360129

Review: July/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER