GIÁRDIA Ag VET FAST

REF VET002

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FINALIDADE

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INSTRUCTIONS FOR USE

PURPOSE

Rapid immunochromatographic test for qualitative determination of *Giardia lamblia* antigen in canine and feline samples of feces.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.
Giardia Ag VET FAST is an *in vitro* diagnostic test for rapid and qualitative detection of *Giardia lamblia* antigens in dog and cat feces. The cassette consists of a nitrocellulose membrane, impregnated in the test line with antibody anti-*Giardia lamblia*. The conjugate consists of antibody anti-*Giardia lamblia*, conjugated with colloidal gold particles. When the sample is added to the cassette, it migrates by capillary action until it reaches the conjugate. In this region, the antigens present in the sample interact with the conjugate, forming the complex: antigen - conjugate.
The sample continues to migrate through the immunochromatographic strip until reaching the test line, where the Giardia antigens also interact with the antibodies immobilized in the test line, forming the complex: immobilized antibody - antigen - conjugate. This shows the formation of a pink-red line on the test line; indicating positive result. Sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed, indicating that the test yielded a valid result.

REAGENTS

Number 1 (R1): Cassette – Store between 15 and 30°C. **Do not freeze.** It contains:

A) 01 plastic cassette with 01 test strip containing
- Control region (C): Antibody anti-mouse antibody immobilized.
- Test region (T): Immobilized monoclonal antibody anti-*Giardia lamblia*.
- Conjugate: Monoclonal mouse antibody against *Giardia lamblia*, colloidal gold conjugate.

B)01 sachet of silica

Number 2 (R2): Diluent – Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Number 3 (R3): Sterile flexible cotton swab.

Accessory:

- Disposable plastic pipette.

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET002-1	10	10 x 1.0 mL
	VET002-2	25	25 x 1.0 mL
	VET002-3*	10	10 x 1.0 mL
Diluent	VET002-4*	25	25 x 1.0 mL
	VET002-5	5	5 x 1.0 mL

* Individualized presentations

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. They are in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.** The cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

- Product for *in vitro* diagnosis, for veterinary use only.**
- Strictly follow the methodology proposed for obtaining results.
- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
- Use clean, dry containers to package biological samples.
- Open the aluminum sachet to remove the cassette (test devices) only at the time of use.
- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.
- The reagents, as well as the samples, should reach room temperature before starting the test.
- Do not use reagents from different lots or kits.
- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
- Kit components are disposable and must not be reused.
- The kit materials that had contact with the sample of the patient must have a special destination. The disposal of sterile flexible cotton swab, the plastic pipettes, the cassettes (test devices) and diluent vials should be disposed of in organic waste.
- Secondary packaging (box). Instruction for use, envelope aluminum and silica sachet can be disposed of in trash common.
- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES

Perform the test with a sample of dog or cat feces.
Test the sample immediately after collection.
If necessary, samples can be stored between 2 to 8°C for up to 2 days or at -20°C for longer periods.

DESCRIPTION OF PROCESS


STABILITY AFTER OPEN
The stability of Giardia Ag VET FAST cassettes (test devices) after opening the aluminized sachet is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE

- The sample should be at room temperature before starting the test.
- Using the swab, collect a sample from 4 different points of feces.
Note: In samples with liquid or pasty consistency, ensure that excess feces are not absorbed by the swab. Sampling directly from the animal's rectum must be carried out only in the absence of a sample in a collection bottle and attention must be paid to the amount of feces sampled, which when insufficient can impair the accuracy of the results.
- Insert the swab into the diluent tube, dissolving the sample.

Note: Use rotating swab movements inside the tube and rubbing it against the wall for 10 seconds to maximize sample dissolution in the diluent.

- Wait for the larger particles to settle at the bottom of the tube (for approximately 5 minutes).
- Note:** Excess particulate material in the sample diluent supernatant can compromise the run.
- Remove the cassette (test device) from the protective packaging and place it on a clean, dry and level surface.
- Use the pipette that comes with the kit to collect the dissolved sample (supernatant).
- Add 4 drops of sample (80 µL), slowly, into the sample application well of the cassette.
- Wait for the formation of the lines. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF RESULTS

Reactive Test: Formation of two lines: one in the control region (C) and one in the test region (T). Any color intensity (light or dark) means the result is reactive. Do not interpret after 20 minutes.	
Not Reactive Test: Line formation only in control region (C). Do not interpret after 20 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

PROCEDURE LIMITATIONS

This product is intended to detect *Giardia lamblia* infections in dogs and cats. A negative result can be observed if the amount of analyte present in the sample is below the detection limit of the test. This can happen, depending on disease period in which the sample was collected, the type of sample used in the test and its storage conditions. A negative result may also occur in case of inefficiency in the sample extraction step in the diluent.

Therefore, for the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE
CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 109 samples characterized by canine giardiasis were analyzed, which 46 samples were diagnosed as positive and 63 were diagnosed as negative. The results show that the clinical sensitivity of the product is 97.83% and the clinical specificity is 96.83%.

Method		Reference		Total
Giardia Ag VET FAST	Result	Reactive	Not reactive	
	Reactive	45	2	
	Not reactive	1	61	
Total result		46	63	109

Clinical Sensitivity: (45/46): 97.83% (IC 95%: 88.47% - 99.94%)
Clinical Specificity: (61/63): 96.83% (IC 95%: 89.00% - 99.61%)
Precision: (106/109): 97.25%

This product has been tested comparing with other methods. A total of 52 samples characterized by feline giardiasis were analyzed, which 41 samples were diagnosed as positive and 11 were diagnosed as negative. The results show that the clinical sensitivity of the product is > 99,9% and the clinical specificity is 97,56%.

Method		Reference		Total
Giardia Ag VET FAST	Result	Reactive	Not reactive	
	Reactive	11	1	
	Not reactive	0	40	
Total result		11	41	52

Clinical Sensivity: (11/11): > 99,9% (IC 95%: 99,9%-100%)
Clinical Specifciy: (40/41): 97,6% (IC 95%: 92,9% – 100%)
Precision: (51/52): 98,07%

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Result	Found Result
Sample 1	10	Reactive	Reactive
Sample 2	10	Reactive	Reactive

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Result	Found Result
Sample 1	03	10	Reactive	Reactive
Sample 2	03	10	Reactive	Reactive

DIAGNOSTIC SIGNIFICANCE

Giardia spp. is a flagellated protozoan that affects humans and animals. Trophozoites affect the small intestine causing acute, chronic, intermittent or persistent diarrhea. Asymptomatic dogs and cats are common, considering that the appearance of clinical signs is related to the condition of the immune system and the individual susceptibility of the animal. Other clinical signs that are also described include abdominal pain, anorexia, flatulence and bloating. *Giardia* infection occurs via the fecal-oral route, through the ingestion of food or water contaminated with parasite cysts, which can survive for long periods in the environment. Infection by *Giardia* can occur simultoncoulsy with other gastrointestinal parasites, leading to worsening of the clinical picture and favoring infections by other pathogens

BIBLIOGRAPHIC REFERENCES

- 1- Hewlett EL, Andrews JS, Ruffier J, Schaefer FW. Experimental Infection of Mongrel Dogs with *Giardia lamblia* Cysts and Cultured Trophozoites. Journal of Infectious Diseases, v.145, n.1, p. 89–93, 1982.
- 2- Hoehler AV, Jex AR, Haydon SR, Stevens MA, Gasser RB. Giardia/giardiasis — A perspective on diagnostic and analytical tools. Biotechnologies Advances, v.32, n.2, p.280-28, 2014.
- 3- Meireles PW. Giardia sp. Giardiase em animais de companhia. [Dissertação] Universidade federal do Paraná. 64p. 2007.
- 4- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product Licensed in the Ministry of Agriculture, Livestock and Supply since 25/03/2020 with the number 10.317/2020.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: March/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER