Bioclin

GLICOSE MONORREAGENTE

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of glucose in biological samples of Plasma (fluoride), serum, cerebrospinal fluid and fluid (ascilitic, pleural and synovial). Colorimetric enzymatic test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Enzymatic - GOD - PAP (Trinder)
The Glucose is enzymatically oxidized by Oxidase Glucose (GOD) according to the following reaction:

The Peroxide of Hydrogen in the presence of Peroxidase (POD) reacts with 4 - Aminoantipyrine and Phenol, forming a cherry red cromogen which color intensity is proportional to the glucose concentration.

REAGENTS

Number 1 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Tampon < 36 mmol/L, Phenol < 20 mmol/L, 4 - Aminoantipirin < 5 mmol/L, Oxidase Glucose > 10.000 U/L, Peroxidase > 700 U/L, stabilizers, surfactant and preservative.
Number 2 - Standard* - Store between 2 and 8°C, well sealed. Contains: Glucose 100,0 mg/dL and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2*
1	250 mL	3 mL
2	2 x 250 mL	3 mL
3	4 x 250 mL	3 mL
4	5 x 20 mL	3 mL
5	5 x 40 mL	3 mL
6	10 x 40 mL	3 mL
7	2 x 60 mL	3 mL
8	4 x 60 mL	3 mL
9	6 x 60 mL	3 mL
10	8 x 60 mL	3 mL
11	10 x 60 mL	3 mL

* The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watches and stopwatches, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. These itens are found at markets in the specialized market of clinical analysis laboratories.

TRANSPORTAION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C. The transport temperature up to 30°C, must not exceed 5 days. Protect from light and avoid moisture. **Do not**

freeze.

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.
- Water level in water bath should be above the level of reagents in test tubes.
- Reagents N° 1 and N°2 contains Sodium Azide which can be irritating to skin and mucous membranes. Handle with care.
- The development of a pinkish clor in reagent N° 1 do not interfere in reagents quality and stability, once the correspondent Blank be used and dosages in regular bases.
- Determinate the factor periodically and in every batch.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Plasma (fluoride), serum, cerebrospinal fluid, fluid (ascilitic, pleural and synovial).

The use of Fluoride Bioclin anticoagulant is recommended as a glycolysis inhibitor. Use 1 drop for every 3 mL of blood. Serum can only be used if it is centrifuged, separated from the cells and dosed immediately after collect. In other biological fluids add a glycolysis inhibitor in the same proportion as described for the blood by centrifuging the sample before starting the dosing. Plasma or Serum: stable for 7 days between 2 and 8°C⁵. Cerebrospinal Fluid: stable for 3 days between 2 and 8°C and 1 month at -20°C⁵.

PROCESS DESCRIPTION

The calibration stability of the Glucose Monoreagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is

suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 tests tubes with the letters B (Blank), A (Sample), P (Standard) and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Reagent N° 2	--	10 µL	--
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL

Homogenize well and place in water-bath at 37°C for 10 minutes.

Read the sample and standard absorbance at 505 nm (490 - 550 nm), matching the zero with the blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Glucose (mg/dL)} = \frac{\text{Sample Absorbance} \times 100}{\text{Standard Absorbance}}$$

Example:

$$\begin{aligned} \text{Sample Absorbance} &= 0.347 \\ \text{Standard Absorbance} &= 0.350 \end{aligned}$$

$$\text{Glucose (mg/dL)} = \frac{0.347 \times 100}{0.350} = 99$$

The reaction follows the Lambert-Beer Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (100 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Glucose (mg/dL)} = \text{Absorbance Sample} \times \text{Calibration Factor}$$

Example:

$$\text{Calibration Factor} = \frac{100}{0.350} = 286$$

$$\text{Glucose (mg/dL)} = 0.347 \times 286 = 99$$

The results are express in mg/dL.

PROCEDURE LIMITATIONS

The proposed method is not indicated for glucose dosage in urine.

INTERFERENCES

Samples with concentration up to 20 mg/dL of Bilirubin, 750 mg/dL of triglycerides and 160 mg/dL of Hemoglobin do not produce significant interference.

In case of interferences produced by the samples, also carry

out the sample of Blank, in order to minimize the action of interferences. Proceed as follows:

Mark 1 tube with Sample of Blank and place 1.0 mL of Sodium Chloride 0.85% with 10 µL of sample. Determine the absorbance at 490 - 550 nm, matching the zero with distilled or deionized water. Subtract the obtained absorbance, of the Sample tube absorbance. Calculate the concentration multiplying the result of Calibration Factor.

The use of highly reductive medicine as such as Ascorbic Acid (Vitamin C) interfere in the reaction, because they compete with the consumption of H₂O₂, supplying values falsely lowered. For this reason, must suspend its use per less than 12 hours before the sample collection.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 917.

REFERENCE VALUES

The reference values in mg/dL for this method were obtained through the determination of glucose in healthy populations of male and female.

Plasma	Premature	20 - 60 mg/dL
	0 to 1 day	40 - 60 mg/dL
	> 1 day	50 - 80 mg/dL
	Children and adults	65 - 99 mg/dL
Liquor		50 - 70 mg/dL

To convert the values of mg/dL into mmol/L (SI), multiply the obtained results by 0.0556.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or

treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Recovery
The recovery analysis was made with 5 sample determinations. The accurate was calculated and was found in a good parity with the reference values, obtaining a recovery between 94% and 104%.

Comparison of Methods
Glucose Monoreagent kit was compared with other commercially available methods for measurement of glucose. For this test, 42 analyses are performed and the results were evaluated. The linear equation obtained was $Y = 1.0053X + 0.4653$ and the correlation coefficient 0.9982. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION
Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	240.75	115.01	73.59
Standard Deviation (mg/dL)	3.72	4.15	1.92
Coefficient of variation (%)	1.55	3.61	2.61

Reproducibility
The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	242.29	114.58	74.06
Standard Deviation (mg/dL)	2.56	0.60	0.56
Coefficient of variation (%)	1.06	0.52	0.76

SENSITIVITY
The sensitivity was calculated from 40 determinations free of a sample free of glucose. An average, which indicates the detection limit, value of 0.55 mg/dL, with standard deviation of 0.25 mg/dL. The sensitivity, which indicates the method quantification limit, corresponds to the average plus three times the standard deviation and is equal to 1.31 mg/dL.

LINEARITY
The reaction is linear up to 500 mg/dL. For samples with values above 500 mg/dL, dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the results obtained by the dilution factor.

DIAGNOSTIC SIGNIFICANCE
The glycemc homeostasis is controlled by the action of several hormones, especially insulin, who keeps the equilibrium of the glucose concentration. Hormonal alterations and other factors generate a variation in the homeostasis, leading to a hyper and hypoglycemic. The Hyperglycemia occurs in several types of Mellitus Diabetes, where frequently are retinopathy, kidney injuries, neuropathies and arteriosclerosis. The Mellitus Diabetes is classified in: Mellitus diabetes insulin dependent (Type I), Mellitus diabetes non insulin dependent (Type II), Mellitus diabetes associated at certain conditions and syndromes (classified previously as secondary diabetes) and Pregnancy diabetes. The Hypoglycemic (HG) levels of glycemia that leads to its manifestations are extremely variable. The manifestations can occur during fasting or postprandial. Occurs fasting HG in insulinome, non pancreatic tumors, hepatic diseases, hypoadrenalin (Addison's disease), hypopituitarism (Simmond's disease), storage of infirmity of glycogenic storaged (Von Gierke's diseases). The postprandial hypoglycemia occurs due to the reactive cause (symptoms of HG 1 to 3 hours after the meals), it can also be of alimental origin or in consequence of mellitus diabetes Type II and the intolerance Glucose.

BIBLIOGRAPHIC REFERENCES
1 - TRINDER, P., Ann. Clin. Biochem., 6-24, 1969.
2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
3 - LOTT, J. A.; TURNER, K., Clin. Chem., 1987, 21:1754, 1970.
4 - CARL, A. B. and EDWARD R. A., Tietz Textbook of Clinical Chem. 2nd ed., 928-997,
5 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE
Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Glucose Monoreagent kit:
10269360384

Review: April/2025

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 LOT NUMBER

 MANUFACTURING DATE

 VALIDITY DATE
(last day of the month)

 TEMPERATURE LIMIT
(store)

 CONTENT IS SUFFICIENT
FOR <N> TEST

 SEE INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC
PRODUCT

 KEEP AWAY
FROM SUNLIGHT

 DO NOT REUSE

 CAUTION

 MADE BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 FLAMMABLE

 CORROSIVE

 TOXIC

 DO NOT USE IF
PACKAGE IS
DAMAGED

 PRODUCT
STERILIZED

 DANGER