Bioclin

HBsAg

REF K145

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Immunochromatographic method for rapid and qualitative determination of Hepatitis B virus surface antigen (HBsAg / subtypes ad and ay) in biological samples of serum, plasma or whole blood. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The HBsAg kit is an immunochromatographic test that contains IgG (anti-HBs) - colloidal gold conjugate that binds to the HBsAg present in the sample and forms an antigen / antibody / colloidal gold complex. The complex migrates through the membrane by capillary action towards the anti-HBsAg antibodies in the test region (T). After the addition of the sample, when binding to the antibodies present in the test region, the appearance of a line that determines the reagent reaction occurs. In the absence of HBsAg surface antigens, the line will not appear. The complex continues to migrate until it reaches the control region (C). The conjugate not bound to the antigen binds to the C region anti-IgG antibody producing a colored line indicating that the test is working and that the result is valid.

REAGENTS

Number 1 (R1): Cassette - Store between 2 and 30 °C. **Do not freeze.** Contains:

- A) 01 plastic cassette (device) with 01 test strip containing:
- Control Region (C): Anti-IgG Antibody;
 - Test Region (T): Goat Anti-HBsAg Antibody;
 - Conjugate: Conjugate IgG (Anti HBs) - colloidal gold.

B) 01 silica sachet

Number 2 (R2): Diluent - Store between 2 and 30°C. **Do not freeze.** Contains: Buffer Solution < 200 mmol/L and preservative.

Accessories for Sample Collection:

- Disposable plastic pipette
- Disposable lancet (Sterile).

PRESENTATIONS

Reagent	Number of Cassettes of Packaging	Diluent
Cassette + Diluent (Pipette included) *	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
	100	2 x 6.0 mL
Cassette + Diluent (Lancet and Pipette included) *	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
	100	2 x 6.0 mL

* The number of accessories is proportional to the number of cassettes.

INSTRUMENTS AND OPERACIONAL INPUTS

Pipettes, tips, watch or stopwatch, 70% alcohol (w/w), cotton and garbage for disposal. These items are found in the specialized market for articles for Clinical Analysis Laboratories. The kit allows visible reading without the aid of equipment.

CONDITIONS OF STORAGE AND TRANSPORTATION

Transport can be carried out at temperatures up to 45°C for a maximum of 45 days. The storage temperature should be 2 to 30°C. Avoid moisture. **Do not freeze.**

WARNINGS

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Use clean and dry containers to store the samples.
- 5- Open the aluminum envelope to remove the cassette (test device), only at the time of use.
- 6- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration.
- 7- The reagents, as well as the sample, must reach room temperature before starting the test.
- 8- Do not mix the reagents (components) of the kits with the same batch or with different batch numbers.
- 9- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.
- 10- The cassette (test device), the pipette and the lancet are disposable and cannot be reused.

11- The kit materials that came into contact with the patient's biological material must have a special destination. Dispose of lancets and pipettes in sharps waste. Cassettes (test devices) and diluent vials must be disposed of in biowaste.

12- Boxes, Instructions for Use, aluminum envelope and silica sachet can be disposed of in regular trash.

13- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the FISPQ (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Advice) by Quibasa.

15- Do not use the product if the packaging is damaged.

16- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma: Samples can be stored under refrigeration, between 2 and 8°C, for a maximum period of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 30 days at -20°C (freezer). To obtain plasma, perform collection using EDTA, Heparin or Citrate as an anticoagulant.


Whole Blood: Collect whole blood by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or by digital puncture with the aid of a lancet. Whole blood collected by lancet should be used for immediate testing. Whole blood collected by venipuncture with anticoagulant can be kept under refrigeration, between 2 and 8°C, for a maximum period of 3 days. **Do not freeze.**

PROCESS DESCRIPTION

The stability of the HBsAg cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE - VENOPUNCTURE

- 1- The sample must be at room temperature before starting the test.
- 2- Remove the cassette (test device) from its protective packaging, place it on a clean and level surface, and identify it properly.
- 3- **For Whole Blood:** Transfer 75 µL of whole blood into the sample well. Hold the vial vertically and apply 2 drops of Diluent (Reagent 2) into the sample well.





- 4- **For Serum or Plasma:** Transfer 80 - 100 µL of serum or plasma into the sample well.

- 5- Wait for the lines to form. Interpret results between 15 and 30 minutes. Do not interpret after 30 minutes.


TECHNIQUE - DIGITAL PUNCH / LANCET

- 1- Remove the cassette (test device) from its protective packaging, place it on a clean and level surface and identify it properly.
- 2- Wash and dry the patient's hands. Apply 70% alcohol (w/w) on the tip of the finger to aseptize the area used.
- 3- Press the tip of the finger that will be pierced by the lancet to collect blood in this region.
- 4- Remove the protective cap from the lancet.
- 5- Position and press the lancet firmly on the area to be lanced. Then the blood will come out through the pierced area.
- 6- Collect the blood using the disposable plastic pipette that comes with the kit. Place the pipette over the drop of blood and let the blood flow by capillarity, without pressing the bulb, up to the line marked on the disposable plastic pipette (equivalent to 75 µL of blood).
- 7- With one hand, cover the pipette hole. With the other hand, dispense the blood into the sample well by pressing the disposable plastic pipette.
- 8- Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (50 µL) of Diluent in the sample well.
- 9- Wait for the lines to form. Interpret results between 15 and 30 minutes. Do not interpret after 30 minutes.


INTERPRETATION OF THE RESULT


Reagent Test: Formation of two red lines, one in the control region (C) and one in the test region (T) after 15 minutes. **Do not interpret after 30 minutes.**

Non-Reagent Test: Formation of a red line in the control region (C) and complete absence of a line in the test region (T) after 15 minutes. **Do not interpret after 30 minutes.**

Invalid Result: The absence of line formation in the control region (C) indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

LIMITATIONS OF THE PROCESS

- 1- The line in the control region will appear before the 15 minutes of the test run, and this does not mean that the results can be interpreted before the expected time.
- 2- It is important to use the correct sample volume, as lower or higher volumes may determine erroneous results.
- 3- The results provided by this kit must be interpreted by the responsible professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.
- 4- If the suspicion of HBV infection persists, a new sample must be collected 30 days after the date of collection of the first sample to carry out a new test.
- 5- This kit has been validated for serum, plasma and whole blood samples only. Do not use other biological liquids or sample pool.
- 6- Do not use samples from the *post-mortem* phase.
- 7- The kit allows the detection of HBsAg in all stages of the infection, acute or chronic, as long as it is at detectable levels.
- 8- Isolated test results cannot be used for definitive diagnosis. To confirm reagent results, follow the guidelines of the Technical Manual for the Diagnosis of Viral Hepatitis approved by Ordinance No. 25, of December 1, 2015.
- 9- A negative result does not exclude the possibility of infection with the Hepatitis B virus (HBV). The concentration of HBsAg in the sample may vary according to the stage of the disease, when the test is carried out, with the presence of anti-HBs antibodies and in patients treated with antivirals, which cause a reduction in the viral load. These conditions can lower the patient's HBsAg levels below the kit's analytical detection limit.

INTERFERENT

Hemolysis with hemoglobin values up to 10 g/L do not interfere with the HBsAg dosage. Significant interferences were not found up to 30 g/L of lipids, 0.2 g/L of Bilirubin and 60 g/L of proteins. Interferences above these values may cause false reagent or false non-reactive results. There was no evidence of interference with Rheumatoid Factor up to 1083 IU/mL, C-Reactive Protein up to 41.5 mg/dL and Anti-Streptolysin O up to 1046 IU/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY Comparison of Methods, Methodological Specificity and Sensitivity

The HBsAg kit was compared with another commercially available HBsAg kit. Eighty-six analyzes were performed and the results were evaluated. The concordance found was > 99.9%. With these results, it can be concluded that the kit has good methodological specificity.

	Kit Reference	HBsAg
Reactive Samples	32	32
Non-Reactive Samples	54	54
Total	86	

Total agreement found: > 99.9% (86/86)

PRECISION Repeatability

The repeatability was calculated from 10 successive determinations, using 3 different samples, obtaining the following results.

Samples	Nº of Repetitions	Expected result	Found result
01 (Serum)	10	Reagent	Reagent
02 (Plasma)	10	Reagent	Reagent
03 (Whole Blood)	10	Reagent	Reagent

The results showed 100% correlation between the expected and found results.

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 different samples, obtaining the following results.

Samples	Nº of Days	Nº of Repetitions	Expected result	Found result
01 (Serum)	03	10	Reagent	Reagent
02 (Plasma)	03	10	Reagent	Reagent
03 (Whole Blood)	03	10	Reagent	Reagent

The results showed 100% correlation between the expected and found results.

SENSITIVITY AND CLINICAL SPECIFICITY

The HBsAg kit has been tested using clinically known samples. For this test, 1365 samples were analyzed, 349 reagent for HBsAg and 1016 non-reagent for HBsAg. The following results were found:

	Expected result	HBsAg
Reactive Samples	349	349
Non-Reactive Samples	1016	1014
Total		1365

The results found suggest the following values:

Clinical Sensitivity: > 99.9% (349/349).

Clinical Specificity: 99.8% (1014/1016).

ANALYTICAL SENSITIVITY

The analytical sensitivity is 2 ng/mL.

HIGH DOSE PROZONE EFFECT

No prozone effect was verified with a sample of high concentration of HBsAg up to 800 µg/mL.

DIAGNOSTIC SIGNIFICANCE

Hepatitis B is an infectious disease caused by the HBV virus that affects the liver. The HBV virus is transmitted through sexual contact, blood exposure (including sharing sharps), and mother-to-child transmission via the placenta or during childbirth. HBV infection results in an apparent number of serological markers and one of the first of these markers is Hepatitis B Surface Antigen (HBsAg). The main ones HBsAg subtypes include ad and ay, all sharing the common determinant 'a'. HBsAg appears between 1 and 10 weeks after infection and before biochemical evidence of liver disease or jaundice. However, for most people, due to its high production, HBsAg can be detected after 30 days of infection. This period is known as the diagnostic window. In the chronic carrier state, HBsAg persists for more than 6 months. Antibodies against HBV antigens appear between 30 and 60 days after infection, varying between individuals. This is the period known as the immunological window. HBV can cause self-limiting acute infection, fulminant hepatitis, chronic hepatitis with progression to cirrhosis and liver failure, besides asymptomatic chronic carrier status. The HBV virus in infected people persists for life and can be transmitted to other people. Thus, hepatitis B has become a public health problem. Therefore, screening for HBsAg is recommended for all donors, pregnant women and people in high-risk groups.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the packaging, as long as they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

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