Bioclin

HIV TRI LINE

REF **K087**

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USAGE INSTRUCTIONS

FUNCTION

Immunochromatographic method for the rapid and qualitative determination of total anti-HIV 1 antibodies (IgM, IgG and IgA), including group M and O and anti-HIV 2 in biological samples of serum, plasma or whole blood. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The 3rd generation HIV Tri Line kit, manufactured with recombinant antigens in both the solid phase and the conjugate, is capable of detecting total antibodies (IgM, IgG and IgA) against HIV 1, including groups M and O and anti-HIV 2.

Recombinant HIV 1 and HIV 2 antigens react with antibodies present in serum, plasma, and whole blood samples. The samples move through a chromatographic membrane by capillary action. The HIV 1 reagent samples will form a red line in the region where the recombinant HIV 1 antigen is immobilized. The HIV 2 reactive samples will form a second red line in the region corresponding to the immobilized HIV 2 antigen. The samples continue to be absorbed through the membrane to the control region, with the formation of another red line, confirming the correct processing of the test.

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 30 °C. **Do not freeze.** Contains:

A) 01 plastic cassette (device) with 01 test strip containing:
 - Control region (C): Polyclonal goat anti-IgG antibody;
 - Test region (T): Recombinant antigens of HIV-1 and HIV-2;
 - Conjugate: Recombinant HIV-1 and HIV-2 antigens conjugated with colloidal gold;

B) 01 sachet of silica;

Number 2 (R2) Diluent - Store between 2 and 30 °C. **Do not freeze.** Contains: Tris buffer (pH 8.0).

Collection Accessories:

- Disposable plastic pipette.
- Disposable lancet

PRESENTATION

The kit has the following presentations:

Reagent	N° of Cassettes per Package	Diluent
Cassette + Diluent (Includes pipettes)*	5	5,0 mL
	10	5,0 mL
	20	5,0 mL
	25	5,0 mL
	30	5,0 mL
	40	2 X 5,0 mL
	50	2 X 5,0 mL
Cassette + Diluent (Includes lancetes and pipettes)*	100	3 X 5,0 mL
	5	5,0 mL
	10	5,0 mL
	20	5,0 mL
	25	5,0 mL
	30	5,0 mL
	40	2 X 5,0 mL
	50	2 X 5,0 mL
	100	3 X 5,0 mL

*These presentations include a lancet and / or disposable plastic pipette. The number of accessories is proportional to the number of cassettes.

OPERATIONAL EQUIPMENT AND SUPPLIES

Pipettes, pointers, watch or stopwatch, alcohol 70% (w/w), cotton and the garbage for disposal. These items are found in the specialized market of articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORTATION CONDITIONS

Transport can be carried out at temperatures up to 45 °C for a maximum of 3 weeks. The storage temperature should be between 2 and 30 °C. Avoid humidity. **Do not freeze.**

SPECIAL CARES

1- Only for *in vitro* diagnostic use.

- 2- Rigorously follow the proposed methodology to obtain exact results.
- 3- The cassettes (device) do not contain contaminating substances, but as the samples are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Use clean and dry containers to settle the samples.
- 5- Open the aluminum envelope to remove the cassette (device) only at the time of use.
- 6- The kit should not be frozen or exposed to high temperatures (above 45 °C), as this will cause it to deteriorate.
- 7- The reagents, as well as the sample, must reach room temperature before starting the test.
- 8- It is recommended the application of local, state and federal environmental protection law for the elimination of reagents and biological material is done in accordance with current legislation.
- 9- To obtain information related to biological safety or in case of accidents with the product, consult the FISPQ (Chemical Product Safety Information Sheet) available on the site www.bioclin.com.br or requesting through the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damage to its package.

11- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum/Plasma: Serum and plasma samples are stable for 3 days when kept between 2 and 8 °C and for 03 months at -20 °C, it can be thawed for use only once. Hemolyzed and lipemic samples should be avoided, as they can cause a false reactive result. Serums and plasmas containing precipitates may give inconsistent results. To obtain the plasma, collect EDTA, Heparin or Citrate as anticoagulant.


Whole Blood: Collect whole blood by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or with the help of a lancet (fingerstick). Whole blood obtained by fingerstick should be used for immediate tests. Whole blood obtained by venipuncture with anticoagulant can be stored at 2-8 °C for up to 2 days. **Do not freeze whole blood.**

PROCESS DESCRIPTION

The stability of the HIV Tri Line cassette (after opening the aluminized envelope) is 2 hours. This stability may vary depending on the environmental conditions after opening the aluminized envelope.


TECHNIQUE - VENOPUNCIÓN

- 1- The sample must be at room temperature before starting the test.
- 2- Remove the cassette (device) from the protective package, place it on a clean and level surface and identify it properly.
- 3- For Serum, Plasma or Whole Blood: Transfer 10 µL of Serum, Plasma or Whole Blood into well A.
- 4- Keep the Diluent (Reagent N° 2) vertically, and apply 2 drops (70 µL) in well B.
- 5- Wait for the lines to form. Interpret the results between 15 and 30 minutes. Do not interpret after 30 minutes.

TECHNIQUE - DIGITAL PUNCTURE / LANCET

- 1- Remove the cassette (device) from the protective package, place it on a clean and level surface and identify it properly.
- 2- Wash and dry the patient's hands. Pass the alcohol 70% (w/w) on the fingertip to asepsis the area used.
- 3- Press the tip of the finger that will be pierced by the lancet for the accumulation of blood in this region.
- 4- Remove the protective cap from the lancet.
- 5- Place and press the lancet firmly on the area to be punctured. Then blood will come out of the pierced area.
- 6- Collect the blood with the help of the disposable plastic pipette that comes with the kit. Aspirate the blood to the marked trace on the disposable plastic pipette.
- 7- In the well A, dispense the blood by pressing the pipette.
- 8- Hold the diluent bottle (Reagent No. 2) vertically and apply 2 drops (70 µL) of diluent to well B.
- 9- Wait for the lines to form. Interpret the results between 15 and 30 minutes. Do not interpret after 30 minutes.

INTERPRETATION OF RESULTS

Reactive Test: Control line formation plus one or more two red lines after 15 minutes in the test region. Do not interpret after 30 minutes.


HIV 1 Reactive: In addition to the red line in the control region (C), another line is formed in the test region 1 (T1).

HIV 2 Reactive: In addition to the red line in the control region (C), another line is formed in the test region 2 (T2).

HIV 1 and 2 Reagents: In addition to the red line in the (C) region, two other lines are formed in the test regions 1 and 2 (T1 and T2). If the color intensity of the T1 line is greater than the T2 line, the result can be considered as reactive HIV-1. If the color intensity of the T2 line is greater than the T1 line, the result can be considered as reactive HIV-2.

Non-Reactive Test: Formation of only one red line, after 15 minutes, in the Control region. Do not interpret after 30 minutes. Note: See notes.

Invalid Test: The absence of line formation in the Control (C) region indicates error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

Notes:

- 1- The line in the control region will appear before 15 minutes of incubation for the test, and this does not mean that the results can be interpreted before the expected time.
- 2- It is important to use the correct sample volume, as lower or higher volumes can lead to erroneous results.
- 3- HIV-2 has a 60% homology in relation to HIV-1. Considering that HIV-1 and HIV-2 coinfection is rare, if the test shows the HIV-1 and HIV-2 reagent lines, this result is likely due to cross-reactivity between them. Therefore, if the HIV Tri Line kit shows a reactive result for HIV-2, follow the guidelines in the Technical Manual for the diagnosis of HIV infection.
- 4- The results of isolated tests cannot be used for a definitive diagnosis. To confirm the results of the reagents, follow the guidelines of Ordinance No. 29, of December 17, 2013, which approves the Technical Manual for the diagnosis of HIV infection.
- 5- An individual infected with the HIV virus may present an immunological window that can vary from 10 days to 6 months. This period includes cases of rapid seroconversion with detection of IgM antibodies and very rare cases of individuals who have a prolonged course of seroconversion. This period may vary depending on the immune response of each individual. The literature reports that 95% of those infected seroconvert within 30 days.
- 6- Considering that the HIV Tri Line kit is of the 3rd generation, that is, manufactured with recombinant antigens both in solid phase and in conjugate, it is capable of detecting total antibodies (IgM, IgG and IgA) against HIV 1, including the group M and O and HIV 2- The literature describes that third-generation kits are capable of detecting antibodies on average 22 to 25 days after infection, since they are capable of detecting IgM antibodies.

- 7- If the suspicion of HIV infection persists, a new sample must be collected 30 days after the date of collection of the first sample to perform a new test.
- 8- False non-reactive results may occur if the test is used on samples from patients known to be seropositive, receiving or have received antiretroviral therapy.
- 9- The results provided by this kit must be interpreted by the responsible professional, and are not the only criteria to determine the diagnosis and / or treatment of the patient.

INTERFERENTS

Hemolysis with hemoglobin values up to 10 g/L does not interfere with the measurement of HIV.

No significant interferences were found up to 30g/L of lipids, 0.2 g/L of bilirubin and 60 g/L of proteins, as well as no interferences were observed with the anticoagulants EDTA, Citrate and Heparin. No cross-reactions were found with HBsAg, HCV, and Syphilis.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. Therefore, the use of controls is recommended, which allow the evaluation, precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods, Methodology and Specificity

To establish the sensitivity and specificity of the HIV Tri Line kit, another commercially available qualitative kit was used as a comparison.

A 100% concordance was evidenced. Comparative studies using 132 truly reactive samples and 93 truly non-reactive samples have shown that the HIV Tri Line kit has >99.9% (95% CI: 98.5 - 100%) sensitivity (132/132) and >99.9% (95 CI: 98.5 - 100%) of specificity (93/93). Cross-reactivity tests were performed using reactive anti-HCV, anti-TP (Syphilis) and anti-HBc sera, and no reactivity was observed.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 different samples, obtaining the following results:

Sample	Nº of Repetitions	Result Expected	Result Found
Serum	10	100% Reactive	100% Reactive
Plasma	10	100% Reactive	100% Reactive
Whole Blood	10	100% Reactive	100% Reactive

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 different samples, obtaining the following results:

Sample	Nº of Days	Nº of Repetitions	Result Expected	Result Found
Serum	03	10	100% Reactive	100% Reactive
Plasma	03	10	100% Reactive	100% Reactive
Whole Blood	03	10	100% Reactive	100% Reactive

DIAGNOSTIC MEANING

The Human Immunodeficiency Virus (HIV) is a retrovirus, identified as the etiological agent of Acquired Immune Deficiency Syndrome (AIDS). AIDS is characterized by changes in the T lymphocyte population, which play a key role in the immune system. In the infected individual, the virus causes a reduction in the subpopulation of T cells, called "helper" T cells, which make these patients susceptible to opportunistic infections and certain malignancies. The main routes of transmission are: unprotected sex, contamination by blood or blood products, and vertical transmission.

HIV is made up of an RNA molecule, protected by a capsid and an envelope. There are two types of HIV virus: HIV 1 and HIV 2. The HIV envelope is the main target of the immune response. The presence of the virus causes the patient's immune system to produce anti-HIV antibodies. Detection of these antibodies can be used as a diagnostic tool.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE

Customer Advisory Service
Tel: 0800 031 5454 | E-mail: sac@bioclin.com.br

Registration number of the HIV Tri Line kit in ANVISA:
10269360148

Revision: July/2022

UNIVERSAL SYMBOLOGY