Bioclin

IMUNOGLOBULINA G

REF K062

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REF: **K062**

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method of determination of IgG in biological samples of serum or plasma (EDTA or Heparin). Immunoturbidimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry.

The reagent allows the quantitative determination of IgG in human serum by immunoprecipitation, in the presence of an activator polymer that increases the sensitivity and speed of the immunoturbidimetric assay. The Immunoglobulins form with the specific antiserum an insoluble complex, producing turbidity, whose intensity increases the absorbance intensity proportionally to the concentration of IgG in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6,000 < 15%, Buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for Immunoglobulin G < 8 g/L, Buffer and preservative.

PRESENTATION

Presentation	Reagent	
	Nº 1	Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
3	3 x 25 mL	3 x 2 mL
4	4 x 25 mL	4 x 2 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport temperature up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Reagents n° 1 and 2 contain Sodium Azide, irritant to skin and mucous. Handle carefully.
- 5- Do not freeze the reagents, the freezing of Antiserum or Buffer can affect their performance.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma obtained with heparin or EDTA. This sample is stable for 2 days if kept in temperatures between 2 and 8°C or 30 days at -20°C. Do not use samples highly hemolyzed or lipemic.

PROCESS DESCRIPTION

The calibration stability of the Immunoglobulin G kit installed on refrigerated equipment is up to 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF WORK REAGENT

Reagents are ready to use.

TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

INTERFERENT

Hemolysis, icterus and lipemia (low or moderate) do not interfere in test performance, using the recommended dilution for the sample. Gross lipemia and turbid samples should be processed with a greater dilution

PROCEDURE LIMITATIONS

Turbidimetric assays of Immunoglobulin cannot discriminate monoclonal, polyclonal and oligoclonal elevations.

However, monotypic Immunoglobulins reach excess Antigen at concentrations much lower than in the case of immunoglobulin polyclonal normal. It is recommended that testing of these proteins be held in serum dilutions which allow readings in lower limits of the calibration curve. It is suggested that such cases results are verified by Immunoelectrophoresis. Usually the phenomenon of excess Antigen is characterized by reduced reaction coefficient of Immunoprecipitation.

It is recommended that tests are performed in duplicate when is not possible to control by Immunoelectrophoresis. Perform proposed dilution of 1:10 and also 1:20. Calculate the concentration of 1:20 dilution found by multiplying the results by 1.961. The result cannot make a significant difference ($\pm 15\%$) of dilution 1:10.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Age	Unit mg/dL
Newborns	750 a 1,500
1 a 3 meses	190 a 780
3 a 5 meses	60 a 860
6 a 12 meses	250 a 1,180
1 a 2 years	350 a 1,200
2 a 4 years	500 a 1,440
4 a 6 years	610 a 1,420
6 a 9 years	570 a 1,410
9 a 11 years	630 a 1,400
11 a 15 years	630 a 1,510
Adults	650 a 1,500

These values should be used as a guideline, and each laboratory should establish its range of reference values, according with the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 96 and 103%.

Comparison of Methods

The Immunoglobulin G kit was compared with other commercially available methods for measurement of IgG. 40 analyzes were performed and the results were evaluated. The linear equation obtained was $y = 1.069X - 4.894$ and the correlation coefficient 0.992. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION**Repeatability**

The repeatability was calculated from 10 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2
Concentration (mg/dL)	1353.20	1772.79
Standard Deviation (mg/dL)	13.20	20.41
Coefficient of Variation (%)	0.98	1.15

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2
Concentration (mg/dL)	1356.99	1777.25
Standard Deviation (mg/dL)	10.50	21.32
Coefficient of Variation (%)	0.77	1.20

SENSITIVITY

The sensitivity was calculated from 15 determinations of a sample free the IgG. The average found, that indicates the method detection limit, was 25.820 mg/dL and the standard deviation of 0.391 mg/dL. The sensitivity, that indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 26.994 mg/dL.

LINEARITY

The reaction is linear up to the level of the highest calibrator concentration.

DIAGNOSTIC SIGNIFICANCE

The Immunoglobulin G is a simple monomeric immunoglobulin which makes 80% of immunoglobulins in the body. It is also distributed in extracellular compartments and is the only one that normally crosses the placenta, giving a large degree of passive immunity to the newborn. It is the main antibody in secondary immune responses and single class antitoxins. FC held their region of complement activation (when attached to the antigen), and assists in phagocytosis by binding to macrophages. With the activation of complement, no generation of chemotaxis neutrophils, increase vascular permeability and the enhanced inflammatory response. The level of IgG is increased in chronic granulomatous infections (example tuberculosis, syphilis, paracoccidioidomycosis), infections in general, hyper immunizations, severe malnutrition, disproteinemia, diseases associated with hematological disorders or granulomas of hypersensitivity and rheumatoid arthritis.

It's level will decreased in agammaglobulinemia aplasia lymphoid, selective IgA deficiency or IgC, IgA myeloma, proteinemia, Bence Jones syndrome, lymphoblastic leukemia and chronic infections by HIV virus.

BIBLIOGRAPHIC REFERENCES

- 1 - KILLINGSWORTH, L. M.; SAVORY, J. :Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem.; 1971, 17, 936.
- 2 - PESCE, A. J.; KAPLAN, L. A.: Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
- 3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
- 4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
- 5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Immunoglobulin G kit: 10269360155

Review: March/2026

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER