ALFA 1 GLICOPROTEÍNA ÁCIDA

REF K060

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for the determination of Alpha 1 Acid Glycoprotein in human samples of serum or plasma (heparin). Immunoturbidimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

The reagent allows the quantitative determination of Alpha 1 Acid Glycoprotein in human serum by immuno-precipitation, in the presence of an activator polymer that increases the sensitivity and immunoturbidimetric test speed. The Alpha1 Glycoprotein Acid to forms along with the specific antiserum an insoluble complex, producing turbidity, whose absorbance is proportional to the concentration of Alpha 1 Acid Glycoprotein in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for Alpha 1 Acid Glycoprotein < 2 g/L, buffer and preservative.

PRESENTATION

Presentation	Reagent	
	Nº 1	Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
3	3 x 25 mL	3 x 2 mL
4	4 x 25 mL	4 x 2 mL
5	1 x 40 mL	1 x 10 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

Calibrator not included in kit.

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Handle with care as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use fresh serum or plasma (heparin). The analyte is stable for 2 days at 2 to 8°C and 30 days at -20°C. Avoid repeated freezing and thawing. For samples thawed, homogenized prior to use. Do not use samples with signs of contamination.

PROCESS DESCRIPTION

WORKING REAGENT PREPARATION

Reagents are ready to use.

TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia (low or moderate) do not not interfere in performance of the tests, using the recommended dilution for sample. Gross lipemia and turbid samples should be processed with a higher dilution.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Serum/Plasma (heparin)	Adults	30 - 120 mg/dL
------------------------	--------	----------------

These values should be used as a guideline, and each laboratory should establish its range of reference values, according with the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 98 and 106%.

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

Alpha 1 Acid Glycoprotein kit was compared with other method for determination of Alpha 1 Acid Glycoprotein dosage commercially available. 40 tests were performed and the results were evaluated. The linear equation obtained was $Y = 0,946X + 5,980$ and correlation coefficient 0,989. With these results we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (mg/dL)	62,87	40,52
Standard Deviation (mg/dL)	0,24	0,42
Coefficient of Variation (%)	0,38	1,05

REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (mg/dL)	62,98	40,80
Standard Deviation (mg/dL)	0,25	0,38
Coefficient of Variation (%)	0,39	0,93

Sensitivity

Sensitivity was calculated from 10 determinations of a sample free of the Alpha 1 Acid Glycoprotein. The average found was 0,53 mg/dL, with standard deviation of 0,05 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,67 mg/dL.

Linearity

The reaction is linear up to the level of the highest calibrator concentration.

DIAGNOSTIC SIGNIFICANCE

Main component of the Winzler Mucoproteins, rises considerably in acute inflammatory processes. It is an important index of rheumatic activity, because it remains high when other tests normalize.

The increase in values also happens in erythematosus lupus, neoplasm metastasis and administration of estrogen. Decrease of values occurs in liver failure, adrenal insufficiency, malnutrition and pituitary insufficiency .

BIBLIOGRAPHIC REFERENCES

- 1 - KILLINGSWORTH, L. M; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17, 936.
- 2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
- 3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
- 4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
- 5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
Email: sac@bioclin.com.br

ANVISA registration for Alpha 1 Acid Glycoprotein kit:
10269360061

Review: February/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MANUFACTURED BY

BATCH CODE

CONTROL

DATE OF MANUFACTURE

POSITIVE CONTROL

USED BY
(last day of month)

NEGATIVE CONTROL

TEMPERATURE LIMITATION
(store at)

BIOLOGICAL RISK

CONTAINS SUFFICIENT
FOR <N> TESTS

INFLAMMABLE

CONSULT INSTRUCTIONS
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC DEVICE

POISON

EUROPEAN AUTHORIZED
REPRESENTATIVE

CE MARK

KEEP AWAY
FROM SUNLIGHT

DO NOT USE IF
PACKAGE IS
DAMAGED