AMILASE

REF K003

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Amylase activity. Colorimetric test for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Caraway Modified - Fixed time Kinetic

The sample is incubated in a Starch substrate. By Iodine addition, the Starch still not hydrolyzed acquires a blue coloration that lowers proportionally to the enzymatic activity, being compared to a control.

REAGENTS

Number 1 - Substrate - Store between 2 and 8°C. Contains: Starch < 5 mmol/L, Phosphate Buffer < 500 mmol/L, stabilizer and preservative.

The reagent in lower temperatures can crystallize. In this case, place the reagent (in its original packing well sealed) in water bath 37°C mixing until crystals dissolution is complete.

Number 2 - Color Reagent - Iodine Stock Solution - Store between 2 and 8°C. Contains: Potassium Iodate < 50 mmol/L, Potassium Iodide < 600 mmol/L, Hydrochloric Acid < 200 mmol/L and stabilizer.

PRESENTATION

Reagent	Volume
Nº 1	50 mL
Nº 2	5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, water-bath 37°C, watch and stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must to be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- The water level in water-bath must be higher than the reagent levels in the test tubes.

6- How the Reagent Nº 1 is easily contaminated by slobber, we recommend the use of masks, pump and automatic pipettes.

7- It is important for the well development of the test, a strict control over the time and temperature.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma with Heparin and Urine. Other anticoagulants (EDTA, Citrate, Oxalate) inhibit the activity of Amylase.

Serum or Plasma: stable for 7 days at 15 to 30°C, 2 months at 2 to 8°C and 1 year at -20°C °.

Urine: stable for 2 days at 15 to 30°C, 10 days at 2 to 8°C and 3 weeks at -20°C °. Gather in the range of 2 to 24 hours.

PROCESS DESCRIPTION

WORK REAGENT PREPARATION

Transfer the ampoule contend (Reagent Nº 2) to empty vessels contained in the kit. Next, add 45 mL of distilled or deionized water and homogenize it. Keep it sealed. Stable per 06 months between 2 and 8°C.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 2 test tubes: C (control), A (sample) and proceed as following:

	Control	Sample
Reagent Nº 1	0,5 mL	0,5 mL
Put in water bath at 37°C per 2 minutes.		
Sample	---	10 µL
Homogenize and incubate at 37°C for exactly 7 minutes and 30 seconds (chronometered).		
Work Reagent	0,5 mL	0,5 mL
Distilled or Deionized Water	4,0 mL	4,0 mL

Homogenize well and determinate the absorbencies of Control and Sample at 660 nm (620 to 700 nm), initializing the instrument with distilled or deionized water. The color is stable per 30 minutes. For urine dosage, measure urine volume collected in an interval of 2 to 24 hours. Adjust the pH between 7,0 and 7,4 with Sodium Carbonate (solid) for acid urine and with Biacid Potassium Phosphate (KH₂PO₄) for alkaline urine. Proceed in the same way for serum.

CALCULATIONS

Ac = Absorbance Control

Aa = Absorbance of Sample

Amylase = $\frac{Ac - Aa}{Ac}$ x 800 (U/dL)

An enzyme unit is the quantity that hydrolyzes totally 10 mg of Starch in 30 minutes at 37°C.

Calculation for urine dosage:

Amylase = $\frac{Ac - Aa}{Ac}$ x $\frac{\text{Urine Volume (mL)}}{\text{Collection Time (hours)}}$ x 8

Example:

Serum Control Absorbance = 0,368

Serum Sample Absorbance = 0,320

Urine Control Absorbance = 0,380

Urine Sample Absorbance = 0,300

Collection Time : 12 hours

Volume: 900 mL

$\frac{0,368 - 0,320}{0,368}$ x 800 = 104 U/dL of Amylase

$\frac{0,380 - 0,300}{0,380}$ x $\frac{900 \times 8}{12}$ = 126 U/h of Amylase

When the value $\frac{Ac - Aa}{Ac}$ is higher than 0,5, dilute the urine with the Sodium Chloride solution 0,85%, repeat the essay and multiply the result by the dilution factor.

The results are expressed in U/dL (serum) or U/h (urine).

PROCEDURE LIMITATIONS

As in every enzymatic reaction, it is of great importance the control of time and temperature of the incubation. The difference of 1 minute in the incubation time carry an error of 13% in the results.

INTERFERENCES

Samples with hemolysis produce results wrongly reduced.

Anticoagulants such as EDTA, Citrate and Oxalate inhibit the amylase activity.

Bilirubin values up to 20 mg/dL and hyperlipemia (triglycerides up to 1500 mg/dL) do not interfere.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

The reference values in U/dL, for this method were obtained through the determination of Amylase in healthy populations of male and female.

Of all ages:

Serum: 60 to 160 U/dL

Urine: 50 to 140 U/h

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

QUALITY CONTROL

Accuracy

COMPARISON OF METHODS AND METHODOLOGICAL SPECIFICITY

The Amylase kit was compared with other commercially available methods for measurement of Amylase. 7 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0,993X + 1,126$ and the correlation coefficient 0,999. With these results, we can conclude that the kit shows good methodological specificity.

Precision

REPEATABILITY

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/dL)	57,95	71,90	121,10
Standard Deviation (U/dL)	0,76	0,79	0,91
Coefficient of Variation (%)	1,31	1,10	0,75

REPRODUCIBILITY

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Average Concentration (U/dL)	58,08	71,90	121,02
Standard Deviation (U/dL)	0,13	0,05	0,08
Coefficient of Variation (%)	0,22	0,07	0,06

Sensitivity

The sensitivity was calculated from twenty (20) determinations of a sample free of Amylase. The average found was 0,0039 U/dL with standard deviation of 0,0001 U/dL. The sensitivity, indicated by the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0,0042 U/dL.

Linearity

The reaction is linear until 400 U/dL. For higher values, dilute the sample with Sodium Chloride 0,85% and proceed with the new determination. Multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Amylase, predominantly of pancreatic and salivary origin, is normally present in blood and in urine in smaller quantities. Raises quickly in plasma after the beginning of symptom of acute pancreatitis, where, parallel to it can be observed the increase of the Urinary Amylase.

Increased values are also, observed in mesentery infarction, perforated gastric ulcers, carcinoma of the pancreas head, mumps, kidney insufficiency and diabetes acidosis.

Plasmatic levels reduced are observed in hepatitis, hepatic cirrhosis, pregnancy toxemia, eclamps and pancreatic carcinoma. Some drugs can interfere in the laboratorial results, elevating the serum levels of Amylase (Morphine, Meperidine, Codeine, Thiazide Diuretics) or reducing, as in the Barbituric poisoning cases.

NUMBER OF TESTS

100 Tests/10 µL of Sample/ 1 mL of Reagent

200 Tests/ 5 µL of Sample/ 500 µL of Reagent

BIBLIOGRAPHIC REFERENCES

1 - CARAWAY, W. T.: Amer J. Clin. Pathol, 1959, 32, 97.

2 - TONKS, D. B.: Clin. Chem., 1983, 9:217.

3 - HENRY, R. J.: Chemistry Technics 4° Ed. New York - Haper e Row, 1964, 498.

4 - HENRY, R. J. and CHIAMORI, N.: Clin. Chem., 1960, 6, 434.

5 - CARL, A. B. and EDWARD, R. A.: Tietz Textbook of Clinical Chem 2nd ed., 1994, 852-863.

6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

 **OBELIS S.A.**
Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE
Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for Amylase kit: 10269360085

Review: January/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED