ANTI ESTREPTOLISINA O

REF K057

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method of quantitative determination of Anti-Streptolysin O (ASO) quantity. Immunoturbidimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

The polystyrene particles coated with Streptolysin O coalesce with samples containing Anti-Streptolysin O. The intensity of scattered light is proportional to the concentration of Anti-Streptolysin O, so that, by comparison with the calibrator concentration known, one can determine the amount of ASO in the sample.

Reagents

- Number 1 - Buffer** - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.
- Number 2 - ASO Latex** - Store between 2 and 8°C. Contains: Suspension of Latex particles sensitized with Streptolysin O < 50 mg/dL, stabilizers and preservative.
- Number 3 - Calibrator** - Store between 2 and 8°C. Contains: Anti Streptolysin O, stabilizers and preservative. **Potentially infectious.**
- Warning:** The concentration of Anti Streptolysin O varies according to the lot - see the bottle's label.

PRESENTATION

Presentation	Reagent		
	Nº 1	Nº 2	Nº 3
1	18 mL	2 mL	1,0 mL
2	45 mL	5 mL	1,0 mL
3	2 x 45 mL	2 x 5 mL	2 x 1,0 mL
4	20 mL	5 mL	1,0 mL
5	2 x 20 mL	2 X 5 mL	1,0 mL
6	3 x 20 mL	3 x 5 mL	2 x 1,0 mL
7	4 x 20 mL	4 x 5 mL	2 x 1,0 mL
8	40 mL	10 mL	1,0 mL
9	2 x 40 mL	2 x 10 mL	2 x 1,0 mL
10	3 x 40 mL	3 x 10 mL	3 x 1,0 mL
11	4 x 40 mL	4 x 10 mL	4 x 1,0 mL
12	60 mL	15 mL	2 x 1,0 mL
13	2 x 60 mL	2 x 15 mL	3 x 1,0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. They can be found at markets specialized in Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain exact results.
- 3-** Water used in material cleaning must to be recent and free of contaminants.
- 4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5-** Reagent Nº 3 should be handled cautiously since it is subject to biological contamination.
- 6-** Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- 7-** The reagent has been tested for HIV, HCV and HBsAg using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**
- 8-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 9-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 10-** Do not use the product in case of damaged packaging.
- 11-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis and intense lipemia. Samples are stable in temperatures between 2 and 8°C for up to 2 days and 3 months at -20°C.

PROCESS DESCRIPTION

The calibration stability of the Anti Streptolysin O kit installed on refrigerated equipment is at least 22 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

The reagents are ready to use.
Thoroughly mix the Reagent Nº 2 before use.

For calibration of the reaction, utilize preferentially the kit's internal calibrator or the Bioclin Multical kit. For control use the Bioclin Multicontrol Kit.

This kit is for use in automated biochemical analyzers. **Check the application sheet for this equipment on www.bioclin.com.br or through SAC.**

PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia (low or moderate) do not interfere with assay performance. Gross lipemia and turbid samples should be processed with a higher dilution and read withing a dynamic range. Values below 20 UI/mL lead to poor reproductive results.

INTERFERENT

No interference was observed by Bilirubin until 60 mg/dL, Hemoglobin until 750 mg/dL, Triglycerides until 900 mg/dL and Rheumatoid Factor 600 UI/L.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 97/662 (Anti-Streptolysin-O, Human, NIBSC Reagent).

REFERENCE VALUES

The reference values for this method were obtained by determining anti-streptolysin O in healthy populations of male and female.

Less than 200 UI/mL

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
QUALITY CONTROL

Accuracy
RECOVERY

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 96 and 105%.

COMPARISON OF METHODS AND METHODOLOGY SPECIFICITY

The Anti Streptolysin O kit was compared with other commercially available methods for determination of Anti-Streptolysin O dosage. 41 tests were performed and the results were evaluated. The linear equation obtained was $Y = 0,947X + 17,127$ and correlation coefficient 0,9779. With these results we can conclude that the kit shows good methodological specificity.

Precision REPEATABILITY

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	96,28	201,59	390,12
Standard Deviation (UI/mL)	1,97	4,35	6,53
Coefficient of Variation (%)	2,04	2,16	1,67

REPRODUCIBILITY

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	101,44	206,18	399,72
Standard Deviation (UI/mL)	4,42	6,87	11,39
Coefficient of Variation (%)	4,35	3,33	2,85

Sensitivity

Sensitivity was calculated from 40 determinations of a sample free of Anti-Streptolysin O. The average found was 2,276 UI/mL, with standard deviation of 1,584 UI/mL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 7,027 UI/mL. Values below 20 UI/mL lead to poor reproductive results.

Linearity

The reaction is linear up to 1200 UL/mL. For samples with higher values than 1200 UI/mL, it is recommended to dilute the sample with Sodium Chloride 0,85%, repeat the dosage and multiply the result by the factor dilution.

High Dose Hook Effect

Hook effect wasn't observed with a high concentration of Antistreptolysin O up to 1800 IU/mL.

DIAGNOSTIC SIGNIFICANCE

The Streptolysin O is a toxic immunogenic exoenzymes produced by group A β -hemolytic streptococci. Given the large amount of exoenzymes released *in vivo*, determination of antibody response against the presence of Streptolysin O has become a routine procedure for the diagnosis and treatment of rheumatic fever, acute glomerulonephritis, scarlet, tonsillitis, erysipelas, sepsius puerperal and other streptococcal infections in group A.

BIBLIOGRAPHIC REFERENCES

- 1 - NASAHICO, Oukumi, Antistreptolysin O Japan Clinical, 1985, 43, 387.
- 2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
- 3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
- 4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
- 5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980, 26, 1459.
- 6 - GALVIN, J. P. et., al Clin. Lab. Assays, 73, 4th, 1983.
- 7 - PASSING, H.; BABLOCK, W., J. Clin. Chem. Clin. Biochem 21, 709; 1983.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - Fax: +55 (31) 3439.5455
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

OBELIS S.A.

Bd. Général Wahis, 53
1030 Brussels, Belgium

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration number for Anti Streptolysin O: 10269360067

Review: January/2020

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED